


Protocol

Isolation of mouse brain-infiltrating leukocytes for single cell profiling of epitopes and transcriptomes



High dimensional compositional and transcriptional profiling of heterogeneous brain-infiltrating leukocytes can lead to novel biological and therapeutic discoveries. High-quality single-cell leukocyte preparations are a prerequisite for optimal single cell profiling. Here, we describe a protocol for epitope and RNA-preserving dissociation of adult mouse brains and subsequent leukocyte purification and staining, which is adaptable to homeostatic and pathogenic brains. Leukocyte preparation following this protocol permits exquisite single-cell surface protein and RNA profiling in applications including CyTOF and CITE-seq.

Ian H. Guldner,
Samantha M.
Golomb, Qingfei
Wang, Emilia Wang,
Siyuan Zhang

iguldner@stanford.edu
(I.H.G.)
szhang8@nd.edu (S.Z.)

Highlights

Mouse brain
dissociation into
single-cell
suspension for
downstream single-
cell profiling

Enrichment for brain-
resident and low-
abundant brain-
infiltrating leukocytes

Protocol tested to
preserve >40 immune
cell surface epitopes

Compatible with
scRNA-seq, flow
cytometry, and
CyTOF

Guldner et al., STAR Protocols
2, 100537
June 18, 2021 © 2021 The
Author(s).
[https://doi.org/10.1016/
j.xpro.2021.100537](https://doi.org/10.1016/j.xpro.2021.100537)



Protocol

Isolation of mouse brain-infiltrating leukocytes for single cell profiling of epitopes and transcriptomes

Ian H. Guldner,^{1,2,3,5,*} Samantha M. Golomb,^{1,2,5} Qingfei Wang,^{1,2} Emilia Wang,^{1,2} and Siyuan Zhang^{1,2,4,5,6,*}

¹Department of Biological Sciences, College of Science, University of Notre Dame, Notre Dame, IN, 46556, USA

²Mike and Josie Harper Cancer Research Institute, University of Notre Dame, South Bend, IN, 46617, USA

³Stanford University, Stanford, CA, 94305, USA

⁴Indiana University Melvin and Bren Simon Cancer Center, Indianapolis, IN, 46202, USA

⁵Technical contact

⁶Lead contact

*Correspondence: iguldner@stanford.edu (I.H.G.), szhang8@nd.edu (S.Z.)
<https://doi.org/10.1016/j.xpro.2021.100537>

SUMMARY

High dimensional compositional and transcriptional profiling of heterogeneous brain-infiltrating leukocytes can lead to novel biological and therapeutic discoveries. High-quality single-cell leukocyte preparations are a prerequisite for optimal single cell profiling. Here, we describe a protocol for epitope and RNA-preserving dissociation of adult mouse brains and subsequent leukocyte purification and staining, which is adaptable to homeostatic and pathogenic brains. Leukocyte preparation following this protocol permits exquisite single-cell surface protein and RNA profiling in applications including CyTOF and CITE-seq.

For complete details on the use and execution of this protocol, please refer to Guldner et al. (2020) and Golomb et al. (2020).

BEFORE YOU BEGIN

The following protocol (Figure 1A) is for one adult mouse brain (male and female, ages 6 weeks to 70 weeks) with a primary emphasis on preparation for Cellular Indexing of Transcriptome and Epitopes by Sequencing (CITE-seq)/RNA Expression and Protein Sequencing (REAP-seq) and a secondary emphasis on preparation for mass cytometry by time of flight spectrometry (CyTOF) (Bandura et al., 2009; Peterson et al., 2017; Stoeckius et al., 2017). CyTOF and CITE-seq are powerful high dimensional single cell profiling strategies. CyTOF profiles surface epitopes (and in some cases, intracellular epitopes) by detection of metal-conjugated antibodies via time of flight mass spectrometry. CyTOF is predominantly used for immune cell population compositional analyses, and limited functional status may be elucidated by use of antibodies targeting functional markers (for example, PD-L1 or CTLA-4). CITE-seq also profiles surface epitopes, but by means of sequencing oligo-tagged antibodies bound to individual cells, and additionally profiles cellular mRNA by standard 3-prime poly-A capture of mRNA and subsequent sequencing by synthesis. The multi-modal nature of CITE-seq allows both compositional analysis and more detailed functional status determination. Importantly, this protocol can be used for brain dissociation and leukocyte enrichment upstream of several other single cell profiling methods including, but not limited to, standard single cell RNA-sequencing (scRNA-seq) and flow cytometry. Reagents should be scaled to the total number of brain preparations accordingly. Steps listed in 'Before You Begin' become increasingly necessary with increasing sample numbers, which is typical for most applications.






Figure 1. Isolation of brain-infiltrating leukocytes for single cell profiling of epitopes and transcriptomes

(A) Overall protocol schematic.

(B) Images of indicated organs without PBS perfusion (right) and with PBS perfusion (left). Scale bar: 2mm.

(C) Photograph of dissociated brain at the end of enzymatic dissociation (note there are no tissue chunks present).

(D) Photographs of Percoll layering sequences (1-4) and layered Percoll gradient before (lower left) and after (lower middle) centrifugation, with enlarged image pointing out the buffy leukocyte layer post-centrifugation (lower right).

(E) Light microscope image of washed leukocytes post-Percoll gradient centrifugation on a hemocytometer (note the lack of debris and single cell nature of suspension). Scale bar: 100µm.

(F) Bar chart showing the percent of viable cells following brain dissociation, leukocyte enrichment, and cell staining, as determined by the percent of cisplatin negative/low cells of all single cells in CyTOF experiments.

See also [Methods video S1](#).

Figure 1. Continued

(G) Bar chart derived from CITE-seq data showing percentage of cells sequenced that were leukocytes (CD45+) from the naive and brain metastasis-burdened brain. Error bars represent SEM.

(H) Stacked bar chart derived from CITE-seq data showing proportions of all identifiable cells sequenced that were leukocytes (CD45+), endothelial cells, pericytes, or metastatic cancer cells from the naive and brain metastasis-burdened brain.

Prepare percoll

⌚ Timing: 5 min

1. Prepare 'working stock' (WS) Percoll by adding 0.71 mL 10× HBSS (without CaCl₂, MgCl₂, or CaSO₄) to 6.42 mL Percoll. Mix well by inverting or vortexing.
2. Prepare 70% Percoll by adding 1 mL 1× HBSS (without Ca²⁺, Mg²⁺, or Phenol Red) to 3.5 mL WS Percoll. Mix well by inverting or vortexing.
3. Prepare '37% Red Percoll' by adding 2.65 mL 'red 1× HBSS' (without CaCl₂, MgCl₂, or CaSO₄; with Phenol Red) to 1.85 mL WS Percoll. Mix well by inverting or vortexing.
4. Prepare 30% Percoll by adding 3 mL 1× HBSS (without Ca²⁺, Mg²⁺, or Phenol Red) to 1.5 mL WS Percoll. Mix well by inverting or vortexing.

⚠ **CRITICAL:** Percoll solutions should only be prepared the day of use as Percoll solutions can become discolored 24h after preparation, hindering the identification of the buffy leukocyte layer. To increase the speed of preparation on the day of the experiment, we recommended aliquoting pre-measured individual components in separate containers so that they can be quickly mixed when needed on the day of the experiment.

⚠ **CRITICAL:** Percoll should be kept at room temperature (22°C–25°C), otherwise Percoll density changes and the gradient centrifugation success will be compromised.

Note: 'red 1× HBSS' is used in step 3 to enhance the visibility and thus identification of the buffy leukocyte layer at the interface of the red and clear Percoll layers. 1× HBSS without Phenol Red may be substituted for Phenol Red 1× HBSS, but the visibility of the leukocyte buffy layer will be less obvious.

Note: The Percoll recipes provided here include slightly more volume than necessary for one brain to account for pipetting error.

Prepare blocking buffer, antibody cocktails and staining solutions, and staining/washing buffers

⌚ Timing: 30 min

5. Prepare CITE-seq staining and washing buffers according to recipes provided in the Materials and Equipment section. Keep all buffers at 4°C or on ice until use.
6. To prepare blocking buffer, dilute stock mouse FcR Blocking Reagent 1:10 in CITE-seq staining buffer. Keep at 4°C or on ice until use.
7. Prepare antibody cocktails. Keep all antibody cocktails at 4°C or on ice until use.
 - a. Prepare CITE-seq antibody cocktail for ~3–4 × 10⁵ cells by adding 0.25 μL (0.125 μg) of each CITE-seq (oligo-tagged) antibody to 50 μL CITE-seq staining buffer. If multiple samples will be prepared and subsequently pooled and 'hashed' to enable multiplexed loading onto the 10× Chromium, first scale up the antibody amounts (minus HTO/"hashing" antibodies) and CITE-seq staining buffer accordingly. Then, aliquot antibody cocktail in 50 μL volumes into individual 1.5 mL tubes. Finally, add 0.3 μL (0.150 μg) unique hashing antibodies, which are oligo-conjugated antibodies that bind ubiquitously expressed proteins MHCI and CD45 and act as sample barcodes, to the appropriate sample tube. For example, add 0.3 μL HTO#1 antibody to




Figure 2. CITE-seq quality control

(A) Schematics of antibody staining process.

(B) Representative bioanalyzer electropherograms of CITE-seq antibody (ADT) library showing ADT product peak ~180 bp, which should be the predominant product peak. Smaller peak at ~140 bp is TSO-RT-Oligo product, which could be removed by further SPRI purification.

(C) Representative bioanalyzer electropherograms of RNA library with expected library size of ~300-700 bp.

(D) Representative Cell Ranger QC summary plot of mean reads per cell versus median genes per cell.

(E) Representative Cell Ranger QC summary plot of mean reads per cell versus sequencing saturation, with the far right data point indicating full sequencing depth of run.

sample tube A, 0.3 μ L HTO#2 antibody to sample tube B, etcetera (Figure 2A). Antibody pools should be mixed by pipetting up and down.

- Prepare CyTOF antibody cocktail by adding volumes of antibodies indicated in Materials and Equipment section to 50 μ L Maxpar Cell Staining Buffer in a 1.5 mL tube (or larger, if needed). Antibody pools should be mixed by pipetting up and down.
- Prepare CyTOF 0.75 μ M Cell-ID Cisplatin work solution for dead/live cell identification by diluting 0.5 μ L Cell-ID Cisplatin in 2 mL Maxpar PBS.

- d. Prepare CyTOF iridium intercalator solution for all cell nuclear stain by diluting 0.5 μ L Cell-ID Intercalator-Ir in 2 mL Maxpar Fix and Perm Buffer.

△ **CRITICAL:** All CyTOF reagents must be kept away from autoclaved and washed materials and be minimally exposed to air to avoid ambient metal contamination (especially Barium contamination). Use only non-autoclaved pipette tips, tubes, etc., and aliquot CyTOF-specific reagents from original reagent bottles by pouring on an as needed basis. Any aliquoted reagents, including antibody cocktails, remaining after experiment should be discarded and never used for subsequent preparations.

△ **CRITICAL:** CITE-seq and CyTOF staining buffers are not interchangeable.

△ **CRITICAL:** If you are preparing multiple samples for CITE-seq and performing hashing, make sure each sample receives a unique HTO antibody. If two samples are accidentally stained with the same HTO antibody, they will not be able to be distinguished from each other during data processing and analysis.

Note: The recommended antibody amounts in this protocol serve as an excellent starting point that has worked well in our hands, but antibodies should ideally be titrated prior to use. Ideally, titration should be performed with CITE-seq and CyTOF antibodies and optimized by sequencing and CyTOF runs. For CITE-seq, CITE-seq antibody manufacturer BioLegend suggests a less expensive and more convenient option is to perform titration experiments by flow cytometry using fluorescently conjugated antibodies of the same clonal identity, the titration results of which translate well to CITE-seq antibodies.

Prepare items for brain perfusion and extraction

⌚ **Timing:** 2 min

8. Fill a 10 mL syringe with a 25 gauge needle with ice cold 1 \times PBS (without Ca²⁺ or Mg²⁺).
9. Gather and clean bone cutting scissors and regular forceps.
10. Prepare 2 mL collection tube for extracted brain by filling with 1 mL 1 \times HBSS (without Ca²⁺, Mg²⁺, or Phenol Red) and place on ice.

Prepare enzyme dissociation mix

⌚ **Timing:** 15 min

11. Prepare stock components of the Multi-Tissue Dissociation Kit I as outlined in manufacturer's datasheet (<https://www.miltenyibiotec.com/upload/assets/IM0018263.PDF>) and reiterated below.
 - a. Prepare Enzyme D by adding 3 mL serum-free RPMI 1640 to lyophilized Enzyme D vial, inverting continuously until fully dissolved.
 - b. Prepare Enzyme R by adding 2.7 mL serum-free RPMI 1640 to lyophilized Enzyme R vial, inverting continuously until fully dissolved.
 - c. Prepare Enzyme A by adding 1 mL Buffer A to lyophilized Enzyme A vial, inverting continuously until fully dissolved.
12. Prepare working enzyme mix by adding 100 μ L Enzyme D, 50 μ L Enzyme R, 12.5 μ L Enzyme A to 2.338 mL 1 \times D-PBS (with Ca²⁺ and Mg²⁺). Prepare immediately prior to beginning brain collection and keep on ice until use.

△ **CRITICAL:** Caution should be used if using other enzymatic tissue dissociation methods, especially if the downstream application includes surface protein staining. We have

observed epitope degradation in using papain-based enzymatic dissociation methods (specifically using the papain-based Neural Tissue Dissociation Kit and Brain Tumor Dissociation Kit from Miltenyi Biotec). We and the manufacturer (Miltenyi Biotec) have experimentally verified that the Multi-Tissue Dissociation Kit I preserves many common immune cell surface epitopes, including all listed in the Key Resources Table.

△ **CRITICAL:** Enzymatic tissue dissociation methods have been shown to reduce the astrocyte and neuron viability and activate and/or alter the transcriptome of microglia and other cell types (Mattei et al., 2020). Users should evaluate their cells to examine whether the digestion has artificially activated or changed their transcriptional status. If enzyme digestion is deemed to alter cellular status, other cold enzymes preparation options could be explored.

Note: We have not observed any detrimental effects as a result of not scaling down the components of the working enzyme mix if less than one brain is used. We have found this enzyme mix can be stretched to dissociate as many as 1.5 brains at the expense of slightly longer dissociation time (~5 min longer).

Note: Unused enzyme aliquots should be immediately stored at -20°C . Enzymes stored at -20°C will remain good for at least 6 months.

KEY RESOURCES TABLE

REAGENT or RESOURCE	SOURCE	IDENTIFIER
Antibodies		
CytoF: anti-CD45-089Y (30-F11)	Fluidigm	Cat#3089005B RRID: N/A
CytoF: rat anti-mouse Ly-6G-141Pr (1A8)	Fluidigm	Cat#3141008B; RRID:AB_2814678
CytoF: anti-mouse CD11c-142Nd (N418)	Fluidigm	Cat#3142003B; RRID:AB_281473
CytoF: anti-mouse CD45R-144Nd (RA3-6B2)	Fluidigm	Cat#3144011B; RRID:AB_2811239
CytoF: anti-mouse CD4-145Nd (RM4-5)	Fluidigm	Cat#3145002B; RRID:AB_2687832
CytoF: rat anti-mouse CD11b-148Nd (M1/70)	Fluidigm	Cat#3148003B; RRID:AB_2814738
CytoF: rat anti-mouse CD44-150Nd (IM7)	Fluidigm	Cat#3150018B; RRID:AB_2827882
CytoF: rat anti-CD25-151Eu (3C7)	Fluidigm	Cat#3151007B; RRID:AB_2827880
CytoF: anti-mouse CD3e-152Sm (145-2C11)	Fluidigm	Cat#3152004B; RRID:AB_2687836
CytoF: anti-mouse PD-L1-153Eu (10F.9G2)	Fluidigm	Cat#3153016B; RRID:AB_2687837
CytoF: anti-CTLA-4-154Sm (UC10-4B9)	Fluidigm	Cat#3154008B RRID: N/A
CytoF: anti-mouse PD-1-159Tb (29F.1A12)	Fluidigm	Cat#3159024B; RRID:AB_2687839
CytoF: anti-Ly-6C-162Dy (HK1.4)	Fluidigm	Cat#3162014B RRID: N/A
CytoF: anti-mouse CX3CR1-164Dy (SA011F11)	Fluidigm	Cat#3164023B; RRID:AB_2832247
CytoF: anti-NK1.1-165Ho (PK136)	Fluidigm	Cat#3165018B RRID: N/A
CytoF: rat anti-mouse ckit-166Er (2B8)	Fluidigm	Cat#3166004B; RRID:AB_2801435
CytoF: rat anti-mouse CD8a-168Er (53-6.7)	Fluidigm	Cat#3168003B; RRID:AB_2811241
CytoF: anti-CD86-172Yb (GL1)	Fluidigm	Cat#3172016B RRID: N/A
CytoF: anti-I-A/I-E-209Bi (M5/114.15.2)	Fluidigm	Cat#3209006B RRID: N/A
CITE-seq: rat anti-CCR2/CD192 (SA203G11)	BioLegend	Cat#150625; RRID:AB_2783122
CITE-seq: rat anti-CD117/c-kit (2B8)	BioLegend	Cat#105843; RRID:AB_2749960
CITE-seq: rat anti-CD11b (M1/70)	BioLegend	Cat#101265; RRID:AB_2734152
CITE-seq: armenian hamster anti-CD11c (N418)	BioLegend	Cat#117355; RRID:AB_2750352
CITE-seq: rat anti-CD172a/SIRP α (P84)	BioLegend	Cat#144033; RRID:AB_2800670
CITE-seq: rat anti-CD25 (PC61)	BioLegend	Cat#102055; RRID:AB_2749982
CITE-seq: rat anti-CD3 (17A2)	BioLegend	Cat#100251; RRID:AB_2750533
CITE-seq: rat anti-CD38 (90)	BioLegend	Cat#102733; RRID:AB_2750556
CITE-seq: rat anti-CD4 (RM4-5)	BioLegend	Cat#100569; RRID:AB_2749956

(Continued on next page)

<i>Continued</i>		
REAGENT or RESOURCE	SOURCE	IDENTIFIER
CITE-seq: rat anti-CD44 (IM7)	BioLegend	Cat#103045; RRID:AB_2734154
CITE-seq: rat anti-CD45 (30-F11)	BioLegend	Cat#103159; RRID:AB_2734156
CITE-seq: rat anti-CD45R/B220 (RA3-6B2)	BioLegend	Cat#103263; RRID:AB_2734158
CITE-seq: rat anti-CD86 (GL-1)	BioLegend	Cat#105047; RRID:AB_2750348
CITE-seq: rat anti-CD8a (53-6.7)	BioLegend	Cat#100773; RRID:AB_2734151
CITE-seq: mouse anti-CD90.1 (OX-7)	BioLegend	Cat#202547; RRID:AB_2783141
CITE-seq: mouse anti-Cx3cr1 (SA011F11)	BioLegend	Cat#149041; RRID:AB_2783121
CITE-seq: rat anti-F4/80 (BM8)	BioLegend	Cat#123153; RRID:AB_2749986
CITE-seq: rat anti-I-A/I-E (M5/114.15.2)	BioLegend	Cat#107653; RRID:AB_2750505
CITE-seq: rat anti-Ly6C (HK1.4)	BioLegend	Cat#128047; RRID:AB_2749961
CITE-seq: rat anti-Ly6G (1A8)	BioLegend	Cat#127655; RRID:AB_2749962
CITE-seq: mouse anti-NK1.1 (PK136)	BioLegend	Cat#108755; RRID:AB_2750536
CITE-seq: rat anti-PD-1 (RMP1-30)	BioLegend	Cat#109123; RRID:AB_2734169
CITE-seq: rat anti-PD-L1 (MIH6)	BioLegend	Cat#153604; RRID:AB_2783125
CITE-seq: rat anti-CD169/Siglec-1 (3D6.112)	BioLegend	Cat#142425; RRID:AB_2783106
CITE-seq: rat anti-Siglec-H (551)	BioLegend	Cat#129615; RRID:AB_275053
CITE-seq: mouse anti-TMEM119 (A16075D)	BioLegend	Cat#853303; RRID:AB_2801201
CITE-seq: mouse anti-XCR1 (Zet)	BioLegend	Cat#148227; RRID:AB_2783120
CITE-seq: rat anti-CD24 (M1/69)	BioLegend	Cat#101841; RRID:AB_2750380
CITE-seq: rat armenian hamster anti-CD103 (2e7)	BioLegend	Cat#121437; RRID:AB_2750349
CITE-seq: rat anti-CD49d (R1-2)	BioLegend	Cat#103623; RRID:AB_2734159
CITE-seq: rat anti-mouse Hashtag 1 antibody (M1/42; 30-F11)	BioLegend	Cat#155801; RRID:AB_2750032
CITE-seq: rat anti-mouse Hashtag 2 antibody (M1/42; 30-F11)	BioLegend	Cat#155803; RRID:AB_2750033
CITE-seq: rat anti-mouse Hashtag 3 antibody (M1/42; 30-F11)	BioLegend	Cat#155805; RRID:AB_2750034
CITE-seq: rat anti-mouse Hashtag 4 antibody (M1/42; 30-F11)	BioLegend	Cat#155807; RRID:AB_2750035
CITE-seq: rat anti-mouse Hashtag 5 antibody (M1/42; 30-F11)	BioLegend	Cat#155809; RRID:AB_2750036
CITE-seq: rat anti-mouse Hashtag 6 antibody (M1/42; 30-F11)	BioLegend	Cat#155811; RRID:AB_2750037
CITE-seq: rat anti-mouse Hashtag 7 antibody (M1/42; 30-F11)	BioLegend	Cat#155813; RRID:AB_2750039
CITE-seq: rat anti-mouse Hashtag 8 antibody (M1/42; 30-F11)	BioLegend	Cat#155815; RRID:AB_2750040
Chemicals, peptides, and recombinant proteins		
Bovine serum albumin	VWR Chemicals	Cat#97061-422 RRID: N/A
Cell-ID Cisplatin	Fluidigm	Cat#201064 RRID: N/A
Cell-ID Intercalator	Fluidigm	Cat#201192B RRID: N/A
Cell Staining Buffer	BioLegend	Cat#20201 RRID: N/A
Dulbecco's PBS (D-PBS), 1x, +Ca, +Mg	Caisson Labs	Cat#PBL02-500M RRID: N/A
Ethylenediaminetetra-acetic acid disodium salt dihydrate (EDTA)	VWR Chemicals	Cat#97061-022 RRID: N/A
FcR Blocking Reagent, mouse	Miltenyi Biotec	Cat#130-092-575 RRID: N/A
Hanks Balanced Salt Solution (HBSS), 1x, -CaCl ₂ , -MgCl ₂ , -MgSO ₄ , +Phenol Red	Thermo Fisher Scientific	Cat#14170161 RRID: N/A
Hanks Balanced Salt Solution (HBSS), 1x, -Ca, -Mg, -Phenol Red	Cytiva Life Sciences	Cat#SH30588.02 RRID: N/A
Hanks Balanced Salt Solution (HBSS), 10x, -CaCl ₂ , -MgCl ₂ , -MgSO ₄ , -Phenol Red	Thermo Fisher Scientific	Cat#14185052 RRID: N/A
Phosphate Buffered Saline (PBS), 1x, -Ca, -Mg, -Phenol Red	Cytiva Life Sciences	Cat#SH3025602 RRID: N/A

(Continued on next page)

Continued

REAGENT or RESOURCE	SOURCE	IDENTIFIER
Maxpar Cell Staining Buffer	Fluidigm	Cat#201068 RRID: N/A
MaxPar Fix and Perm Buffer	Fluidigm	Cat#201067 RRID: N/A
Maxpar PBS	Fluidigm	Cat#201058 RRID: N/A
MaxPar Water	Fluidigm	Cat#201069 RRID: N/A
Percoll	GE Healthcare	Cat#17-0891-02 RRID: N/A
RPMI 1640	Thermo Fisher Scientific	Cat#11875093 RRID: N/A
Tween 20	VWR Chemicals	Cat#97062-332 RRID: N/A
1x EQ beads	Fluidigm	Cat#201078 RRID: N/A
Critical commercial assays		
Chromium Single Cell 3' Library and Gel Bead Kit v3	10x Genomics	PN-1000092 RRID: N/A
Single Cell 3' GEM Library and Gel Bead Kit v3.1	10x Genomics	PN-1000092 RRID: N/A
Chromium Chip B Single Cell Kit	10x Genomics	PN-1000074 RRID: N/A
Chromium i7 Multiplex Kit	10x Genomics	PN-120262 RRID: N/A
Multi-tissue Dissociation Kit I	Miltenyi Biotec	130-110-201 RRID: N/A
Deposited data		
Mouse CITE-seq data	(Guldner et al., 2020)	GSE134285
Experimental models: cell lines		
E0771	CH3 BioSystems	SKU:940001-Vial RRID: N/A
Experimental models: organisms/strains		
C57BL/6	Jackson Laboratories	Stock#:000664 RRID: N/A
Oligonucleotides		
Primer: ADT_Truseq_i7_UDI01 - CAAGCAGAAGACGGCATAACGAGAACC GCGGGTACTGGAGTTCCTTGGCACCCGAGAATTCCA	This paper	N/A
Primer: ADT_Truseq_i7_UDI04 - CAAGCAGAAGACGGCATAACGAGTTGGACTTTGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA	This paper	N/A
Primer: ADT-RPI-1 - CAAGCAGAAGA CGGCATAACGAGATCGTGATGTGACTGGAGTTCCTTGGCACCCGAGAATTCCA	This paper	N/A
Primer: ADT-RPI-2 - CAAGCAGAAGACGGCATAACGAGATACATCGGTGACTGAGTTCCTTGGCACCCGAGAATTCCA	This paper	N/A
Primer: HTO_Nextera_i7_UDP01 - CAAGCAGAAGACGGCATAACGAGATCGCTCATGTGACTGGAGTTCAGACGTGTGC	This paper	N/A
Primer: HTO_Nextera_i7_UDP02 - CAAGCAGAAGACGGCATAACGAGATTATCTGACGTGACTGGAGTTCAGACGTGTGC	This paper	N/A
Primer: HTO_Nextera_i7_UDP03 - CAAGCAGAAGACGGCATAACGAGATATATGAGTGACTGGAGTTCAGACGTGTGC	This paper	N/A
Primer: HTO_Nextera_i7_UDP04 - CAAGCAGAAGACGGCATAACGAGATCTTAGTGACTGGAGTTCAGACGTGTGC	This paper	N/A
Primer: HTO-N701 - CAAGCAGAAGACGGCATAACGAGATTCCGCTTAGTGACTGGAGTTCAGACGTGTGC	This paper	N/A
Primer: HTO-N702 - CAAGCAGAAGACGGCATAACGAGATCTAGTACGTGACTGGAGTTCAGACGTGTGC	This paper	N/A
Software and algorithms		
FlowJo	BD	https://www.flowjo.com

(Continued on next page)

Continued

REAGENT or RESOURCE	SOURCE	IDENTIFIER
R studio	(Ihaka and Gentleman, 1996)	https://rstudio.com
CellRanger 3.1.0	10x Genomics	https://support.10xgenomics.com/single-cell-gene-expression/software/downloads/latest
Cytobank Premium	Beckman Coulter	https://cytobank.org
GraphPad Prism	GraphPad Software	https://www.graphpad.com
Seurat_3.1.1	(Satija et al., 2015)	https://satijalab.org/seurat/
Other		
CyTOF2	Fluidigm	N/A
Illumina NovaSeq 6000	Illumina	N/A
10x Chromium Controller and Accessory Kit	10x Genomics	PN-120223

MATERIALS AND EQUIPMENT

CITE-seq Staining Buffer

Component	Final Concentration or Percentage	Amount
PBS (1 ×, without Mg ²⁺ or Ca ²⁺)	1 ×	10 mL
BSA	2.0%	0.2 g
Tween 20	0.02%	2 μL
Total	n/a	10 mL

Store at 4°C for up to one month.

CITE-seq Wash Buffer A

Component	Final Concentration or Percentage	Amount
PBS (1 ×, without Mg ²⁺ or Ca ²⁺)	1 ×	10 mL
EDTA (0.5 M)	2.0mM	40 μL
BSA	2.0%	0.2 g
Tween 20	0.02%	2 μL
Total	n/a	10 mL

Store at 4°C for up to one month.

CITE-seq Wash Buffer B

Component	Final Concentration or Percentage	Amount
PBS (1 ×, without Mg ²⁺ or Ca ²⁺)	1 ×	10 mL
EDTA (0.5 M)	1.0mM	20 μL
BSA	2.0%	0.2 g
Tween 20	0.02%	2 μL
Total	n/a	10 mL

Store at 4°C for up to one month.

CITE-seq Wash Buffer C

Component	Final Concentration or Percentage	Amount
PBS (1 ×, without Mg ²⁺ or Ca ²⁺)	1 ×	10 mL
EDTA (0.5 M)	0.1mM	2.0 μL
BSA	1.0%	0.1g
Tween 20	0.02%	2 μL
Total	n/a	10 mL

Store at 4°C for up to one month.

Working Stock Percoll

Component	Final Concentration or Percentage	Amount
Percoll (aka original stock)	90%	6.42 mL
10× HBSS	1×	0.71 mL
Total	n/a	7.13 mL

Store at room temperature (22°C–25°C) for 3 h or less.

70% Percoll

Component	Final Concentration or Percentage	Amount
Working Stock Percoll	70% (relative to original stock)	3.5 mL
1× HBSS	22%	1.0 mL
Total	n/a	4.5 mL

Store at room temperature for 3 h or less.

37% Percoll

Component	Final Concentration or Percentage	Amount
Working Stock Percoll	37% (relative to original stock)	1.85 mL
1× HBSS (with Phenol Red)	58%	2.65 mL
Total	n/a	4.5 mL

Store at room temperature for 3 h or less.

30% Percoll

Component	Final Concentration or Percentage	Amount
Working Stock Percoll	30% (relative to original stock)	1.5 mL
1× HBSS (with Phenol Red)	66.66%	3.0 mL
Total	n/a	4.5 mL

Store at room temperature for 3 h or less.

Cell-ID Cisplatin Working Stock

Component	Final Concentration or Percentage	Amount
Cell-ID Cisplatin	0.75 μM	0.5 μL
Maxpar PBS	n/a	2 mL
Total	n/a	2.5 mL

Make up fresh each use.

Cell-ID Intercalator-IR Working Stock

Component	Final Concentration or Percentage	Amount
Cell-ID Intercalator IR	0.02%	0.5 μL
Maxpar Fix and Perm Buffer	n/a	2 mL
Total	n/a	2.5 mL

Make up fresh each use.

CyTOF Antibody Volumes (per 50 μL Maxpar Cell Staining Buffer/sample)

Component	Working dilution	Amount
anti-CD45-089Y (30-F11)	1:100	0.5 μL
rat anti-CD25-151Eu (3C7)	1:100	0.5 μL
anti-CTLA-4-154Sm (UC10-4B9)	1:100	0.5 μL

(Continued on next page)

Continued

Component	Working dilution	Amount
anti-mouse PD-L1-153Eu (10F.9G2)	1:100	0.5 µL
anti-mouse CD11c-142Nd (N418)	1:100	0.5 µL
anti-mouse CD4-145Nd (RM4-5)	1:100	0.5 µL
anti-mouse CD3e-152Sm (145-2C11)	1:100	0.5 µL
anti-mouse PD-1-159Tb (29F.1A12)	1:100	0.5 µL
CytoF: rat anti-mouse ckit-166Er (2B8)	1:100	0.5 µL
rat anti-mouse CD11b-148Nd (M1/70)	1:100	0.5 µL
rat anti-mouse CD44-150Nd (IM7)	1:50	1.0 µL
anti-NK1.1-165Ho (PK136)	1:50	1.0 µL
CytoF: anti-mouse CX3CR1-164Dy (SA011F11)	1:50	1.0 µL
anti-CD86-172Yb (GL1)	1:50	1.0 µL
anti-Ly-6C-162Dy (HK1.4)	1:200	0.25 µL
rat anti-mouse Ly-6G-141Pr (1A8)	1:200	0.25 µL
anti-I-A/I-E-209Bi (M5/114.15.2)	1:200	0.25 µL
rat anti-mouse CD8a-168Er (53-6.7)	1:200	0.25 µL
anti-mouse CD45R-144Nd (RA3-6B2)	1:200	0.25 µL

Make up fresh each use.

CITE-seq Antibody Volumes (per 50 µL CITE-seq Staining Buffer/sample)

Component	Working dilution	Amount
rat anti-CCR2/CD192 (SA203G11)	1:200	0.25 µL
rat anti-CD117/c-kit (2B8)	1:200	0.25 µL
rat anti-CD11b (M1/70)	1:200	0.25 µL
armenian hamster anti-CD11c (N418)	1:200	0.25 µL
rat anti-CD172a/SIRPα (P84)	1:200	0.25 µL
rat anti-CD25 (PC61)	1:200	0.25 µL
rat anti-CD3 (17A2)	1:200	0.25 µL
rat anti-CD38 (90)	1:200	0.25 µL
rat anti-CD4 (RM4-5)	1:200	0.25 µL
rat anti-CD44 (IM7)	1:200	0.25 µL
rat anti-CD45 (30-F11)	1:200	0.25 µL
rat anti-CD45R/B220 (RA3-6B2)	1:200	0.25 µL
rat anti-CD86 (GL-1)	1:200	0.25 µL
rat anti-CD8a (53-6.7)	1:200	0.25 µL
mouse anti-CD90.1 (OX-7)	1:200	0.25 µL
mouse anti-Cx3cr1 (SA011F11)	1:200	0.25 µL
rat anti-F4/80 (BM8)	1:200	0.25 µL
rat anti-I-A/I-E (M5/114.15.2)	1:200	0.25 µL
rat anti-Ly6C (HK1.4)	1:200	0.25 µL
rat anti-Ly6G (1A8)	1:200	0.25 µL
mouse anti-NK1.1 (PK136)	1:200	0.25 µL
rat anti-PD-1 (RMP1-30)	1:200	0.25 µL
rat anti-PD-L1 (MIH6)	1:200	0.25 µL
rat anti-CD169/Siglec-1 (3D6.112)	1:200	0.25 µL
rat anti-Siglec-H (551)	1:200	0.25 µL
mouse anti-TMEM119 (A16075D)	1:200	0.25 µL
mouse anti-XCR1 (Zet)	1:200	0.25 µL
rat anti-CD24 (M1/69)	1:200	0.25 µL
rat armenian hamster anti-CD103 (2e7)	1:200	0.25 µL
rat anti-CD49d (R1-2)	1:200	0.25 µL
rat anti-mouse Hashtag 1 antibody (M1/42; 30-F11)	1:200	0.30 µL

(Continued on next page)

Continued

Component	Working dilution	Amount
rat anti-mouse Hashtag 2 antibody (M1/42; 30-F11)	1:200	0.30 μ L
rat anti-mouse Hashtag 3 antibody (M1/42; 30-F11)	1:200	0.30 μ L
rat anti-mouse Hashtag 4 antibody (M1/42; 30-F11)	1:200	0.30 μ L
rat anti-mouse Hashtag 5 antibody (M1/42; 30-F11)	1:200	0.30 μ L
rat anti-mouse Hashtag 6 antibody (M1/42; 30-F11)	1:200	0.30 μ L
rat anti-mouse Hashtag 7 antibody (M1/42; 30-F11)	1:200	0.30 μ L
rat anti-mouse Hashtag 8 antibody (M1/42; 30-F11)	1:200	0.30 μ L

Make up fresh each use.

Stock Enzyme D

Component	Final Concentration or Percentage	Amount
Lyophilized Enzyme D	n/a	1 vial
serum-free RPMI 1640	n/a	3 mL
Total	n/a	3 mL

Store at -20°C for up to one year.

Stock Enzyme R

Component	Final Concentration or Percentage	Amount
Lyophilized Enzyme R	n/a	1 vial
serum-free RPMI 1640	n/a	2.7 mL
Total	n/a	2.7 mL

Store at -20°C for up to one year.

Stock Enzyme A

Component	Final Concentration or Percentage	Amount
Lyophilized Enzyme A	n/a	1 vial
Buffer A	n/a	1 mL
Total	n/a	1 mL

Store at -20°C for up to one year.

Enzyme Mix

Component	Final Concentration or Percentage	Amount
Enzyme D	n/a	100 μ L
Enzyme R	n/a	50 μ L
Enzyme A	n/a	12.5 μ L
D-PBS	n/a	2.338 mL
Total	n/a	2.5 mL

Make up fresh each use.

Alternative Enzyme Dissociation Kits*

Reagent or Resource	Source	Identifier
Neural Tissue Dissociation Kit (P)	Miltenyi Biotec	130-092-628
Brain Tumor Dissociation Kit	Miltenyi Biotec	130-095-942

*Demonstrated to destroy some epitopes

STEP-BY-STEP METHOD DETAILS

Animal perfusion and brain extraction

⌚ Timing: 3 min

This step will perfuse the mouse to eliminate intra-vascular blood cells and results in brain extraction.

1. Anesthetize the mouse using methods approved by the institutional regulatory committee.
2. Spray chest of the mouse with water to mat down fur.
3. Perform a lateral incision down the center of the rib cage, cutting the skin.
4. Perform a lateral incision through the rib cage just to the left or right of the sternum up the entire rib cage to fully expose the chest cavity. Take care to not puncture heart, lungs, or other organs while making this excision. Remove as much of the rib cage as necessary to be able to see and hold the heart with forceps.
5. Holding the heart steady with forceps, make a small incision in the right atrium and immediately begin to perfuse the mouse with 10 mL ice-cold 1× PBS (without Mg²⁺ or Ca²⁺) by inserting 25 gauge needle into the apex of the left ventricle, slowly depressing the plunger, perfusing the entire 10 mL volume in no less than 1 min.
6. To extract the brain, first cut the vertebral column/spinal cord at the base of the skull (approximately C1-C3 region) with large bone scissors to detach the skull from the body. This cut should expose the brain stem. Next, insert one blade of fine scissors into the brain stem opening, bracing the sharp part of the scissor blade against the left inner part of the skull. Snip the skull laterally from the opening of the brain stem until you reach the left eye of the mouse, resulting in the left side of the skull as an 'opening point' and the right side of the skull as a 'hinge'. Next, using fine scissors or forceps, push the skull cap upwards, using the previous cut as an insertion/lever point. This should completely remove the skull cap, exposing the top of the brain and entire left side of the brain. Finally, slide forceps between the bottom of the brain and skull to pop out the brain. Keep the brain in a 2 mL tube with cold 1× HBSS on ice until the next step.

△ CRITICAL: Perfusion performed in step 5 must be successful to ensure leukocytes isolated in later steps were brain parenchyma-infiltrating rather than those within the vasculature. If perfusion is successful, the liver should turn from dark brown/pink to pale brown/beige (Figure 1B). If these color changes do not occur, attempt to perfuse the mouse again with at least 10 mL ice cold PBS.

Optional: Mouse can be perfused with >10mL of ice-cold PBS in step 5, which is recommended if time allows to further ensure removal of intra-vascular leukocytes.

Brain dissociation

⌚ Timing: ~25 min

This step will dissociate the mouse brain into a single cell suspension.

7. Mechanical Dissociation
 - a. Using sharp scissors (iris scissors or similar), thoroughly mince mouse brain in 2 mL tube filled with 1 mL 1× HBSS (without Ca²⁺, Mg²⁺, or Phenol Red) until chunks remain that are generally small enough to be pipetted using a p1000 micropipette.
 - b. Triturate (mechanically dissociate) minced brain using a p1000 micropipette to pipette up and down until there is absolutely no resistance upon pipetting upward. Ideally, this should be done with the sample tube sitting in ice or a pre-chilled tube rack.

- c. Transfer brain slurry contents to 15 mL conical tube and use fresh 1 × HBSS to wash out remaining brain slurry or fine chunks from 2 mL tube into 15 mL tube. As it is important to recover as much of the brain slurry as possible, several washes of 1-2 mL 1 × HBSS of the 2 mL tube may be performed. Because this will be centrifuged and pelleted in the next step, there should be no concern about the total volume of 1 × HBSS in the 15 mL transfer tube.
 - d. Add at least 2 mL 1 × HBSS to 15 mL conical tube containing brain slurry.
 - e. Centrifuge for 3 min at 300RCF and 4°C with acceleration and brake set to low settings.
8. Enzymatic Dissociation
- a. Aspirate supernatant from the brain cell pellet.
 - b. Resuspend each brain cell pellet in 2.5 mL pre-prepared Multi-Tissue Dissociation Kit I enzyme mix.
 - c. Rotate or rock sample on a rocker or tube rotator at a medium speed for 10 min at 37°C.
 - d. Triturate sample for 30 seconds to 1 min using a p1000 micropipette, eliminating most chunks of tissue. Any remaining chunks will be further dissociated in subsequent steps.
 - e. Rotate or rock sample at a medium speed for 10 min at 37°C.
 - f. Triturate sample for 30 seconds to 1 min using a p1000 micropipette. No tissue chunks should be present at this point (Figure 1C), but if there are, rotate or rock sample at a medium speed for an additional 2–3 min at 37°C, triturate again, and repeat process until no tissue chunks remain.
 - g. Add at least 8 mL ice cold 1 × HBSS (without Ca²⁺, Mg²⁺, or Phenol Red) to the sample to help slow the enzymatic reaction. If any tissue chunks remain, remove them by filtering the sample through a 40 μm cell strainer and into a new tube.
 - h. Centrifuge for 10 min at 500RCF at 4°C with acceleration and brake set to low settings.

⚠ **CRITICAL:** Mechanical dissociation in steps 7a and 7b should not be replaced with other methods, such as sonication or a Dounce homogenization. Sonication could damage cells, and we observed that cell preparations following Percoll gradient centrifugation contained significant debris with dounce homogenization that would be unsuitable for downstream CyTOF or CITE-seq experiments.

⚠ **CRITICAL:** Mechanical dissociation in steps 7a and 7b must be well done to ensure thorough dissociation and reduce overall enzymatic dissociation time. We have found inadequate mechanical dissociation significantly reduces the effectiveness of enzymatic digestion and final leukocyte yield. An example of the effect of insufficient mechanical dissociation on the final dissociated product (post-enzymatic dissociated) compared to ideal mechanical dissociation is shown in Figure 1C.

⚠ **CRITICAL:** During enzymatic dissociation in step 8e, ensure that the sample is continuously agitated. If the sample becomes stagnant due to either rotator/rocker failure, becoming lodged in the tip of the conical tube, etc., the dissociation rate is significantly reduced, potentially compromising cell viability and cell yield.

Note: We use a large volume of ice cold 1 × HBSS in step 8g to help stop the enzymatic reaction. We theorize this large volume should dilute the reaction and the cold temperature will reduce enzymatic activity. We caution against using serum-containing reagents to help stop the reaction as the serum could artificially stimulate leukocytes in the preparation.

Leukocyte purification by percoll gradient centrifugation

⌚ **Timing:** 35 min

This step will purify brain-infiltrating leukocytes, both enriching these cells in a condensed, extractable layer in the Percoll gradient while removing debris such as myelin (Figure 1D and Methods video S1).

- Aspirate supernatant from cell pellet.
- Resuspend cell pellet in 3.5 mL 70% Percoll by pipetting up and down. Bring total volume to 4 mL with 70% Percoll.
- Carefully layer a total of 4 mL of 37% Red Percoll on top of the 70% layer without allowing intermixing of the two Percoll densities. This can be achieved by using a p1000 micropipette and slowly dripping the 37% Percoll onto the sidewall of the conical tube while the conical tube is held at a slight (~45 degree) angle ([Methods video S1](#)).
- Carefully layer a total of 4 mL of 30% Percoll on top of the 37% Red Percoll layer without allowing intermixing of the two Percoll densities. This can be achieved as described above. The final layered solution should consist of three distinct layers as shown in [Figure 1D](#).
- Centrifuge sample in a swinging bucket centrifuge for 25 min at 800 RCF at room temperature (23°C–25°C) with the acceleration at the lowest setting and no break. Following centrifugation, there will be a buffy leukocyte layer at the interface of the 37% Red Percoll and clear 70% Percoll layers and a thick myelin layer at the top of the sample ([Figure 1D](#)).
- Using a vacuum aspirator connected to a p1000 micropipette tip, remove layers above buffy leukocyte layer at the interface of the 37% Red Percoll and 70% clear Percoll layers. Perform this step slowly, keeping the micropipette tip at the surface of the liquid (as opposed to submerged) to ensure the myelin layer is thoroughly removed and minimizing disturbance of the buffy leukocyte layer. Once within 1–2 mL of the buffy leukocyte layer, use a p1000 micropipette with a wide-bore tip to gently remove the remaining Percoll above the buffy leukocyte layer.
- Using a p1000 micropipette with a fresh tip, extract the buffy leukocyte layer and place it in a new 15 mL conical tube. No more than 2 mL of Percoll on the top or bottom of the buffy leukocyte layer should be extracted.
- Bring the volume of extracted cells up to 15 mL with cold 1× HBSS (without Ca²⁺, Mg²⁺, or Phenol Red) and centrifuge for 7 min at 350–400RCF at 4°C with the acceleration and brake set to medium settings.
- Aspirate supernatant from cell pellet, leaving as minimal of a volume of residual liquid as possible without disturbing the cell pellet. A small aliquot of cells must be taken out to verify the preparation largely consists of debris-free single cells and to obtain a cell count and examine viability ([Figure 1E](#) and [1F](#)). Proceed with the subsequent steps only if the cell preparation is a clean, single cell suspension with >80% viability.

△ **CRITICAL:** Centrifuge used for Percoll gradient centrifugation in step 15 must be set according to the indicated speed and temperature settings. A very low acceleration and no break will prevent early and late disturbances of the Percoll gradient, which could significantly reduce cell yield. Also, centrifugation must occur at room temperature as the density of Percoll would be changed at lower or higher temperatures, compromising cell yield.

△ **CRITICAL:** Once the Percoll gradient is fully layered in step 12, the sample should be centrifuged within 10 min as the 70% layer containing cells can begin to partition, slightly compromising cell yield. Prepared gradients should be maintained at room temperature before centrifugation.

△ **CRITICAL:** Minimize the amount of Percoll taken up when extracting the buffy leukocyte layer in step 15. Extracting excessive amounts of Percoll will make it impossible to pellet the cells, which is necessary to optimally remove Percoll.

Optional: Centrifugation time in step 15 can be reduced to 20 min with no appreciable effect on cell yield.

Optional: The 30% Percoll layer placed in step 12 can be omitted without compromising cell yield or purity. This layer is mostly used to extend the buffer space between the buffy leukocyte layer and myelin debris, reducing the chance of intermixing during leukocyte extraction.

Note: Related to steps 10 and 11, cells can be resuspended in 37% Percoll instead of 70% Percoll and layered as described above with minimal effect on cell yield. This method has been described and evaluated elsewhere ([Lee and Tansey, 2013](#)).

Leukocyte blocking, staining, and washing

⌚ Timing: 50 min

This step will block non-specific antibody binding and stain leukocytes with selected antibodies.

18. CITE-seq Applications

- a. Resuspend cell pellet in 50 μ L pre-prepared Fc Receptor block. Incubate on ice for 20 min.
- b. Prepare CITE-seq antibody cocktail as shown in [Figure 2A](#). Add 50 μ L antibody cocktail to sample, pipetting up and down to mix with the sample. Incubate on ice for 25 min.
- c. Add at least 5 mL CITE-seq staining buffer to sample and centrifuge for 5 min at 350–400RCF at 4°C with the acceleration and brake set to medium settings.
- d. Aspirate supernatant from the cell pellet. Resuspend in at least 5 mL CITE-seq Wash Buffer A and centrifuge for 5 min at 350–400RCF at 4°C with the acceleration and brake set to medium settings.
- e. Aspirate supernatant from the cell pellet. Resuspend in at least 5 mL CITE-seq Wash Buffer B and centrifuge for 5 min at 350–400RCF at 4°C with the acceleration and brake set to medium settings.
- f. Aspirate supernatant from the cell pellet. If multiple samples were prepared and hashed, combine individual samples in one tube. Resuspend individual sample or pooled samples in 500 μ L CITE-seq Wash Buffer C and aliquot \sim 40 μ L of the sample to perform a cell count. Add at least 4.5 mL CITE-seq Wash Buffer C and centrifuge sample for 5 min at 350–400RCF at 4°C with the acceleration and brake set to medium settings. While cells are pelleting, count the aliquoted cells.
- g. Resuspend cells in CITE-seq staining buffer at a concentration of 1,600 cells/ μ L for 10X superloading.

19. CyTOF applications

- a. Add at least 5 mL Maxpar PBS to cell pellet and centrifuge for 5 min at 350–400RCF at 4°C with the acceleration and brake set to medium settings.
- b. Aspirate supernatant from the cell pellet. Resuspend in 1 mL Maxpar PBS. Take a small aliquot of cells to count.
- c. Transfer an aliquot of 1 million cells or fewer in 100 μ L or less Maxpar PBS to a new 15 mL tube. Bring the total volume of cell suspension to 100 μ L with Maxpar Cell Staining Buffer.
- d. Add 100 μ L of pre-prepared 0.75 μ M Cell-ID cisplatin to the cell suspension and mix by pipetting. Incubate cell suspension in Cell-ID cisplatin for 5 min at room temperature.
- e. Add 2 mL Maxpar Cell Staining Buffer to cell suspension. Centrifuge for 5 min at 350–400RCF at 4°C with the acceleration and brake set to medium settings.
- f. Aspirate supernatant from the cell pellet. Resuspend cell pellet in 50 μ L pre-prepared Fc Receptor block. Incubate on ice for 20 min.
- g. Add 50 μ L or pre-prepared antibody cocktail to sample, pipetting up and down to mix with the sample. Incubate on ice for 25 min.
- h. Add at least 1 mL Maxpar Cell Staining Buffer to sample and centrifuge for 5 min at 350–400RCF at 4°C with the acceleration and brake set to medium settings.
- i. Aspirate supernatant from the cell pellet. Resuspend in 150 μ L Maxpar PBS and then add 150 μ L 3.2% PFA centrifuge for 5 min at 350–400RCF at 4°C with the acceleration and brake set to medium settings. Incubate at room temperature for 20 min.
- j. Add at least 1 mL Maxpar Cell Staining Buffer to sample and centrifuge for 5 min at 800RCF at room temperature with the acceleration and brake set to medium settings.

- k. Aspirate supernatant from the cell pellet. Resuspend in 300 μ L pre-prepared Cell-ID Cell-ID Intercalator-Ir to label all nucleated cells. Incubate at room temperature for 1 h or 4°C for 16 h.
- l. Add 1 mL Maxpar Cell Staining Buffer to sample and centrifuge for 5 min at 800RCF at room temperature.
- m. Aspirate supernatant from the cell pellet. Resuspend in 1 mL Maxpar Water and centrifuge for 5 min at 800RCF at room temperature.
- n. Aspirate supernatant from the cell pellet. Resuspend in 200 μ L Maxpar Water and aliquot a small volume of cell suspension to count. Bring total cell suspension volume up to 1 mL with Maxpar Water and centrifuge for 5 min at 800RCF at room temperature.
- o. Aspirate supernatant from the cell pellet. If running sample on a CyTOF instrument immediately, resuspend cells at a concentration of 500 cells/ μ L in Milli-Q water containing 0.1 \times EQ Beads.

Note: Samples prepared for CyTOF can be stored for days to weeks following the first wash post-Cell-ID Intercalator-Ir staining (step 19l). To store, aspirate supernatant following the first wash post-Cell-ID Intercalator-Ir staining, overlay the cell pellet with \sim 50 μ L Maxpar Water, seal the tube cap with parafilm, and store at 4°C.

Note: For multiplexed analysis, user should consider to barcode each sample using Cell-ID™ 20-Plex Pd Barcoding Kit available from Fluidigm.

10 \times genomics-based cell partitioning, library preparation, sequencing, and analysis

⌚ Timing: 1–2 days

This step will partition single cells, create a library for sequencing, and sequence the single cell library.

20. For 10 \times Genomics-based cell partitioning and library preparation, follow current versions of 10 \times Genomics protocols and CITE-seq protocols (10 \times Genomics cell partitioning and library preparation protocol: https://assets.ctfassets.net/an68im79xiti/4tjk4KvXzTWgTs8f3tvUjq/2259891d68c53693e753e1b45e42de2d/CG000183_ChromiumSingleCell3_v3_UG_Rev_C.pdf; CITE-seq without hashing library preparation protocol: https://citeseq.files.wordpress.com/2019/02/cite-seq_190213.pdf; CITE-seq with hashing library preparation protocol: https://citeseq.files.wordpress.com/2019/02/cite-seq_and_hashing_protocol_190213.pdf). The qualities of cDNA, ADT and HTO libraries should be verified using bioanalyzer before proceeding for sequencing. Examples of high quality cDNA, ADT and HTO libraries are presented in Figures 2B and 2C. Specifically, product peaks should be \sim 180 bp for ADT libraries and \sim 300–700 bp for mRNA libraries, and these should be the predominant peaks within each library (Figures 2B and 2C). If other product peaks are present, libraries may be re-cleaned with SPRI beads to eliminate unwanted fragments and retain the desired fragments.
21. Sequence the CITE-seq library to the desired sequence depth. In the case of CITE-seq, a critical consideration is the ratio of mRNA to antibody libraries sequenced to obtain optimal sequencing depth for both libraries. Although the amount of each library sequenced will vary based on several factors (for example: the flow cell used and the number of individual libraries if samples were hashed). To this end, we have routinely sequenced hashed samples on a NovaSeq platform. Sequencing a CITE-seq library consisting of 9 HTO-hashed/multiplexed individual samples (68% cDNA library, 25% ADT library, 7%, HTO library) on an NovaSeq S1 flow cell routinely provides the following reads: \sim 1000 \times 10⁶ reads for cDNA, \sim 400 \times 10⁶ reads for ADT, and \sim 100 \times 10⁶ reads for HTO. In our experience, with an expected recovery of \sim 16,000 single cells per NovaSeq S1 sequencing lane, such sequencing depth is sufficient to reach >4000 genes/cell and >50% cDNA sequencing saturation of each singlet CD45+ immune cell




Figure 3. RNA and antibody-based epitope expression in brain infiltrating leukocytes

(A) Example of gating leukocyte populations using CITE antibodies. Briefly, CNS native myeloid cells (microglia and BAM) were gated based on both populations lowly expressing CD45 and then differentiated from each other based on higher I-A-I-E and CD38 expression in BAM. Peripheral leukocytes were gated based on high CD45 expression. Within the peripheral leukocyte population: B cells were defined as having CD19 and B220 expression. T cells were defined based on high CD3 expression and low NK1.1 expression, then subpopulations further distinguished on the basis of CD4 and CD8 expression. Myeloid cells were distinguished from adaptive peripheral immune cells by high CD11b expression, and then further subdivided on the basis of Ly6C expression levels and Ly6G expression.

(B) Example of gating leukocyte populations using CyTOF antibodies.

Figure 3. Continued

(C) Biaxial plots of CITE epitope expression on the x axis and its corresponding mRNA expression on the y axis. Discordant expression of mRNA and protein can be observed.

(D) UMAPs derived from CITE-seq data clustered by RNA transcription cluster color coded by transcriptional cluster (left) or canonical cell ID determined by CITE-antibody gating (right). While canonically identified cells generally group within the same cluster, there is appreciable dispersion of canonical cell types among clusters, indicating the utility in using CITE antibody tags rather than mRNA alone to identify canonical cell types.

recovered (Figures 2C and 2D). Multiple CITE-seq libraries can be multiplexed and sequenced on Nova-seq S2 flow cell to further reduce the per cell sequencing cost.

22. To analyze data (analyze differential gene expression, generate dimensional reduction plots, etc.) a number of excellent options are available, including Seurat (Satija et al., 2015) (R coding (Ihaka and Gentleman, 1996)), Scanpy (Wolf et al., 2018) (Python coding), and Cytobank (cloud-based software program). A gating strategy for a brain leukocyte CITE antibody panel is displayed in Figure 3A.

Pause point: Amplified cDNA may be stored at -20°C until ready for final sequencing library preparation. Prepared sequencing library may be stored at -20°C until ready for sequencing.

Note: The upstream work of this protocol should be adaptable to most, if not all, forms of cell partitioning and library preparation.

Note: 10x Genomics regularly updates their reagent kits and protocols related to step 20. Be sure to check 10x Genomics reagent kits that will be used for your experiment and check the associated 10x Genomics protocol before proceeding with droplet generation or library preparation. Likewise, check for updates to the CITE-seq library preparation protocol, which is continuously being adjusted.

Note: When performing cell hashing/sample multiplexing using HTO antibodies to enable sample pooling, we routinely use the HTODemux function in the R package Seurat (Satija et al., 2015) to demultiplex HTO hashed samples as part of the analysis performed in step 22. However, we have observed a low frequency of HTO cross-contamination. Assignment of HTO to respective samples needs to be interpreted with caution and should be cross-validated with other computational methods.

CytoF data acquisition and analysis

⌚ **Timing:** ~2 h

This step will run cell suspensions on a CyTOF instrument to ultimately obtain a count matrix of the number of antibodies bound to each cell.

23. To run cell suspensions on a CyTOF instrument, follow the most up to date user manual of the specific CyTOF instrument. See <https://www.fluidigm.com/binaries/content/documents/fluidigm/consumables/pages/mass-cytometry/mass-cytometry/fluidigm%3Aresources/cytof-2-mass-cytometer/fluidigm%3Afile> in the case of operating a CyTOF2.
24. To analyze data (perform gating to analyze population differences and surface marker expression differences, generate dimensional reduction plots, etc.) a number of excellent options are available, including FlowJo and Cytobank. A gating strategy for a brain leukocyte CyTOF panel is displayed in Figure 3B.

EXPECTED OUTCOMES

Single cell profiling requires sample preparations of high quality (a single cell suspension with high cell viability and little debris) that are ideally enriched for the cell population(s) of interest. Adopting

this protocol, a debris-free, high viability single cell preparation enriched for leukocytes will be obtained (Figures 1E–1H), meeting the pre-requisite for successful single cell profiling. The number of leukocytes obtained will depend on the pathological state of the mouse, but should yield a minimum of $\sim 3 \times 10^5$ leukocytes regardless of state, with a viability >90% (Figure 1F). This cell yield and viability meet or exceed the minimum requirements of applications such as CyTOF and CITE-seq. Ultimately, leukocyte preparations derived from following this protocol will lead to exquisite single cell profiling in applications such as CITE-seq/REAP-seq and CyTOF, as demonstrated in our recent publications (Golomb et al., 2020; Guldner et al., 2020). Critically, we observed our brain dissociation protocol preserved many common immune cells identifying or functional surface protein epitopes, as evidenced by positive staining of >40 cell surface epitopes for which we stained in CITE-seq experiments, a subset of which are shown in Figures 3A–3C. When comparing mRNA expression and protein expression of a particular gene, we note a level or discordance between mRNA and protein expression levels (Figure 3C). We also observe that cells identified by use of surface marker expression via CITE antibodies show an appreciable level of dispersion across multiple RNA-based clusters on a UMAP (Figure 3D). Together, the data in Figures 3C and 3D demonstrate the importance of incorporating CITE antibodies into scRNA-seq experiments for proper identification of leukocytes based on the field standard of surface proteins, which may have expression levels discordant with their respective mRNAs.

QUANTIFICATION AND STATISTICAL ANALYSIS

Pearson correlation between the two features displayed on the biaxial plots in Figure 3C is presented in the top left corner of each plot. This analysis was performed in the R package, Seurat (Satija et al., 2015). This protocol is primarily aimed at describing a strategy for isolating brain infiltrating leukocytes for single cell profiling of transcriptomes and epitopes. For further information on statistical analyses for single cell RNA seq and CyTOF datasets, please refer to application specific guidelines.

LIMITATIONS

High-dimensional profiling strategies require cell preparations of high quality (a single cell suspension with high cell viability and little debris) that are ideally enriched for the cell population(s) of interest. Further, the single cell suspensions must be stained with optimized concentrations of antibodies to obtain interpretable results and reduce unnecessary overuse of expensive reagents. The protocol described herein has allowed us to prepare high quality single cell brain leukocyte preparations that have resulted in optimal CyTOF and CITE-seq profiling results. The key broad steps in our protocol include brain dissociation, leukocyte enrichment, and cell staining. Although this protocol has yielded results that meet or exceed current standards, like any methodology, this protocol is not without limitations.

The isolation of mouse brain-infiltrating leukocytes described herein has proven successful in our hands in the contexts of homeostasis, aging, and brain metastasis (Golomb et al., 2020; Guldner et al., 2020). Importantly, these represent conditions in which the brain has an abundance of myelin (aging) or large epithelial lesions (brain metastasis) that can present challenges to brain dissociation and the generation of single cell suspensions. Notably, despite the varying brain conditions we have studied, our protocol was adaptable to each scenario without modification. We have not tested this protocol in the context of protein aggregation pathologies or primary brain tumors, which may present unique challenges to dissociation and subsequent generation of single cell suspensions.

Brain dissociation is a critical step in this protocol that lays the foundation for successful downstream leukocyte enrichment, staining, and profiling. We have tested several mechanical (for example: sonication and Dounce homogenization) and enzymatic (for example: papain-based and collagenase-based) methods for brain dissociation, finding that a combination of mechanical (mincing and triturating) and collagenase-based enzymatic (Multi-Tissue Dissociation Kit I) methods provided optimal results. Importantly, we have optimized this dissociation method to not require any

specialized lab equipment. However, this method does have limitations in that (1) it is time consuming and (2) it could destroy some surface protein epitopes. In relation to time, 12 mouse brains can be processed using this protocol with two persons within an appropriate amount of time that does not substantially affect cell yield or viability. Unfortunately, besides the antibodies/surface epitopes we validated here, the resistance of epitopes to enzymatic dissociation must be determined empirically, but can easily be accomplished by less expensive methods such as flow cytometry using identical antibody clones. While manual dissociation alone can be accomplished faster and completely on ice, it could damage and kill cells and may not dissociate tissue thoroughly enough, impacting cell yield and sample cleanliness, both of which are important factors for applications such as CyTOF or scRNA-seq. Lastly, it is worth mentioning that enzyme-based dissociation that occurs at elevated temperatures, such as performed in this protocol at 37°C, could potentially artificially activate/alter the transcriptome of leukocytes (Mattei et al., 2020). If this is observed, we recommend exploring other dissociation options that can be accomplished at lower temperatures.

Even in the inflamed brain, leukocytes represent a small proportion of cells relative to neuronal and glial cells, thus making leukocyte enrichment a critical step. Percoll gradient centrifugation as performed in this protocol can remove myelin debris while confining leukocytes to a small, visible band within the Percoll gradient. This method is a relatively fast and efficient way to achieve a debris-free enriched leukocyte population, but could include some contaminating cells. In fact, we observed that brain metastasis cancer cells can be found within the leukocyte buffy layer and represent approximately 15% of the cells in this layer, and endothelial cells and pericytes can compose approximately 7% and 2% of cells in the buffy layer, respectively (Figures 1G and 1H). Still, neurons, astrocytes, and oligodendrocytes are largely absent from the leukocyte buffy layer, as evidenced by not being identified in sequencing data. Alternative enrichment methods include fluorescent activated cell sorting (FACS) or magnetic-based sorting strategies which have the potential to be more specific, but are at least as time consuming as Percoll gradient centrifugation, require the use of more specialized equipment, and not as adaptable to performing in a parallelized manner, which is important when several samples are being processed.

Single cell profiling is rapidly and continuously evolving, with existing methods being refined and new techniques being invented. Regardless of the changes that single cell profiling undergo, the requirement for high quality single cell preparation upstream of the profiling will remain constant. This protocol has been demonstrated to reliably dissociate the murine brain, whether homeostatic, aged, or tumor-burdened, into a single cell suspension, and highly enrich brain-infiltrating leukocytes. We have subsequently performed CyTOF and CITE-seq on leukocytes prepared based on the protocol described herein, achieving high quality results, as shown in the figures of this protocol, and our previous publications (Golomb et al., 2020; Guldner et al., 2020).

TROUBLESHOOTING

Problem 1

Chunks remain after 20–25 min of enzymatic dissociation (step 8).

Potential solution

Obtaining a tissue suspension free of tissue chunks after enzymatic dissociation requires sufficient upstream preparation, namely, very fine mincing of the brain and thorough trituration after mincing and in the middle of enzymatic dissociation. The brain should initially be minced and triturated upon extraction so that there is no resistance when pipetting the dissociated brain slurry using a P1000 micropipette. Additionally, trituration in the middle of enzymatic dissociation should be performed until nearly all brain chunks are broken up, and titration after enzymatic dissociation should be performed until no more than ~6 tiny brain chunks remain.

Problem 2

Percoll layers blend together during layering (steps 11 and 12).

Potential solution

Percoll layers tend to blend together due to intermixing during the process of Percoll layering. Slow, gentle placement of less dense Percoll layers on top of the denser Percoll layers is the key to prevent the blending of layers. Additionally, during layering, it can be helpful to (1) tilt tube at a 45 degree angle and (2) drip Percoll dropwise down the sidewall of the tube ~5 mm away from the surface of the Percoll already in the tube.

Problem 3

Accidental myelin contamination (step 15).

Potential solution

When aspirating myelin, the thick milky-colored layer at the top of all the Percoll layers, and the Percoll above the buffy leukocyte layer, if not all of the myelin was aspirated from the side wall of the conical tube, it may slide down the sidewall of the tube and mix with the buffy leukocyte layer. This issue primarily lies within insufficient removal of the myelin layer. To aid in complete myelin removal, either (1) take extra care to aspirate all of the myelin from the sidewall of the tube and/or (2) add 2 mL 1× HBSS above the 30% Percoll layer before centrifugation. This HBSS layer will remain on top of the myelin layer after centrifugation and aid in cleaner aspiration of myelin.

Problem 4

Cells do not pellet after first percoll wash (step 16).

Potential solution

If too much extra Percoll was extracted with the buffy leukocyte layer (step 14), there may be insufficient dilution of the Percoll in the subsequent washing step, resulting in a cushion of dense Percoll between the bottom of the conical tube and the cells after centrifugation. More specifically, there will be a buffy layer suspended a little above the tip of the conical tube rather than a cell pellet. To avoid this issue, take care not to extract more than ~3 mL of Percoll total. If this Percoll cushion forms, carefully extract the suspended cells (still look like a buffy layer) in as minimal of a volume as possible and either (1) re-wash or (2) proceed with the next step. As Percoll may interfere with subsequent steps, we advise to re-wash and ensure Percoll is largely removed from the cells.

Problem 5

Unexpected leukocyte activation (step 22).

Potential solution

If leukocytes, especially microglia from a naive brain, display an unexpectedly activated transcriptome, it may be due to any part of the procedure. It is important to work on ice when possible throughout this procedure, as heat could artificially activate leukocytes. If working on ice to the maximum possible does not resolve the artificial activation issue, other dissociation options that avoid incubation at higher temperature should be explored.

RESOURCE AVAILABILITY

Lead contact

Further information and requests for resources and reagents should be directed to and will be fulfilled by the lead contact, Siyuan Zhang (szhang8@nd.edu).

Materials availability

This study did not generate new materials.

Data and code availability

This protocol did not generate new data or employ custom code.

SUPPLEMENTAL INFORMATION

Supplemental information can be found online at <https://doi.org/10.1016/j.xpro.2021.100537>.

ACKNOWLEDGMENTS

This work was funded by NIH R01 CA194697-01 (S.Z.), NIH R01 CA222405 - 01A1 (S.Z.), Notre Dame CRND Catalyst Award (S.Z. and I.H.G.), and NIH CTSI core facility pilot grants (S.Z.). We would like to acknowledge and thank the Dee Family Endowment (S.Z.). We are grateful for the use of the following core facilities and their members' technical insight: Notre Dame Genomics and Bioinformatics Core Facility, Notre Dame Freimann Life Sciences Center, and Indiana University Simon Cancer Center Core Facility.

AUTHOR CONTRIBUTIONS

I.H.G., Q.W., and S.Z. designed experimental protocols. I.H.G., S.M.G., Q.W., E.W., and S.Z. performed experiments. I.H.G., S.M.G., E.W., and S.Z. analyzed data. I.H.G., S.M.G., E.W., and S.Z. wrote and revised the protocol. S.Z. supervised the study. All authors reviewed the manuscript.

DECLARATION OF INTERESTS

The authors declare no competing interests.

REFERENCES

- Bandura, D.R., Baranov, V.I., Ornatsky, O.I., Antonov, A., Kinach, R., Lou, X., Pavlov, S., Vorobiev, S., Dick, J.E., and Tanner, S.D. (2009). Mass cytometry: technique for real time single cell multitarget immunoassay based on inductively coupled plasma time-of-flight mass spectrometry. *Anal. Chem.* *81*, 6813–6822.
- Golomb, S.M., Guldner, I.H., Zhao, A., Wang, Q., Palakurthi, B., Aleksandrovic, E.A., Lopez, J.A., Lee, S.W., Yang, K., and Zhang, S. (2020). Multi-modal single-cell analysis reveals brain immune landscape plasticity during aging and gut microbiota dysbiosis. *Cell Rep.* *33*, 108438.
- Guldner, I.H., Wang, Q., Yang, L., Golomb, S.M., Zhao, Z., Lopez, J.A., Brunory, A., Howe, E.N., Zhang, Y., Palakurthi, B., et al. (2020). CNS-native myeloid cells drive immune suppression in the brain metastatic niche through Cxcl10. *Cell* *183*, 1234–1248.e25.
- Ihaka, R., and Gentleman, R. (1996). R: a language for data analysis and graphics. *J. Comput. Graph. Stat.* *5*, 299–314.
- Lee, J.-K., and Tansey, M.G. (2013). Microglia isolation from adult mouse brain. *Methods Mol. Biol.* *1041*, 17–23.
- Mattei, D., Ivanov, A., van Oostrum, M., Pantelyushin, S., Richetto, J., Mueller, F., Beffinger, M., Schellhammer, L., Vom Berg, J., Wollscheid, B., et al. (2020). Enzymatic dissociation induces transcriptional and proteotype bias in brain cell populations. *Int. J. Mol. Sci.* *21*.
- Peterson, V.M., Zhang, K.X., Kumar, N., Wong, J., Li, L., Wilson, D.C., Moore, R., McClanahan, T.K., Sadekova, S., and Klappenbach, J.A. (2017). Multiplexed quantification of proteins and transcripts in single cells. *Nat. Biotechnol.* *35*, 936–939.
- Satija, R., Farrell, J.A., Gennert, D., Schier, A.F., and Regev, A. (2015). Spatial reconstruction of single-cell gene expression data. *Nat. Biotechnol.* *33*, 495–502.
- Stoeckius, M., Hafemeister, C., Stephenson, W., Houck-Loomis, B., Chattopadhyay, P.K., Swerdlow, H., Satija, R., and Smibert, P. (2017). Simultaneous epitope and transcriptome measurement in single cells. *Nat. Methods* *14*, 865–868.
- Wolf, F.A., Angerer, P., and Theis, F.J. (2018). SCANPY: large-scale single-cell gene expression data analysis. *Genome Biol.* *19*, 15.