



Carrier-assisted One-pot Sample Preparation for Targeted Proteomics Analysis of Small Numbers of Human Cells

Kendall Martin¹, Tong Zhang¹, Pengfei Zhang², William B. Chrisler¹, Fillmore L. Thomas³, Fen Liu⁴, Tao Liu¹, Wei-Jun Qian¹, Richard D. Smith¹, Tujin Shi¹

¹Biological Sciences Division, Pacific Northwest National Laboratory, Richland, Washington 99354, USA

²NHC Key Laboratory of Cancer Proteomics, Department of Oncology, Xiangya Hospital, Central South University, Changsha, Hunan, 410008, P.R. China

³Environmental Molecular Sciences Laboratory, Pacific Northwest National Laboratory, Richland, Washington 99354, USA

⁴Department of Critical Care Medicine, The First Affiliated Hospital of Nanchang University, Nanchang, Jiangxi, 330006, P.R. China

Abstract

Protein analysis of small numbers of human cells and single cells is primarily achieved by targeted proteomics with antibody-based immunoassays, whereas they have inherent limitations (e.g., low multiplex and unavailability of antibodies for new proteins). Mass spectrometry (MS)-based targeted proteomics has emerged as an alternative in terms of being antibody-free, high multiplex, high specificity, and quantitation accuracy. Recent advances in MS instrumentation makes MS-based targeted single-cell proteomics possible for multiplexed quantification of highly abundant proteins in single cells. However, there is a technical issue for effective processing of single cells with minimal sample loss for MS analysis. To address this issue we have recently developed a simple protein carrier-assisted one-pot sample preparation coupled with liquid chromatography (LC) – selected reaction monitoring (SRM) termed cLC-SRM for targeted proteomics analysis of small numbers of human cells including single cells. This method capitalizes on using the combined excessive exogenous protein as a carrier and low-volume one-pot processing to greatly reduce surface adsorption losses and high-specificity LC-SRM to effectively address the increased dynamic concentration range resulted from the addition of exogeneous carrier protein. Herein we have provided a detailed protocol for cLC-SRM analysis of small numbers of human cells including cell sorting, cell lysis and digestion, LC-SRM analysis, and data analysis. cLC-SRM was demonstrated to enable multiplexed accurate quantification of most moderately abundant proteins in small numbers of cells (e.g., 10–100 cells) and highly abundant proteins in single cells. Most importantly, this method can be easily implemented in any MS and proteomics laboratories at no additional cost for instrument or reagents. Further improvements in detection sensitivity and sample throughput are needed towards multiplexed targeted single-cell proteomics analysis. We

Corresponding author: Tujin Shi (tujin.shi@pnl.gov).

Disclosures

The authors have nothing to disclose.

anticipate that eventually it can be broadly applied to biomedical research and systems biology with the potential of facilitating precision medicine.

Summary

We recently developed a protein carrier-assisted one-pot sample preparation coupled to Liquid Chromatography (LC) – Selected Reaction Monitoring (SRM) termed cLC-SRM for multiplexed targeted proteomics analysis of small numbers of cells including single cells. This method capitalizes on using excessive exogenous protein as a carrier and high-specificity LC-SRM for targeted quantification.

Keywords

Protein Carrier; Carrier-assisted One-pot Sample Preparation; cLC-SRM; Small Numbers of Cells; Single Cells; Mass-limited Sample; PCR Tube; Targeted Single-cell Proteomics

Introduction

Recent technological advances in genomics (transcriptomics) allow for precise measurements of whole genome (transcriptome) in single cells^{1–3}. However, single-cell proteomics technologies are lagging far behind but as equally important as genomics (transcriptomics) technologies^{4–8}. Furthermore, protein abundance cannot necessarily be inferred from mRNA abundance⁹, and proteome is more complex and dynamic than transcriptome¹⁰. A large number of mixed populations of cells (i.e., bulk cells) are generally used to generate comprehensive proteome data^{11–13}. However, such bulk measurements average out stochastic variations of individual cells, thus obscure important cell-to-cell variability (i.e., cell heterogeneity)^{4,14}. Limitations of such bulk measurements become even more severe when the cells of interest only account for a small portion of the total populations of cells (e.g., cancer stem cells within tumors at an early-stage cancer). Therefore, there is a huge knowledge gap between single-cell proteomics and genomics (transcriptomics).

Antibody-based immunoassays (e.g., flow or mass cytometry) are predominantly used for targeted proteomics analysis of single cells^{6,7,15–18}, whereas they share common shortcomings with other antibody-based methods (e.g., low multiplex and unavailability of antibodies for new proteins). Mass spectrometry (MS)-based targeted proteomics has emerged as an alternative for accurate protein quantification in terms of being antibody-free, high multiplex (200 proteins in a single analysis¹⁹), high quantitation accuracy (absolute amounts or concentrations), high specificity and reproducibility (10% CV)^{20–23}. However, MS-based single-cell proteomics is still at the early infancy stage. Recent significant progress in sample preparation has made MS-based single-cell proteomics possible for quantitative analysis of highly abundant proteins from single human cells. For example, the most advanced MS platform coupled with ultralow-flow RPLC flow rates can only allow label-free MS detection and quantification of ~670–870 proteins out of the total 12,000 proteins in single HeLa cells^{24,25}.

Currently, there are five MS-based single-cell proteomics approaches available, in which three are for global proteomics (nanoPOTS: nanowell-based Preparation in One pot for Trace Samples²⁶; iPAD-1: integrated proteome analysis device for single-cell analysis²⁷; SCoPE-MS: single cell proteomics by mass spectrometry²⁸) and the other two are for targeted proteomics (cLC-SRM: carrier-assisted liquid chromatography (LC) – selected reaction monitoring (SRM)²⁹; cPRISM-SRM: carrier-assisted high pressure, high-resolution separations with intelligent selection and multiplexing coupled to SRM³⁰). However, all these approaches have inherent technical drawbacks. nanoPOTS and iPAD-1, which downscale sample processing volume to 2–200 nL, are not easily adoptable for broad benchtop applications^{26,27}. For SCoPE-MS, a TMT (tandem mass tag) carrier is added after single-cell processing so it cannot effectively prevent surface adsorption losses during sample processing²⁸, resulting in low reproducibility with a correlation coefficient of only ~0.2–0.5 between replicates³¹. Furthermore, due to the inability to fractionate ultrasmall single-cell samples TMT-based single-cell analysis lacks quantitation accuracy from ratio compression or distortion caused by coeluting interferences³². For cLC-SRM and cPRISM-SRM, using exogenous proteins as a carrier is more suitable for targeted proteomics because peptides from excessive exogenous proteins are frequently sequenced by MS/MS, which greatly reduces the chance for sequencing low abundant endogenous peptides^{29,30}. Unlike global single-cell proteomics with relative quantification, the two targeted proteomics approaches can provide accurate or absolute protein analysis of single cells. Most importantly, they can be easily implemented in any MS and proteomics laboratories at no additional cost for instrument or reagents. In contrast to cPRISM-SRM that requires using high-resolution RPLC for PRISM fractionation, cLC-SRM has a significant advantage in multiplexing (100s of proteins in a single analysis) but with slightly lower detection sensitivity²⁹. Therefore, cLC-SRM is more readily to access and should have more broad utilities for accurate multiplexed protein analysis of small numbers of cells as well as mass-limited precious clinical specimens.

Herein we describe a detailed protocol to perform cLC-SRM for convenient targeted proteomics analysis of small numbers of human cells including single cells. The protocol consists of the following major steps: cell sorting by FACS (fluorescence activated cell sorting), cell lysis and digestion processed in low-volume single PCR tube, LC-SRM data collection, and SRM data analysis using publicly available Skyline software (Figure 1). Its broad utility was demonstrated by absolute targeted quantification of EGFR/MAPK pathway proteins in 1–100 MCF7 or MCF10A cells along with our previously well-established SRM assays and determined pathway protein copies per cell at a wide dynamic range of concentrations²⁹. We anticipate that with the detailed protocol most proteomics researchers can readily implement cLC-SRM in their laboratories for accurate protein analysis of ultrasmall samples (e.g., rare tumor cells) to meet their project needs.

Protocol

The step-by-step cLC-SRM analysis is shown in Figure 1. Below is the detailed protocol:

1. Pretreatment of PCR tubes

1. Add 100 μL of nonhuman (e.g., *S. oneidensis*) cell lysate digests at 0.2 $\mu\text{g}/\mu\text{L}$ to PCR tubes. Incubate at room temperature for overnight to coat PCR tube surface.
2. Remove the cell lysate digests, rinse PCR tubes with HPLC-grade water for 3 times, and then airdry PCR tubes in a fume hood.
3. Store coated PCR tubes into 4 $^{\circ}\text{C}$ freezer till further use.

2. FACS sorting

1. Align a fluorescence-activated cell sorter (FACS) (e.g., BD Influx flow cytometer) into a Hard-Shell 96-Well PCR Plates using fluorescent beads.
2. Verify that the FACS machine is well aimed for sorting cells into the bottom of PCR plates: 1) Set the angle of the sort stream to the minimum possible, which increases the chance that the cells will be sorted into the bottoms of wells rather than hit the sides; 2) To aim the machine, use the test stream to ‘sort’ droplets of PBS buffer (e.g., 100 droplets) onto the well of empty plate to make sure that the collected droplets are on the bottom. If the droplets do not land onto the right position, recalibrate the cell sorter and repeat until the droplets in the sort stream are deposited correctly.
3. Place the coated PCR tubes into 96-well PCR tube rack with caps open.
4. Sort of desired numbers of cells into precoated PCR tubes.
5. Immediately centrifuge the sorted cells at 850 rpm for 10 min at 4 $^{\circ}\text{C}$ to keep them at the bottom of the tube to avoid potential cell loss.
6. Store the sorted cells into -80°C freezer until further analysis.

3. Addition of protein carrier, heavy internal standards, and TFE

1. Add 4 μL of 25 mM NH_4HCO_3 into PCR tubes containing collected cells.
2. Add 1 μL of 10 ng/ μL BSA in 25 mM NH_4HCO_3 into sample tubes.
3. Add 0.3 μL of 100 fmol/ μL crude heavy peptide standards (total 30 fmol) into sample tubes.

Note: Dispense of the small volume of crude heavy peptide standards into the middle of the solution and make sure there is no left-over liquid in the pipet tip. High-purity light peptides (>95%) are used to calibrate crude heavy peptide concentrations for absolute quantification.

4. Add 9 μL of 100% TFE into sample tubes with the final ~60% TFE.
5. Centrifuge at 5000 rpm for 5 mins and then gently mixed at 850 rpm for 3 mins.

4. Cell lysis and protein denaturation

1. Sonicate samples on ice for 1 min with 1-min interval for a total of 5 cycles

2. Incubate samples at incubated at 75–90 °C for 1 h for protein denaturation using the RT-PCR thermocycler with the heated-lid option.
3. Cool samples to room temperature with centrifugation at 5000 rpm for 3 mins.

5. Reduction and alkylation (optional)

1. Add 0.6 µL of 50 mM DTT with the final concentration of 2 mM.
Note: 10 µL of 500mM DTT is added into 90 µL of 25 mM NH₄HCO₃ to make 50 mM DTT.
2. Centrifuge at 5000 rpm for 5 mins and then gently mixed at 850 rpm for 3 mins.
3. Incubate at 56°C for 1 h with gentle shaking at 850 rpm.
4. Cool samples to room temperature and centrifuge at 5000 rpm for 3 mins.
5. Add 0.5 µL of 60 mM IAA was added to the RT-PCR tube with the final concentration of ~2 mM.
Note: 15 µL of 400 mM IAA is added into 85 µL of 25 mM NH₄HCO₃ to make 60 mM IAA.
6. Incubate in the dark at room temperature for 30 mins with gentle shaking at 850 rpm.

6. Trypsin digestion

1. Reduce the sample volume down to ~4 µL for greatly reducing TFE volume using a Speed Vac concentrator.
Note: TFE is a volatile solvent with boiling point of 74°C. High percentage TFE inhabits trypsin digestion.
2. Add 9 µL of 25 mM NH₄HCO₃ and 1–3 µL of 15 ng/µL trypsin for digestion with the final trypsin concentration 1–3 ng/µL.
3. Note: When 10 human cells (~1 ng) are digested, the ratio of trypsin enzyme over protein is 15:1, 750-fold higher than that in standard trypsin digestion (the ratio of 1:50).
4. Mix gently at 850 rpm for 3 mins and incubate for overnight (~16 hours) at 37 °C.
5. Add 0.5 µL of 5% formic acid to stop enzymatic reaction.
6. Centrifuge for 1 h at 5000 rpm.
7. Store in –80 °C freezer until further LC-SRM analysis or in 4 °C freezer for immediate analysis.

7. Preparation for direct LC-SRM analysis

1. Remove the cap of PCR tube.
2. Insert PCR tube into LC vial.

3. Close the LC vial for LC-SRM analysis.

Note: For small numbers of cells processed in 96-well PCR plate, the mat cover will be added on the top of plate for direct LC-SRM analysis.

8. LC-SRM analysis

1. Analyze samples using a nanoACQUITY UPLC® (Waters Corporation, Milford, MA) coupled to a TSQ Vantage triple quadrupole mass spectrometer (Thermo Scientific, San Jose, CA).
2. The Waters BEH 1.7 μm C18 columns (75 μm i.d. \times 20 cm for standard gradient and 100 μm i.d. \times 10 cm for short gradient) are connected to a chemically etched 20 μm i.d. fused silica electrospray emitter via a stainless metal union.
3. Replace 5 μL sample loop with 20 μL sample loop for directly loading all the samples (~15 μL in total) into LC column to maximize detection sensitivity.
4. Use 0.1% formic acid in water and 0.1% formic acid in 90% acetonitrile as mobile phases A and B for capillary RPLC separation.
5. The binary LC gradient is used at flow rates of 300 nL/min for standard gradient and 400 nL/min for short gradient. Standard gradient: 5–20% B in 26 mins, 20–25% B in 10 mins, 25–40% B in 8 mins, 40–95% B in 1 min and at 95% B for 7 mins for a total of 52 mins and the analytical column re-equilibrated at 99.5% A for 8 mins. Short gradient for fast separation: 5–95% B in 5 mins and the analytical column re-equilibrated at 99.5% A for 5 mins for a total of 10 mins.
6. The TSQ Vantage mass spectrometer is operated with ion spray voltages of 2,400 \pm 100 V, a capillary offset voltage of 35 V, a skimmer offset voltage of –5 V and a capillary inlet temperature of 220° C.
7. The tube lens voltages are obtained from automatic tuning and calibration without further optimization.
8. Before LC-SRM analysis of samples, different concentrations of crude heavy peptide standards are spiked into 20–50 ng/ μL of nonhuman cell lysate digests to determine peptide LC retention time (RT) for building scheduled SRM method and potential interference transitions for accurate quantification.
9. The RT scheduled SRM mode is used for multiplexed quantification with the scan window of 6 mins.

Note: The scan window will be adjusted accordingly depending on the number of transitions in a certain time window.

10. The cycle time is set to 1 s, and the dwell time for each transition is automatically adjusted depending on the number of transitions scanned at different retention time windows. A minimal dwell time 10 ms is used for each SRM transition. All target proteins can be simultaneously monitored in a single LC-SRM analysis

9. Data Analysis

1. Import raw data files into publicly available Skyline software³³ for data visualization of each target peptide to determine their detectability.
2. Use the best transition without matrix interference for SRM quantification.
3. Use two criteria to determine SRM peak detection and integration: 1) same retention time; 2) approximately the same relative SRM peak intensity ratios across multiple transitions between endogenous (light) peptides and heavy peptide internal standards.
4. Use multiple technical or biological replicates to obtain the standard deviation (SD) and coefficient of variation (CV).
5. Inspect all the SRM data manually to ensure correct SRM peak assignment and SRM peak boundary for reliable quantification.
6. The signal-to-noise ratio (S/N) is calculated by the peak apex intensity over the highest background noise within a retention time region of ± 15 s for standard gradient and ± 6 s for short gradient. The limit of detection (LOD) and the limit of quantification (LOQ) are defined as the lowest concentration points of each target protein at which the S/N ratio of surrogate peptides was at least 3 and 7, respectively.
7. For conservatively determining the LOQ values, in addition to the requirement of S/N ≥ 7 , one additional criteria is applied: surrogate peptide response over different cell numbers must be within the linear dynamic range.
8. All calibration curves are plotted using Microsoft Excel 2010. The RAW data from TSQ Vantage are loaded into Skyline software to display graphs of extracted ion chromatograms (XICs) of multiple transitions of target proteins monitored.

Representative results

The performance of cLC-SRM was first evaluated by targeted quantification of EGFR/MAPK pathway proteins in 0.5–20 ng of MCF7 cell lysates (equivalent to 5–200 cells). As shown in Figure 2A, extracted ion chromatograms (XICs) clearly shows transitions for ATADDELSFK derived from GRB2 present at $\sim 220,000$ copies per MCF7 cell³⁴. cLC-SRM was able to reproducibly quantify endogenous ATADDELSFK down to 5 MCF7 cell equivalents with an S/N ratio of 14 and $\sim 1,800$ zmol of quantification sensitivity. The resultant calibration curves displayed excellent linearity over a wide dynamic range of MCF7 cell equivalents with LOQs of ~ 5 –10 cells for high-abundance proteins and ~ 10 –20 cells for moderate-abundance proteins (Figure 2B). The median SRM technical CV for all target peptides across all the data points was $\sim 9\%$, consistent with the technical reproducibility of well-characterized LC-SRM with CV below 10%.^{20,23,35–38}

We next applied cLC-SRM to measure EGFR pathway proteins in 10 and 100 intact MCF7 cells sorted by FACS. cLC-SRM enabled detection of high- and moderate-abundance

proteins in 10 intact MCF7 cells (Figure 2C). Even low-abundance AKT1 at ~3700 copies per cell was also detected with the average S/N ratio of 7, suggesting ~60 zmol of absolute sensitivity of cLC-SRM. With the cell number increased to 100, a majority of previously identified important EGFR pathway proteins (22 out of the total 32 proteins)³⁴ were reliably detected and quantified by cLC-SRM. Furthermore, we have recently tested whether a short gradient time (e.g., 5 mins vs standard 45 mins) is sufficient for rapid cLC-SRM analysis. As shown in Figure 3, XICs clearly shows endogenous detection of VLTPTQVK peptide derived from PEBP1 at ~744,000 copies per cell³⁴ in single MCF10A cells sorted by FACS with an S/N ratio of 5 and ~1,240 zmol of quantification sensitivity. As expected, with the cell number increased to 50 and 75 SRM signal became much stronger with detection of all three transitions which have the same pattern as its corresponding heavy internal standard. All these results have shown that cLC-SRM can be used for multiplexed, sensitive, absolute quantification of target proteins in small numbers of human cells including single cells. A short gradient is feasible for cLC-SRM presumably because sample complexity for 20 ng of tryptic peptides (10 ng BSA and 100 human cells \approx 10 ng proteins) from carrier-assisted small numbers of cells can be effectively address by high-resolution capillary RPLC separation with high loading capacity of 200 ng.

Discussion

cLC-SRM is a convenient targeted proteomics method that enables accurate multiplexed protein analysis of small numbers of cells including single cells. This method capitalizes on protein carrier-assisted one-pot sample preparation that combines all steps into one pot (e.g., single tube or single well) including cell collection, multistep cell lysis and digestion, transfer of peptide digests to capillary LC column for MS analysis (Figure 1). This ‘all-in-one’ low-volume one-pot processing presumably maximizes recovery of small numbers of cells for quantitative targeted proteomics analysis by greatly reducing possible surface absorption losses.

There are two critical steps for cLC-SRM analysis: 1) after FACS sorting immediate centrifugation is required to ensure collected cells at the bottom of tubes or wells because low volume (~15 μ L) is used to process small numbers of cells; 2) prior to the addition of trypsin, it was necessary to reduce the TFE levels to 10% for effective trypsin digestion. To avoid drying out samples with unrecoverable loss, it was necessary to frequently check the sample volume (~4 μ L of the remaining volume) during Speed Vac concentrating. One-pot sample preparation can be further simplified by removal of reduction and alkylation steps with neglectable effect on digestion efficiency (Figure 3). One drawback is that cysteine containing peptides cannot be selected as surrogate peptides for target proteins because of instability of free cysteine. To avoid Speed Vac concentrating step for TFE removal we consider replacing TFE with MS-compatible surfactants (e.g., DDM: n-Dodecyl β -D-maltoside) for cell lysis because frequent checking sample volume severely prevents automation of sample processing with low throughput. In addition, single carrier protein (e.g., BSA) may not be sufficient to unbiasedly reduce surface adsorption losses for every protein. This was evidenced by a few undetected proteins at >30,000 copies per cell (e.g., MAP2K1) in 10 MCF7 cells that should have been detected²⁹. Nonhuman cell lysates presumably can be used as a more effective proteome carrier in contrast to our current BSA

carrier because there are many different types of nonhuman proteins with less interference for SRM detection of human proteins.

When compared to other available MS-based single-cell proteomics methods that require specific devices and only provide relative quantification^{26,27}, cLC-SRM is readily adopted and implemented in any MS proteomics laboratory without additional costs and can provide accurate or absolute quantification of target proteins in small numbers of cells including single cells. But unlike other global single-cell proteomics methods, cLC-SRM is not suitable for global proteome profiling. Furthermore, cLC-SRM has significant advantages over conventional targeted single-cell proteomics using antibody-based immunoassays in terms of multiplexing (200 proteins in a single analysis) and quantitation specificity and accuracy. Future developments will focus on significant improvements in detection sensitivity and sample throughput for enabling rapid absolute quantification of most proteins (e.g., important negative feedback regulators) in single human cells. Automation of cLC-SRM is another direction to improve its robustness, reproducibility, and the overall sample throughput. We anticipate that cLC-SRM can be broadly used for multiplexed accurate quantification of target proteins in small subpopulations of cells and rare tumor cells as well as ultrasmall clinical specimens. In turn, it will be greatly beneficial to current biomedical research and systems biology.

Acknowledgements

This work was supported by NIH Grant R21CA223715 (TS). Dr. Pengfei Zhang acknowledges the fellowship support from China Scholarship Council. The experimental work described herein was performed in the Environmental Molecular Sciences Laboratory, Pacific Northwest National Laboratory, a national scientific user facility sponsored by the United States of America Department of Energy under Contract DE-AC05-76RL0 1830.

References

1. Wang Z, Gerstein M & Snyder M RNA-Seq: a revolutionary tool for transcriptomics. *Nat Rev Genet* 10 (1), 57–63, (2009). [PubMed: 19015660]
2. Miyamoto DT et al. RNA-Seq of single prostate CTCs implicates noncanonical Wnt signaling in antiandrogen resistance. *Science* 349 (6254), 1351–1356, (2015). [PubMed: 26383955]
3. Wang Y & Navin NE Advances and applications of single-cell sequencing technologies. *Mol Cell* 58 (4), 598–609, (2015). [PubMed: 26000845]
4. Macaulay IC, Ponting CP & Voet T Single-Cell Multiomics: Multiple Measurements from Single Cells. *Trends Genet* 33 (2), 155–168, (2017). [PubMed: 28089370]
5. Gavasso S, Gullaksen SE, Skavland J & Gjertsen BT Single-cell proteomics: potential implications for cancer diagnostics. *Expert Rev Mol Diagn* 16 (5), 579–589, (2016). [PubMed: 26895397]
6. Wu M & Singh AK Single-cell protein analysis. *Curr Opin Biotechnol* 23 (1), 83–88, (2012). [PubMed: 22189001]
7. Hughes AJ et al. Single-cell western blotting. *Nature Methods* 11 (7), 749–U794, (2014). [PubMed: 24880876]
8. Peterson VM et al. Multiplexed quantification of proteins and transcripts in single cells. *Nat Biotechnol* 10.1038/nbt.3973, (2017).
9. Liu Y, Beyer A & Aebersold R On the Dependency of Cellular Protein Levels on mRNA Abundance. *Cell* 165 (3), 535–550, (2016). [PubMed: 27104977]
10. Choudhary C & Mann M Decoding signalling networks by mass spectrometry-based proteomics. *Nat Rev Mol Cell Biol* 11 (6), 427–439, (2010). [PubMed: 20461098]

11. Mertins P et al. Proteogenomics connects somatic mutations to signalling in breast cancer. *Nature* 534 (7605), 55–62, (2016). [PubMed: 27251275]
12. Zhang H et al. Integrated Proteogenomic Characterization of Human High-Grade Serous Ovarian Cancer. *Cell* 166 (3), 755–765, (2016). [PubMed: 27372738]
13. Zhang B et al. Proteogenomic characterization of human colon and rectal cancer. *Nature* 513 (7518), 382–387, (2014). [PubMed: 25043054]
14. Gaudet S & Miller-Jensen K Redefining Signaling Pathways with an Expanding Single-Cell Toolbox. *Trends Biotechnol* 34 (6), 458–469, (2016). [PubMed: 26968612]
15. Giesen C et al. Highly multiplexed imaging of tumor tissues with subcellular resolution by mass cytometry. *Nat Methods* 11 (4), 417–422, (2014). [PubMed: 24584193]
16. Di Palma S & Bodenmiller B Unraveling cell populations in tumors by single-cell mass cytometry. *Curr Opin Biotechnol* 31 122–129, (2015). [PubMed: 25123841]
17. Schulz KR, Danna EA, Krutzik PO & Nolan GP Single-cell phospho-protein analysis by flow cytometry. *Curr Protoc Immunol Chapter 8 Unit 8 17 11–20*, (2012). [PubMed: 22855358]
18. Willison KR & Klug DR Quantitative single cell and single molecule proteomics for clinical studies. *Curr Opin Biotechnol* 24 (4), 745–751, (2013). [PubMed: 23810371]
19. Lee JY et al. Detection of Head and Neck Cancer Based on Longitudinal Changes in Serum Protein Abundance. *Cancer Epidemiol Biomarkers Prev* 10.1038/msb4100014 doi: 10.1158/1055-9965.EPI-1120-0192., (2020).
20. Shi T et al. Advancing the sensitivity of selected reaction monitoring-based targeted quantitative proteomics. *Proteomics* 12 (8), 1074–1092, (2012). [PubMed: 22577010]
21. Shi T et al. Advances in targeted proteomics and applications to biomedical research. *Proteomics* 16 (15–16), 2160–2182, (2016). [PubMed: 27302376]
22. Lange V, Picotti P, Domon B & Aebersold R Selected reaction monitoring for quantitative proteomics: a tutorial. *Mol Syst Biol* 4 222, (2008). [PubMed: 18854821]
23. Picotti P & Aebersold R Selected reaction monitoring-based proteomics: workflows, potential, pitfalls and future directions. *Nat Meth* 9 (6), 555–566, (2012).
24. Zhu Y et al. Proteomic Analysis of Single Mammalian Cells Enabled by Microfluidic Nanodroplet Sample Preparation and Ultrasensitive NanoLC-MS. *Angew Chem Int Ed Engl* 57 (38), 12370–12374, (2018). [PubMed: 29797682]
25. Cong Y et al. Improved Single-Cell Proteome Coverage Using Narrow-Bore Packed NanoLC Columns and Ultrasensitive Mass Spectrometry. *Anal Chem* 92 (3), 2665–2671, (2020). [PubMed: 31913019]
26. Zhu Y et al. Nanodroplet processing platform for deep and quantitative proteome profiling of 10–100 mammalian cells. *Nat Commun* 9 (1), 882, (2018). [PubMed: 29491378]
27. Shao X et al. Integrated Proteome Analysis Device for Fast Single-Cell Protein Profiling. *Anal Chem* 90 (23), 14003–14010, (2018). [PubMed: 30375851]
28. Budnik B, Levy E, Harmange G & Slavov N SCoPE-MS: mass spectrometry of single mammalian cells quantifies proteome heterogeneity during cell differentiation. *Genome Biol* 19 (1), 161, (2018). [PubMed: 30343672]
29. Zhang P et al. Carrier-Assisted Single-Tube Processing Approach for Targeted Proteomics Analysis of Low Numbers of Mammalian Cells. *Anal Chem* 91 (2), 1441–1451, (2019). [PubMed: 30557009]
30. Shi T et al. Facile carrier-assisted targeted mass spectrometric approach for proteomic analysis of low numbers of mammalian cells. *Commun Biol* 1 103, (2018). [PubMed: 30271983]
31. Vitrinel B, Iannitelli DE, Mazzoni EO, Christiaen L & Vogel C A simple method to quantify proteins from one thousand cells. *bioRxiv* doi: 10.1101/753582 (2019).
32. Rauniyar N & Yates JR 3rd. Isobaric labeling-based relative quantification in shotgun proteomics. *J Proteome Res* 13 (12), 5293–5309, (2014). [PubMed: 25337643]
33. MacLean B et al. Skyline: an open source document editor for creating and analyzing targeted proteomics experiments. *Bioinformatics* 26 (7), 966–968, (2010). [PubMed: 20147306]
34. Shi T et al. Conservation of protein abundance patterns reveals the regulatory architecture of the EGFR-MAPK pathway. *Sci Signal* 9 (436), rs6, (2016). [PubMed: 27405981]

35. Addona TA et al. Multi-site assessment of the precision and reproducibility of multiple reaction monitoring-based measurements of proteins in plasma. *Nat Biotechnol* 27 (7), 633–641, (2009). [PubMed: 19561596]
36. He J et al. Analytical platform evaluation for quantification of ERG in prostate cancer using protein and mRNA detection methods. *J Transl Med* 13 54, (2015). [PubMed: 25889691]
37. Shi T et al. Antibody-free, targeted mass-spectrometric approach for quantification of proteins at low picogram per milliliter levels in human plasma/serum. *Proc Natl Acad Sci U S A* 109 (38), 15395–15400, (2012). [PubMed: 22949669]
38. Shi T et al. Targeted quantification of low ng/mL level proteins in human serum without immunoaffinity depletion. *J Proteome Res* 12 (7), 3353–3361, (2013). [PubMed: 23763644]




Figure 1: Overview of all the steps in cLC-SRM (carried-assisted one-pot sample preparation coupled with liquid chromatography – selected reaction monitoring).

Nonhuman cell lysate digests (e.g., *S. oneidensis*) is used to pretreat PCR tubes for coating tube surface. Small numbers of human cells or single cells sorted by FACS are collected into pretreated PCR tubes. BSA protein (or nonhuman cell lysate proteome) carrier, heavy internal standards, and TFE (or DDM) are added into sample tubes sequentially for facilitating cell lysis and reducing surface adsorption losses. Conceptually the combined DDM and nonhuman cell lysate proteome carrier will work well for cLC-SRM. Cell lysis is conducted by sonication and protein denaturation is achieved by heating at high temperature. DDT and IAA reagents are used for reduction and alkylation, respectively (the step is optional). Trypsin is added for digestion with much higher ratios of trypsin: protein than that for standard trypsin digestion. The cap of sample tube is removed and then PCR tube is inserted into LC vial for direct LC-SRM analysis. Collected SRM data are analysis using publicly available Skyline software.




Figure 2: Sensitivity and accuracy of cLC-SRM for multiplexed quantification of EGFR/MAPK pathway proteins.

(A) Extracted ion chromatograms (XICs) of transitions monitored for ATADDELSFK derived from GRB2 at different numbers of MCF7 cell equivalents: 548.8/924.4 (*blue*), 548.8/853.4 (*purple*), 548.8/738.4 (*chestnut*). (B) Calibration curves for quantifying high abundance GRB2 and moderate-abundance PTPN11 with the use of the best responsive interference-free transitions, ATADDELSFK (548.8/924.4) for GRB2 and SNPGDFTLSVR (596.9/496.3) for PTPN11. Three and two SRM replicates were performed for 0–10 and 20–200 MCF7 cell equivalents, respectively. (C) Comparison of SRM signal between 10 and 100 MCF7 cells sorted by FACS. Each sample consists of two biological replicates with the addition of ~30 fmol of heavy peptide standards per replicate. XICs of transitions monitored for ATADDELSFK derived from GRB2: 548.8/924.4 (*blue*), 548.8/853.4 (*purple*), 548.8/738.4 (*chestnut*); XICs of transitions monitored for FYGAEIVSALDYLHSEK derived from AKT1: 648.0/897.9 (*blue*), 648.0/816.4 (*purple*), 648.0/283.1 (*chestnut*); XICs

of transitions monitored for LPSADVYR derived from SOS1: 460.7/807.4 (*blue*), 460.7/710.3 (*purple*), 460.7/404.2 (*chestnut*).

Author Manuscript

Author Manuscript

Author Manuscript

Author Manuscript

VLTPQTQVK (~744,000)




Figure 3: An example of targeted proteomics analysis of small numbers of MCF10A cells including single cells with cLC-SRM at short LC gradient.

Comparison of SRM signal among 1, 50 and 75 MCF10A cells sorted by FACS. ~30 fmol of internal standard was added to each sample. XICs of transitions monitored for VLTPQTQVK derived from PEBP1: 443.3/673.4 (*blue*), 443.3/572.3 (*purple*), 443.3/213.2 (*chestnut*). For single cells, the transition 443.3/213.2 was removed due to severe matrix interference.

Materials

Name	Company	Catalog Number	Comments
2 mL glass LC vial	Microsolv	9502S-WCV	Vessel to hold PCR tube for autosample injection
BSA	Sigma-Aldrich	P0834-10×1mL	Carrier protein for greatly reducing surface adsorption losses
DTT	Thermo Scientific	A39255	Reagent for reduction
Formic acid	Thermo Scientific	28905	Reagent for stopping enzyme reaction
IAA	Thermo Scientific	A39271	Reagent for alkylation
Peptide internal standards	New England peptide		Targeted quantification of EGFR/ MAPK pathway proteins
RT-PCR tube	GeneMate Bioexpress	T-3035-1	0.2 mL PCR tube for one-pot sample preparation
Skyline software	University of Washington		Publicly available for SRM data analysis
Sonicator	Hielscher Ultrasound Technology	UTR200	Sonication on ice for cell lysis
Speed Vac concentrator	Thermo Scientific		Reduction of the percentage of TFE for effective trypsin digestion
TFE	Sigma-Aldrich	18370-10×1mL	60% TFE for cell lysis
Thermocycler w/ heated lid	Peltier Thermal Cycler	PTC-200	Heating for protein denaturation
Trypsin Gold	Promega	V528A	Enzyme for protein digestion
Waters BEH C18 column	Waters		C18 column for peptide separation