Technical Assessment of the Affymetrix Yeast Expression GeneChip YE6100 Platform in a Heterologous Model of Genes That Confer Resistance to Antimalarial Drugs in Yeast

MARTIN E. NAU,¹ LYNDAL R. EMERSON,² RODGER K. MARTIN,³ DENNIS E. KYLE,³ DYANN F. WIRTH,² and MARYANNE VAHEY⁴*

Henry M. Jackson Foundation for the Advancement of Military Medicine, Rockville, Maryland¹; Department of Immunology and Infectious Diseases, Harvard School of Public Health, Boston, Massachusetts²; and Division of Experimental Therapeutics³ and Division of Retrovirology,⁴ Walter Reed Army Institute of Research, Washington, D.C.

Received 4 January 2000/Returned for modification 17 February 2000/Accepted 24 February 2000

The advent of high-density gene array technology has revolutionized approaches to drug design, development, and characterization. At the laboratory level, the efficient, consistent, and dependable exploitation of this complex technology requires the stringent standardization of protocols and data analysis platforms. The Affymetrix YE6100 expression GeneChip platform was evaluated for its performance in the analysis of both global (6,000 yeast genes) and targeted (three pleiotropic multidrug resistance genes of the ATP binding cassette transporter family) gene expression in a heterologous yeast model system in the presence and absence of the antimalarial drug chloroquine. Critical to the generation of consistent data from this platform are issues involving the preparation of the specimen, use of appropriate controls, accurate assessment of experiment variance, strict adherence to optimized enzymatic and hybridization protocols, and use of sophisticated bioinformatics tools for data analysis.

A universal challenge to drug therapy is the development of drug resistance. Efforts to understand the molecular mechanisms of the emergence of resistance to drugs span the fields of infectious disease, cancer, and toxicology. The eventuality of drug resistance necessitates the ongoing development of new drugs and interventions. A decade of research has identified a class of genes associated with multidrug resistance (8, 9).

The multidrug resistance genes (mdr genes) are part of the ATP binding cassette (ABC) transporter genes in mammalian cells (4, 7, 10). To facilitate the detection of drug resistance and to expedite the development of new drugs, several in vitro model systems have been developed that examine the activity of mdr and ABC transporters. One such system is the heterologous yeast model in which the genes PDR5, PDR10, and SNQ2, members of the pleiotropic drug resistance (pdr) family in yeast, have been associated with drug resistance (2, 9, 10, 15, 16, 17, 18). Observations that there may be 30 or more genes in yeast that are related by sequence homology to the ABC transporter gene family complicate the association of drug resistance with a particular gene (3). The Saccharomyces cerevisiae genome sequencing project revealed 31 ABC genes, which have been classified into six distinct subfamilies based on phylogenetic analysis (3, 7, 14, 19, 20). The pdr family is the largest of these subgroups, with 10 members. In total there are 12 ABC genes that have been associated with modulation of resistance to xenobiotics to date. The PDR5 gene has been linked to resistance to cycloheximide, mycotoxins, and cerulenin, and its product has been found to transport glucocorticoids (2, 3, 4, 10, 13). A second member of the *pdr* group, SNQ2, has been found to be linked to resistance to 4-nitrosoquinoline-N-oxide, methyl-nitro-nitrosoguanidine, and metal

ions such as Na⁺, Li⁺, and Mn⁺ (3, 16, 18). The $\Delta snq2 \Delta pdr5$ deletion strain exhibits a more pronounced sensitivity to metal ions and other drug substrates (3). *PDR10* is closely related to *PDR5* (65% sequence identity); however, the functional relatedness of these genes remains to be determined. Interestingly, *PDR10* has been found to localize to the cell surface like *PDR5* and *SNQ2* (3, 9).

With the introduction of the Affymetrix yeast expression GeneChip YE6100 platform (YE6100 platform), it has become feasible to plan experiments to simultaneously assess the changes in the expression patterns of not only the pleiotropic drug resistance gene family but also 6,000 yeast genes (5). Previously, Wodicka et al., at Affymetrix, characterized the basic performance characteristics of a prototype for the YE6100 platform to generate a global survey of 6,000 yeast genes (22). This platform was refined and exploited by Cho et al. to survey the complete yeast genome (6). Holstege et al., using an elegant battery of controls, exploited the commercially available YE6100 platform to assess the transcriptional control of yeast cell division (11). Winzeler et al. used a customized gene array platform for direct allelic scanning of the entire yeast genome (21).

To test the practical potential of the commercially available YE6100 platform to address drug resistance, a well-defined heterologous yeast model system was chosen. The expression profiles of two strains of *S. cerevisiae* were evaluated in the presence and absence of the antimalarial drug chloroquine. Strain YPH 499 (499) is wild type and refractory to the drug chloroquine. Strain YHW 1052 (1052) is a mutant with deletions in the *PDR5*, *PDR10*, and *SNQ2* genes and is thus more susceptible to chloroquine. The aim of this paper is to detail the technical aspects of the utilization of the YE6100 platform that are critical to the generation of consistent and reliable gene expression data in the study of drug resistance. The implementation of the methods and protocols presented in this paper will facilitate more intensive efforts to elucidate the

^{*} Corresponding author. Mailing address: Gene Array Laboratory, Walter Reed Army Institute of Research, 1600 East Gude Dr., Rockville, MD 20850. Phone: (301) 251-5058. Fax: (301) 762-7460. E-mail: mvahey@pasteur.hjf.org.

TABLE 1. Cell densities and mRNA yields

| Time point ^a | Strain | Cell de | mRNA (µg) | | | | |
|----------------------------|-------------|---|--|--|-------------|------------------|--|
| | | Introduction of drug | Har | vest | Cantual | Treated cultures | |
| | | | Control cultures | Treated cultures | cultures | | |
| Early | 499 1052 | $4.3 \times 10^{6} \\ 4.5 \times 10^{6}$ | 1.7×10^{7} 1.5×10^{7} | $\begin{array}{c} 1.1\times10^7\\ 8.2\times10^6\end{array}$ | 9.9 25.5 | 14.7 13.3 | |
| Middle | 499 1052 | $\begin{array}{c} 4.2\times10^6\\ 4.1\times10^6\end{array}$ | $\begin{array}{c} 3.2\times10^7\\ 4.2\times10^7\end{array}$ | $\begin{array}{c} 1.4\times10^7\\ 1.1\times10^7\end{array}$ | 6.6 8.0 | 14.4 15.4 | |
| Late | 499 1052 | $\begin{array}{c} 3.0\times10^6\\ 3.2\times10^6\end{array}$ | $\begin{array}{c} 3.5\times10^7\\ 3.2\times10^7 \end{array}$ | $\begin{array}{c} 2.6\times10^7\\ 1.8\times10^7 \end{array}$ | 13.1 9.4 | 6.4 7.5 | |

^{*a*} Early, 2 h; middle, 3 h; late, 4.5 h. Treated cultures received 1.5, 2.5, and 2.5 mg/ml for the early, middle and late time points, respectively.

details of the molecular interactions involved in the emergence of drug resistance. Two levels of data analysis, the global assessment of functional gene families and the targeted assessment of particular genes, will be addressed to demonstrate the type of information gleaned from each.

MATERIALS AND METHODS

Strains and media. The strains, 1052 and 499, used in this study were the kind gifts of Karl Kuchler of The University and Biocenter of Vienna, Vienna, Austria. The yeast strain 1052 (Δpdr5::TR1 Δsnq2::hisG Δpdr10::hisG) was utilized for this study along with its isogenic parental strain 499 (MATa ade2-101cc his $3\Delta 200 \ leu 2-\Delta 1 \ lys 2-801$ am trp $1-\Delta 1 \ ura 3-52$). Strain 1052 is deficient in three ABC transporters encoded in the pdr pathway (PDR5, PDR10, and SNQ2). In strain 1052, the deletion in PDR5 is from nucleotide (nt) +399 through nt +4456. The deletion in *PDR10* is from nt -90 through nt +4307. The deletion in *SNQ2* is from nt -6 through nt +3899. The 50% inhibitory concentrations of the drug chloroquine are 127 mg/ml for 499 and 50.00 mg/ml for 1052 as determined in nonaerated liquid medium and in solid medium culture. In liquid culture the 50% inhibitory concentrations of the drug chloroquine are 4.75 \pm 0.75 mg/ml for 499 and 1.38 \pm 0.13 mg/ml for 1052. Starter cultures were taken from colonies lifted from freshly streaked agar plates and grown overnight (to confluence at 2×10^8 cells/ml) at 30°C and 300 rpm in 5 to 10 ml of yeast-peptone-dextrose medium. The 5- to 10-ml starter cultures were diluted into 1,200 ml of prewarmed and aerated yeast-peptone-dextrose medium in a 4-liter flask to a density of 1.5×10^6 cells/ml. Cultures were grown at 30°C and 300 rpm for 2 h or until the cell density reached 3.0×10^6 cells/ml. At this juncture the culture was split into two 600-ml aliquots in two prewarmed 2-liter flasks. Chloroquine was added to the treatment flask to a concentration of 1.5 or 2.5 mg/ml from a 200-mg/ml concentrated stock of chloroquine diphosphate salt (Sigma, St. Louis, Mo.) dissolved in sterile double-distilled water. This solution had a pH of approximately 4.0. An exact volume of sterile double-distilled water, adjusted to the pH of the chloroquine solution, was added to the control flask. Table 1 shows the cell densities from critical points in the growth and treatment of the cultures used in the study. The assay points in the study are defined as early (2 h with or without 1.5 mg of drug per ml), middle (3 h, with or without 2.5 mg of drug per ml), and late (4.5 h, with or without 2.5 mg of drug per ml).

Cell harvesting and preparation of poly(A) RNA. Cultures were harvested identically at three time points: 2, 3, and 4.5 h. It is imperative that all cultures be treated exactly the same during the harvesting procedure. The overnight yeast culture was dispensed into 12 50-ml polypropylene conical tubes (Falcon/Becton Dickinson Labware, Franklin Lakes, N.J.) and centrifuged in a clinical centrifuge for 5 min at 4°C and at 2,000 × g. The pellet was resuspended in 5 ml of Tri-Reagent (Molecular Research Center, Woodlands, Tex.), and an equal volume of 400- μ m-diameter acid-washed glass beads was added. The mixture, and the manufacturer's instructions for the preparation of total RNA were followed. Poly(A) RNA (mRNA) was prepared from total RNA using the Oligotex (Qiagen, Valencia, Calif.) method according to the manufacturer's instructions.

cDNA synthesis. Double-stranded cDNA was synthesized in two steps using the Superscript Choice System (GibcoBRL, Rockville, Md.) and the reverse transcription primer T7-(dt)₂₄ [5'GGCCAGTGAATTGTAATACGACTCACT ATAGGGAGGCGG(T)₂₄ 3'] (GENSET Corp., LaJolla, Calif.). First-strand synthesis was carried out in a 20-µl reaction mixture. Approximately 3.0 µg of mRNA was annealed to 7 µg of T7-(dt)₂₄ primer at 70°C for 10 min. Reverse transcription was carried out at 37°C for 1 h in a mixture with final concentrations of 50 mM Tris-HCl (pH 8.3), 75 mM KCl, 3 mM MgCl₂, 10 mM dithiothreitol,

500 μ M each dATP, dCTP, dGTP, and dTTP, and 20,000 to 30,000 U of Superscript II reverse transcriptase per ml, and the reaction was terminated by placing the tube on ice. Second-strand synthesis was carried out in 150 μ l, incorporating the entire 20- μ l first-strand reaction mixture and a 130- μ l secondstrand reaction mixture for final concentrations of 25 mM Tris-HCl (pH 7.5), 100 mM KCl, 5 mM MgCl₂, 10 mM (NH₄)₂SO₄, 0.15 mM β -NAD⁺, 250 μ M each dATP, dCTP, dGTP, and dTTP, 1.2 mM dithiothreitol, 65 U of DNA ligase per ml, 250 U of DNA polymerase I per ml, and 13 U of RNase H per ml. The mixture was incubated at 16°C for 2 h, whereupon 2 μ l of T4 DNA polymerase at 5 U/ μ l was added and the incubation was continued at 16°C for 5 min. To terminate the reaction, 10 μ l of 0.5 M EDTA was added. The cDNA was purified using phenol-chloroform-isoamyl alcohol (24:23:1) saturated with 10 mM Tris-HCl (pH 8.0)–1 mM EDTA (AMBION, Inc., Austin, Tex.). The purified cDNA was precipitated with 5 M ammonium acetate and absolute ethanol at -20° C for 20 min. The pellet was resuspended in 7 to 9 μ l of RNase-free water to achieve a final concentration of between 0.25 and 0.65 μ g/ μ l.

In vitro transcription and fluorescent labeling. Synthesis of biotin-labeled cRNA was carried out by in vitro transcription using the MEGAscript T7 In Vitro Transcription Kit (AMBION, Inc.). According to the manufacturer's instructions, 0.4 to 1.0 µg of double-stranded cDNA was placed in a 20-µl reaction mix, at room temperature, containing Ambion $1 \times$ reaction buffer and enzyme mix (proprietary). The labeling mix consisted of 7.5 mM ATP, 7.5 mM GTP, 5.6 mM UTP, 1.9 mM biotinylated UTP (ENZO Diagnostics, Farmingdale, N.Y.), 5.6 mM CTP, and 1.9 mM biotinylated CTP (ENZO). The reaction mixture was incubated at 37°C for 5 h. The biotin-labeled cRNA was purified using RNeasy spin columns (Qiagen) according to the manufacturer's protocol. The biotinlabeled cRNA was fragmented in a 40-µl reaction mixture containing 40 mM Tris-acetate (pH 8.1), 100 mM potassium acetate, and 30 mM magnesium acetate, incubated at 94°C for 35 min, and then put on ice. One microliter of the intact biotin-labeled cRNA and 2 μl of the fragmented sample were run on a 1%agarose gel to evaluate both the yield and size distribution of the intact and fragmented products.

Hybridization, staining, and scanning of the GeneChip. The biotin-labeled and fragmented cRNA was hybridized to the YE6100 Yeast GeneChip array (Affymetrix, Santa Clara, Calif.) according to the manufacturer's instructions. Briefly, a 220-µl hybridization solution of 1 M NaCl, 10 mM Tris (pH 7.6), 0.005% Triton X-100, 50 pM control oligonucleotide B2 (5' bioGTCAAGATG CTACCGTTCAG 3') (Affymetrix), control cRNA (Bio B [150 pM], Bio C [500 pM], Bio D [2.5 nM], and Cre X [10 nM]) (American Type Tissue Collection, Manassas, Va., and Lofstrand Labs, Gaithersburg, Md.), 0.1 mg of herring sperm DNA per ml, and 0.05 µg of the fragmented labeled sample cRNA per µl was heated to 95°C, cooled to 40°C, and clarified by centrifugation before being applied to each of the four subarrays (A, B, C, and D) that comprise the YE6100 Yeast GeneChip platform. Hybridization was at 40°C in a rotisserie hybridization oven (model 320; Affymetrix) at 60 rpm for 16 h. Following hybridization, the GeneChip arrays were washed 10 times at 25°C with 6× SSPE-T buffer (1 M NaCl, 0.006 M EDTA, 0.06 M Na₃PO₄, 0.005% Triton X-100, pH 7.6) using the automated fluidics station protocol. GeneChip arrays were incubated at 50°C in $0.5 \times$ SSPE-T for 20 min at 60 rpm in the rotisserie oven and then stained for 15 min room temperature and 60 rpm with streptavidin phycoerythrin (Molecular Probes, Inc., Eugene, Oreg.) stain solution at a final concentration of 10 µg/ml in 6× SSPE-T buffer and 1.0 mg of acetylated bovine serum albumin (Sigma) per ml. The GeneChip arrays were washed twice at room temperature with $6 \times$ SSPE-T buffer and then scanned with a GeneArray Scanner (Hewlett-Packard, Santa Clara, Calif.), controlled by GeneChip 3.1 software (Affymetrix).

Assay monitoring and controls. The TEST 1 GeneChip (Affymetrix) was used according to the manufacturer's instructions to assess critical features of the mRNA preparations and the cDNA generated from the yeast strains and to evaluate the stringency of staining and hybridization. In addition, a battery of three types of GeneChip controls present on the TEST 1 GeneChip and on each of the four arrays in the YE6100 GeneChip set were employed according to the manufacturer's instructions. Details of the use and performance of these critical controls are given in Results. A method of mathematical scaling was employed by the GeneChip 3.1 software (Affymetrix) to normalize the fluorescence signal from each probe cell on each GeneChip and thus facilitate the reliable comparison of data from independent experiments.

Data analysis algorithm for the assessment of variance. The Affymetrix raw data set was scrutinized to eliminate any transcripts with fewer than 50% of probe cells contributing to the data. Subsequently, the first step in raw data mining for the assessment of variance captured all gene transcripts that were present on both GeneChips being compared (PP data set). The second step required that a decision be made to define what degree of change would be considered significant. We chose to approach this issue objectively, using a distribution analysis of the complete PP data set which defined a mean for the population of values and subsequently determined quartile percentages of 25, 50, and 75% above and below that mean. For the assessment of variance, outliers were defined as values exceeding the mean by 10-fold and were eliminated from the data set. When the PP data set was examined in this way, a value of 3.0-fold was applied to all subsequent analyses. Variances between GeneChips (intraexperimental variance) and between independent mRNA targets (interex-



FIG. 1. Performance of the battery of GeneChip controls with two independent preparations of templates from strains 1052 and 499. The ordinate shows the relative fluorescence values for each of the control markers listed on the abcissa. The linear regression r^2 value for the standard curve generated by the Bio B, Bio C, Bio D, and Cre X markers is 0.86.

perimental variance) were assessed by scoring the percentage of PP transcripts that exhibit no change relative to the total number of PP transcripts.

Data analysis algorithm for interrogation. For experiments in which differences in expression profiles between the drug-treated and untreated yeast strains were examined, the data analysis captured data from genes that were present in both cases (PP data set), as well as genes present in one case and absent in the other (PA or AP data set). All values above the 3.0-fold cutoff were included in the analysis of experimental expression profiles. The experimental design employed the analysis of data from the untreated control as a baseline for comparison to the treated strain in all cases. The cumulative fold change for the expression of all genes in a particular functional family was the sum of the levels of change of gene expression, using the values for the untreated strain as the control.

Bioinformatics analyses. GeneSpring version 2.1 (Silicon Genetics, San Carlos, Calif.) was used to derive global trends in the expression profiles and to specifically assess the expression patterns of the *pdr* gene targets. We used the temporal analysis of all of the raw data from the Affymetrix platform normalized to a single mean by the GeneSpring software.

RESULTS

Consistent cell harvests and mRNA yields. Table 1 shows the yields of cells and of mRNA across the three time points of the experiment and at the two concentrations of chloroquine used in the study. The amounts of cells harvested were comparable and equivalent at all time points.

Assessment of GeneChip performance. A battery of controls was used for all experiments. Three types of GeneChip controls are present on the TEST 1 GeneChip and on each of the four GeneChips in the YE6100 set. The first set of controls consists of four synthetically generated plasmid templates that are subjected to reverse transcription to incorporate fluorescent label according to the manufacturer's instructions (Affymetrix). These four cRNA templates, Bio B, Bio C, Bio D, and Cre X, are mixed in a cocktail to generate final concentrations of 150 pM, 500 pM, 2.5 nM, and 10 nM, respectively. These concentrations generate a standard curve and can thus be used to standardize interexperimental variation and efficiency of cDNA synthesis and labeling and to provide the dynamic range of the assay. Ultimately, the standard curve generated by these templates can be used to quantitate the level of RNA expression for a given gene on a per-cell basis. The second set of controls used on the GeneChip assesses the efficiency of cDNA synthesis by quantitating the amounts of 3' and 5' portions of target sequences generated during cDNA synthesis by assessing the expression of the yeast actin gene. Optimal synthesis reactions will generate equivalent amounts of signal in the 3' and 5' prime targets. The third set of GeneChip controls involves the evaluation of the integrity of the mRNA preparation used in the analysis and reports the GeneChip-based determination of equivalent amounts of mRNA used in the test. This is achieved by the assessment of the 18S rRNA gene expression profile, which is divided on the GeneChip into five sets of probe cells or segments (a through e).

The results of the analysis of TEST 1 GeneChip controls for two independent evaluations of strains 499 and 1052 are shown in Fig. 1. The ordinate indicates the relative fluorescence intensity reported by the GeneArray Scanner. The data from the 18S rRNA series show less than a twofold range in segment a and no significant difference in segment b, c, or d, except for the 1052 data point, which is less than onefold lower in segment c. This data set supports the hypothesis that equivalent amounts of mRNA were used in the cDNA reaction in preparation for GeneChip analysis. Also shown in Fig. 1 are the

| | a | Relative fluorescence intensity ^b with strain: | | | |
|---|----------|---|---|--|--|
| Control type | n" | 1052 | 499 | | |
| Bio B, Bio C, Bio D, and Cre X 18S ribosomal genes (a-e) | 20 40 | $4.13 \pm 0.53 \ (r^2 = 0.83) \\ 0.98 \pm 0.38$ | $3.27 \pm 0.56 \ (r^2 = 0.81)$ 1.04 + 0.43 | | |
| Yeast actin, 3' and 5' | 36 | 2.74 ± 0.55 | 3.04 ± 0.55 | | |

TABLE 2. Descriptive statistical analysis of the GeneChip control platform

^{*a*} Total number of data points in the analysis.

^b Data are expressed as the median ± the standard error. Regression values for the standard curve are for linear regression calculations.

results of the assessment of 3' and 5' segments of the actin gene expression. There is no significant difference between the fluorescence values for the 3' end of the yeast actin gene and for the 5' end of the yeast actin gene in this experiment. This data set indicates an optimal yield from the cDNA synthesis reaction. The manufacturer (Affymetrix) suggests that the yield of 3' product may vary by as much as fourfold. In our hands, optimization of the cDNA synthesis step routinely yielded less than a 0.5-fold difference between 3' and 5' segments.

The standard curve generated by the synthetic templates Bio B, Bio C, Bio D, and Cre X is shown in Fig. 1. The curve has an r^2 value of 0.86 and was remarkably consistent between strains, between GeneChips, and for two independent template preparations. Table 2 summarizes data on the performance of the battery of the three sets of controls that were generated by between 20 and 40 independent GeneChip assessments. A descriptive statistical analysis of the data set shows stringent inter- and intraexperimental consistency.

Assessment of assay variance. Table 3 presents data on the results of two independent expression profiles for each strain, 499 and 1052, in the absence of drug. These data were generated using one of the four GeneChips that comprise the complete YE6100 GeneChip platform (GeneChips A through D). In each case an independent growth and harvest of yeast cells followed by an independent preparation of GeneChip-ready template was carried out. Genes were scored as being present in both sets of data (PP), exhibiting no change in expression between the two sets of data, having increased or decreased, and, finally, having increased or decreased by threefold. For strain 1052, the total number of PP genes was 1,450, of which 32 increased by threefold, 116 decreased by threefold, and 1,302 (89%) remained unchanged, thus generating a variance between the two runs of 10.2%. For strain 499, the total number of PP genes was 1,439, of which 72 increased by threefold, 153 decreased by threefold, and 1,214 (84%) remained unchanged, thus generating a variance between the two runs of 15.6%. To further reduce these levels of interexperimental variance, the original culture was split into two cultures and reassessed for percentage of variance. As a result of splitting the original culture in this way, rather than growing two sideby-side cultures, the variance was reduced to zero for both strains, since there were no genes that changed greater than threefold between the two runs.

Global expression profiles of strains 1052 and 499 in the presence and absence of chloroquine. Shown in Fig. 2 and 3 are the results of a global survey of the 6,000 genes on the YE6100 GeneChip platform as assessed in strains 1052 and 499, respectively, in the presence and absence of the drug chloroquine and at each of the three time points and two drug concentrations used in the study. The control in each case was the value from the strain in the absence of the drug. Cumulative fold change values for the functional families are arrived at by simple summation of the levels of change from the control for each gene in a functional family.

As compared with the middle and late time points, the early time points for both 1052 and 499 exhibit a lower level of expression, with some increase in genes associated with membranes in strain 499. At the middle time point, however, both strains exhibit an increase in gene expression, with few genes showing a decrease. Genes associated with membranes, metabolism, and ribosomes showed the most increase in strain 1052 at the middle time point. The levels of the cumulative increase in expression were 2- to 10-fold higher in strain 499 at the middle time point. Increases in the expression of genes associated with membranes, metabolism, and ribosomes were similar in pattern but greater in magnitude to the changes at this time point in strain 1052. The most dramatic change occurred in strain 499 at the middle time point in the increase in expression of genes associated with synthetic pathways. In strain 1052, the late time point data set was dominated by a large decrease in the expression of genes associated with membranes. In contrast to the case for the two earlier time points, most expression levels were reduced in strain 1052 at the late time point. The expression of genes in strain 499 was also decreased at the late time point compared with the two earlier time points. The largest decline in expression was in the genes associated with translation and transcription.

Targeted expression profiles of the *pdr* genes *PDR5*, *PDR10*, and *SNQ2* in strains 1052 and 499 in the presence and absence of chloroquine. Figure 4 shows the expression profiles at three time points and in the absence or the presence of two different concentrations of the antimalarial drug chloroquine. The expression of the gene *PDR5* was decreased in the 1052 mutant strain in the presence and absence of the drug. In contrast, the expression of the gene *PDR10* was increased in strain 1052 in the presence and the absence of chloroquine. The expression

TABLE 3. Calculation of variance for two independently grown and tested samples of each strain

| Strain | No. of genes present in both profiles | No. (%) of genes: | | | | | Variance | | |
|-------------|--|--------------------------|--------------------------|--------------------------|------------------------------|--------------------------|--------------------------|------------|--------------|
| | | With no change | That increase | That decrease | With <3-fold or no change | That increase ≥3-fold | That decrease ≥3-fold | No. | % |
| 1052 499 | 1,450 1,439 | 786 (54.2) 690 (47.9) | 394 (27.2) 421 (29.3) | 270 (18.6) 328 (22.8) | 1,302 (89.8) 1,214 (84.4) | 32 (2.2) 72 (5.0) | 116 (8.0) 153 (10.6) | 148 225 | 10.2 15.6 |



FIG. 2. Cumulative change of gene expression levels in the mutant strain 1052 in the presence of chloroquine. The ordinate shows the cumulative fold changes for the expression levels of genes categorized by the functional family designation shown on the abcissa. The functional families are cell cycle and division proteins (CCD), drug resistance membrane proteins (DRM), kinases (KIN), membrane proteins (MEM), metabolic pathway proteins (MET), ribosomal proteins (RIBO), respiratory chain proteins (RSP), synthetic metabolic pathways (SMP), transcription and translation proteins (TRAN), pathology-related proteins (PATH), and stress-related proteins (SR). The expression level of genes in the untreated sample is used to determine the baseline for the degree of change of gene expression. The profiles for the early time point, the middle time point, and the late time point are shown.

of the gene *SNQ2* was moderate but level in strain 1052 in the presence of drug and moderate with a minor increase in slope in the absence of the drug. The wild-type strain 499 exhibited an increase in the expression levels of *PDR5* in the presence of drug but not in the absence of drug. In the absence of drug, the expression of the gene *PDR5* was moderate and level across all time points. The expression levels of *PDR10* and *SNQ2* in strain 499 remained low and level in both the presence and absence of the drug.

DISCUSSION

Template preparation. Several approaches to the extraction of total RNA and the subsequent preparation of mRNA are currently available. We found that the combination of two commercially available kits, the Tri-Reagent and Qiagen Oligotex kits, gave the most dependable results with yeast. The most critical aspects of the preparation of template for the Affymetrix GeneChip YE6100 platform are the quality of the mRNA and the degree to which it is representative of the biological nature of the sample. To ensure a representative sample, it is imperative to standardize the growth and handling of the yeast cultures. Holstege et al. first suggested that the attention to detail involved in the growth and harvest of yeast cultures for expression profiling was critical to the dependability of the data generated (11). We confirm and extend that observation by emphasizing the added importance of standardizing the treatment of these strains with the drug chloroquine

and minimizing experimental variance by splitting single cultures for drug treatment. It is imperative to ascertain the phenotypes of the wild-type and mutant strains in the presence of a drug prior to the characterization of the expression profiles generated as a result of treatment with that drug.

Quality control and assessment. The Affymetrix GeneChip YE6100 is exquisitely sensitive and necessitates the use of powerful controls to assure that all aspects of the procedure are consistent and reliable. Of the four types of controls available for expression profiling, using the Affymetrix GeneChip YE6100, we chose to apply three. The only control that we did not utilize involved the addition of synthetic total RNA template to the RNA samples extracted from the yeast strains. Instead, we chose to use data from the 3' and 5' ends of the yeast actin gene as a more accurate and less intrusive measure of the yield, quality, and representative nature of the mRNA. The data generated by these controls result directly from the sample tested and are not enhanced or quenched by the presence of artificial template.

We have determined that the battery of three controls that we routinely employ are essential to the interpretation, consistency, and reliability of expression profiling experiments. Perhaps the most powerful of the sets of controls is the standard curve generated by the synthetic templates Bio B, Bio C, Bio D, and Cre X. These data points offer the investigator the power to express GeneChip data on a semiquantitative level. The 18S ribosomal protein series and the yeast actin 3' and 5'



FIG. 3. Cumulative change of gene expression levels in the wild-type strain 499 in the presence of chloroquine. The ordinate shows the cumulative fold change for the expression levels of genes categorized by the functional family designation shown on the abcissa. The functional families are cell cycle and division proteins (CCD), drug resistance membrane proteins (DRM), kinases (KIN), membrane proteins (MEM), metabolic pathway proteins (MET), ribosomal proteins (RIBO), respiratory chain proteins (RSP), synthetic metabolic pathways (SMP), transcription and translation proteins (TRAN), pathology-related proteins (PATH), and stress-related proteins (SR). The expression level of genes in the untreated sample is used to determine the baseline for the degree of change of gene expression. The profiles for the early time point, the middle time point, and the late time point are shown.

end targets provide critical information on the preparation of the RNA and on the representative quality of the cDNA subsequently produced. The fact that all of these controls reside on each GeneChip further supports and ensures the generation of dependable data both within and between experiments. Most importantly, remarkably low levels of intraexperimental variance can be achieved, despite the enormous number of complex steps involved in generating an expression profile, by faithful attention to optimized laboratory protocols and by the vigilant use of the battery of GeneChip controls.

Interpretation of GeneChip expression profiles. We employed a well-characterized heterologous yeast model to assess the impact of the drug chloroquine on the yeast *pdr* genes *PDR5*, *PDR10*, and *SNQ2*. We assessed the expression profile data on two levels: (i) the global analysis of cumulative changes in expression of genes classified into broad functional families and (ii) the targeted expression analysis of the three *pdr* genes across the three time points and two drug concentrations used in the study. Jelinsky and colleagues used the global assessment of expression profiles to assess changes in gene expression in yeast in response to alkylating agents (12). Alon and colleagues employed targeted expression and cluster analysis to define expression patterns in colon tumors (1).

The assessment of the global alterations in expression profiles of broadly defined functional families in each of the strains in the presence of drug clearly identifies that in the mutant, the functional family most significantly affected by the drug is the membrane protein group. Strain 1052 exhibits a 250-fold reduction in the cumulative gene expression in the membrane protein group. The functional family of drug resistance-related membrane proteins is also reduced in cumulative gene expression by 75-fold. In contrast to these observations, the wild-type strain exhibits an increase in the expression of membraneassociated proteins and, most significantly, in proteins involved with the processes of transcription and translation. By the late time point, the wild-type strain exhibits a 100-fold decrease in the expression of proteins related to transcription and translation. Clearly these two strains respond with distinct strategies to the presence of drug. The assessment of the degree of cumulative change in the expression profiles of broadly defined functional families of genes can be readily made from the data reported by the Affymetrix GeneChip YE6100 platform. This information is most useful in suggesting the focus of further data mining to elucidate the specifics of a biological pathway affected by the drug.

The GeneSpring bioinformatics platform commercially available from Silicon Genetics interrogates the Affymetrix GeneChip YE6100 data in a significantly more powerful way. This tool allows for the identification of the patterns and magnitude of expression of any single gene assessed by the Affymetrix GeneChip YE6100 over the course of the study. The expression profile of individual targeted genes as well as the patterns or clusters of related genes can also be elucidated by the analysis. In the model system employed in this study, the promoter region of the *PDR10* gene was disrupted. An unchanged or reduced expression of this gene might be predicted as a result of this deletion. The expression profiles derived by GeneSpring analysis of the *PDR10* gene in the mutant strain



FIG. 4. Expression profiles of the yeast *pdr* genes *PRD5*, *PDR10*, and *SNQ2*. The ordinate shows the relative fluorescence intensity for (i) each of the study time points (early [E], middle [M], and late [L]), (ii) the two experimental treatments (drug treated [T] and untreated [U], and (iii) each of the two strains (1052 and 499), as shown on the abcissa.

1052 exhibit an unexpectedly high level of expression in both the presence and absence of chloroquine. Several explanations for this observation can be proposed.

The elevated levels of the mutant PDR10 gene expression may reflect the bias of the GeneChip to assess the 3' region of a gene. It is important to take into account that the Affymetrix GeneChip YE6100 platform interrogates 25-mer regions that cover the last 600 bp of the 3' end of the gene (5). This region is distal to the deletion made at the 5' promoter region of the PDR10 gene. Alternatively, there may be a difference in the efficiency of the promoter region, or in the stability or rate of turnover of the gene product, in the mutant as compared to that of the intact gene in the wild-type strain. In the wild-type strain, there is an increase in the production of PDR5 in response to drug treatment, while the PDR10 and SNQ2 expression levels remain moderate and unchanged, respectively. This pattern may reflect the specificity of the PDR5 response to the drug chloroquine in this strain (9). In contrast, expression levels of PDR5 and SNQ2 in the mutant strain show little or no response to the presence of the drug. Mechanistic explanations of the biological function of the gene products of PDR5, PDR10, and SNQ2 in the mutant and wild-type strains warrant further investigation. These observations show the complexity of the interpretation of expression profile data and underscore the necessity of ascertaining, by an independent assessment, information on the functional status of a gene target.

In summary, the utilization of optimized laboratory protocols, monitored by stringent controls, generates a powerful data set from the Affymetrix Expression GeneChip platform. The interpretation of the patterns and magnitudes of expression profiles represented in the data set requires the application of bioinformatics tools and a fundamental knowledge of the model being examined. The power of the method resides in the sensitivity, accuracy, and speed with which the expression of over 6,000 genes in response to experimental conditions can be simultaneously assessed. Confirmation of the trends observed in the data generated by expression profiling serves as a point of departure for further analysis of gene function and thus of the molecular mechanisms of drug action.

ACKNOWLEDGMENTS

We thank Karl Kuchler for the yeast strains YHW 1052 and YPH 499 used in this study, Anthony Lailin of Silicon Genetics and Brian Shimada and Mark Hurt of Affymetrix for their technical advice, and Deborah L. Birx and Nelson L. Michael for helpful discussions.

This work was supported in part by Cooperative Agreement no. DAMD17-93-V-3004 between the U.S. Army Medical Research and Materiel Command and the Henry M. Jackson Foundation for the Advancement of Military Medicine.

REFERENCES

- Alon, U., N. Barkai, D. A. Notterman, K. Gish, S. Ybarra, D. Mack, and S. J. Levine. 1999. Broad patterns of gene expression revealed by clustering analysis of tumor and normal colon tissues probed by oligonucleotide arrays. Proc. Natl. Acad. Sci. USA 96:6745–6750.
- Balzi, E., M. Wang, S. Leterme, L. Van Dyck, and A. Goffeau. 1994. *PDR5*, a novel yeast multidrug resistance conferring transporter controlled by the transcription regulator PDR1. J. Biol. Chem. 269:2206–2214.
- Bauer, B. E., H. Wolfger, and K. Kuchler. 1999. Inventory and function of yeast ABC proteins: sex, stress, pleiotropic drug and heavy metal resistance. Biochim. Biophys. Acta 1461:217–236.
- Bissinger, P. H., and K. Kuchler. 1994. Molecular cloning and expression of the Saccharomyces cerevisiae STS1 gene product. A yeast ABC transporter conferring mycotoxin resistance. J. Biol. Chem. 269:4180–4186.

- Chee, M., R. Yang, E. Hubbell, A. Berno, X. C. Huang, D. Stern, J. Winkler, D. J. Lockhart, M. S. Morris, and S. P. A. Fodor. 1996. Accessing genetic information with high density DNA arrays. Science 274:610–614.
- Cho, R. J., M. Fromont-Racine, L. Wodicka, B. Feierbach, T. Strearns, P. Legrain, D. J. Lockhart, and R. W. Davis. 1998. Parallel analysis of genetic selections using whole genome oligonucleotide arrays. Proc. Natl. Acad. Sci. USA 95:3752–3757.
- Decottignies, A., and A. Goffeau. 1997. Complete inventory of the yeast ABC proteins. Nat. Genet. 15:137–145.
- Delling, U., M. Raymond, and E. Schurr. 1998. Identification of Saccharomyces cerevisiae genes conferring resistance to quinoline ring containing antimalarial drugs. Antimicrob. Agents Chemother. 42:1034–1041.
- Egner, R., F. E. Rosenthal, A. Kralli, D. Sanglard, and K. Kuchler. 1998. Genetic separation of FK506 susceptibility and drug transport in the yeast *PDR5* ATP binding cassette multidrug resistance transporter. Mol. Biol. Cell 9:523–543.
- Hirata, D., K. Yano, K. Miyahara, and T. Miyakawa. 1994. Saccharomyces cerevisiae YDR1, which encodes a member of the ATP-binding cassette (ABC) superfamily, is required for multidrug resistance. Curr. Genet. 26: 285–294.
- Holstege, F. C. P., E. G. Jennings, J. J. Wyrick, T. I. Lee, C. J. Hengartner, M. R. Green, T. R. Golub, E. S. Lander, and R. S. Young. 1998. Dissecting the regulatory circuitry of a eukaryotic genome. Cell 95:717–728.
- Jelinsky, S. A., and L. Samson. 1999. Global response of Saccharomyces cerevisiae to an alkylating agent. Proc. Natl. Acad. Sci. USA 95:1486–1491.
- Kralli, A., S. P. Bohen, and K. R. Yamamoto. 1995. LEM1, an ATP-bindingcassette transporter, selectively modulates the biological potency of steroid hormones. Proc. Natl. Acad. Sci. USA 92:4701–4705.
- Kuchler, K., and R. Egner. 1997. Unusual protein secretion and translocation pathways in yeast: implication of ABC transporters, p. 49–85. In K.

Kuchler, A. Rubartelli, and B. Holland (ed.), Unusual secretory pathways: from bacteria to man. Landes Bioscience, Austin, Tex.

- Mahe, Y., Y. Lemoine, and K. Kuchler. 1996. The ATP binding cassette transporters *Pdr5* and *Snq2* of *Saccharomyces cerevisiae* can mediate transport of steroids in vivo. J. Biol. Chem. 271:25167–25172.
- Miyahara, K., M. Mizunuma, D. Hirata, E. Tsuchiya, and T. Miyakawa. 1996. The involvement of the *Saccharomyces cerevisiae* multidrug resistance transporters Pdr5p and Snq2p in cation resistance. FEBS Lett. 399:317–320.
- Ruetz, S., U. Delling, M. Brault, E. Schurr, and P. Gros. 1996. The pmdr1 gene of *Plasmodium falciparum* confers cellular resistance to antimalarial drugs in yeast cells. Proc. Natl. Acad. Sci. USA 93:9942–9947.
- Servos, J., E. Haase, and M. Brendel. 1993. Gene SNQ2 of Saccharomyces cerevisiae, which confers resistance to 4-nitroquinoline-N-oxide and other chemicals, encodes a 169 kDa protein homologous to ATP-dependent permeases. Mol. Gen. Genet. 236:214–218.
- Taglicht, D., and S. Michaelis. 1998. Saccharomyces cerevisiae ABC proteins and their relevance to human health and disease. Methods Enzymol. 292: 130–162.
- Wendler, F., H. Bergler, K. Prutej, H. Jungwith, G. Zisser, K. Kuchler, and G. Hogenauer. 1997. Diazaborine resistance in the yeast *Saccharomyces cerevisiae* reveals a link between YAP1 and the pleiotropic drug resistance genes *PDR1* and *PDR3*. J. Biol. Chem. 272:27091–27098.
- Winzeler, E. A., D. R. Richards, A. R. Conway, A. L. Goldstein, S. Kalman, M. J. McCullough, J. H. McCusker, D. A. Stevens, L. Wodicka, D. J. Lockhart, and R. W. Davis. 1998. Direct allelic variation scanning of the yeast genome. Science 281:1194–1197.
- Wodicka, L. H. Dong, M. Mittman, M. H. Ho, and D. J. Lockhart. 1997. Genome-wide expression monitoring in *Saccharomyces cerevisiae*. Nat. Biotechnol. 15:1359–1367.