## VIDAS<sup>®</sup> TB-IGRA reagents induce a CD4+ and CD8+ T-cell IFN- $\gamma$ response for both TB infection and active TB

Dear Editor,

TB is one of the leading causes of death and the WHO has estimated an additional 9.9 million new TB cases in 2020.1 Global control and elimination of TB can be improved by identifying latent TB infection (LTBI) and treating those who are at highest risk of developing active TB.2-5 Many diagnostic assays are available for LTBI, including the tuberculin skin test (TST) and interferon-gamma (IFN-γ) release assays (IGRAs). QuantiFERON-TB Gold Plus (QTF-Plus; Qiagen, Hilden, Germany) and T-SPOT.TB (Oxford Immunotec, Abingdon, UK) are the most routinely used IGRAs.<sup>2-4</sup> The design and development of IGRAs are based on the knowledge that both CD4+ and CD8+ T-cells are important in TB immune pathogenesis.<sup>6,7</sup> CD4+ Tcell response plays an important role, as shown by the high risk of progression to TB disease in subjects with CD4+ T-cell impairment, such as HIV-infected individuals.1 CD8-specific T-cells are involved in the control of Mycobacterium tuberculosis (MTB) bacterial load, and are increased at TB disease diagnosis and reduced at the end of therapy.<sup>8,9</sup> A new assay has recently been developed (VIDAS® TB-IGRA, bioMérieux, Marcy-l'Étoile, France) based on early secretory antigenic target-6 and culture filtrate protein-10 peptides<sup>10</sup> for the diagnosis of individuals infected with MTB, both disease and infection. In a recent multicentre study, this test showed a higher accuracy for active TB than QFT-Plus, used as surrogate for TB infection, with 97.5% sensitivity.<sup>11,12</sup> The assay is completely automated; the collected whole blood is dispensed from a single tube into three different wells: NIL (negative control), AG (MTB antigen) and Mitogen (positive control). After on-board incubation, IFN- $\gamma$  is measured using VIDAS 3 instrument.<sup>12</sup> Here, we aimed to characterise the T-cell-specific response induced by VIDAS TB-IGRA stimulation reagents using flow cytometry to analyse the CD4+ and CD8+ T-cell response in relation to a subject's TB status.

Between November 2020 and January 2021, peripheral blood mononuclear cells (PBMCs) were taken from patients with active TB and LTBI. Active TB patients were enrolled within 7 days from starting treatment, and those with LTBI before starting

EP and CF contributed equally to this study.

preventive treatment. Active TB was microbiologically confirmed using culture testing and/or nucleic acid tests. All enrolled individuals were tested using QTF-Plus assay. PBMCs were isolated, frozen and stored in liquid nitrogen until use. For flow cytometry assay, 1x10<sup>6</sup> PBMCs were stimulated with VIDAS TB-IGRA reagents (bioMérieux): NIL (the inert compound contained in the buffer employed to dilute the antigen peptides that was used as negative control, 10 µg/ml), AG (MTB antigen, 100 µg/ml), Mitogen (positive control, 20 µg/ml). Anti-CD28 and anti-CD49d (1 µg/ml each) monoclonal antibodies and Golgi  $Plug^{{}^{\mathrm{\scriptscriptstyle T\!M}}}$  (1  $\mu g/ml)$  were added to PBMCs, as previously described.<sup>7</sup> After 16–24 h, PBMCs were stained using CD3 PE-Cy7, CD4 BV711 and CD8 APC-H7. Cytofix/Cytoperm<sup>™</sup> was used for IFN-γ-APC intracellular staining according to the manufacturer's recommendations (Figure 1) (all antibodies and reagents from BD Biosciences, San Jose, CA, USA). Samples were acquired using a BD Lyric cytometer (BD Biosciences). Dead cells were excluded by both side/forward scatter using BD Horizon<sup>™</sup> Fixable Viability Stain 700 (BD Bioscience). Data were analysed using FlowJo software v10 (Tree Star, Ashland, OR, USA). Cytokine background was subtracted from the stimulated tubes. Detection of at least 10 events in the CD4+ IFN- $\gamma^+$  or CD8+ IFN- $\gamma^+$ gate was considered a positive response.<sup>14</sup> Statistical analysis was performed using Prism software (Graph-Pad Prism 8, La Jolla, CA, USA). Median and interquartile ranges (IQRs) were calculated for continuous measures,  $\chi^2$  test for proportions, and Mann-Whitney U-test for pairwise comparison. Data were considered statistically significant when P <0.05. All participants provided written informed consent; the study was approved by the Ethical Committee of Istituto Nazionale Malattie Infettive Lazzaro Spallanzani-IRCCS, Rome, Italy (approval 126/2020 and 58/2018) and complied with the principles of the Declaration of Helsinki.<sup>13</sup>

We enrolled 10 LTBI and 10 active TB individuals. Among the LTBI subjects, 9 were recent contacts of pulmonary TB patients and 1 was a remote infection.<sup>2</sup> We did not find significant differences for age (active TB: median age 33 years, IQR 24–49; LTBI: median age 44 years, IQR 28–52), sex (active TB: female 20%, LTBI: female 30%), bacilli Calmette-Guérin vaccination (active TB: 60%; LTBI: 60%) and

Article submitted 8 August 2021. Final version accepted 4 October 2021.



**Figure 1** Gating strategy for IFN- $\gamma$  analysis in T-cell subsets and IFN- $\gamma$  production by CD4<sup>+</sup> and CD8<sup>+</sup> T-cells in patients with active TB and LTBI. Plots show the distribution and frequency of IFN- $\gamma$  in PBMCs of a representative active TB and LTBI subject. Lymphocytes were gated based on their SSC-A and FSC-A parameters, following which only live cells were considered according to the negativity of the Viability Stain antibody and single cells were gated as described. CD3<sup>+</sup> T-cells were gated according to the expression of their surface marker. SSC-A = side scatter; FSC-A = forward scatter; IFN- $\gamma$  = interferon-gamma; LTBI = latent TB infection; PBMC = peripheral blood mononuclear cell.

origin (West Europe active TB: 40%, LTBI: 40%; East Europe active TB: 60%; LTBI: 50%) between the two groups. No significant differences in QTF-Plus scores between active TB and LTBI were reported (active TB: 70%, LTBI: 100%). We performed an intracellular staining to evaluate the frequency of CD4<sup>+</sup> IFN- $\gamma^+$  and CD8<sup>+</sup> IFN- $\gamma^+$  T-cells (Figure 2). Regardless of the CD4<sup>+</sup> or CD8<sup>+</sup> T-cell subset involved, we scored 8/10 active TB responders (80%), and 4/10 LTBI responders (40%) (Figure 3). In particular, 87.5% of the active TB patients and 50% of the LTBI subjects had a CD4+ T-cell-positive response to the AG (Figure 3). Similarly, 62.5% of active TB and 50% of LTBI samples showed a CD8+ T-cell response (Figure 3). Although we did not find any significant difference, the frequency of CD4+ IFN- $\gamma^+$  and CD8+ IFN- $\gamma^+$  T-cells was higher in LTBI than in active TB. All samples scored positive for the mitogen-control, indicating that participants were immunocompetent (data not shown), and had a cell viability between 89.4% and 99.7% (data not shown).

Here, for the first time, using a flow cytometry approach we have demonstrated that the AG reagent of VIDAS TB-IGRA induces both CD4+ and CD8+ T-cell responses in a cohort of subjects with active TB and LTBI enrolled in a low TB endemic country such as Italy. Previous studies have reported the ability of other IGRAs, such as the QFT-Plus, to induce both a CD4+ and CD8+ T-cell response in active TB and LTBI subjects with or without HIVinfection.<sup>2,3</sup> In the present study, we confirmed this by using a different test. We also observed a higher number of CD8+ T-cell responses to VIDAS TB-IGRA AG in active TB than in LTBI; this has been previously shown using the QFT-Plus assay in both HIV-uninfected and -infected subjects.<sup>3,8</sup> Interestingly, we also found CD8+T-cell responders in LTBI. This is due to the fact that TB comprises a spectrum of different stages characterised by host-pathogen



**Figure 2** Percentage of IFN- $\gamma^+$  CD4<sup>+</sup> T-cells and IFN- $\gamma^+$  CD8<sup>+</sup> T-cells in untreated cells and in response to AG stimulation in LTBI subjects and active TB patients. LTBI = latent TB infection; AG = MTB antigen; IFN- $\gamma$  = interferon-gamma; FSC-A = forward scatter.



Figure 3 Left: Box showing both numbers and percentage of participants who showed a CD4+ or/and CD8+ T-cell response. Middle and right graphs: percentage of CD4<sup>+</sup>, CD8<sup>+</sup>/IFN-γ<sup>+</sup> T-cells in the responders to VIDAS TB-IGRA AG. Each dot/square represents a single active TB or LTBI patient, black lines represent median values. Individuals who showed both CD4+ and CD8+ responses are represented with the same colour in all graphs. Individuals who showed a single CD4<sup>+</sup> or CD8<sup>+</sup> T-cell response are indicated in white. No differences between groups were observed using the Mann–Whitney test. AG = (VIDAS TB-IGRA) antigen tube reagent; LTBI = latent TB infection; IFN- $\gamma$  = interferon-gamma.

equilibrium;<sup>15</sup> it is therefore not uncommon to detect a CD8<sup>+</sup> T-cell response during LTBI as well. In this study, only 4/10 of the enrolled LTBI subjects scored positive on VIDAS TB-IGRA AG using flow cytometry, although all LTBI subjects enrolled in our study had a positive response to the QFT-Plus test. This is likely due to the use of frozen PBMCs that may limit the detection of specific responses compared to fresh samples.<sup>16</sup> In conclusion, we have shown that VIDAS TB-IGRA reagents induce both CD4<sup>+</sup> and CD8<sup>+</sup> T-cell IFN-γ response in both LTBI and active TB patients. Our findings provide further details on the immunological features of the responses associated with the test and on the contribution of T-cell subpopulations. These results are important in shedding light on VIDAS TB-IGRA characteristics, a new generation of completely automated IGRA.12

E. PETRUCCIOLI,<sup>1</sup> C. FARRONI,<sup>1</sup> G. CUZZI,<sup>1</sup> V. VANINI,<sup>1,2</sup> F. PALMIERI,<sup>3</sup> P. VITTOZZI,<sup>3</sup> D. GOLETTI<sup>1</sup> <sup>1</sup>Translational Research Unit, and <sup>2</sup>UOS Professioni Sanitarie Tecniche, and <sup>3</sup>Respiratory Infectious Diseases Unit, National Institute for Infectious Diseases Lazzaro Spallanzani, Istituto di Ricovero e Cura a Carattere Scientifico, Rome, Italy Correspondence to: D Goletti, Translational Research Unit of the Research Department, National Institute for Infectious Diseases, Padiglione del Vecchio, Room 13, Via Portuense 292, Rome 00149, Italy. e-mail: delia.goletti@inmi.it; delia.goletti@gmail.com

## Acknowledgements

LTB

The authors would like to thank the patients for their generosity and nurses for their contribution regarding sample collection from patients: L Dell'Arciprete, F Parisotto, L Tinari, C Fassina and E Ercoli; and A Aiello and S N Fard for their help in revising the manuscript.

Conflicts of interest: DG is a consultant of bioMérieux (Marcy l'Etoile, France), Qiagen (Hilden, Germany), and Diasorin and received a grant from Biomerieux. This study was supported by bioMérieux, and the Italian Ministry of Health (Ricerca Corrente-Linea 4).

## References

- 1 World Health Organization. Global Tuberculosis Report 2021: WHO, 2021 https://www.who.int/publications/i/item/ 9789240037021.
- 2 Petruccioli E, et al. Effect of HIV-infection on QuantiFERONplus accuracy in patients with active tuberculosis and latent infection. J Infect 2020; 80: 536-546.
- 3 Petruccioli E, et al. First characterization of the CD4 and CD8 Tcell responses to QuantiFERON-TB plus. J Infect 2016; 73: 588-597.
- 4 Chiacchio T, et al. Characterization of QuantiFERON-TB-plus results in latent tuberculosis infected patients with or without immune-mediated inflammatory diseases. J Infect 2019; 79: 15-23.
- 5 Behr MA, et al. Latent Tuberculosis: Two Centuries of Confusion. Am J Respir Crit Care Med. 2021; 204: 142-148.
- 6 Chiacchio T, et al. Impact of antiretroviral and tuberculosis therapies on CD4(+) and CD8(+) HIV/M. tuberculosis-specific T-cell in co-infected subjects. Immunol Lett 2018; 198: 33-43.
- 7 Petruccioli E, et al. IFN-y/TNF-a specific-cells and effector memory phenotype associate with active tuberculosis. J Infect 2013; 66: 475-486.
- 8 Rozot V, et al. Combined use of Mycobacterium tuberculosisspecific CD4 and CD8 T-cell responses is a powerful diagnostic tool of active tuberculosis. Clin Infect Dis 2015; 60: 432-437.
- 9 Cho S, et al. Antimicrobial activity of MHC class I-restricted CD8+ T cells in human tuberculosis. Proc Natl Acad Sci USA 2000; 97: 12210-12215.
- 10 Hamada Y, et al. Tests for tuberculosis infection: landscape analysis. Eur Respir J 2021; https://doi.org/10.1183/13993003. 00167-2021. [In press]
- 11 Dheda K, et al. Evaluation of Vidas® TB-IGRA for tuberculosis infection detection, bioMérieux white paper. Marcy l'Etoile, France: bioMérieux, 2021 https://go.biomerieux.com/tb-igrawhite-paper.
- 12 VIDAS® TB-IGRA: Reliable and fully automated solution for Tuberculosis (TB) infection detection. Marcy l'Etoile, France: bioMérieux, 2020. https://www.biomerieux-diagnostics.com/ vidasr-tb-igra.

- 13 World Medical Association. Declaration of Helsinki ethical principles for medical research involving human subjects. Helsinki, Finland: WMA, 2018. https://www.wma.net/ policies-post/wma-declaration-of-helsinki-ethical-principlesfor-medical-research-involving-human-subjects/2018.
- 14 Roederer M. How many events is enough? Are you positive? Cytometry A 2008; 73: 384–385.
- 15 Barry CE, 3rd, et al. The spectrum of latent tuberculosis: Rethinking the biology and intervention strategies. Nat Rev Microbiol 2009; 7: 845–855.
- 16 Ford T, et al. Cryopreservation-related loss of antigen-specific IFN-γ producing CD4(+) T-cells can skew immunogenicity data in vaccine trials: lessons from a malaria vaccine trial substudy. Vaccine 2017; 35: 1898–1906.