

## Long Term Accuracy of SARS-CoV-2 Interferon- $\gamma$ Release Assay and its Application in Household Investigation

Kanagavel Murugesan<sup>1</sup>, Prasanna Jagannathan<sup>2</sup>, Jonathan Altamirano<sup>3</sup>, Yvonne A Maldonado<sup>3</sup>, Hector F. Bonilla<sup>2</sup>, Karen B. Jacobson<sup>2</sup>, Julie Parsonnet<sup>2</sup>, Jason R. Andrews<sup>2</sup>, Run-Zhang Shi<sup>1</sup>, Scott Boyd<sup>1</sup>, Benjamin A. Pinsky<sup>1,2,4</sup>, Upinder Singh<sup>2,5</sup>, Niaz Banaei<sup>1,2,6</sup>

<sup>1</sup> Department of Pathology, Stanford University School of Medicine, Stanford, CA, USA

<sup>2</sup> Division of Infectious Diseases and Geographic Medicine, Department of Medicine, Stanford University School of Medicine, Stanford, CA, USA

<sup>3</sup> Division of Infectious Diseases, Department of Pediatrics, Stanford University School of Medicine, Stanford, CA, USA

<sup>4</sup> Clinical Virology Laboratory, Stanford Health Care, Stanford, CA, USA

<sup>5</sup> Department of Microbiology and Immunology

<sup>6</sup> Clinical Microbiology Laboratory, Stanford Health Care, Stanford, CA, USA

### **\*Corresponding Author:**

Niaz Banaei MD

3375 Hillview Ave, Rm. 1602,

Palo Alto, Ca 94304 USA

Phone 650-736-8052

[nbanaei@stanford.edu](mailto:nbanaei@stanford.edu)

**Summary:**

SARS-CoV-2 immunodiagnostics are needed to identify infected individuals in order to understand the transmission dynamics of emerging variants and to assess vaccine response. Interferon-gamma release assay maintains sensitivity 10 months post-infection in convalescents and detects more household contacts than IgG.

Accepted Manuscript

## Abstract

**Background.** An immunodiagnostic assay that sensitively detects a cell-mediated immune response to SARS-CoV-2 is needed for epidemiological investigation and for clinical assessment of T cell-mediated immune response to vaccines, particularly in the context of emerging variants that might escape antibody responses.

**Methods.** The performance of a whole blood interferon-gamma (IFN- $\gamma$ ) release assay (IGRA) for the detection of SARS-CoV-2 antigen-specific T cells was evaluated in COVID-19 convalescents tested serially up to 10 months post-infection and in healthy blood donors. SARS-CoV-2 IGRA was applied in contacts of households with index cases. Freshly collected blood in the lithium heparin tube was left unstimulated, stimulated with a SARS-CoV-2 peptide pool, and stimulated with mitogen.

**Results.** The overall sensitivity and specificity of IGRA were 84.5% (153/181; 95% confidence interval [CI] 79.0-89.0) and 86.6% (123/142; 95% CI;80.0-91.2), respectively. The sensitivity declined from 100% (16/16; 95% CI 80.6-100) at 0.5-month post-infection to 79.5% (31/39; 95% CI 64.4-89.2) at 10 months post-infection ( $P < 0.01$ ). The IFN- $\gamma$  response remained relatively robust at 10 months post-infection (3.8 vs. 1.3 IU/mL, respectively). In 14 households, IGRA showed a positivity rate of 100% (12/12) and 65.2% (15/23), and IgG of 50.0% (6/12) and 43.5% (10/23) in index cases and contacts, respectively, exhibiting a difference of +50% (95% CI +25.4-+74.6) and +21.7% (95% CI, +9.23-+42.3), respectively. Either IGRA or IgG was positive in 100% (12/12) of index cases and 73.9% (17/23) of contacts.

**Conclusions.** The SARS-CoV-2 IGRA is a useful clinical diagnostic tool for assessing cell-mediated immune response to SARS-CoV-2.

**Keywords:** IGRA, COVID-19, T-cell response, Immunocompromised patients, whole blood assay

Accepted Manuscript

## Introduction

The emergence of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) has demanded the development of sensitive laboratory diagnostics to detect active and remote infections to control the pandemic [1]. In the arena of immunodiagnosics, serologic assays that probe antibody responses to SARS-CoV-2 are far more available than assays that measure T cell responses [1]. However, serological tests for anti-SARS-CoV-2 antibodies may not accurately predict the magnitude, functionality, and durability of T cell-mediated immune response to SARS-CoV-2 [2, 3], particularly in immunocompromised patients with impaired B cell function [3-6]. Recent studies have shown that T cell responses are more sensitive markers of past SARS-CoV-2 infection compared with antibody responses and postulated to represent a correlate of protective immunity [5, 7-11]. Furthermore, T cell responses are more robust than antibody responses in convalescents with mild or asymptomatic COVID-19 infection [10-13]. Thus, clinical T cell assays for SARS-CoV-2 are needed to evaluate individuals for the cell-mediated immune response against SARS-CoV-2 as evidence for past infection and immune response to vaccination [14].

Interferon gamma release assay (IGRA) is an in vitro blood diagnostic used clinically to measure IFN- $\gamma$  released by antigen-specific T cells after stimulation with pathogen-specific peptides. IGRA is best known for its role in diagnosing latent *Mycobacterium tuberculosis* infection [15, 16] in which either purified or whole blood mononuclear cells are stimulated overnight followed by IFN- $\gamma$  enzyme-linked immunosorbent spot or ELISA, respectively, to measure IFN- $\gamma$  response from sensitized T cells [9, 17]. Recently, we and another group reported on the accuracy of a laboratory-developed whole blood IGRA for the detection of the SARS-CoV-2 specific T cell responses in convalescents weeks after infection [18, 19]. Since then, our assay has been further optimized with commercial peptide pools and

simplified with in-tube stimulation, and offered clinically in a CLIA laboratory at Stanford Health Care. In this study, we report the longitudinal accuracy of SARS-CoV-2 IGRA up to 10 months after infection and evaluated its utility in household contacts exposed to an index case. We show that IGRA response in convalescents is durable and it identifies more infected household contacts compared with antibody testing.

## Methods

**Ethics.** This study was conducted under two Institutional Review Board-approved protocols (IRB-55619 and IRB-55479) named “A Phase 2 Randomized, Single-Blind Study of a Single Dose of Peginterferon Lambda-1a Compared with Placebo in Outpatients with Mild COVID-19” and “Natural History of Shedding and Household Transmission of COVID-19: Constructing Patterns of Viral Spread and Evolution” which collected blood from infected and exposed individuals, respectively. Informed consent was obtained from participants before blood collection.

**Study design.** This was a case-control study to evaluate the longitudinal sensitivity and specificity of SARS-CoV-2 IGRA (Figure 1). Additionally, the utility of SARS-CoV-2 IGRA was assessed in household contacts of COVID-19 index cases. Blood samples from outpatient COVID-19 patients with positive SARS-CoV-2 reverse transcriptase (RT)-PCR enrolled in an IFN- $\lambda$  vs. placebo therapeutic clinical trial (referred to as cases) and healthy blood donors with no COVID-19 symptoms (referred to as controls) were tested with SARS-CoV-2 IGRA and IgG ELISA. Cases were tested at 0.5, 1, 4, 7 and 10 months post-infection between April 6, 2020 and April 30, 2021. Controls were tested once between May 11, 2020 and November 11, 2021. None of the samples were tested in a previous study [18]. Based on

prior studies, IFN- $\lambda$  has been shown to not negatively impact adaptive response to SARS-CoV-2 [20]. Household contacts of RT-PCR-positive COVID-19 index cases with mild COVID-19 from 14 households were tested with SARS-CoV-2 IGRA and IgG ELISA days after the index case was diagnosed. Members of each household were tested on the same day.

**Specimens.** Whole blood was collected in lithium heparin blood collection tube and transported to the Stanford Health Care Clinical Microbiology Laboratory at room temperature for IGRA testing. Blood from healthy blood donors was purchased from Stanford Blood Center.

**Antigen and Mitogen.** The SARS-CoV2 peptide pools consisting of PepTivator SARS-CoV-2 Prot S, S1, N and M were purchased from Miltenyi Biotec (Bergisch Gladbach, Germany). Per the manufacturer, S1 pool consists of 15-mer peptides overlapping by 11-residues covering the S1, while the remaining pools contain 15-mer peptides corresponding to immunodominant epitopes. Peptides were developed for stimulation of both CD4 and CD8 T cells. Each 60nmol vial, which contains 60 peptides, was reconstituted aseptically in 2mL of sterile de-ionized water and all 4 vials were combined into a single mega pool. For in-tube stimulation 80 $\mu$ L of the mega pool was dispensed into a BD vacutainer no additive tube (Beckton Dickenson; city, state). Optimization studies were performed with an individual S1 peptide pool preparation and in combination with S, and S, M, N (Supplementary Figure 1A). Phytohemagglutinin PHA-P Mitogen (Sigma) was dissolved in sterile de-ionized water at 1mg/mL. For in-tube stimulation, 50 $\mu$ L was dispensed into a BD Vacutainer No Additive Tube. Prepared tubes and peptide reagents were stored at -80 °C.

**Interferon Gamma Release Assay (IGRA).** SARS-CoV-2 IGRA was performed as described previously [18]. One mL of freshly collected blood was transferred to a 24 well tissue culture plate or BD vacutainer with no additive tubes at 1mL per well/tube. One well or tube was left unstimulated (nil), one well or tube was stimulated with SARS-CoV-2 antigen mega pool at 7.5 nmol/mL (approximately 1 µg/mL), and one well or tube was stimulated with mitogen at 50 µL/mL. Antigen concentration was based on optimization studies stimulating with 3.75 and 7.5 nmol/mL (Supplementary Figure 1B). In-tube stimulation was validated using in-plate stimulation as the reference method (Supplementary Figure 2) and used only for clinical testing after implementation in the clinical laboratory [21]. The blood samples were mixed gently and incubated at 37 °C with a relative humidity of 95% for 20 to 24 h. Plasma was separated and stored at 4 °C. The IFN-γ concentration was measured with an automated enzyme-linked immunosorbent assay (ELISA) instrument (DSX; Dynex Technologies, Chantilly, VA) using the QuantiFERON-TB ELISA kit (Qiagen, Germantown, MD)). A four-point standard curve was used to calculate IFN-γ concentration in international units (IU)/mL. Linear range was 0 to 10 IU/mL. IFN-γ response was defined as positive if antigen-nil  $\geq 0.35$  IU/mL; negative if antigen-nil  $< 0.35$  and mitogen-nil  $\geq 0.5$  IU/mL, and indeterminate if nil  $> 8$  IU/mL or antigen-nil  $< 0.35$  and mitogen-nil  $< 0.5$  IU/mL.

**Enzyme-linked immunosorbent assay (ELISA).** Anti-SARS-CoV-2 Spike S1 domain IgG ELISA was performed on lithium heparin plasma using the EUROIMMUN instrument and reagents (Lübeck, Germany) per the manufacturer's instructions. The anti-SARS-CoV-2 IgG antibody level was reported as the ratio of optical density (OD) of the sample over the OD of the calibrator. Linear range was 0 to 12 OD ratio. The ratio was defined as follows:  $< 0.8$  negative;  $\geq 0.8$  to  $< 1.1$  borderline;  $\geq 1.1$  positive.



**Statistical Analysis.** The Mann-Whitney U test was used to compare median IFN- $\gamma$  response between groups. Adjusted P values were determined as per the Bonferroni's multiple comparison test. Fisher's exact test was used to analyze differences between proportions. Statistical analysis was done with GraphPad Prism 8.0.1 software (San Diego, CA, USA).

## Results

### Patient Cohorts

To assess the accuracy of SARS-CoV-2 IGRA longitudinally after natural infection, a total of 94 unique adult COVID-19 convalescents were tested at 0.5 (n=16), 1 (n=8), 4 (n=48), 7 (n=70) and 10 (n=39) months post infection for total of 181 tests. Eighteen convalescents participated in one timepoint, 67 participated in 2 timepoints, 7 participated in 3 timepoints, and 2 participated in 4 timepoints. A total of 142 uninfected healthy blood donors were tested one time. COVID-19 convalescents had either mild (91%) or asymptomatic (9%) COVID-19. The median age of COVID-19 convalescents was 38 years (interquartile range [IQR], 29-55 years) and 45% were female. To evaluate the utility of SARS-CoV-2 IGRA in household contacts on index cases with COVID-19, 12 index cases (index cases from 2 households did not participate) and 23 contacts belonging to 14 households were tested. Demographic and clinical characteristics of the household participants are summarized in Table 1 and Supplementary Table 1A and B. The median age of household participants was 46 years (IQR, 20-62 years) and 49% were female. The median time from index case RT-PCR positivity to blood collection for this study was 17 days (IQR, 5-96 days) for index cases. The median time to blood collection for 17 contacts was 10 days (IQR, 6-23

days); for the remaining six contacts, median time to blood collection was 187 days (IQR, 149-202 days).

Six contacts (26.1%) were RT-PCR positive after exposure to index case, but four of them were RT-PCR positive on the day of sample collection for this study.

### **Accuracy of SARS-CoV-2 IGRA**

In COVID-19 convalescents tested with SARS-CoV-2 IGRA, the overall sensitivity over 10 months post-infection was 84.5% (153/181; 95% confidence interval [CI] 79.0%–89.0%) (Table 2 and Figure 2A). When only their first test result is considered, the overall sensitivity was 85.1% (80/94; 95% CI 76.4–91.0). In healthy blood donors the specificity was 86.6% (123/142; 95% CI 80.0-91.2) (Table 2 and Figure 2A). No indeterminate results were obtained in both groups. The sensitivity at 0.5 and 1-month post infection was 100% and dropped to 79.5% ( $P < 0.01$ ) at 10-month post infection (Table 2 and Figure 2A). Similar trends were observed when convalescents with three or more timepoints were analyzed (Supplementary Figure 3A). Compared with median IFN- $\gamma$  response in healthy blood donors (0.02 IU/mL, IQR 0.01-0.06), the IFN- $\gamma$  response was significantly higher in convalescents at 0.5-month (3.77 IU/mL, IQR 1.38-6.75,  $P < 0.001$ ), 1-month (2.96 IU/mL, IQR 1.31-6.10,  $P < 0.001$ ), 4-month (2.15 IU/mL, IQR 0.54-3.53,  $P < 0.001$ ), 7-month (1.16 IU/mL, IQR 0.53-3.86,  $P < 0.001$ ), and 10-month (1.26 IU/mL, IQR 0.38-2.96,  $P < 0.001$ ) post-infection (Figure 2A). The percent change in median responses between successive timepoints was -21.5% (0.5- to 1-month), -27.4% (1- to 4-month), -46.0% (4- to 7-month), and +8.6% (7- to 10-month). Among convalescents, compared with 0.5-month post-infection, IFN- $\gamma$  response was not significantly different at 1 and 4-month post-infection (Figure 2A). However, there was a

significant decline in IFN- $\gamma$  response at 7-month ( $P<0.05$ ), and 10-month ( $P<0.01$ ) post-infection (Figure 2A).

The overall sensitivity of IgG over 10 months post-infection was 83.9% (152/181; 95% CI 77.9-88.6) in convalescents. The specificity was 97.8% (3/181; 95% CI 93.9-99.4) in healthy blood donors (Table 2 and Figure 2B). The sensitivity was high at 0.5-month (87.5% [14/16]; 95% CI 63.9-97.8) and 1-month (100% [8/8]; 95% CI 67.6-100) post infection compared with 71.8% (28/39; 95% CI 56.2-83.5;  $P<0.001$ ) at 10-month (Figure 2B). Compared with median IgG OD ratio in healthy blood donors 0.24 (IQR 0.18-0.32) the median IgG OD ratio was 6.97 (IQR 4.25-10.36,  $P<0.001$ ) at 0.5-month, 9.04 (IQR 8.14-10.06,  $P<0.001$ ) at 1-month, 3.27 (IQR 1.94-4.52,  $P<0.001$ ) at 4-month, 3.11 (IQR 2.2-5.42,  $P<0.001$ ) at 7-month, and 2.07 (IQR 0.94-4.46,  $P<0.001$ ) at 10-month time points (Figure 2B). The percent change in median IgG OD ratios between successive timepoints was +29.7% (0.5- to 1-month), -63.8% (1- to 4-month), -4.9% (4- to 7-month), and -33.4% (7- to 10-month). Compared with 0.5-month post-infection, IgG OD ratio was significantly lower at 4 ( $P<0.01$ ), 7 ( $P<0.001$ ), and 10-month ( $P<0.001$ ) post-infection (Figure 2B).

When comparing IGRA with IgG, the overall concordance between qualitative results was 79.6% (144/181) in convalescents and 85.9% (122/142) in healthy blood donors. Overall percent agreement was 82.4% (266/323), positive percent agreement was 78.5% (135/172), and negative percent agreement was 86.8% (131/151). Quantitatively, no significant correlation was found between IGRA IFN- $\gamma$  response and IgG OD ratio ( $r=0.11$ ,  $P>0.05$ ) (Supplementary Figure 4).

### **SARS-CoV-2 IGRA in Household contacts**

In household members tested one time after recruitment, 100% (12/12) of the index cases and 65.2% (15/23) of the contacts were positive with the SARS-CoV-2 IGRA (Figure 3A and Supplementary Figure 5A). No indeterminate results were obtained. In 13 contacts with history of RT-PCR-positive SARS-CoV-2 infection (7 with remote history), 12 (92.3%) were SARS-CoV-2 IGRA positive. The median IFN- $\gamma$  response was 2.81 IU/mL (IQR, 0.67-5.73) in the index cases compared with 1.38 IU/mL (IQR, 0.03-2.7,  $P \geq 0.05$ ) in the contacts (Figure 3C).

SARS-CoV-2 IgG ELISA was positive in 50.0% (6/12) of the index cases and 43.5% (10/23) of the contacts (Figure 3B and Supplementary Figure 5B). In 13 contacts with history of SARS-CoV-2 infection, 6 (46.2%) were SARS-CoV-2 IgG positive. The median IgG antibody ratio was 1.11 (IQR, 0.28-1.91) in the index cases compared with 0.75 (IQR, 0.2-7.87,  $P > 0.05$ ) in the contacts (Figure 3D).

In the household cohort, concordance between IGRA and IgG results was 50.0% (6/12) and 60.9% (14/23) in the index cases and the contacts, respectively. When combining IGRA and IgG results, either test was positive in 100% (12/12) of the index cases and 73.9% (17/23) of the contacts (Figure 3C and 3D). Compared with IgG, IGRA had a differential positivity rate of +50% (95% CI +25.4-+74.6) and +21.7% (95% CI, +9.23-+42.3) in the index cases and the contacts, respectively. The positivity rate of IGRA was significantly higher than IgG in index cases and contacts ( $P < 0.05$  for both).

## Discussion

Accurate immunodiagnosics are needed to assess T cell-mediated immune response to SARS-CoV-2 after natural infection and after vaccination to inform providers on possible protective immunity to SARS-CoV-2 [2] and to accurately understand transmission dynamics of emerging variants [9]. Furthermore, T cell responses may be important for understanding immunity to SARS-CoV-2 infection, particularly in the context of emerging variants that might escape antibody responses [22, 23]. Using a whole blood SARS-CoV-2 IGRA, we show that T cell-mediated immune response to SARS-CoV-2 was sustained longitudinally with 85% and 80% sensitivity at 4- and 10-month, respectively, post-infection in convalescents who had mild or asymptomatic COVID-19. Our findings are consistent with a prior study showing a robust memory T-cell response months after SARS-CoV-2 infection in individuals with mild or asymptomatic infection [10, 11, 13, 24, 25], and similar to prior studies we observed a mild decline in the cellular and humoral immune response in convalescents [6, 26, 27]. Consistent with prior studies we also found a higher proportion of household contacts with positive SARS-CoV-2 IGRA compared with IgG [11, 12]. These findings indicate that the SARS-CoV-2 IGRA used in this study is more sensitive than IgG testing for detection of asymptomatic or mild infection amongst close contacts, at least in the early period, and thus may be an important tool for epidemiological studies aimed at acutely understanding the transmission dynamics of emerging SARS-CoV-2 variant such as the delta variant [28].

SARS-CoV-2 IGRA may also serve an important role in the clinical assessment of protective T-cell-mediated immune response to the SARS-CoV-2 vaccines [9, 29, 30]. Humoral immune response has been shown to decline after SARS-CoV-2 vaccination [26, 27, 31]. Furthermore, it was recently shown that the antibody response rate to the SARS-CoV-2 vaccine in

transplant patients was low after the second dose (40%) and third dose (68%) [32]. Preliminary analysis of SARS-CoV-2 vaccine response in immunocompromised patients at our health system using the in-tube version of our whole blood SARS-CoV-2 IGRA, which allows immediate stimulation and provides higher sensitivity [18], has revealed a significantly higher IGRA positivity rate compared with IgG positivity rate [21]. These findings support the role of SARS-CoV-2 IGRA in vaccine response assessment, particularly in elderly and immunocompromised patients given that immunosenescence and immunosuppression, respectively, may dampen adaptive immune responses and leave the host vulnerable to subsequent infection [4, 33, 34]. Thus, IGRA and IgG serology may serve a role in informing providers on the status of cell-mediated and humoral response to SARS-CoV-2 vaccine and the need for revaccination.

Although the findings are promising, this study has several limitations. First, the whole blood SARS-CoV-2 IGRA did not distinguish between CD4 and CD8 T cell responses. Prior studies in convalescents have shown that IFN- $\gamma$  response in IGRA is predominantly CD4 T cell-derived [18, 35]. However, a more complex assay design that allows measurement of CD4 and CD8 T cells response is possible if clinically indicated. Second, given that healthy blood donors were recruited during the pandemic, we would not know whether the 16% positivity rate with SARS-CoV-2 IGRA in this group was due to true-positive results due to past SARS-CoV-2 infection, false-positive results, or cross-reactivity due to past infection with seasonal CoV [12, 36, 37]. Third, the IgG response in this study was not further assessed for its ability to neutralize the virus. Such characterization was not relevant in the context of investigating IGRA accuracy. Forth, SARS-CoV-2 IGRA was not applied to the investigation of vaccine response in this study. Studies are underway to measure the vaccine response with IGRA in

immunocompromised patients. Lastly, the household study had a small sample size and lacked multiple time points to more accurately assess the performance of IGRA vs IgG, and we did not perform risk assessment in household contacts to correlate IGRA positivity with exposure risk.

In summary, the whole blood SARS-CoV-2 IGRA was shown to maintain sensitivity in convalescents up to 10 months post-infection and was shown to have a higher positivity rate than IgG in household contacts of COVID-19 cases. SARS-CoV-2 IGRA, including commercial options [38], is a simple and robust clinical immunodiagnostic test that can be applied to accurately understand the transmission of emerging variants and to assess vaccine response in vulnerable populations.

Accepted Manuscript

## NOTES

**Funding:** This work was not funded.

**Conflict of interest:** YAM reports grants or contracts from Respiratory Syncytial Virus (RSV) Prefusion F Subunit Vaccine in Infants Born to Women Vaccinated During Pregnancy, *Pfizer, C3671008, 07/15/2020 – 07/15/2025*, Stanford Precision Health for Ethnic and Racial Equity (SPHERE) Transdisciplinary Collaborative Center, *NIH, U54 MD010724, 04/11/16-8/31/2022*, Monitoring COVID19 and Building Capacity with Northern Plains Tribes and the Future of Pandemics, *NIH, U54 MD010724- 05S1, 04/11/16-8/31/2022*, A genomic tool for identifying pathogenic circulating vaccine-derived polioviruses, *NIH R21AI148810, 01/01/2020– 12/31/2021*, Stanford Aging & Ethnogeriatrics Transdisciplinary Collaborative Center (SAGE), *NIH P30AG059307, (Periyakoil), 09/30/2018 – 06/30/2023*, Stanford Human Systems Immunology Center, *Bill & Melinda Gates Foundation, OPP1113682 / (Davis), 11/04/2014 – 09/30/2020*, Congenital and Perinatal Infections Rare Diseases Clinical Research Consortium (RDCRC)-, *University of Alabama, NIH, 000522211-022 Maldonado, PI, Kimberlin, D, lead PI 09/01/2019 – 08/30/2022*, A Phase 1, Open-label dose finding study to evaluate safety, tolerability, and immunogenicity and phase 2/3 Placebo-controlled, observer-blinded safety, Tolerability, and Immunogenicity study of a SARS-Vov-3 RNA vaccine candidate against COVID-19 in healthy children <12 years of age. , *Pfizer, C4591007, Maldonado, site PI, 04/02/2021-04/01/2024*, COVID-19: Understanding When to Return



Safely to Work and School, *Chan Zuckerberg Foundation, 12089sc, Maldonado- consortium PI, Rutherford-lead PI, 04/30/2020 – 11/30/2021* all outside of the submitted work; payment or honoraria from American Academy of Pediatrics, National Conference meeting; and Pfizer, Data Safety Monitoring Board, non-COVID Vaccine trial. All other authors have no conflict of interest.

Accepted Manuscript

## References

1. Safiabadi Tali SH, LeBlanc JJ, Sadiq Z, et al. Tools and techniques for severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2)/COVID-19 detection. *Clin Microbiol Rev*, 2021; 34: 10.
2. Dan JM, Mateus J, Kato Y, et al. Immunological memory to SARS-CoV-2 assessed for up to 8 months after infection. *Science*, 2021; 371: 10.
3. Le Bert N, Clapham HE, Tan AT, et al. Highly functional virus-specific cellular immune response in asymptomatic SARS-CoV-2 infection. *J Exp Med*, 2021; 218: 10.
4. Abbasi J. The promise and peril of antibody testing for COVID-19. *JAMA*, 2020; 323: 1881-3.
5. Bange EM, Han NA, Wileyto P, et al. CD8(+) T cells contribute to survival in patients with COVID-19 and hematologic cancer. *Nat Med*, 2021; 27: 1280-9.
6. Bonifacius A, Tischer-Zimmermann S, Dragon AC, et al. COVID-19 immune signatures reveal stable antiviral t cell function despite declining humoral responses. *Immunity*, 2021; 54: 340,354.e6.
7. Cox RJ, Brokstad KA. Not just antibodies: B cells and T cells mediate immunity to COVID-19. *Nature reviews. Immunology*, 2020; 20: 581-2.
8. Cañete PF, Vinuesa CG. COVID-19 makes B cells forget, but T cells remember. *Cell (Cambridge)*, 2020; 183: 13-5.

9. Bertoletti A, Le Bert N, Qui M, Tan AT. SARS-CoV-2-specific T cells in infection and vaccination. *Cell Mol Immunol*, 2021; 18: 2307-12.
10. Wang Z, Yang X, Zhong J, et al. Exposure to SARS-CoV-2 generates T-cell memory in the absence of a detectable viral infection. *Nat Commun*, 2021; 12: 1724-z.
11. Sekine T, Perez-Potti A, Rivera-Ballesteros O, et al. Robust T cell immunity in convalescent individuals with asymptomatic or mild COVID-19. *Cell*, 2020; 183: 158,168.e14.
12. Gallais F, Velay A, Nazon C, et al. Intrafamilial exposure to SARS-CoV-2 associated with cellular immune response without seroconversion, france. *Emerging infectious diseases*, 2021; 27: 113-21.
13. Rodda LB, Netland J, Shehata L, et al. Functional SARS-CoV-2-specific immune memory persists after mild COVID-19. *Cell*, 2021; 184: 169,183.e17.
14. Ameratunga R, Woon S, Jordan A, et al. Perspective: Diagnostic laboratories should urgently develop T cell assays for SARS-CoV-2 infection. *Expert review of clinical immunology*, 2021; 17: 421-30.
15. Pai M, Denkinger CM, Kik SV, et al. Gamma interferon release assays for detection of mycobacterium tuberculosis infection. *Clin Microbiol Rev*, 2014; 27: 3-20.
16. Shafeque A, Bigio J, Hogan CA, Pai M, Banaei N. Fourth-generation QuantiFERON-TB gold plus: What is the evidence? *J Clin Microbiol*, 2020; 58.
17. Bert NL, Tan AT, Kunasegaran K, et al. SARS-CoV-2-specific T cell immunity in cases of COVID-19 and SARS, and uninfected controls. *Nature*, 2020; 584: 457-62.

18. Murugesan K, Jagannathan P, Pham TD, et al. Interferon- $\gamma$  release assay for accurate detection of severe acute respiratory syndrome coronavirus 2 T-cell response. *Clinical Infectious Diseases*, 2020.
19. Petrone L, Petruccioli E, Vanini V, et al. A whole blood test to measure SARS-CoV-2-specific response in COVID-19 patients. *Clinical microbiology and infection*, 2021; 27: 286.e7,286.e13.
20. Andreakos E, Tsiodras S. COVID -19: Lambda interferon against viral load and hyperinflammation. *EMBO Molecular Medicine*, 2020; 12.
21. Yang LM, Costales C, Ramanathan M, et al. Cell-mediated and humoral immune response to SARS-CoV2 mRNA vaccination and booster dose in immunosuppressed patients. *medRxiv*, 2022.
22. Tarke A, Sidney J, Methot N, et al. Impact of SARS-CoV-2 variants on the total CD4(+) and CD8(+) T cell reactivity in infected or vaccinated individuals. *Cell Rep Med*, 2021; 2: 100355.
23. Geers D, Shamier MC, Bogers S, et al. SARS-CoV-2 variants of concern partially escape humoral but not T-cell responses in COVID-19 convalescent donors and vaccinees. *Science Immunology*, 2021; 6: 10.
24. Stephens DS, McElrath MJ. COVID-19 and the path to immunity. *JAMA*, 2020; 324: 1279-81.

25. Schwarzkopf S, Krawczyk A, Knop D, et al. Cellular immunity in COVID-19 convalescents with PCR-confirmed infection but with undetectable SARS-CoV-2-specific IgG. *Emerging infectious diseases*, 2021; 27: 10.
26. Choe PG, Kang CK, Suh HJ, et al. Waning antibody responses in asymptomatic and symptomatic SARS-CoV-2 infection. *Emerging Infectious Diseases*, 2021; 27.
27. Chia WN, Zhu F, Ong SWX, et al. Dynamics of SARS-CoV-2 neutralising antibody responses and duration of immunity: A longitudinal study. *The Lancet Microbe*, 2021; 2: e240-9.
28. Farinholt T, Doddapaneni H, Qin X, et al. Transmission event of SARS-CoV-2 delta variant reveals multiple vaccine breakthrough infections. *medRxiv*, 2021.
29. Polack FP, Thomas SJ, Kitchin N, et al. Safety and efficacy of the BNT162b2 mRNA covid-19 vaccine. *N Engl J Med*, 2020; 383: 2603-15.
30. Kalimuddin S, Tham CYL, Qui M, et al. Early T cell and binding antibody responses are associated with COVID-19 RNA vaccine efficacy onset. *Med (N Y)*, 2021; 2: 682,688.e4.
31. Meyer B. Waning antibodies to SARS-CoV-2 - don't panic. *The Lancet Regional Health Europe*, 2021; 4: 100115.
32. Kamar N, Abravanel F, Marion O, Couat C, Izopet J, Del Bello A. Three doses of an mRNA covid-19 vaccine in solid-organ transplant recipients. *The New England journal of medicine*, 2021; 385: 661-2.

33. Winter AK, Hegde ST. The important role of serology for COVID-19 control. *The Lancet infectious diseases*, 2020; 20: 758-9.
34. Chen L, Xiong J, Bao L, Shi Y. Convalescent plasma as a potential therapy for COVID-19. *The Lancet infectious diseases*, 2020; 20: 398-400.
35. Grifoni A, Weiskopf D, Ramirez SI, et al. Targets of T cell responses to SARS-CoV-2 coronavirus in humans with COVID-19 disease and unexposed individuals. *Cell*, 2020; 181: 1489,1501.e15.
36. Braun J, Loyal L, Frensch M, et al. SARS-CoV-2-reactive T cells in healthy donors and patients with COVID-19. *Nature*, 2020; 587: 270-4.
37. Mateus J, Grifoni A, Tarke A, et al. Selective and cross-reactive SARS-CoV-2 T cell epitopes in unexposed humans. *Science*, 2020; 370: 89-94.
38. Huzly D, Panning M, Smely F, et al. Validation and performance evaluation of a novel interferon- $\gamma$  release assay for the detection of SARS-CoV-2 specific T-cell response. *medRxiv*, 2021.

**Table 1.** Demographic and clinical characteristics of household members.

Demographic and Clinical Characteristics		Index cases (n=12)		Contacts (n=23)	
		n	%	n	%
<b>Median age (IQR)</b>		56 (24-66)		44 (16-54)	
<b>Male</b>		7	58.3	11	47.8
<b>Race</b>	<b>White</b>	8	67.7	17	73.9
	<b>Hispanic</b>	4	33.3	4	17.4
	<b>Asian</b>	-	-	2	8.7
<b>Symptoms</b>	<b>Fever</b>	4	33.3	2	8.7
	<b>Cough</b>	6	50.0	6	26.1
	<b>Sore throat</b>	3	25.0	3	13.0
	<b>Running nose</b>	6	50.0	3	13.0
	<b>Chest pain</b>	3	25.0	1	4.3
	<b>Joint pain</b>	1	8.3	2	8.7
	<b>Fatigue</b>	4	33.3	5	21.7
	<b>Shortness of breath</b>	3	25.0	2	8.7
	<b>Loss of smell</b>	2	17.7	1	4.3
	<b>Diarrhea</b>	-	-	1	4.3
	<b>Headache</b>	1	8.3	3	13.0

Accepted Manuscript

**Table 2.** Accuracy of SARS-CoV-2 IGRA and IgG in convalescents and healthy blood donors.

	Months Since SARS-CoV-2 Infection (No. of Subjects)					Healthy Blood Donors (No. of Subjects)	
	0.5 (16)	1 (8)	4 (48)	7 (70)	10 (39)		
IGRA	Sensitivity (n/N; 95% CI)	100 (16/16; 80.6-100)	100 (8/8; 67.6-100)	85.4 (41/48; 72.8-92.8)	81.4 (57/70; 70.7-88.8)	79.5 (31/39; 64.5-89.2)	-
	Specificity (n/N; 95% CI)	-	-	-	-	-	86.6 (123/142; 80.0-91.2)
	IgG	87.5 (14/16; 63.9-97.8)	100 (8/8; 67.6-100)	87.5 (42/48; 75.3-94.1)	85.7 (60/70; 75.7-92.1)	71.8 (28/39; 56.2-83.5)	-
	Specificity (n/N; 95% CI)	-	-	-	-	-	97.8 (3/181; 93.9-99.4)

No, number; CI, confidence interval

Accepted Manuscript



**Figure legends:**

**Figure 1. Patient cohorts included in this study to evaluate the performance of whole blood SARS-CoV-2 IGRA.**


**Figure 2. SARS-CoV-2 IGRA and IgG ELISA results in convalescents and healthy blood donors.** IFN- $\gamma$  response with whole blood IGRA (A) and IgG antibody optical density (OD) ratio (B) in healthy blood donors (HBD) (n=142) and COVID-19 convalescents at 0.5 (n=16), 1 (n=8), 4 (n=48), 7 (n=70), and 10 (n=39) months post-infection. Overall sensitivity over 10 months post-infection was 84.5% (153/181; 95% CI 79.0%–89.0%) with IGRA and 83.9% (152/181; 95% CI 77.9-88.6) with IgG. Dotted lines represent the assay cutoffs (0.35 IU/mL for IGRA and 1.1 OD ratio for IgG). Quantitative results >10 IU/mL are shown as 10 IU/mL and >12 OD ratio are shown as 12 OD ratio. Horizontal lines show the median. whiskers show the range of IFN- $\gamma$  response (IU/mL). Response in convalescents was compared to controls (blue). Response in convalescents were compared to 0.5 month time point (red) \*, P < 0.05; \*\*, P < 0.01; \*\*\*, P < 0.001.

**Figure 3. SARS-CoV-2 IGRA and IgG ELISA result in index cases and household contacts.**

Venn diagrams show positive results with IGRA and/or IgG in 12 index cases (A) and 23 contacts (B).

Graphs show IFN- $\gamma$  response with whole blood IGRA (C) and IgG antibody optical density (OD) ratio (D) in index cases and contacts. Dotted lines represent the assay cutoffs (0.35 IU/mL for IGRA and 1.1 OD ratio for IgG). Quantitative results >10 IU/mL are shown as 10 IU/mL and >12 OD ratio are shown as 12 OD ratio. Horizontal lines show the median. Whiskers show the range. Response in index cases was compared to household contacts. ns, not significant.

Figure 1



Accepted

Figure 2




Figure 3

