

ARTICLE



Molecular Diagnostics

Sensitive and robust liquid biopsy-based detection of *PIK3CA* mutations in hormone-receptor-positive metastatic breast cancer patients

Christoph Suppan¹, Ricarda Graf², Stephan Jahn³, Qing Zhou², Eva Valentina Klocker¹, Rupert Bartsch⁴, Angelika Terbuch¹, Karl Kashofer³, Peter Regitnig³, Joerg Lindenmann⁵, Florian Posch¹, Hanno Gerritsmann⁶, Philipp J. Jost^{1,7}, Ellen Heitzer^{2,8,9}, Nadia Dandachi^{1,10} and Marija Balic^{1,11}

© The Author(s), under exclusive licence to Springer Nature Limited 2021

BACKGROUND: The benefit of alpelisib in hormone-receptor-positive (HR+) metastatic breast cancer patients provided clinical evidence for the increasing importance of *PIK3CA* testing. We performed a comparison of liquid biopsy and tissue-based detection of *PIK3CA* mutations.

MATERIALS AND METHODS: *PIK3CA* hotspot mutation analysis using a high-resolution SiMSen-Seq assay was performed in plasma from 93/99 eligible patients with HR+/HER2– breast cancer. Additionally, mFAST-SeqS was used to estimate the tumour fractions in plasma samples. In 72/93 patients, matched tissue was available and analysed using a customised Ion Torrent panel.

RESULTS: *PIK3CA* mutations were detected in 48.6% of tissue samples and 47.3% of plasma samples, with identical *PIK3CA* mutation detected in 24/72 (33.3%) patients both in tissue and plasma. In 10 (13.9%) patients, mutations were only found in plasma, and in 6 (8.3%) patients, *PIK3CA* mutations found in tissue were not detectable in ctDNA. In 49/93 plasma samples without detectable *PIK3CA* mutations, 22 (44.9%) samples had elevated tumour fractions, implying true negative results.

CONCLUSION: SiMSen-Seq-based detection of *PIK3CA* mutations in plasma shows advantageous concordance with the tissue analyses. A combination with an untargeted approach for detecting ctDNA fractions may confirm a negative *PIK3CA* result and enhance the performance of the SiMSen-Seq test.

British Journal of Cancer (2022) 126:456–463; <https://doi.org/10.1038/s41416-021-01601-9>

INTRODUCTION

Breast cancer (BC) remains the leading cause of cancer-related death in the global female population [1]. Despite the substantial progress in systemic treatment, metastatic disease remains incurable, but novel treatment options offer the chance for prolonged disease stabilisation [2]. However, particularly for the treatment of hormone-receptor (HR)-positive disease, substantial progress in systemic treatment has been made. After the introduction of CDK4/6 inhibitors, targeted signal-transduction pathways have been established in the meantime [3]. Phosphatidylinositol 3-kinases (PI3Ks) are among the most important targets in this regard, apparently mediating the conversion of phosphatidylinositol 4,5-bisphosphate (PI2) to phosphatidylinositol 3,4,5-triphosphate (PIP3) [4]. The PI3K/AKT/mTOR pathway's ensuing upregulation plays a critical role in essential processes of

disease progression, including cell growth, proliferation, metabolism and invasion [5].

In breast cancer, the PI3K/AKT/mTOR signaling pathway can be altered by different mechanisms. In approximately a third of early breast cancers activating mutations in the helical or kinase domain of *PIK3CA* can be found (45% in luminal A, 29% in luminal B, 39% in HER2-enriched and 9% in the basal carcinomas, respectively) [6, 7]. Similar rates of mutations can be found in metastatic breast cancer, although fewer data are available [6].

Recent data suggest that *PIK3CA* alterations have an adverse impact on the prognosis of HR-positive disease and may play a role in endocrine resistance [8–10]. Several drugs targeting the PI3 kinase pathway have been tested, but to date, a clinically meaningful benefit of adding an α -specific PI3K inhibitor (alpelisib) with an acceptable tolerability profile was only demonstrated in

¹Division of Oncology, Department of Internal Medicine, Medical University of Graz, Graz, Austria. ²Institute of Human Genetics, Diagnostic and Research Center for Molecular Biomedicine (Austria), Medical University of Graz, Graz, Austria. ³Institute of Pathology, Medical University of Graz, Graz, Austria. ⁴Department of Medicine 1, Clinical Division of Oncology, Medical University of Vienna, Vienna, Austria. ⁵Division of Thoracic Surgery and Hyperbaric Surgery, Department of Surgery, Medical University of Graz, Graz, Austria. ⁶Novartis Pharma GmbH, Jakob-Lind-Straße 5, 1020 Vienna, Austria. ⁷Department of Medicine III, Klinikum rechts der Isar, TUM School of Medicine, Technical University of Munich, Munich, Germany. ⁸BioTechMed-Graz, Graz, Austria. ⁹Christian Doppler Laboratory for Liquid Biopsies for Early Detection of Cancer, Graz, Austria. ¹⁰Research Unit Epigenetic and Genetic Cancer Biomarkers, Medical University of Graz, Graz, Austria. ¹¹Research Unit Translational Breast Cancer, Medical University of Graz, Graz, Austria. ✉email: nadia.dandachi@medunigraz.at; marija.balic@medunigraz.at

Received: 24 February 2021 Revised: 6 October 2021 Accepted: 12 October 2021

Published online: 9 November 2021

the SOLAR-1 phase III trial [11]. In this trial, only a minor proportion of patients was pretreated with CDK4/6 inhibitors, but the benefit of PI3K inhibition in this patient's subset was confirmed in the phase II BYLieve study [12]. With these data, the importance of the *PIK3CA* encoding for the α -isoform of the catalytic subunit of the class IA PI3K kinase as a biomarker was postulated; the best testing approach, however, remains elusive. The stability of mutations throughout the disease was only recently addressed by Mosele et al. [9], and prospective collection of metastatic tissue is ongoing in several large clinical trials such as AURORA (Aiming to Understand the Molecular Aberrations in Metastatic Breast Cancer, NCT02102165) [13, 14]. As tissue from metastatic sites may not always be readily available, interest in liquid biopsy techniques has grown. In the SOLAR-1 trial, in parallel to the tissue-based detection, cell-free DNA was analysed for the presence of *PIK3CA* mutations. Along with the approval of alpelisib, FDA has approved the companion diagnostic test, therascreen® *PIK3CA* RGQ PCR Kit (QIAGEN Manchester, Ltd.), to select patients who have *PIK3CA* mutations in tumour tissue specimens and/or in circulating tumour DNA (ctDNA) isolated from plasma specimens.

Therefore, our study's primary aim was to compare plasma-based detection of *PIK3CA* hotspot mutations with the tissue-based detection from the most recent tumour tissue.

MATERIALS AND METHODS

Patients and sample collections

Between February 2017 and June 2021, 99 locally advanced or metastatic breast cancer patients were recruited for this prospective observational study at the Division of Oncology, Department of Internal Medicine, Medical University of Graz, Austria. One patient who had received a bone marrow transplantation and three patients who had triple-positive or triple-negative tumours were excluded from the study (Fig. 1). The study was approved by the ethics committee of the Medical University Graz (ethical approval number 21–227 ex 09/10), and written informed consent was obtained from all patients. All analyses were carried out considering the guidelines for good scientific practice as officially required from the Medical University of Graz. Clinical data and histopathological tumour characteristics were collected from clinical and pathological records in an electronic database system. Blood samples were obtained either at first diagnosis of metastases, during varying lines of treatment, and/or at a




Fig. 1 CONSORT diagram for the study. The numbers of patients and samples included in the study and reasons for the exclusion are shown.

further moment of progression/development of new metastases, along with a planned introduction of a new line of systemic treatment. The most recent tissue was used for tissue analyses, either metastatic tissue, or if metastatic tissue was not available the primary tumour tissue was analysed.

Preparation of plasma DNA

Isolation of plasma DNA was performed as previously described [15]. Briefly, whole blood (18 mL) was collected in PAXgene Blood ccfDNA Tubes (Qiagen, Hilden, Germany). Plasma was isolated from blood samples by centrifugation at $200 \times g$ for 10 min, followed by $1600 \times g$ for 10 min with brake and acceleration powers set to slow. The supernatant was then removed and centrifuged again at $1600 \times g$ for an additional 10 min. Plasma was then stored in 2 mL tubes at -80°C . DNA was isolated from 2 mL plasma using the QIAamp Circulating Nucleic Acid Kit (Qiagen, Hilden, Germany) according to the manufacturer's recommendations or using the QIASymphony PAXgene Blood ccfDNA Kit (Qiagen, Hilden, Germany). Qubit dsDNA HS Assay Kit (Life Technologies, Vienna, Austria) was used for the quantification of plasma DNA.

High-resolution screening for *PIK3CA* mutations in plasma

To screen for *PIK3CA* mutation, the SiMSen-Seq method (simple, multiplexed, PCR-based barcoding of DNA for sensitive mutation detection using sequencing) was used [16]. This technology is based on molecular barcoding of individual DNA template strands to track all sequencing reads back to a single original template and correct for PCR errors during library preparation. To assess the most common hotspots ($n = 11$), three assays covering exons 7, 9 and 20 were designed and validated as described in ref. [16] (see Table S1, Table S2 and Supplementary methods). After a comprehensive evaluation of the assay performance using SeraCare reference materials and control samples, we set the limit of detection (LOD) to 0.25%. In addition, samples were only considered *PIK3CA* mutated if ≥ 5 consensus reads carried the mutation (for details, see supplementary material, Supplementary Figs. S1 and S2).

For the analysis of plasma samples, an average input of 26.3 ng (range, 8.5–103.8 ng) cDNA per sample was used for library preparation. After purification with a 0.8-fold volume of AMPure XP Beads (Beckman Coulter, Brea, CA, USA), libraries were quantified with qPCR using the QIAseq Library Quant Assay (QIAGEN, Hilden, Germany). Libraries were sequenced on an Illumina MiSeq, or NextSeq using 150 bp single-end reads aiming for 10 million reads per sample. Barcoded SiMSen-Seq data were analysed using Debarcer [16]. After consensus read generation, we obtained an average coverage of 2185x (range 49–17227x) for all interrogated positions. The molecular coverage of the SiMSen-Seq assay significantly correlated with the amount of input DNA used for sequencing (Spearman's $\rho = 0.66$, $p < 0.001$; Supplementary Fig. S3). The presence of mutations was reported as variant allele frequency (VAF, percentage of alternate reads relative to the total number of reads at a given position).

The raw sequencing datasets have been deposited at the European Genome-phenome Archive (EGA; <http://www.ebi.ac.uk/ega/>) under study accession number EGAS00001004940 as following dataset: *PIK3CA* SiMSen-Seq data of plasma samples.

mFAST-SeqS

LINE1-amplicon libraries for mFAST-SeqS were prepared as previously described [17]. Briefly, 1 ng of plasma DNA was amplified with Phusion Hot Start II Polymerase in eight PCR cycles using target-specific L1 primers. The 10 μL of the purified (AMPure Beads, Beckman Coulter, Brea, CA, USA) PCR products were used for the second PCR in which Illumina-specific adaptors and indices were added. L1 amplicon libraries were pooled equimolarly and sequenced on an Illumina MiSeq or NextSeq, generating 150 bp single reads or 76 bp paired-end reads aiming for 100,000 reads per sample. Sequence reads were aligned to the hg19 genome, and reads with a mapping quality > 15 were counted per chromosome arm using an in-house script. Normalised read counts were compared to a control population ($n = 35$) using z-score statistics by subtracting the mean and dividing by the standard deviation to assess over- and under-representation of LINE1-sequences. Because little or no reads aligned to the short arms of acrocentric chromosomes 13p, 14p, 15p, 21p and 22p, these were omitted from the analysis. Finally, all chromosome arm-specific z-scores were squared and summed up, resulting in a genome-wide z-score, which acts as a surrogate for tumour fraction [18].

Tissue-based analyses

DNA was extracted from macrodissected tumour areas of unstained FFPE sections. Adequate tumour cell content was verified from HE stained parallel sections. Extraction was performed using the Maxwell 16 instrument (Promega, Wisconsin, USA) and the Maxwell RSC DNA FFPE Kit (Promega, CatNr: AS1450) or the Ion AmpliSeq Direct FFPE DNA Kit (Thermo Fisher Scientific, Massachusetts, USA). DNA was quantified by Picogreen fluorescence, and typically, 40 ng DNA were used for library preparation. NGS libraries were prepared using a custom-designed AmpliSeqHD panel covering the following mutation hotspots: *EGFR* (L858), *ERBB2* (S310, L755, V762, V777, P780), *ERBB3* (V104, E928), *KRAS* (G12, G13), *PTEN* (T319), *ESR1* (E207, K303, E380, Y537, D538), *PIK3CA* (E103, N345, C420, E542, E545, Q546, E726, H1047, G1049), *AKT1* (E17), *TP53* (R175, H179, R213, Y220, R248, R273). The panel includes all mutations evaluated in the SOLAR-1 trial [11]. Molecular barcoding enables a detection threshold of 0.1% mutant allele. Sequencing was performed on an Ion S5XL benchtop sequencer (Thermo Fisher Scientific) to a length of 200 bp. We used the open-source Ion Torrent software (Thermo Fisher Scientific, open source, GPL, <https://github.com/iontorrent/>) for initial data analysis. Briefly, this included base calling, alignment to the reference genome (hg19) using the TMAP mapper and variant calling by a modified diBayes approach taking the flow space information into account. Called variants were deconvoluted with regards to molecular barcodes and annotated using open-source software ANNOVAR [19] and SnpEff [20]. All coding, nonsynonymous mutations were further evaluated and visually inspected in IGV (<http://www.broadinstitute.org/igv/>), and variant calls resulting from technical read errors or sequence effects were excluded from the result. This data are available on request from the authors.

Statistical analyses

All statistical analyses were performed with Stata 16.1 (Stata Corp., Houston, TX, USA) and Graphpad. Continuous variables were reported as medians [25–75th percentile] and count data as absolute frequencies (%). The distribution of baseline variables between patients with or without a *PIK3CA* mutation (binary variable) was assessed with rank-sum tests (continuous variables) and χ^2 , or with Fisher's exact tests (categorical variables). The agreement between plasma and tissue was evaluated with a kappa score and 95% CIs. The exact McNemar's test for comparing paired sample proportions was used to assess the statistical significance of differences between plasma and tissue results. The correlation between z-scores and mutant *PIK3CA* allele fractions was determined by the Spearman rank correlation coefficient. In an explorative analysis, we associated *PIK3CA* mutation status with progression-free survival using Kaplan–Meier estimators, log-rank tests and Cox proportional hazards regression. Figures were generated with Stata or GraphPad Prism (v9).

RESULTS

Patient characteristics

Ninety-nine patients were enrolled, and 95 advanced or metastatic breast cancer patients with a median age of 65.4 years [25–75th percentile: 54.9–74.7] were analysed in this study. Table 1 shows the main patients' characteristics. The majority of patients were female (95.8%), and all patients had HR+/-HER2– tumours (100%). The most frequent metastatic site was bone (66.3%), followed by lymph nodes (31.6%), lung (24.2%), liver (22.1%) and pleura (10.5%). Fifty patients (52.6%) had multiple metastatic sites, and most patients (80.0%) received first or second-line palliative treatment at study inclusion. Seventy-four tissue samples and ninety-three plasma samples were available for *PIK3CA* mutation analysis (Fig. 1). For the comparison cohort of 72 paired tissue and plasma samples, the median time interval between tissue and plasma sample collection was 5.5 months [25–75th percentile: 0.7–34.7 months].

PIK3CA mutations and z-scores in patients with evaluable plasma samples (N = 93)

Of all 93 patients with evaluable plasma samples, 47.3% (44/93) had detectable *PIK3CA* mutations. In detail, 48 *PIK3CA* mutations were detected, and in four samples, the co-occurrence of two *PIK3CA*

Table 1. Baseline characteristics of the study population (N = 95).

	Total cohort N = 95
Age at inclusion (years)	65.4 (54.9–74.7)
Female gender	91 (95.8%)
Histological type	
IDC/NST	59 (62.1%)
ILC	23 (24.2%)
Other or not reported	13 (13.7%)
Molecular subtype	
HR+/HER2–	95 (100.0%)
Bone metastases	63 (66.3%)
Lung metastases	23 (24.2%)
Liver metastases	21 (22.1%)
Pleura metastases	10 (10.5%)
Lymph node metastases	30 (31.6%)
Number of metastatic sites	
0–1 sites ^a	45 (47.4%)
≥2 sites	50 (52.6%)
Treatment line at inclusion	
≤2 lines	76 (80.0%)
>2 lines	19 (20.0%)
<i>PIK3CA</i> mutation detected in plasma (N = 93) ^b	44 (47.3%)
<i>PIK3CA</i> mutation detected in tissue (N = 74) ^c	36 (48.6%)
Source of tissue (N = 74)	
Metastasis	45 (60.8%)
Primary tumour	29 (39.1%)
Interval between tissue and blood sample collection (months)	7.7 (0.7–40.5)

Data are medians [25–75th percentile] for continuous data and absolute frequencies (%) for count data.

^aIncludes one patient with locally advanced disease.

^bNGS failed in two plasma samples.

^cNo sufficient tumour material was available for NGS testing in 21 cases.

mutations was observed (Table 2). Double or compound *PIK3CA* mutations are described in *cis* on the same allele; they are supposed to increase PI3K activity and enhance downstream signaling, cell proliferation and tumour growth [21]. The VAF of detected *PIK3CA* mutations in plasma ranged from 0.27 to 49.85%, with a median of 1.67%. The most frequent mutation reported in plasma samples was H1047R (32.6%), followed by E542K (9.7%) and E545K (7.5%).

Metastatic sites, including bone, liver, lung, pleura and lymph nodes, were not significantly associated with *PIK3CA* mutated ctDNA. Multiple metastases were numerically (56.3% versus 37.8%) but not statistically associated with *PIK3CA* mutations (χ^2 test $p = 0.075$). Finally, there was no difference between the line of treatment and the presence of *PIK3CA* mutations in plasma.

The median z-score resulting from mFAST-SeqS analyses was 2.13 [25–75th percentile: 1.53–4.18]. Based on a previously established cut-off at ≥ 3 , which depending on the amount and amplitude of copy number alterations, correlates to a tumour fraction of 5–10% [22], 33/93 (35.5%) patients had elevated z-scores. We found no statistically significant correlation between z-scores and mutant *PIK3CA* allele frequencies (Spearman's $\rho = -0.11$, $p = 0.285$). Surprisingly, z-scores were higher in samples where no *PIK3CA* mutation was detected (median 2.43, range 0.85–69.68) than *PIK3CA* mutated samples (median 1.82, range 0.75–16.59; Wilcoxon rank-sum test $p = 0.042$; Supplementary Fig. S4). 22/49 (44.9%) of patients with no

Table 2. Patients with double *PIK3CA* mutations detected in plasma.

Patient ID	First mutation in plasma (AF %)	Second mutation in plasma (AF %)	First mutation in tissue (%MAF)	Second mutation in tissue (%MAF)	Δ time (months)
B211	E542K (12.76)	H1047R (0.33)	E542K (28.57)	None detected	20.3
B235	E542K (1.34)	H1047R (0.77)	None detected	None detected	22.6
B259	E542K (1.71)	H1047R (4.42)	E542K (13.52)	None detected	60.3
B387	E545K (16.36)	E542K (9.87)	E545K (62.09)	E542K (26.66)	17.6

detectable plasma *PIK3CA* mutation presented z-scores ≥ 3 , indicating a high tumour content. Since ctDNA but no *PIK3CA* mutation could be detected in these samples, these results can be considered true negatives for the assessed *PIK3CA* mutations. In contrast, for 27/49 (55.1%) patients with low (<3) z-scores and, therefore, low tumour fractions in plasma, the presence of *PIK3CA* mutations below the LOD cannot be entirely excluded. In these cases, tissue-based *PIK3CA* analysis might be indicated. Taken together, in 27 out of 93 patients (29.0%), the *PIK3CA* status was negative with a z-score < 3 based on ctDNA analyses. For these 27 patients, 20 tissue samples were evaluable for *PIK3CA* testing, and *PIK3CA* mutations were detected only in 4/20 (20.0%) of these patients.

***PIK3CA* mutations in patients with evaluable tissue samples (N = 74)**

PIK3CA mutations were detected in 48.6% (36/74) of archival tumour tissue samples (29 primary, 45 metastatic). In one tumour sample, two *PIK3CA* mutations were observed (Table 2). VAFs ranged from 4.36 to 72.93% (median 28.64%). The most frequent *PIK3CA* mutation detected in tissue was H1047R (23.0%), followed by E545K (13.5%), E542K (9.5%) and H1047L (2.7%). The location and number of metastases, as well as treatment lines, were not associated with *PIK3CA* mutations in tissue samples.

Concordance of *PIK3CA* mutation status in plasma and tissue

Comparative analyses of plasma and tissue were performed on samples from 72 patients for the 11 *PIK3CA* hotspot mutations covered by both assays (Supplementary Fig. S5). We found a moderate agreement between tissue and plasma samples, with a Cohen kappa = 0.58 (95% CI, 0.426–0.651) and an overall concordance rate of 72.2% (52/72, Table 3). Overall concordance was similar if only metastatic tissue samples (32/45, 71.1%) or only primary tissue samples (20/27, 74.1%) were considered. In detail, 24/72 patients (33.3%) had the same *PIK3CA* alteration in plasma and tissue, while no *PIK3CA* mutation in either sample type was found in 28/72 patients (38.9%). Discordant results (20/72, 27.8%; Table 4) included 10 (13.9%) patients with a *PIK3CA* mutation in plasma only and six patients (8.3%) with a *PIK3CA* mutation in tissue only. In addition, in four patients (5.6%), a different *PIK3CA* mutation was detected in the plasma sample compared to the tissue sample. The overall comparison showed no statistical difference between the results of tissue and plasma (McNemar's test $p = 0.317$). For concordant samples, the *PIK3CA* VAF in ctDNA were generally higher (median 2.78, range 0.57–17.70) than for the discordant samples (median 0.97, range 0.27–2.97; Wilcoxon rank-sum test $p = 0.001$), where mutations were only detected in plasma (Supplementary Fig. S6). Z-scores were similar in discordant samples (median 2.09, range 0.77–7.90) compared to concordant samples (median 2.12, range 0.85–69.68; Wilcoxon rank-sum test $p = 0.475$). Within the 20 discordant samples, 6 had z-scores ≥ 3 . Finally, the time interval between plasma and tissue collection did not differ significantly and ranged from –9.9 to 139.2 months (median 1.9 months) for discordant samples and from –14.8 to 229.3 months (median 9.6 months) for concordant samples (Wilcoxon rank-sum test $p = 0.119$).

Overall in the comparison cohort (N = 72), 38 patients had *PIK3CA* mutations (52.8%) detectable in plasma samples,

Table 3. Summary of *PIK3CA* mutation concordance data comparing plasma with tissue samples (N = 72).

	No	%
Total concordance (no alteration or identical alteration)	52	72.2
No alterations found in plasma and tissue	28	38.9
Identical alterations found in plasma and tissue	24	33.3
Total discordant samples	20	27.8
Different alterations found in plasma and tissue	4	5.6
Plasma only	10	13.9
Tissue only	6	8.3
Total number of alterations found in tissue	35	–
Total number of alterations found in plasma	42	–
Samples with double mutations in tissue	1	1.4
Samples with double mutations in plasma	4	5.6

compared to 34 patients with *PIK3CA* mutations in tissue samples (47.2%). In four plasma samples, double *PIK3CA* mutations were detected, while only one double mutation occurred in tumour tissue samples. H1047R was the most frequent mutation found in both plasma and tissue samples (34.7% in plasma, 22.2% in tissue), followed by E545K (8.3% in plasma, 12.5% in tissue), E542K (8.3% in both plasma and tissue) and C420R (1.4% in both plasma and tissue). The H1047L variant was detected in two tissue samples (2.8%) but not in plasma samples.

Focusing only on the question of whether we were able to identify patients qualifying for alpelisib treatment based on the detection of any of the 11 *PIK3CA* mutations included in the SOLAR-1 trial [11], we detected at least one mutation in 28/72 patients both in tissue and plasma. In 10 patients, such mutations were found only in plasma. Among six patients with tissue *PIK3CA* mutations not detectable in ctDNA, four patients had low tumour fractions (z-score < 3). Comparing our results to the liquid biopsy-based results in the SOLAR-1 study, we found that the proportions of discordant plasma samples compared to wild-type tissue samples were significantly higher in our study ($p < 0.001$). In contrast, the proportion of discordant plasma samples compared to mutant tissue samples was considerably lower in our study ($p = 0.002$, Fig. 2). (https://www.accessdata.fda.gov/cdrh_docs/pdf19/P190004B.pdf) These results indicate that our SIMSen-Seq assay was more sensitive than the theascreen assay used in the SOLAR-1 study [11].

Finally, we performed an exploratory analysis assessing the association between *PIK3CA* mutation and progression-free survival in a subset of 47 HR+/HER2– breast cancer patients receiving CDK inhibitor treatment. In this subset of patients, the presence of a *PIK3CA* mutation in plasma was not significantly associated with PFS (HR 1.79, 95% CI 0.73–4.37, log-rank $p = 0.197$; Supplementary Fig. S7)

DISCUSSION

This analysis with patient-paired plasma and tissue samples shows that a liquid biopsy-based approach holds promise to identify

Table 4. Discordant *PIK3CA* mutation results ($N = 20$) between plasma and tissue samples.

Patient ID	Mutation(s) in plasma (AF %)	z-score	Mutation in tissue (MAF %)	Source of tissue	Time interval between plasma and tissue collection ^a
B176	H1047R (0.27)	7.39	ND	Metastasis	–11 days
B195	H1047R (2.97)	4.18	ND	Metastasis	6 days
B200	H1047R (0.98)	1.75	ND	Metastasis	4.4 months
B201	E542K (0.46)	2.11	H1047R (21.1)	Metastasis	22 days
B205	H1047R (0.72)	7.90	ND	Metastasis	14 days
B207	H1047R (0.97)	3.30	ND	Metastasis	–3.8 months
B209	H1047R (1.17)	1.22	ND	Primary	6.6 months
B212	ND	3.19	E545K (30.7)	Primary	–9.9 months
B218	ND	3.40	E542K (23.7)	Metastasis	16.1 months
B235	E542K (1.34), H1047R (0.77)	1.69	ND	Primary	22.6 months
B242	ND	2.07	E545K (55.8)	Metastasis	10.1 months
B275	ND	2.51	H1047R (25.6)	Primary	2.2 months
B285	H1047R (0.31)	1.81	ND	Metastasis	–4.6 months
B296	H1047R (0.58)	0.77	H1047L (30.5)	Metastasis	2.3 months
B303	H1047R (0.42)	0.81	E542K (29.5)	Metastasis	7 days
B304	H1047R (1.08)	2.20	ND	Primary	139.2 months
B307	H1047R (0.41)	1.51	H1047L (25.9)	Primary	88.5 months
B341	ND	2.43	E545K (30.9)	Primary	20 days
B350	H1047R (0.47)	1.29	ND	Metastasis	1.8 months
B354	ND	0.94	H1047R (15.9)	Metastasis	119.0 months

ND no mutation detected.

^aNegative numbers mean that plasma samples were analysed at an earlier date than tissue samples.

patients with the most common activating *PIK3CA* mutations. Using the high-resolution NGS-based SiMSen-Seq approach, *PIK3CA* mutations were detected non-invasively in 47% of our patients with a moderate concordance to tissue data (overall concordance rate 72%). Although digital PCR may currently be regarded as the gold standard for mutation detection at hotspots, it can be particularly challenging when analyzing multiple targets since either multiple reactions are needed, or no specific sequence information can be obtained. Although SiMSen-Seq is more expensive than digital PCR, all targets can be analysed within one reaction, and sequence information of the mutations is provided. A combination with an untargeted, mutation-independent approach for detecting ctDNA fractions might confirm a negative *PIK3CA* result, thereby enhancing the performance of the SiMSen-Seq test. This approach allowed the optimisation of the selection of candidate patients for *PIK3CA* targeted treatment.

The concept of ctDNA based detection of *PIK3CA* mutations as predictive biomarkers for the PI3K inhibitors was already evaluated in the phase III trials BELLE-2 and BELLE-3 with the Pan-PI3K inhibitor buparlisib [8, 23], and in the SANDPIPER trial with the alpha-specific PI3K inhibitor taselisib [24]. Although these trials were positive, based on the toxicity and the magnitude of benefit, neither buparlisib nor taselisib have been clinically developed further in HR+/HER2– breast cancer. However, *PIK3CA* mutation testing based on liquid biopsy has shown promising results as a predictive biomarker for PI3K inhibitors [25], and these data were more recently confirmed in the SOLAR-1 trial. In contrast, *PIK3CA* mutations in ctDNA did not predict the clinical benefit of everolimus in the BOLERO2 study. Among 550 patients included in the ctDNA analyses, median PFS was comparable in *PIK3CA* wild-type and mutant patients receiving everolimus [10].

The FDA approval of alpelisib includes the companion diagnostic test theascreen® *PIK3CA* RGQ PCR Kit [26] and liquid

biopsy testing. Details on the detection of *PIK3CA* mutations in ctDNA in the SOLAR-1 trial have not yet been published, and only conference papers refer to the sensitivity of the liquid biopsy-based testing to detect *PIK3CA* mutant patients [27]. Out of 328 patients identified as positive in the tissue, in 54.9% of patients, *PIK3CA* mutations could be confirmed in ctDNA via the qPCR assay. These data indicate that the sensitivity of the test performed in patients included in the SOLAR-1 study was substantially lower than the sensitivity of the NGS-based SiMSen-Seq assay used in our study. Furthermore, by adding the mFAST-Seq5 z-score as a surrogate marker for tumour fraction in plasma to the results of SiMSen-Seq, we could narrow down the need for tissue testing to only 27/93 (29%) patients. In 4 out of these 20 patients with evaluable tissue samples, we identified *PIK3CA* mutation in the tissue. This finding may be of great importance for ctDNA testing in general since the combination of mutation-specific analyses with an untargeted, fast and cost-effective assessment of tumour fractions may help data interpretation. In addition to the relevance of its one-time measurement, longitudinal assessment of tumour fractions using mFAST-Seq5 sequencing harbors important prognostic information [28].

PIK3CA mutations occur on average in ~40% of HR+ patients, with a slightly higher proportion in Luminal A patients [6]. In a recent analysis of metastatic breast cancer tissue samples from the SAPHIR study, the positivity of *PIK3CA* mutated samples was lower in the HR+ HER2– patients, with 28%. However, recent data show a highly heterogenous mutation frequency between distinct ethnicities in the first line setting in the HR+ HER2– subgroup of patients (46% in Chinese vs 22% in US patients) [29].

Preliminary data from the AURORA demonstrated a similar overall range of *PIK3CA* mutations in metastatic breast cancer without providing information on the distribution with intrinsic subtypes [13]. Nevertheless, the activating *PIK3CA* mutations were




Fig. 2 Concordance of *PIK3CA* mutation status between ctDNA and tissue samples in our study and in the SOLAR-1 study. Total numbers of concordant (turquoise) and discordant (red) cases; ctDNA, circulating tumour DNA; wt, wild type; mt, mutant.

associated with diminished sensitivity to chemotherapy and worse OS [9]. Tzanikou et al. have shown that the frequency of detection of *PIK3CA* mutations in liquid biopsy is higher in metastatic patients than in early stages, possibly either due to a higher proportion of mutated patients but also to higher tumour load in peripheral blood [30]. In the recently published analyses of the first 381 matched pairs of primary tumours and metastases from the AURORA study, an enrichment of *PIK3CA* mutations in metastasis was observed [14, 31]. Thus, it is reasonable to expect an increase in *PIK3CA* mutations in progressive or recurrent disease [31]. In our present study, we did not observe this increase since the most recent tissue available was used for the NGS, and discordant results were observed independently, as summarised in Table 4. Yet, it may contribute to a higher proportion of positive patients overall. Another study demonstrated that the proportion of multiple *PIK3CA* mutated tumours increases with metastatic disease, thereby further increasing the importance of continuous testing at the later stage [21]. Our results demonstrate that co-occurrence of mutations (e.g. double mutations) were primarily found in ctDNA, thereby substantiating the assumption that DNA may be more representative of the heterogeneity of the disease. However, with our approach, we can neither technically nor biologically confirm that these mutations occurred on the same allele. Additionally, the VAF of the second mutation was mainly low, except for one patient where the detected VAF of the second mutation not found in tissue was higher in plasma. One possible reason for the discrepancy, in this case, could be that the time difference between the primary tumour and plasma analysis was significantly longer compared to the other samples. Peripheral blood for ctDNA was also partly collected at a later advanced and/or metastatic stage, whereas a significant proportion of the analysed tissues (61%) were from the most recently biopsied single metastatic lesion. Since compound mutations were shown to be associated with improved benefit from *PIK3CA* inhibitors, this fact may also be in favour of the definition of *PIK3CA* mutation status based on the ctDNA analyses [21].

The detection of *PIK3CA* mutations in plasma implies the presence of tumour-specific mutations since germline *PIK3CA* mutations lead to early lethality and are not expected in humans

[32]. Analytical sensitivity of the assays may allow detection of very low amounts of mutated DNA, with a detection limit of 0.25% VAF in our study. However, the clinical relevance of detecting very low VAF is not well understood. There are several possibilities for low VAF, including a low tumour DNA fraction and a low proportion of *PIK3CA* mutated cells in the overall tumour mass. Concerning clonality, data suggest that most *PIK3CA* mutations are clonal or truncal mutations [33], leaving the possibility that low tumour ctDNA burden is the primary cause of low-level detection. In our study, we observed that discordant samples with low-level VAF were mainly retrieved from patients with low ctDNA burden as measured by z-scores. The correlation between *PIK3CA* VAF levels and clinical benefits from the alpelisib and the best time point to start the treatment is currently unknown and remains a topic for future studies.

Our study has several limitations, including the small sample size, a higher percentage of *PIK3CA* mutated cancers, and the focus on *PIK3CA* mutations only. A broader approach providing information on various important genes would be beneficial, but with a potential of decreased sensitivity for *PIK3CA* mutation detection compared to our present analyses. Our future studies will also include a comparison of SiMSen-Seq based detection of *PIK3CA* mutations with NGS technologies analyzing a larger number of genes important for treating metastatic breast cancer, as was recently performed in the PlasmaMATCH study [34]. Finally, for the clinical implementation of our SiMSen-Seq assays additional in-depth technical validation experiments are necessary to assess sensitivity and specificity. These experiments are already ongoing at our institution.

Nevertheless, the PlasmaMATCH study demonstrated groundbreaking clinical progress for selected metastatic breast cancer patients, and we believe that this study paved the way for future large clinical trials. On the other hand, we currently need to establish precise and low-cost technologies to select patients for currently available PI3K treatments, and our SiMSen-Seq based assay performed well in the current study.

In conclusion, within our study, SiMSen-Seq-based detection of *PIK3CA* hotspot mutations from plasma has shown promising results. In combination with an untargeted, mutation-independent approach for detecting ctDNA fractions, we believe that for samples with z-score >3 and no *PIK3CA* hotspot mutations, no additional tissue analyses are needed. Adding ctDNA fractions thus further increased the validity of negative plasma results. SiMSen-Seq alone and this combinatory approach for plasma analyses substantially improved the robustness of liquid biopsy-based *PIK3CA* mutations testing, particularly compared to the FDA-approved companion diagnostic test theascreen® *PIK3CA* used in the SOLAR-1 study. This approach should be validated in a larger cohort of patients and compared with larger panel NGS ctDNA sequencing methods to establish the most promising approach for selecting individual treatments for our patients.

DATA AVAILABILITY

Data are available at request from authors and were partially (SiMSen-Seq and mFAST-SeqS) deposited at the European Genome-phenome Archive (EGA; <http://www.ebi.ac.uk/ega/>) under study accession number EGAS00001004940 as following dataset: *PIK3CA* SiMSen-Seq data of plasma samples.

REFERENCES

- Torre LA, Bray F, Siegel RL, Ferlay J, Lortet-Tieulent J, Jemal A. Global cancer statistics, 2012. *CA Cancer J Clin.* 2015;65:87–108.
- Gennari A, Conte P, Rosso R, Orlandini C, Bruzzi P. Survival of metastatic breast carcinoma patients over a 20-year period: a retrospective analysis based on individual patient data from six consecutive studies. *Cancer.* 2005;104:1742–50.
- Ribnikar D, Volovat SR, Cardoso F. Targeting CDK4/6 pathways and beyond in breast cancer. *Breast.* 2019;43:8–17.

4. Karakas B, Bachman KE, Park BH. Mutation of the PIK3CA oncogene in human cancers. *Br J Cancer*. 2006;94:455–9.
5. Fruman DA, Chiu H, Hopkins BD, Bagrodia S, Cantley LC, Abraham RT. The PI3K pathway in human disease. *Cell*. 2017;170:605–35.
6. Koboldt DC, Fulton RS, McLellan MD, Schmidt H, Kalicki-Verizer J, McMichael JF, et al. Comprehensive molecular portraits of human breast tumours. *Nature*. 2012;490:61–70.
7. Saal LH, Holm K, Maurer M, Memeo L, Su T, Wang X, et al. PIK3CA mutations correlate with hormone receptors, node metastasis, and ERBB2, and are mutually exclusive with PTEN loss in human breast carcinoma. *Cancer Res*. 2005;65:2554–9.
8. Di Leo A, Johnston S, Lee KS, Ciruelos E, Lonning PE, Janni W, et al. Buparlisib plus fulvestrant in postmenopausal women with hormone-receptor-positive, HER2-negative, advanced breast cancer progressing on or after mTOR inhibition (BELLE-3): a randomised, double-blind, placebo-controlled, phase 3 trial. *Lancet Oncol*. 2018;19:87–100.
9. Mosele F, Stefanovska B, Lusque A, Tran Dien A, Garberis I, Droin N, et al. Outcome and molecular landscape of patients with PIK3CA-mutated metastatic breast cancer. *Ann Oncol*. 2020;31:377–86.
10. Moynahan ME, Chen D, He W, Sung P, Samoila A, You D, et al. Correlation between PIK3CA mutations in cell-free DNA and everolimus efficacy in HR(+), HER2(–) advanced breast cancer: results from BOLERO-2. *Br J Cancer*. 2017;116:726–30.
11. Andre F, Ciruelos E, Rubovszky G, Campone M, Loibl S, Rugo HS, et al. Alpelisib for PIK3CA-mutated, hormone receptor-positive advanced breast cancer. *N Engl J Med*. 2019;380:1929–40.
12. Rugo HS, Lerebours F, Ciruelos E, Drullinsky P, Ruiz-Borrego M, Neven P, et al. Alpelisib plus fulvestrant in PIK3CA-mutated, hormone receptor-positive advanced breast cancer after a CDK4/6 inhibitor (BYLieve): one cohort of a phase 2, multicentre, open-label, non-comparative study. *Lancet Oncol*. 2021;22:489–98.
13. Aftimos PG, Antunes De Melo e Oliveira AM, Hilbers F, Venet D, Vingiani A, Nili Gal Yam E, et al. First report of AURORA, the breast international group (BIG) molecular screening initiative for metastatic breast cancer (MBC) patients (pts). *Ann Oncol*. 2019;30:iii48.
14. Aftimos P, Oliveira M, Irrthum A, Fumagalli D, Sotiriou C, Nili Gal-Yam E, et al. Genomic and transcriptomic analyses of breast cancer primaries and matched metastases in AURORA, the Breast International Group (BIG) molecular screening initiative. *Cancer Discov*. 2021; 11:2796–2811.
15. Heitzer E, Ulz P, Belic J, Gutsch S, Quehenberger F, Fischereder K, et al. Tumor-associated copy number changes in the circulation of patients with prostate cancer identified through whole-genome sequencing. *Genome Med*. 2013;5:30.
16. Stahlberg A, Krzyzanowski PM, Egyud M, Filges S, Stein L, Godfrey TE. Simple multiplexed PCR-based barcoding of DNA for ultrasensitive mutation detection by next-generation sequencing. *Nat Protoc*. 2017;12:664–82.
17. Belic J, Koch M, Ulz P, Auer M, Gerhalter T, Mohan S, et al. Rapid identification of plasma DNA samples with increased ctDNA levels by a modified FAST-Seq approach. *Clin Chem*. 2015;61:838–49.
18. Belic J, Graf R, Bauernhofer T, Cherkas Y, Ulz P, Waldispuehl-Geigl J, et al. Genomic alterations in plasma DNA from patients with metastasized prostate cancer receiving abiraterone or enzalutamide. *Int J Cancer*. 2018;143:1236–48.
19. Wang K, Li M, Hakonarson H. ANNOVAR: functional annotation of genetic variants from high-throughput sequencing data. *Nucleic Acids Res*. 2010;38:e164.
20. Cingolani P, Platts A, Wang le L, Coon M, Nguyen T, Wang L, et al. A program for annotating and predicting the effects of single nucleotide polymorphisms, SnpEff: SNPs in the genome of *Drosophila melanogaster* strain w1118; iso-2; iso-3. *Fly (Austin)*. 2012;6:80–92.
21. Vasan N, Razavi P, Johnson JL, Shao H, Shah H, Antoine A, et al. Double PIK3CA mutations in cis increase oncogenicity and sensitivity to PI3Kalpha inhibitors. *Science*. 2019;366:714–23.
22. Suppan C, Brcic I, Tiran V, Mueller HD, Posch F, Auer M, et al. Untargeted assessment of tumor fractions in plasma for monitoring and prognostication from metastatic breast cancer patients undergoing systemic treatment. *Cancers (Basel)*. 2019;11:1171.
23. Baselga J, Im SA, Iwata H, Cortes J, De Laurentiis M, Jiang Z, et al. Buparlisib plus fulvestrant versus placebo plus fulvestrant in postmenopausal, hormone receptor-positive, HER2-negative, advanced breast cancer (BELLE-2): a randomised, double-blind, placebo-controlled, phase 3 trial. *Lancet Oncol*. 2017;18:904–16.
24. Baselga J, Dent SF, Cortés J, Im Y-H, Diéras V, Harbeck N, et al. Phase III study of taselelisib (GDC-0032) + fulvestrant (FULV) v FULV in patients (pts) with estrogen receptor (ER)-positive, PIK3CA-mutant (MUT), locally advanced or metastatic breast cancer (MBC): primary analysis from SANDPIPER. *J Clin Oncol*. 2018;36:LBA1006–LBA1006.
25. Baselga J, Sellami D, El-Hashimy M, Dharan B, Wang A, Scheuer N, et al. Abstract A050: PIK3CA mutation status in tumor tissue and ctDNA as a biomarker for PFS in patients with HR+, HER2- ABC treated with buparlisib or placebo plus fulvestrant: results from the BELLE-2 and BELLE-3 randomized studies. *Mol Cancer Ther*. 2018;17:A050–A050.
26. Martinez-Saez O, Chic N, Pascual T, Adamo B, Vidal M, Gonzalez-Farre B, et al. Frequency and spectrum of PIK3CA somatic mutations in breast cancer. *Breast Cancer Res*. 2020;22:45.
27. Juric D, Ciruelos E, Rubovszky G, Campone M, Loibl S, Rugo H, et al. Abstract GS3-08: Alpelisib + fulvestrant for advanced breast cancer: subgroup analyses from the phase III SOLAR-1 trial. *Cancer Res*. 2019;79:GS3-08-GS3.
28. Dandachi N, Posch F, Graf R, Suppan C, Klocker EV, Müller HD, et al. Longitudinal tumor fraction trajectories predict risk of progression in metastatic HR + breast cancer patients undergoing CDK4/6 treatment. *Mol Oncol*. 2020;15:2390–2400.
29. Li H, Davis AA, Liu X-R, Xie F, Gui X-Y, Chen Y-F, et al. Cell-free DNA comparative analysis of hormone receptor-positive, first-line metastatic breast cancer genomic landscape in the United States and China. *J Clin Oncol*. 2020;38:1059.
30. Tzanikou E, Markou A, Politaki E, Koutsopoulos S, Psyrri A, Mavroudis D, et al. PIK3CA hotspot mutations in circulating tumor cells and paired circulating tumor DNA in breast cancer: a direct comparison study. *Mol Oncol*. 2019;13:2515–30.
31. Markou A, Farkona S, Schiza C, Efstathiou T, Kounelis S, Malamos N, et al. PIK3CA mutational status in circulating tumor cells can change during disease recurrence or progression in patients with breast cancer. *Clin Cancer Res*. 2014;20:5823–34.
32. Bi L, Okabe I, Bernard DJ, Wynshaw-Boris A, Nussbaum RL. Proliferative defect and embryonic lethality in mice homozygous for a deletion in the p110alpha subunit of phosphoinositide 3-kinase. *J Biol Chem*. 1999;274:10963–8.
33. O'Leary B, Cutts RJ, Liu Y, Hrebien S, Huang X, Fenwick K, et al. The genetic landscape and clonal evolution of breast cancer resistance to palbociclib plus fulvestrant in the PALOMA-3 trial. *Cancer Discov*. 2018;8:1390–403.
34. Turner NC, Kingston B, Kilburn LS, Kernaghan S, Wardley AM, Macpherson IR, et al. Circulating tumour DNA analysis to direct therapy in advanced breast cancer (plasmaMATCH): a multicentre, multicohort, phase 2a, platform trial. *Lancet Oncol*. 2020;21:1296–308.

ACKNOWLEDGEMENTS

We thank Klara Balic for English proofreading. In addition, we want to acknowledge that part of this manuscript was published as a doctoral thesis.

AUTHOR CONTRIBUTIONS

CS, AT, ND and MB designed the study. RG and KK performed sequencing experiments, SJ and PR performed pathological evaluation, CS, EVK, AT, FP and MB contributed to patient care and data acquisition. ND and FP performed the statistical analysis, CS, AT, ND, MB, EH and PJJ performed interpretation of results. CS, ND, MB and FP drafted and revised the manuscript. ND, EH and MB supervised the study. All authors reviewed the manuscript and approved its final version.

FUNDING

This work was supported by the Christian Doppler Research Fund for Liquid Biopsies for Early Detection of Cancer led by EH supported by the Federal Ministry of Digital and Economic Affairs and by Novartis and Pfizer.

ETHICS APPROVAL AND CONSENT TO PARTICIPATE

The study was approved by the ethics committee of the Medical University Graz (ethical approval number 21–227 ex 09/10), and written informed consent was obtained from all patients.

COMPETING INTERESTS

CS received travel expenses, consulting fees and honoraria from Amgen, Astra Zeneca, Eli Lilly, Novartis, Pfizer, Roche and Samsung, SJ received honoraria from Novartis and Roche, EVK has received travel expenses and honoraria from Astra Zeneca and Pierre-Fabre, RB had a consulting role, received honoraria, research funding and travel expenses from Astra Zeneca, Daiichi, Eisai, Eli Lilly, MSD, Novartis, Pfizer, Pierre-Fabre, Puma, Roche, AT received honoraria and research funding from Astra Zeneca, BMS, Merck and Roche. KK received honoraria and advisory from Astra Zeneca, Thermo Fischer, Novartis and BMS. PR received consulting fees, honoraria and travel expenses from Diaceutics, Novartis, and Roche, PJJ has had a consulting or advisory role, received honoraria, research funding and/or travel/accommodation expenses from Abbvie, Bayer, Boehringer, Novartis, Pfizer, Servier, Roche, BMS and Celgene, HG is a full time/part-time employee by Novartis, EH Servier, has received funding from Servier, Freenome, CA and PreAnalytiX and had an advisory role for Roche. MB has received honoraria, consulting fees research funds and/or travel expenses from Amgen, Astra Zeneca, Daiichi, Eli Lilly, MSD, Novartis, Pierre-Fabre, Pfizer, Roche and Samsung.

ADDITIONAL INFORMATION

Supplementary information The online version contains supplementary material available at <https://doi.org/10.1038/s41416-021-01601-9>.

Correspondence and requests for materials should be addressed to Nadia Dandachi or Marija Balic.

Reprints and permission information is available at <http://www.nature.com/reprints>

Publisher's note Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional affiliations.