

1 Testing Frequency Matters | An Evaluation of the Diagnostic Performance of a SARS-
2 CoV-2 Rapid Antigen Test in United States Correctional Facilities

3 Margaret L. Lind, MPH, Ph.D.¹

4 Olivia L. Schultes, B.A.²

5 Alexander J. Robertson, B.S.¹

6 Amy J. Houde, M.S.W, L.C.S.W³

7 Derek A.T. Cummings, Ph.D.^{4,5}

8 Albert I. Ko, M.D.^{2,6*}

9 Byron S. Kennedy, M.D., Ph.D.^{3*}

10 Robert P. Richeson, D.C.^{3*}

11 *These authors contributed equally

12

13 1. Department of Epidemiology of Microbial Diseases, Yale School of Public Health,
14 New Haven, CT, USA

15 2. Department of Epidemiology, University of Washington's School of Public Health,
16 Seattle, WA, USA

17 3. Connecticut Department of Correction, Wethersfield, CT, USA

18 4. Department of Biology, University of Florida, Gainesville, FL, USA

19 5. Emerging Pathogens Institute, University of Florida, Gainesville, FL, USA

20 6. Instituto Gonçalo Moniz, Fundação Oswaldo Cruz, Salvador, BA, Brazil

21 **Corresponding Author:**

22 Margaret L. Lind, PhD

23 Margaret.Lind@yale.edu

1 Yale University School of Public Health
2 60 College Street
3 New Haven, CT 06510
4 United States of America
5

6 **Post Publication Corresponding Author:**

7 Robert P. Richeson, DC
8 Robert.Richeson@ct.gov
9 Connecticut Department of Corrections
10 24 Wolcott Hill Rd.
11 Wethersfield, CT 06109
12 United States of America
13
14
15

16 **Short Title:** BinaxNOW Accuracy in Jails and Prisons

17 **Summary:** Through the evaluation of the BinaxNOW rapid antigen test under varying
18 collection frequency strategies, we provide evidence of the utility of serial rapid antigen
19 test collection within congregate facility settings for outbreak investigation, screening,
20 and when rapid detection is required.

Abstract

Background: The CDC recommends serial rapid antigen assay collection within congregate facilities. Though modeling and observational studies from communities and long-term care facilities have shown serial collection provides adequate sensitivity and specificity, the accuracy within correctional facilities remains unknown.

Methods: Using Connecticut Department of Corrections (DOC) data from November 21st 2020 to June 15th 2021, we estimated the accuracy of a rapid assay, BinaxNOW, under three collection strategies, single test collection and serial collection of two and three tests separated by 1-4 days. The sensitivity and specificity of the first (including single), second, and third serially collected BinaxNOW tests were estimated relative to RT-PCRs collected within one day of the BinaxNOW test. The accuracy metrics of the testing strategies were then estimated as the sum (sensitivity) and product (specificity) of tests in each strategy.

Results: Of the 13,112 residents who contributed ≥ 1 BinaxNOW test during the study period, 3,825 contributed ≥ 1 RT-PCR paired BinaxNOW test. In relation to RT-PCR, the three-rapid antigen test strategy had a sensitivity of 95.9% (95% confidence intervals (CI): 93.6-97.5%) and specificity of 98.3% (CI: 96.7-99.1%). The sensitivity of the two- and one-rapid antigen test strategies were 88.8% and 66.8%, respectively, and the specificities were 98.5% and 99.4%, respectively. The sensitivity was higher among symptomatic residents and when RT-PCRs were collected before BinaxNOW tests.

Conclusions: We found serial antigen test collection resulted in high diagnostic accuracy. These findings support serial collection for outbreak investigation, screening, and when rapid detection is required (such as intakes or transfers).

Keywords: Diagnostic Accuracy, COVID-19, Rapid Antigen Test, Correctional Facility

Introduction

Within the United States, state and federal run correctional facilities have experienced high transmission rates of SARS-CoV-2 and remain high-risk settings for COVID-19.[1,2] In fact, data from September through November 2020 show that residents of Federal Bureau of Prisons were 4.7 times more likely to become infected with SARS-CoV-2 and 2.6 times more likely to die from COVID-19 than general US residents.[2] Despite the development of COVID-19 vaccines and vaccination programs for incarcerated populations, vaccine coverage remains below that needed for population level protection.[3–9] Rapid and accurate SAR-CoV-2 testing will therefore remain a key component of infection prevention within correctional facilities.

In August 2020, the Food and Drug Administration (FDA) issued an emergency use authorization of Abbott's BinaxNOW, a COVID-19 rapid antigen test.[10,11] Compared with reverse transcription polymerase chain reaction (RT-PCR), the rapid turnaround time and low cost of rapid antigen tests makes them a cost-effective strategy for congregate settings, where transmission risk is high and implementation of serial mass screening is feasible.[12] Unfortunately, prior studies from community, educational, and long-term care facility settings found the sensitivity of rapid antigen tests to be poor to moderate (53-77%) and to be lower among asymptomatic individuals and individuals early or late in their course of infection.[13–18] These findings, thus, call into question the use of rapid antigen tests as single point of care tests.

Single test collection is not, however, the intended testing strategy for rapid antigen tests outside of healthcare settings.[11,19,20] Instead, both the manufacturer and the FDA advice serial collection of at least two tests.[11,19] In alignment with the intended use, the Centers for Disease Control and Prevention (CDC) recommends

serial collection of rapid antigen testing within congregate settings for screening and during outbreaks regardless of symptom presentation.[21] This guidance is supported by modeled evidence that outbreak control depends largely on frequency and speed of testing and observational data showing that serial testing improves the sensitivity of rapid antigen tests in both nursing home and community settings.[12,14,22,23] However, the value of serial testing within correctional facilities remains unknown.

Herein, we present findings of a study that evaluated the accuracy of serial BinaxNOW rapid antigen testing during a mass screening and testing program implemented by the Connecticut Department of Corrections (DOC) during the COVID-19 pandemic wave in late 2020 and early 2021. We specifically compared the accuracy of the rapid antigen test relative to RT-PCR under three collection strategies: serial testing of up to three negative rapid antigen tests, serial testing of up to two negative rapid antigen tests, and rapid antigen collection in isolation.

Methods

Setting and Specimen Collection

On November 21, 2020, the DOC initiated the implementation of the Abbott BinaxNOW rapid antigen test (headquarters: Chicago, IL) collection into their SARS-CoV-2 testing program for symptomatic testing, contact tracing, testing of residents undergoing admission to the correctional facilities and inter-facility transfer and mass voluntary asymptomatic screening. For each instance of rapid antigen use, the DOC guidelines recommend the serial collection of up to three negative rapid antigen tests taken on day one, four, and seven (where day one is the day of exposure, if exposed). Due to concerns around the sensitivity of the rapid antigen test, phased implementation

of rapid antigen testing with confirmatory RT-PCR was performed. While undergoing serial test collection, residents of DOC facilities were placed in quarantine or isolation.

Trained medical staff collected anterior nasal and nasopharyngeal swab specimens for rapid antigen and RT-PCR testing. Quest Diagnostic facilities performed SARS-CoV-2 RT-PCR testing of nasopharyngeal swab specimens using the SARS-CoV-2 RNA (COVID-19), Qualitative NAAT assay and defined positive tests as a cycle threshold value of less than 40.[24] At the time rapid antigen test collection was implemented, the DOC oversaw 17 facilities with a resident census of around 10,000 residents (9,945 residents on July 1st, 2020).[25,26] The Yale University Institutional Review Board classified this study as public health surveillance.

Data Collection and Cohort Development

Using resident and testing data queried from internally maintained DOC databases, we retrospectively identified all rapid antigen (BinaxNOW) and RT-PCR testing records from November 21st, 2020 to June 15th, 2021. We included tests collected among residents with prior documented SARS-CoV-2 infections but excluded duplicated records (same day, assay, and assay results) and reporting errors (negative recorded on the same day as a positive).

Following the DOC recommendations for serial rapid antigen test collection, we identified series of rapid antigen tests (rapid antigen tests collected within one and four days of each other or tests collected in the absence of any test in the prior or following four days; Figure 1.A). We then paired the first, second and third rapid antigen tests of the identified series with RT-PCR tests collected between one day prior to or following the rapid antigen test (Figure 1.B). In the event of multiple RT-PCR matches per rapid

antigen test, ordered preference was given to positive RT-PCRs (Figure 1.B; Option 1), RT-PCRs collected on the same day as the rapid antigen test (Figure 1.B; Option 2), and RT-PCRs collected before the rapid antigen test (Figure 1.B; Option 3). We defined symptoms as the presence or absence of COVID-19 related symptoms reported at the time of rapid antigen testing.[27]

Statistical Approach: Diagnostic Accuracy

The characteristics of residents with at least one RT-PCR paired rapid antigen test were summarized by presence of RT-PCR positive SARS-CoV-2 events using counts and percentages for categorical factors and means and standard errors for continuous factors. We estimated the sensitivity and specificity for each testing strategy using the following equations[28]:

One-Test Sensitivity Equation: Weighted Average($Test1_{sen}$, $Test2_{sen}$, $Test3_{sen}$)

*Two-Test Sensitivity Equation: $Test1_{sen} + (1 - Test1_{sen}) * Test2_{sen}$*

*Three-Test Sensitivity Equation: $Test1_{sen} + (1 - Test1_{sen}) * Test2_{sen} + (1 - (Test1_{sen} + Test2_{sen})) * Test3_{sen}$*

One-Test Specificity Equation: Weighted Average($Test1_{spec}$, $Test2_{spec}$, $Test3_{spec}$)

*Two-Test Specificity Equation: $Test1_{spec} * Test2_{spec}$*

*Three-Test Specificity Equation: $Test1_{spec} * Test2_{spec} * Test3_{spec}$*

Where $Test1_{sen/spec}$ was the sensitivity or specificity of RT-PCR paired first rapid antigen tests, $Test2_{sen/spec}$ was the sensitivity or specificity of RT-PCR paired second rapid antigen tests, and so on (eFigure1). The sensitivity and specificity of RT-PCR paired first, second, and third rapid antigen tests were estimated using generalized estimating equations (GEE) with robust standard errors and a logit link. We propagated the uncertainty through the serial testing equations using posterior simulation of 1000 random draws of the GEE estimate.[29,30] Additionally, we estimated the diagnostic accuracy stratified by the following *a priori* selected factors: age, symptom presence,

and test order. We tested for additive differences in the diagnostic accuracy of the stratified samples by subtracting the draws of each sample from a selected reference category and defined significance as confidence intervals above or below the null (zero).

For the PPV and NPV, we simulated 1000 average daily prevalence estimates using a Poisson regression with an outcome of positive SARS-CoV-2 test (either rapid antigen or RT-PCR test). To reduce the risk of including multiple positive tests from the same testing event, we excluded positive events within 5 days of each other. With the estimated prevalence, sensitivity, and specificity, we estimated the PPV and NPV for each of the sensitivity, specificity, and prevalence drawn estimates (n =1000).[31] For each accuracy metric, we calculated 95% confidence intervals as the 2.5 and 97.5 percentiles.[30] All analyses were conducted in R 4.1.0 using the geepack and multcomp packages.[32,33]

Sensitivity Analyses

To test the robustness of our findings, we performed sensitivity analyses where we selected first, second and third tests at random instead of the observed order, limited to rapid antigen test series where tests were collected exactly three days apart, and invoked different selection approaches in the event of multiple RT-PCRs linked to a rapid antigen test. Finally, in a post-hock analysis, we performed an age stratified analysis among residents whose RT-PCR and rapid antigen tests were collected on the same day.

Results

Between November 21st 2020 and June 15th 2021, 128,986 RT-PCR and 30,112 rapid antigen tests were collected among 17,669 DOC residents (eFigure2). Of these residents, 3,825 contributed at least one RT-PCR paired rapid antigen test. The majority of residents were male (76.5%), and the most frequently observed race was Black (40.9%). Residents who experienced a paired RT-PCR positive SARS-CoV-2 event were demographically similar to residents who did not experience a paired RT-PCR positive SARS-CoV-2 event (Table 1).

Among the 18,186 identified rapid antigen test series, 4,919 consisted of two or more tests and 2,911 consisted of three or more test. A total of 3,844 first, 677 second, and 314 third rapid antigen tests from identified series were collected within one day of a RT-PCR and considered RT-PCR paired (Figure 2). The test collection order and proportion of tests collected among symptomatic and asymptomatic residents were similar between RT-PCR paired first, second, and third rapid antigen tests (eTable1).

Most RT-PCR paired rapid antigen tests were negative (first test: 91.5%, second test: 88.2%, third test: 89.8%). Of the negative rapid antigen tests, most were paired with negative RT-PCRs (first test: 91.5%; second test: 89.2%; third test: 89.8%; eTable2). Relative to RT-PCR, the sensitivity of the first, second and third rapid antigen tests were 67.7% (95% confidence interval (CI): 62.9-72.0%), 65.4% (CI: 55.9-74.1%) and 62.7% (CI: 49.1-75.2%), respectively. This resulted in sensitivities of 66.8% (CI: 62.8-70.6%), 88.8% (CI: 85.1-91.9%) and 95.9% (CI: 93.6-97.5%) for the one-, two-, and three-rapid antigen test collection strategies, respectively. The specificity of the first, second, and third rapid antigen tests were above 99%. The specificity of two and three serially collected rapid antigen tests was 98.5% (CI: 97.3-99.1%) and 98.3% (CI: 96.7-

99.1%), respectively (Table2). The PPV, based on an observed prevalence of 57 (CI: 53-62) cases per 100,000 residents, was highest for the one-rapid antigen test collection strategy (76.7%; CI: 69-83.1%) and lowest for the three-rapid antigen test collection strategy (62.9%; CI: 46.6-75.5%). The NPV for each rapid antigen test collection strategy was high (one-rapid antigen test: 99.0%, two-rapid antigen tests: 99.0%, three-rapid antigen tests: 98.9%, eTable3).

The sensitivity of each collection strategy was significantly higher among symptomatic residents than among asymptomatic residents (difference for one-rapid antigen test: 30.8% (CI: 23.5-37.4%), two-rapid antigen tests: 16.2% (CI: 11.0-22.0%), three-rapid antigen tests: 8.6% (CI: 5.2-13.0%). The sensitivity was higher when rapid antigen tests were collected on the same day as RT-PCRs than when rapid antigen tests were collected on the day before RT-PCRs (difference for one-rapid antigen test: 28.3% (18.3-38.0%), two-rapid antigen tests: 25.4% (CI: 12.1-38.2%), three-rapid antigen tests: 27.0% (CI: 11.0-41.2%)). The sensitivity of the two-rapid antigen test collection strategy was significantly higher for residents less than 37 years of age than for residents more than 37 years of age (difference: 6.8% (CI: 0.2-13.9%), Table3).

The specificity for each rapid antigen test collection strategy was significantly higher among asymptomatic residents than among symptomatic residents (difference for one-rapid antigen test: 4.2% (0.9-11.2%), two-rapid antigen tests: 16.1% (CI: 3.5-42.9%), three-rapid antigen tests: no sample). Interestingly, the specificity of RT-PCR paired second rapid antigen tests collected among symptomatic residents was low (85.8%; CI: 58.4-96.4%) but the specificity of RT-PCR paired first rapid antigen tests

was high (97.5; CI: 90.5-99.3%). The specificities did not vary significantly by age or timing of tests (Table3).

Sensitivity Analyses

We performed multiple sensitivity analyses to test the robustness of our findings to varying rapid antigen test collection orders, rapid antigen test collection windows, and RT-PCR pair selection. The sensitivity of the one-rapid antigen test, two-rapid antigen test, and three-rapid antigen test collection strategies ranged from 57.2-74.6%, 84.7-97.6%, and 92.8-98.7%, respectively (eTables4-6/eFigure3). The highest and lowest sensitivities were observed when we selected the first, second, and third tests at random (eFigure3). As with the primary analysis, we observed high specificities for each testing strategy (one-rapid antigen test: 98.9-99.8%, two-rapid antigen test: 98.0-99.5%, three-rapid antigen test: 96.8-99.3%; eTables4-6/eFigure3). In a post-hock analysis, we found that the accuracy of the rapid antigen test did not differ significantly by age group among residents whose RT-PCR and rapid antigen tests were collected on the same day (eTable7).

Discussion

We evaluated the predictive accuracy of the rapid COVID-19 antigen test, BinaxNOW, in relation to RT-PCR under three different collection strategies (a single test in isolation and two and three serial tests separated by 1–4-day intervals) among residents of Connecticut state prisons and jails between November 21st, 2020 and June 15th, 2021. The three-rapid antigen test collection strategy is currently recommended by Connecticut DOC.[34] In alignment with diagnostic accuracy estimates from other congregate settings, we found that rapid antigen tests had a moderate sensitivity when

collected in isolation, but the sensitivity increased significantly when rapid antigen tests were collected in pairs and triplets.[13,14,22]

In relation to RT-PCR, we found the current DOC rapid antigen test collection strategy, three-rapid antigen test, had a high sensitivity (96%) and specificity (98%). While the specificity of the two and one-rapid antigen test collection strategies were higher (99% and 99%, respectively), the sensitivity for these less intensive collection strategies were significantly lower (89% and 67%, respectively). These findings suggest that, among 100 residents infected with SARS-CoV-2, 96 would be captured by the three-rapid antigen test collection strategy. Compared with this strategy, the two-rapid antigen test and the one-rapid antigen test collection strategies would miss an additional seven and 29 infected residents, respectively. Conversely, the three-rapid antigen test collection strategy would only misdiagnose two out of 100 uninfected residents. Though this is one more than the one-rapid antigen test collection strategy, the cost of false negatives (missed infections) far outweighs the cost of false positives (excess isolation) under scenarios of highly transmissible infectious diseases, such as SARS-CoV-2.[35]

In alignment with prior studies, we found the sensitivity of rapid antigen testing to be significantly higher later in the course of infection (rapid antigen tests collected on the same day or after RT-PCR) and among residents exhibiting COVID-19 symptoms.[36,37] However, contrary to previous findings, we found the specificity was lower among symptomatic residents than among asymptomatic residents.[37–39] While surprising, prior studies found the specificity of the Access Bio CareState test was significantly and the BinaxNOW was non-significantly lower among residents who had experienced symptoms for a prolonged period of time (>7 days).[39,40] Though

symptom duration data were unavailable, the observed difference was driven by the specificity of second, serially collected antigen tests. This suggests that, like other rapid antigen tests, the specificity of the BinaxNOW may decline with increased symptom duration. Despite the observed differences, the specificity among symptomatic residents remained high for each testing strategy with available data (one-rapid antigen test: 95.3%, two-rapid antigen tests: 82.9%). Thus, this difference in specificity does not invalidate serial tests collection among symptomatic residents of correctional facilities.

Taken collectively, our findings support the use of serial rapid antigen testing under scenarios when PCR turnaround time is long, rapid detection is required, or when isolating/quarantining all exposed residents is unfeasible. Such scenarios include times of unknown exposure (intake and transfers), contact tracing, asymptomatic and symptomatic screening, and during outbreak investigations. These recommendations stem from the rapid turnaround time and low cost of a single test, the collectively high diagnostic accuracy of serial collection, and the ability of serial collection strategies to detect events under continuous exposure scenarios. Though our findings speak predominately to the value of serial collection resulting from a single exposure, serial collection provides additional benefit through capture of infections from exposures that occurred immediately prior to or following the first collected test.[12] The combination of these benefits, thus, may result in more rapid isolation of infected residents, and in turn, a reduction in facility wide transmission.

Limitations and Strengths

Our study was subject to limitations typical of retrospective diagnostic validation analyses. First, race/ethnicity was missing for a large portion of the population and we

were unable to test for differences in test accuracy by race. However, it is unlikely that race would impact the diagnostic accuracy of the rapid antigen test. Second, we relied on a reference outcome of RT-PCR positivity, which is an imperfect indicator of infection. Third, we did not have access to cycle threshold values and were unable to test the impact of viral load on rapid antigen test performance. Fourth, our accuracy estimates relied on collected tests that may be biased towards department specific testing practices. However, we observed similar results when we selected the first, second, and third test of each series at random (eFigure3; eTable5). Finally, our study was conducted prior to the large Delta and Omicron waves of 2021/2022. While the diagnostic accuracy of rapid antigen tests likely varies between the different variants, we believe the benefits of serial collection will hold.

Our study had several strengths including our large sample of paired assays collected among a diverse population of Connecticut State correctional facility residents. This large sample allowed us to estimate and compare the diagnostic accuracy of three different collection strategies, including the three-test collection strategy employed by the Connecticut DOC. Our large sample also allowed us to examine characteristics associated with the accuracy of the rapid antigen test within correctional facilities settings and speak to the use of different collection strategies based on these characteristics. Additionally, through the inclusion of numerous sensitivity analyses, we were able to show that our findings were not the result of the data cleaning or modeling assumptions we employed within this analysis. Finally, we were able to include all unique rapid antigen test sets and account for within person correlation in our uncertainty intervals using generalized estimating equations.

Conclusion

Compared with singularly collected tests, we found that serial collection of BinaxNOW rapid antigen tests resulted in meaningfully higher sensitivities and comparably high specificities among residents in state correctional facilities. We found this held for both asymptomatic and symptomatic residents and regardless of rapid antigen test collection time relative to RT-PCR collection. These findings speak to the utility of serially collected rapid antigen tests within correctional facilities for asymptomatic and symptomatic screening, contact tracing, and during outbreak investigations. If employed under such scenarios, rapid antigen testing may result in faster isolation of infected residents and reduce transmission within facilities.

NOTES

Author Contributions: Margaret L. Lind had full access to all the data in the study and takes responsibility for the integrity of the data and the accuracy of the data analysis.

Concept and design: Ko, Kennedy, Richeson

Acquisition, analysis, or interpretation of data: Schultes, Lind

Drafting of the manuscript: Lind, Schultes, Robertson, Ko, Cummings

Critical revision of the manuscript for important intellectual content: All authors.

Statistical analysis: Lind

Administrative, technical, or material support:

Supervision: Ko, Kennedy, Richeson, Cummings

Additional Contributions: We thank the Connecticut Department of Corrections and the healthcare workers in each correctional facility for administrative support and

administering RT-PCR and rapid antigen tests. None of the individuals received compensation for their contribution.

Data Use Agreement: The data used in this study belongs to the Connecticut Department of Corrections. Qualified researchers may submit a data share request for de-identified patient level data by contacting the corresponding author with a detailed description of the research question.

Funding: Connecticut Department of Public Health Emerging Infections Program (EIP): COVID-19 contract (DPH log # 2021-0071-3)

Conflict of Interest Disclosures: A.I.K serves as an expert panel member for Reckitt Global Hygiene Institute, scientific advisory board member for Revelar Biotherapeutics and a consultant for Tata Medical and Diagnostics and Regeneron Pharmaceuticals, an unpaid leadership or fiduciary role on the Board of Directors for the American Society of Tropical Medicine and Hygiene and has received grants from Merck, Regeneron Pharmaceuticals and Tata Medical and Diagnostics for research related to COVID-19, all of which are outside the scope of the submitted work. D.A.T.C. reports contracts or grants from Merck for research unrelated to this manuscript. Other authors declare no conflict of interest.

References

1. Hagan LM, McCormick DW, Lee C, et al. Outbreak of SARS-CoV-2 B.1.617.2 (Delta) Variant Infections Among Incarcerated Persons in a Federal Prison - Texas, July-August 2021. *MMWR Morb Mortal Wkly Rep* **2021**; 70:1349–1354.
2. Toblin RL, Hagan LM. COVID-19 Case and Mortality Rates in the Federal Bureau of Prisons. *Am J Prev Med* **2021**; 61:120–123.
3. Chin ET, Leidner D, Ryckman T, et al. Covid-19 Vaccine Acceptance in California State Prisons. *New England Journal of Medicine* **2021**; 385:374–376.
4. Goldberg Y, Mandel M, Bar-On YM, et al. Waning Immunity after the BNT162b2 Vaccine in Israel. *New England Journal of Medicine* **2021**; Available at: <https://www.nejm.org/doi/10.1056/NEJMoa2114228>. Accessed 29 October 2021.
5. Garcia-Beltran WF, Lam EC, St. Denis K, et al. Multiple SARS-CoV-2 variants escape neutralization by vaccine-induced humoral immunity. *medRxiv* **2021**; :2021.02.14.21251704.
6. Tartof SY. Effectiveness of mRNA BNT162b2 COVID-19 vaccine up to 6 months in a large integrated health system in the USA: a retrospective cohort study. :10.
7. Fowlkes A. Effectiveness of COVID-19 Vaccines in Preventing SARS-CoV-2 Infection Among Frontline Workers Before and During B.1.617.2 (Delta) Variant Predominance — Eight U.S. Locations, December 2020–August 2021. *MMWR Morb Mortal Wkly Rep* **2021**; 70. Available at: <https://www.cdc.gov/mmwr/volumes/70/wr/mm7034e4.htm>. Accessed 7 October 2021.
8. Kustin T, Harel N, Finkel U, et al. Evidence for increased breakthrough rates of SARS-CoV-2 variants of concern in BNT162b2-mRNA-vaccinated individuals. *Nat Med* **2021**; 27:1379–1384.
9. Hagan LM, Dusseau C, Crockett M, Rodriguez T, Long MJ. COVID-19 vaccination in the Federal Bureau of Prisons, December 2020—April 2021. *Vaccine* **2021**; 39:5883–5890.
10. FDA Media. COVID-19 Update: FDA Authorizes First Diagnostic Test Where Results Can Be Read Directly From Testing Card. FDA, 2020. Available at: <https://www.fda.gov/news-events/press-announcements/covid-19-update-fda-authorizes-first-diagnostic-test-where-results-can-be-read-directly-testing-card>. Accessed 12 August 2021.
11. Abbott. Abbott's BinaxNOW™ Rapid Antigen Self Test Receives FDA Emergency Use Authorization for Asymptomatic, Over-the-Counter, Non-Prescription, Multi-Test Use. Available at: <https://abbott.mediaroom.com/2021-03-31-Abbotts-BinaxNOW-TM-Rapid-Antigen-Self-Test-Receives-FDA-Emergency-Use-Authorization-for-Asymptomatic-Over-the-Counter-Non-Prescription-Multi-Test-Use>. Accessed 8 November 2021.

12. Larremore DB, Wilder B, Lester E, et al. Test sensitivity is secondary to frequency and turnaround time for COVID-19 screening. *Sci Adv* **2021**; 7:eabd5393.
13. Okoye NC, Barker AP, Curtis K, et al. Performance Characteristics of BinaxNOW COVID-19 Antigen Card for Screening Asymptomatic Individuals in a University Setting. *J Clin Microbiol* **2021**; 59:e03282-20.
14. McKay SL, Tobolowsky FA, Moritz ED, et al. Performance Evaluation of Serial SARS-CoV-2 Rapid Antigen Testing During a Nursing Home Outbreak. *Ann Intern Med* **2021**; 174:945–951.
15. Sood N, Shetgiri R, Rodriguez A, et al. Evaluation of the Abbott BinaxNOW rapid antigen test for SARS-CoV-2 infection in children: Implications for screening in a school setting. *PLoS One* **2021**; 16:e0249710.
16. Dinnes J, Deeks JJ, Berhane S, et al. Rapid, point-of-care antigen and molecular-based tests for diagnosis of SARS-CoV-2 infection. *Cochrane Database of Systematic Reviews* **2021**; Available at: <https://www.cochranelibrary.com/cdsr/doi/10.1002/14651858.CD013705.pub2/full>. Accessed 16 August 2021.
17. Brümmer LE, Katzenschlager S, Gaeddert M, et al. Accuracy of novel antigen rapid diagnostics for SARS-CoV-2: A living systematic review and meta-analysis. *PLoS Med* **2021**; 18:e1003735.
18. Pilarowski G, Lebel P, Sunshine S, et al. Performance Characteristics of a Rapid Severe Acute Respiratory Syndrome Coronavirus 2 Antigen Detection Assay at a Public Plaza Testing Site in San Francisco. *J Infect Dis* **2021**; 223:1139–1144.
19. FDA. Fact Sheet For Individuals | Abbott Diagnostics | BinaxNOW COVID-19 Ag Card. 2022; Available at: <chrome-extension://efaidnbmnnnibpcajpcglclefindmkaj/viewer.html?pdfurl=https%3A%2F%2Fwww.fda.gov%2Fmedia%2F141568%2Fdownload&clen=274942&chunk=true>. Accessed 18 February 2022.
20. FDA. Fact Sheet For Healthcare Providers | Abbott Diagnostics | BinaxNOW COVID-19 Ag Card. 2022; Available at: <chrome-extension://efaidnbmnnnibpcajpcglclefindmkaj/viewer.html?pdfurl=https%3A%2F%2Fwww.fda.gov%2Fmedia%2F141568%2Fdownload&clen=274942&chunk=true>. Accessed 18 February 2022.
21. CDC. Interim Guidance for Antigen Testing for SARS-CoV-2. 2022. Available at: <https://www.cdc.gov/coronavirus/2019-ncov/lab/resources/antigen-tests-guidelines.html>. Accessed 18 February 2022.
22. Smith RL, Gibson LL, Martinez PP, et al. Longitudinal Assessment of Diagnostic Test Performance Over the Course of Acute SARS-CoV-2 Infection. *The Journal of Infectious*

- Diseases **2021**; Available at: <https://doi.org/10.1093/infdis/jiab337>. Accessed 12 August 2021.
23. Shah MM, Salvatore PP, Ford L, et al. Performance of Repeat BinaxNOW Severe Acute Respiratory Syndrome Coronavirus 2 Antigen Testing in a Community Setting, Wisconsin, November 2020–December 2020. *Clinical Infectious Diseases* **2021**; 73:S54–S57.
 24. Quest Diagnostics. SARS-CoV-2 RNA, Qualitative Real-Time RT-PCR (Test Code 39433). Available at: <https://www.fda.gov/media/136231/download>. Accessed 29 October 2021.
 25. Connecticut State Department of Correction. Total Supervised Population Count. Available at: <https://portal.ct.gov/DOC/Report/Total-Supervised-Population-Count>. Accessed 18 February 2022.
 26. Connecticut State Department of Correction. Facilities. Available at: <https://portal.ct.gov/DOC/Miscellaneous/Facilities>. Accessed 18 February 2022.
 27. CDC. Coronavirus Disease 2019 (COVID-19) – Symptoms. 2021. Available at: <https://www.cdc.gov/coronavirus/2019-ncov/symptoms-testing/symptoms.html>. Accessed 4 August 2021.
 28. Weinstein S, Obuchowski NA, Lieber ML. Clinical Evaluation of Diagnostic Tests. *American Journal of Roentgenology* **2005**; 184:14–19.
 29. Genders TSS, Spronk S, Stijnen T, Steyerberg EW, Lesaffre E, Hunink MGM. Methods for calculating sensitivity and specificity of clustered data: a tutorial. *Radiology* **2012**; 265:910–916.
 30. Burkart KG, Brauer M, Aravkin AY, et al. Estimating the cause-specific relative risks of non-optimal temperature on daily mortality: a two-part modelling approach applied to the Global Burden of Disease Study. *The Lancet* **2021**; 398:685–697.
 31. Tenny S, Hoffman MR. Prevalence. In: StatPearls. Treasure Island (FL): StatPearls Publishing, 2022. Available at: <http://www.ncbi.nlm.nih.gov/books/NBK430867/>. Accessed 16 February 2022.
 32. Højsgaard S, Halekoh U, Yan J. The R Package geepack for Generalized Estimating Equations. *Journal of Statistical Software* **2006**; 15:1–11.
 33. Hothorn T, Bretz F, Westfall P. Simultaneous inference in general parametric models. *Biom J* **2008**; 50:346–363.
 34. CDC. SARS-CoV-2 Antigen Testing in Long Term Care Facilities Considerations for Use in Nursing Homes and other Long-Term Care Facilities care Workers. 2020. Available at: <https://www.cdc.gov/coronavirus/2019-ncov/hcp/nursing-homes-antigen-testing.html>. Accessed 3 November 2021.

35. Hu B, Guo H, Zhou P, Shi Z-L. Characteristics of SARS-CoV-2 and COVID-19. *Nat Rev Microbiol* **2021**; 19:141–154.
36. Mitchell SL, Orris S, Freeman T, et al. Performance of SARS-CoV-2 antigen testing in symptomatic and asymptomatic adults: a single-center evaluation. *BMC Infectious Diseases* **2021**; 21:1071.
37. Pray IW. Performance of an Antigen-Based Test for Asymptomatic and Symptomatic SARS-CoV-2 Testing at Two University Campuses — Wisconsin, September–October 2020. *MMWR Morb Mortal Wkly Rep* **2021**; 69. Available at: <https://www.cdc.gov/mmwr/volumes/69/wr/mm695152a3.htm>. Accessed 16 February 2022.
38. Prince-Guerra JL. Evaluation of Abbott BinaxNOW Rapid Antigen Test for SARS-CoV-2 Infection at Two Community-Based Testing Sites — Pima County, Arizona, November 3–17, 2020. *MMWR Morb Mortal Wkly Rep* **2021**; 70. Available at: <https://www.cdc.gov/mmwr/volumes/70/wr/mm7003e3.htm>. Accessed 4 August 2021.
39. Pollock NR, Jacobs JR, Tran K, et al. Performance and Implementation Evaluation of the Abbott BinaxNOW Rapid Antigen Test in a High-Throughput Drive-Through Community Testing Site in Massachusetts. *J Clin Microbiol* **2021**; 59:e00083-21.
40. Pollock NR, Tran K, Jacobs JR, et al. Performance and Operational Evaluation of the Access Bio CareStart Rapid Antigen Test in a High-Throughput Drive-Through Community Testing Site in Massachusetts. *Open Forum Infectious Diseases* **2021**; 8:ofab243.

Tables

Table1: Demographic Characteristics of Residents of Connecticut State Correctional Facilities with Time Matched Rapid Antigen and RT-PCR tests by Occurrence of SARS-CoV-2 Infection

Characterisitcs	Full Population	RT-PCR Positive SARS-CoV-2	
		Yes	No
<i>Participants</i>	(N=3825)	(N=522)	(N=3425)
Sex (N, %)			
Female	897 (23.5%)	118 (22.6%)	832 (24.3%)
Male	2928 (76.5%)	404 (77.4%)	2593 (75.7%)
Age (Mean, SD)	37 (12)	38 (11)	37 (12)
Race (N, %)			
American Indian	15 (0.4%)	3 (0.6%)	13 (0.4%)
Asian	19 (0.5%)	2 (0.4%)	18 (0.5%)
Black	1564 (40.9%)	207 (39.7%)	1392 (40.6%)
White	1186 (31.0%)	186 (35.6%)	1053 (30.7%)
Unknown/Missing	1041 (27.2%)	124 (23.8%)	949 (27.7%)
<i>RT-PCR Paired Rapid Antigen Tests</i>	(N=4835)	(N=557)	(N=4278)
Sex (N, %)			
Female	1385 (28.6%)	141 (25.3%)	1244 (29.1%)
Male	3450 (71.4%)	416 (74.7%)	3034 (70.9%)
Age (Mean, SD)	37 (12)	38 (11)	37 (12)
Race (N, %)			
American Indian	21 (0.4%)	3 (0.5%)	18 (0.4%)
Asian	22 (0.5%)	2 (0.4%)	20 (0.5%)
Black	1903 (39.4%)	219 (39.3%)	1684 (39.4%)
White	1584 (32.8%)	202 (36.3%)	1382 (32.3%)
Unknown/Missing	1305 (27.0%)	131 (23.5%)	1174 (27.4%)
Rapid Antigen Tests in Series (N) ^a			
First	3844 (79.5%)	402 (72.2%)	3442 (80.5%)
Second	677 (14.0%)	104 (18.7%)	573 (13.4%)
Third	314 (6.5%)	51 (9.2%)	263 (6.1%)

^a Rapid antigen test series defined as tests collected within 1-4 days of each other or in the absence of a test in the 4 days prior

Table 2: Rapid Antigen Test Accuracy Relative to RT-PCR Among Residents of Connecticut State Correctional Facilities Under Varying Collection Strategies

	No. RT-PCR Positive Pairs	Sensitivity (95% CI)		No. RT-PCR Negative Pairs	Specificity (95% CI)	
		Single Test ^a	Testing Strategy ^b		Single Test ^a	Testing Strategy ^b
First Rapid Antigen Test	402	67.7 (62.9, 72.0%)	66.8 (62.8, 70.6%)	3442	99.4 (99.1, 99.6%)	99.4 (99.1, 99.6%)
Second Rapid Antigen Test	104	65.4 (55.9, 74.1%)	88.8 (85.1, 91.9%)	573	99.1 (97.9, 99.6%)	98.5 (97.3, 99.1%)
Third Rapid Antigen Test	51	62.7 (49.1, 75.2%)	95.9 (93.6, 97.5%)	263	100 (98.6, 100%)	98.3 (96.7, 99.1%)

^a 95% confidence intervals estimated using generalized estimating equations with robust standard errors when >1 test pair per person was present, alternatively, Wald confidence intervals were estimated

^b Serial testing sensitivity was estimated as the additive probability (positive for any rapid antigen test); serial testing specificity was estimated as the multiplicative probability (negative for all rapid antigen test), posterior simulation of 1000 draws was used to propagate uncertainty through the equations

^c Positive predictive value (PPV) and negative predictive value (NPV) estimated using the estimated prevalence, sensitivity, and specificity

Table 3: Rapid Antigen Test Accuracy Relative to RT-PCR Among Residents of Connecticut State Correctional Facilities Under Varying Collection Strategies by Demographic and Clinical Characteristics

Characteristics	No. RT-PCR Positive Pairs	Sensitivity (95% CI)			No. RT-PCR Negative Pairs	Specificity (95% CI)		
		Single Test ^a	Serial Testing ^b	Difference in Serial Testing Accuracy ^c		Single Test ^{a,c}	Serial Testing ^{b,c}	Difference in Serial Testing Accuracy ^b
First Rapid Antigen Test								
Sample collection sequence								
Rapid Antigen Test Collected Before RT-PCR	95	49.6% (39.8, 59.3%)	45.9% (37.3, 54.8%)	-28.3% (-38.0, -18.3)	832	99.6% (98.9, 99.9%)	99.5% (98.7, 99.8%)	0.1% (-0.8, 0.6)
Rapid Antigen Test Collected Same Day as RT-PCR	236	76.3% (70.4, 81.2%)	74.2% (69.4, 78.5%)	-	1702	99.5% (99.0, 99.7%)	99.4% (98.9, 99.7%)	-
Rapid Antigen Test Collected After RT-PCR	71	63.4% (51.7, 73.7%)	65.0% (54.7, 74.1%)	-9.2% (-20.5, 1.0)	908	99.1% (98.2, 99.6%)	99.0% (98.2, 99.5%)	-0.4% (-1.3, 0.3)
Symptom presentation ^d								
Symptomatic	131	88.6% (81.9, 92.9%)	88.1% (82.4, 92.2%)	30.8% (23.5, 37.4)	79	97.5% (90.5, 99.3%)	95.3% (88.5, 98.2%)	-4.2% (-11.2, -0.9)
Asymptomatic	269	58.0% (52.0, 63.7%)	57.2% (52.2, 62.0%)	-	3357	99.5% (99.1, 99.7%)	99.5% (99.2, 99.7%)	-
Age ^e								
≤37 years	199	67.3% (60.6, 73.6%)	68.5% (63.0, 73.6%)	3.7% (-3.8, 11.5)	2020	99.6% (99.1, 99.8%)	99.4% (99.0, 99.6%)	0.1% (-0.4, 0.7)
>37 years	203	67.9% (61.5, 74.0%)	64.7% (59.2, 70.0%)	-	1422	99.2% (98.6, 99.6%)	99.3% (98.7, 99.6%)	-
Second Rapid Antigen Test								
Sample collection sequence								
Rapid Antigen Test Collected Before RT-PCR	17	35.2% (16.9, 59.2%)	67.7% (55.1, 80.4%)	-25.4% (-38.2, -12.1)	130	99.2% (94.8, 99.9%)	98.8% (94.3, 99.6%)	0.1% (-4.3, 2.3)
Rapid Antigen Test Collected Same Day as RT-PCR	69	70.9% (59.4, 80.7%)	93.2% (89.7, 95.6%)	-	267	99.2% (97.1, 99.8%)	98.7% (96.5, 99.4%)	-
Rapid Antigen Test Collected After RT-PCR	18	72.4% (48.6, 88.0%)	90.1% (80.0, 95.8%)	-3.1% (-13.5, 3.7)	176	98.9% (95.5, 99.7%)	97.9% (94.6, 99.0%)	-0.7% (-4.2, 1.6)

Symptom presentation ^d								
Symptomatic	26	88.5% (69.5, 96.2%)	98.7% (96.1, 99.6%)	16.2% (11.0, 22.0)	14	85.8% (58.2, 96.4%)	82.9% (56.4, 94.3%)	-16.0% (-42.9, -3.5)
Asymptomatic	78	57.8% (46.6, 68.1%)	82.3% (76.7, 87.0%)	-	559	99.5% (98.3, 99.8%)	98.9% (97.8, 99.4%)	-
Age ^e								
≤37 years	57	75.3% (62.6, 84.8%)	92.0% (87.4, 95.3%)	6.8% (0.2, 13.9)	353	98.6% (96.6, 99.4%)	98.1% (96.2, 99.0%)	-0.8% (-2.9, 0.9)
>37 years	47	53.0% (38.9, 66.8%)	85.1% (79.2, 89.9%)	-	220	100% (98.4, 100%)	99.0% (97.5, 99.5%)	-

Third Rapid Antigen Test

Sample collection sequence

Rapid Antigen Test Collected Before RT-PCR	5	0.1% (0.0, 51.6%)	70.9% (56.7, 86.8%)	-27.0% (-41.2, -11.0)	73	100% (95.1, 100%)	98.0% (92.1, 99.6%)	-0.2% (-6.0, 3.2)
Rapid Antigen Test Collected Same Day as RT-PCR	46	69.7% (55.0, 81.1%)	97.9% (96.3, 98.9%)	-	132	100% (97.3, 100%)	98.2% (95.3, 99.3%)	-
Rapid Antigen Test Collected After RT-PCR	0	-	-	-	58	99.9% (94.0, 100%)	97.0% (91.0, 98.9%)	-1.2% (-7.3, 2.3)

Symptom presentation^d

Symptomatic	17	88.4% (63.1, 97.1%)	99.9% (99.3, 100%)	8.6% (5.2, 13.0)	0	-	-	-
Asymptomatic	34	49.9% (33.7, 66.0%)	91.2% (86.8, 94.6%)	-	263	100% (98.6, 100%)	98.7% (97.1, 99.3%)	-
Age ^e								
≤37 years	28	64.3% (45.2, 79.6%)	97.2% (94.6, 98.6%)	2.9% (-0.9, 7.6)	169	100% (97.8, 100%)	97.7% (95.2, 98.9%)	-0.7% (-3.4, 2.9)
>37 years	23	60.8% (40.3, 78.2%)	94.2% (89.9, 97.0%)	-	94	100% (96.2, 100%)	98.5% (94.9, 99.5%)	-

^a 95% confidence intervals estimated using Generalized Estimating Equations with robust standard errors when >1 test pair per person was present, alternatively, Wald confidence intervals were estimated

^b One rapid antigen test collection accuracy metrics estimated as the weighted average of the metrics for the first, second, and third rapid antigen test of the included series; serial collection (2 or 3 tests) sensitivity was estimated as the additive probability (positive any rapid antigen test); serial testing specificity was estimated as the multiplicative probability (negative for all rapid antigen test), posterior simulation of 1000 draws was used to propagate uncertainty through the equations

^c Difference in the accuracy metric between the stratified groups calculated on the draw level; uncertainty intervals estimated using percentile method

^d Symptomatic infection was defined as reporting one or more COVID-19-related symptoms (6) at time of rapid antigen test collected closest in time to the RT-PCR

^e Stratified according to mean age of the population (37 years of age)

Figure Legends:

Figure 1: Depiction of the selected rapid antigen tests and the paired RT-PCRs

Legend: (A) The first three tests of rapid antigen test series (rapid antigen tests collected within one and four days of each other or tests collected in the absence of any test in the prior or following four days) were selected. (B) Rapid antigen tests were matched to RT-PCRs collected within one day prior to or following the rapid antigen test. If more than one RT-PCR was paired to the rapid antigen test, we preferentially selected positive RT-PCRs (Option 1) followed by those collected on the same day as the rapid antigen test (Option 2) and those collected prior to the rapid antigen test (Option 3).

Figure 2: Flowchart of residents of Connecticut state run correctional facilities between November 21, 2020, and June 15, 2021, and rapid antigen tests included in the analysis.

Legend: Flowchart of residents of Connecticut run correctional facilities and rapid antigen tests included in our analysis.

The first three tests of rapid antigen test series (rapid antigen tests collected within one and four days of each other or tests collected in the absence of any test in the prior or following four days) were selected. RT-PCR paired rapid antigen tests were defined as rapid antigen tests with RT-PCR tests collected between one day prior to and one day after the rapid antigen test.


A. Rapid Antigen Tests

* Indicates the selected rapid antigen test


First Rapid Antigen Test of a Series *



Second Rapid Antigen Test of Series *




Third Rapid Antigen Test of a Series *




B. Rapid Antigen Test/RT-PCR Pairing

* Indicates the selected RT-PCR test

Option A



Option B



Option C




Figure 1
229x99 mm (1.9 x DPI)




Figure 2
229x53 mm (1.9 x DPI)