



Article

# Development of a Sensitive Screening Method for Simultaneous Determination of Nine Genotoxic Nitrosamines in Active Pharmaceutical Ingredients by GC-MS

Anna B. Witkowska<sup>1,2,\*</sup> , Joanna Giebułtowicz<sup>3</sup> , Magdalena Dąbrowska<sup>1</sup> and Elżbieta U. Stolarczyk<sup>1,4</sup>

<sup>1</sup> Research Analytics Team, Analytical Department, Łukasiewicz Research Network—Industrial Chemistry Institute, 8 Rydygiera Street, 01-793 Warsaw, Poland

<sup>2</sup> Department of Bioanalysis and Drugs Analysis, Doctoral School, Medical University of Warsaw, 61 Żwirki i Wigury, 02-091 Warsaw, Poland

<sup>3</sup> Department of Bioanalysis and Drugs Analysis, Faculty of Pharmacy, Medical University of Warsaw, 1 Banacha, 02-097 Warsaw, Poland

<sup>4</sup> Spectrometric Methods Department, National Medicines Institute, Chelmska 30/34, 00-725 Warsaw, Poland

\* Correspondence: anna.witkowska@wum.edu.pl

**Abstract:** A worldwide crisis with nitrosamine contamination in medical products began in 2018. Therefore, trace-level analysis of nitrosamines is becoming an emerging topic of interest in the field of quality control. A novel GC-MS method with electron ionization and microextraction was developed and validated for simultaneous determination of nine carcinogenic nitrosamines (NDMA, NMEA, NDEA, NDBA, NMOR, NPYR, NPIP, NDPA, and *N*-methyl-npz) in active pharmaceutical ingredients (APIs): cilostazol, sunitinib malate, and olmesartan medoxomil. The method was validated according to the International Council for Harmonisation of Technical Requirements for Pharmaceuticals for Human Use (ICH) guidelines, demonstrating good linearity in the range of LOQ up to 21.6 ng/mL (120% of specification limit). The limits of detection for the nine nitrosamines were determined to be in the range 0.15–1.00 ng/mL. The developed trace level GC-MS method turned out to be specific, accurate, and precise. The accuracy of all the tested APIs ranged from 94.09% to 111.22% and the precision evaluated by repeatability, intermediate precision, and system precision was  $RSD \leq 7.65\%$ . Nitrosamines were not detected in cilostazol and sunitinib, whereas in olmesartan medoxomil NDEA was detected at the level of LOQ. The novel protocol was successfully applied for nitrosamines determination in selected APIs and can be used for the routine quality control of APIs under Good Manufacturing Practices rules, ensuring the safety and effectiveness of pharmaceutical products.

**Keywords:** nitrosamines; gas chromatography–mass spectrometry; active pharmaceutical ingredient; ionization; microextraction; validation



**Citation:** Witkowska, A.B.; Giebułtowicz, J.; Dąbrowska, M.; Stolarczyk, E.U. Development of a Sensitive Screening Method for Simultaneous Determination of Nine Genotoxic Nitrosamines in Active Pharmaceutical Ingredients by GC-MS. *Int. J. Mol. Sci.* **2022**, *23*, 12125. <https://doi.org/10.3390/ijms232012125>

Academic Editors: Geoffrey Brown, Enikő Kallay and Andrzej Kutner

Received: 31 August 2022

Accepted: 4 October 2022

Published: 12 October 2022

**Publisher's Note:** MDPI stays neutral with regard to jurisdictional claims in published maps and institutional affiliations.




**Copyright:** © 2022 by the authors. Licensee MDPI, Basel, Switzerland. This article is an open access article distributed under the terms and conditions of the Creative Commons Attribution (CC BY) license (<https://creativecommons.org/licenses/by/4.0/>).

## 1. Introduction

A worldwide crisis with nitrosamine contamination in medical products began in 2018 [1]. Since then, widespread investigations by regulatory agencies, including the European Medicines Agency (EMA) and United States Food and Drug Administration (US FDA), were undertaken. In drugs such as angiotensin II receptor blockers (ARBs), ranitidine, metformin, rifampin, rifapentine, and, recently, verenicline, *N*-nitrosodimethylamine (NDMA) and other nitrosamines have been detected [2,3]. Nitrosamines are known as probable human carcinogens (IARC 2A group for NDMA), potent genotoxic agents, and “cohort of concern” compounds according to ICH M7 guidance [4,5]. They can form in reaction of amine source (e.g., 2°, 3° amines, amidines, hydrazines) and nitrosating agents (e.g., nitrite in acidic conditions, nitrogen oxide, nitrous acid). According to newly established guidelines, every marketing authorization holder has to identify the risk of nitrosamine formation or cross-contamination in the manufacturing/storage of chemical and biological

medicines. Root causes of nitrosamine impurities in APIs (active pharmaceutical ingredients) and drug products can be manufacturing process related as well as stability of the drug substance/product or excipient compatibility related (Figure 1) [6–8].



**Figure 1.** Root causes of nitrosamine impurities in APIs and drug products.


The scale of the problem is revealed by the fact that, till now, more than one thousand batches of medical products have been recalled from the market due to the detection of harmful nitrosamines [9].

A large number of studies also resulted in the introduction of new, more restrictive legislation that regulated the concentration of nitrosamines in drug products. Therefore, trace-level analysis of nitrosamines is becoming an emerging topic of interest in the field of quality control. Till now, a wide range of quantitative analytical methods to measure nitrosamine traces in pharmaceuticals have been established, including gas chromatography with single quadrupole or tandem mass spectrometers (GC-MS/MS) with a headspace system or direct injection [10–16], liquid chromatography–tandem mass spectrometers (LC-MS/MS) or high-resolution mass spectrometers [17–20], and liquid chromatography with UV detectors and supercritical fluid chromatography (SFC) [21–23]. Solid-phase extraction (SPE), solid-phase microextraction (SPME), and dispersive liquid–liquid microextraction (DLLME) are also commonly used with GC-MS to produce comparable results of the LOQ (range 0.5–1.5 ppb) to LC-MS/MS [24–26]. Many already established methods are complex and can be applied to only one drug product/API, such as ranitidine and metformin. Our newly established, sensitive analytical screening procedure involving microextraction is fast, simple, environmentally friendly, and with limits of quantitation (LOQs) close to multi-step methods involving SPME or SPE extraction. Moreover GC-MS offers low operational costs, and is commonly used for analysis of volatile impurities in APIs. We demonstrate, for the first time, that our analytical procedure can be successfully applied for various APIs selected for quantitative testing after risk assessment.

## 2. Results and Discussion

### 2.1. Risk Assessment

Figure 2 presents the structures of tested APIs and N-nitrosamines examined in this study.



**Figure 2.** Structures of the APIs and nitrosamines presented in this study. Red circles—types of nitrosamine-related substructures present in APIs.

Since 2019, national regulatory agencies (US FDA, EMA) have obligated every marketing authorization holder to review their manufacturing process in order to identify the risk of nitrosamine formation or cross-contamination and, if necessary, to perform confirmatory testing. Kao et al. [27] in their recent study reported a substructure-based approach as a method for screening and investigating nitrosamines in medical products. According to their work, there are eight substructures present in APIs: DBA (dibutylamine), DEA (diethylamine), DIPA (diisopropylamine), DMA (dimethylamine), IPEA (isopropylethylamine), MPA (methylphenylamine), NMP (N-methyl-2-pyrrolidone), and a tetrazole ring can be suspected to be precursors of nitrosamines. A tetrazole ring is present in the structure of olmesartan medoxomil and cilostazol, while DEA in sunitinib malate (Figure 2). Another risk factor for nitrosamine contamination of sunitinib malate is the use of pyrrolidine in the synthesis. Many suppliers of raw materials for the synthesis of APIs do not provide information about the reagents used in the synthesis, which is why it is so important to selectively check batches with risk elements. Regulatory agencies have implemented analytical methods for sartans, antihistamines, and antidiabetic drugs, but they also continue to test other drug products; so, our new, sensitive screening method for determination of the nine nitrosamines in cilostazol and sunitinib malate is a very useful and suitable tool.


## 2.2. Method Development

A key parameter in the development of the microextraction method for the samples was optimization of the sample solvent. The solubility studies of olmesartan medoxomil, sunitinib malate, and cilostazol are presented in Supplementary Tables S1–S3. The choice of solvent involved the selection of the medium in which centrifugal extraction of the analytes under study was most efficient. Finally, methanol was chosen as the solvent in which cilostazol, olmesartan medoxomil, and sunitinib malate were slightly soluble. During method development, the centrifuge speed, centrifugation time, and amount of solvent taken for extraction were optimized. The best sensitivity was obtained for a centrifugation speed of 15,000 rpm, time of 10 min, and 250  $\mu$ L of solvent. At lower centrifugation speeds and shorter times (5000 rpm, 10,000 rpm, 2 and 5 min), there was no complete extraction of the analytes (Table S4). At longer centrifugation times, no increase in method sensitivity was observed. During the extraction, no evaporation of the solvent in the centrifuge chamber was noticed. From the capillary columns commonly used in the analysis of nitrosamines, such as DB-5 ms, DB-624 ms, and VF-WAXms [28–30], the best chromatographic separation was obtained on the column VF-WAXms (30 m  $\times$  0.32 mm, film thickness of 1.0  $\mu$ m). For the quantification analysis on a single quadrupole mass spectrometer, we used electron

ionization in selected ion monitoring mode (SIM) of particular ions, ensuring the highest sensitivity was chosen (see Section 3.2).

### 2.3. Validation Results

During the development of the method, satisfactory separation of nine nitrosamines was obtained in selected APIs. The standard added to the sample did not cause peaks splitting and the matrix did not affect the result (Figures S1–S27). Figure 3 shows a chromatogram of the standard solution.



**Figure 3.** SIM mode chromatogram of the N-nitrosamines standard—18 ng/mL.

Results for linearity, limit of detection (LOD), and LOQ are summarized in Table 1. The determination coefficients ( $R^2$ ) of the linear regression for the nine nitrosamines were over 0.995, which clearly demonstrates that the method was suitable for quantitative analysis.

**Table 1.** Validation results for linearity, LOD, and LOQ.

Analyte	Range (ng/mL)	* $Y = ax + b$ ( $R^2 \geq 0.990$ )	LOD (ng/mL)		LOQ (ng/mL)	
		** $t_a$ ** $t_b$	$S/N \geq 3$	( $3.3 \times \sigma/S$ )	$S/N \geq 10$	( $10 \times \sigma/S$ )
NDMA	0.45–21.6	$Y = 0.0127x + 0.0005$ ( $R^2 = 0.9996$ ) 158.53 0.5695	0.15	0.47	0.45	1.43
NMEA	1.5–21.6	$Y = 0.0104x + 0.0011$ ( $R^2 = 0.9997$ ) 159.85 1.4050	0.50	0.40	1.50	1.20
NDEA	1.0–21.6	$Y = 0.0088x + 0.0009$ ( $R^2 = 0.9997$ ) 153.81 1.3608	0.30	0.45	1.00	1.37
NDPA	3.0–21.6	$Y = 0.0029x + 0.0007$ ( $R^2 = 0.9988$ ) 70.521 1.4061	1.00	0.78	3.00	2.38

Table 1. Cont.

Analyte	Range (ng/mL)	*Y = ax + b (R <sup>2</sup> ≥ 0.990)	LOD (ng/mL)		LOQ (ng/mL)	
		** t <sub>a</sub> ** t <sub>b</sub>	S/N ≥ 3	(3.3 × σ/S)	S/N ≥ 10	(10 × σ/S)
NMOR	3.0–21.6	Y = 0.0039x – 0.0002 (R <sup>2</sup> = 0.9991) 83.655 0.7650	1.00	0.67	3.00	2.03
NPYR	1.5–21.6	Y = 0.0085x + 0.0007 (R <sup>2</sup> = 0.9994) 105.24 0.7193	0.50	0.60	1.50	1.82
NPIP	1.5–21.6	Y = 0.0087x + 0.0005 (R <sup>2</sup> = 0.9992) 91.615 0.4482	0.50	0.69	1.50	2.10
NDBA	1.5–21.6	Y = 0.0032x + 0.0002 (R <sup>2</sup> = 0.9986) 71.346 0.4260	0.50	0.88	1.5	2.66
N-methyl-npz	1.5–21.6	Y = 0.0065x + 0.0013 (R <sup>2</sup> = 0.9996) 137.02 2.2741	0.50	0.46	1.50	1.40

\* a—slope; b—y-intercept. \*\* t<sub>a</sub>, t<sub>b</sub>—parameters of the Student's *t*-test; t<sub>b</sub> = b/S<sub>b</sub>, t<sub>a</sub> = a/S<sub>a</sub>, where S<sub>a</sub>—the standard deviation of a, and S<sub>b</sub>—the standard deviation of b.

The LODs and LOQs calculated both ways from the linearity and ratio peak/noise height were very similar. LODs were determined to be within the range 0.15–1.00 ppb (Table 2); such low LODs are comparable to the results obtained with GC-MS/MS by Liu et al. [11] (Table S5).

The accuracy of the method was verified at three levels—80%, 100%, and 120%—of the specification limit. The results for accuracy for all tested APIs are summarized in Table 2. The results exhibited that recoveries for all nitrosamines in sunitinib malate ranged from 94.09% to 111.22%, for cilostazol from 95.91% to 104.87%, and for olmesartan medoxomil from 92.20% to 108.72%, respectively.

Method precision was evaluated by repeatability, intermediate precision, and system precision. As summarized in Table 3, this GC-MS method exhibited satisfactory RSD% values for repeatability in all APIs in the range of 0.76–5.75%, and for intermediate precision and system precision in the range of 1.13–7.65% and 1.03–3.65% (area), respectively. Since the acceptance criteria for recovery and precision for all APIs were established for 60–115% and RSD < 21%, all the results meet the criteria according to the Association of Official Analytical Chemists [31].

The robustness of the GC-MS method was evaluated on standard solution to determine the impact on the result in case of small variations in the chromatographic parameters: oven initial column temperature ± 5 °C; first temperature ramp ± 1 °C/min; and gas pressure ± 10%. The robustness results (Table S6) showed that the retention time for all nitrosamines change within a range of ± 1.1 min and the separation of the nitrosamines was not affected.

**Table 2.** Validation results of the accuracy for sunitinib, cilostazol, and olmesartan medoxomil.

Analyte	Accuracy Mean Recovery 60–115% RSD ≤ 21% (n = 9)		
	Sunitinib Malate	Cilostazol	Olmesartan Medoxomil
NDMA	94.45 1.98	98.26 2.11	104.66 5.68
NMEA	94.09 2.10	99.14 1.78	95.85 4.24
NDEA	97.97 3.57	97.15 2.02	97.34 4.44
NDPA	97.07 5.05	95.91 2.15	92.20 6.50
NMOR	111.22 2.89	97.10 3.66	103.60 5.16
NPYR	95.87 3.82	100.62 3.00	108.72 5.29
NPIP	99.37 5.13	96.66 6.02	ND
NDBA	99.30 6.10	104.87 2.38	102.17 8.98
<i>N</i> -methyl-npz	102.94 7.53	99.90 5.17	ND

ND—not detected.

**Table 3.** Validation results of precision for sunitinib malate, cilostazol, and olmesartan medoxomil.

Matrix Analyte	Sunitinib		Cilostazol		Olmesartan Medoxomil		System Precision	
	Repeatability RSD (%)	Intermediate Precision RSD (%)	Repeatability RSD (%)	Intermediate Precision RSD (%)	Repeatability RSD (%)	Intermediate Precision RSD (%)	Area Nitrosamine/Area NDMA-d6 RSD (%)	Retention Time RSD (%)
NDMA	2.12	2.05	0.76	1.13	3.07	2.51	1.31	0.03
NMEA	2.25	1.99	1.27	1.42	2.47	2.21	1.08	0.02
NDEA	2.48	3.37	1.10	1.37	3.00	4.04	1.31	0.02
NDPA	3.95	6.26	2.97	2.25	5.75	4.24	1.03	0.01
NMOR	3.18	5.55	2.73	4.49	3.07	3.97	3.52	0.01
NPYR	2.51	5.08	1.77	3.91	2.55	7.50	2.42	0.01
NPIP	3.67	4.06	3.91	4.60	ND	ND	2.01	0.01
NDBA	4.32	6.81	1.68	1.46	4.51	4.71	3.07	0.01
<i>N</i> -methyl-npz	2.22	6.45	5.60	7.65	ND	ND	3.65	0.01

ND—not detected.

Moreover, during validation, the stability of the standard solution and sample solution was investigated. Standard solution stability studies for all nitrosamines in methanol demonstrated that after 24 and 48 h in 25 °C, the RSD% (area) ranged from 0.63 to 4.79%, demonstrating that the standard solution was stable up to 48 h after preparation (Tables S7–S15). Stability of the sample solution of cilostazol, olmesartan medoxomil, and sunitinib malate spiked with nitrosamines at 100% of a specification level was evaluated after 24 h in 25 °C by measuring the peak area quotient of nitrosamines and NDMA-d6. The RSD% for all the nitrosamines in olmesartan medoxomil was in range of 1.36–18.33%, in cilostazol, 0.69–9.14%, and in sunitinib malate, 0.38–10.63%, which indicated the increase

in RSD with time and the necessity of sample solution preparation directly before analysis (Tables S16–S24).

Supplementary Material Table S5 shows the comparison of the results of this publication with various reported GC-MS and GC-MS/MS methods for the detection of N-nitrosamines. The novel GC-MS method with microextraction showed superior results in the validation test compared with other GC-MS and GC-MS/MS methods. The developed method offers several benefits over previously established methodologies, including enhanced sensitivity (LOD below 1 ppb), with no requirements for derivatization and multistep sample preparation (SPE, SPME or DLLEM); a reduction in solvent usage (only 250  $\mu$ L); and application to a wide variety of APIs. We have also demonstrated that the developed method using a single quadrupole is as good as the method with tandem mass spectrometry.

#### 2.4. Application of the Method for Nitrosamine Determination in Selected APIs

The validated GC-MS method was applied to determine nine nitrosamines in three active pharmaceutical ingredients: cilostazol, sunitinib malate, and olmesartan medoxomil. In cilostazol and sunitinib, no nitrosamines were found. NDEA was detected in olmesartan medoxomil at a concentration of 1.66 ng/mL. Therefore, based on the proposed European Pharmacopoeia Commission limit (30 ppb), the acceptance criteria for nitrosamines in APIs have been met. Olmesartan medoxomil has never been recalled from the market due to detection of nitrosamines beyond acceptable intake limits [2].

### 3. Materials and Methods

#### 3.1. Materials

Cilostazol (99.8%), sunitinib malate (99.7%), and olmesartan medoxomil (99.7%) were purchased from Pharmaceutical Research Institute (Poland). The nitrosamines mix (containing NDMA, NMEA, NDEA, NDPA, NMOR, NPYR, NPIP, and NDBA) EPA 8270 (Appendix IX), certified reference material; 2000  $\mu$ g/mL in methanol) was obtained from Sigma-Aldrich (USA). 1-Methyl-4-nitrosopiperazine was obtained from Angene International Limited (London, England). N-Nitrosodimethylamine-d6 (1000  $\mu$ g/mL in methanol) was purchased from LGC (Kielcin, Poland). Methanol LC-MS reagent was obtained from J.T Baker (Oslo, Norway).

#### 3.2. Equipment and Method

The GC-MS analysis was performed using a Shimadzu gas chromatograph GCMS-QP2010 ULTRA with MS detector. The separation of nitrosamines was performed on a VF-WAXms polyethylene glycol column (30 m  $\times$  0.32 mm, film thickness of 1.0  $\mu$ m) from Agilent Technologies. In the study, centrifuge tubes with a 0.2  $\mu$ m PVDF filter from Anchem and an M-Science centrifuge from MPW Med. Instruments were used. All the GC-MS method's experimental conditions are presented in Table 4. The ions selected for each nitrosamine for quantification and confirmation are shown in Table 5.

#### 3.3. Standard Solution Preparation

The internal standard (IS) solution (70 ng/mL in methanol) was prepared from the internal stock solution NDMA-d6 (1000  $\mu$ g/mL in methanol) and diluted with methanol. The nitrosamine intermediate dilution was prepared by dissolving the appropriate amount of the nitrosamines mix (certified reference material) to reach 10  $\mu$ g/mL of the nitrosamine mix in MeOH. A standard solution of nine nitrosamines at the concentration level of 18.0 ng/mL was prepared from the nitrosamine intermediate dilution in a 10 mL volumetric flask with methanol.

**Table 4.** GC-MS method experimental conditions.

GC Parameters	
Carrier gas	Helium
Column flow rate	2.1 mL/min
Oven initial temperature	75 °C (hold time: 2 min)
First temperature ramp	9.0 °C/min
Second temperature	200 °C/min (hold time: 2 min)
Second temperature ramp	15.0 °C/min
Final temperature	240 °C (hold time: 3 min)
Pressure	45.5 kPa
Injection mode	Splitless, high-pressure injection 130 kPa
Injection port temperature	150 °C
MS Parameters	
Ion Source Temperature	240
Interface Temperature	240
Acquisition type	SIM
Solvent delay	6.5 min

**Table 5.** GC-MS SIM for the detected nitrosamines.

Compound	Quantification Ion ( <i>m/z</i> )	Confirmation Ion ( <i>m/z</i> )
NDMA	74	44
NMEA	88	56
NDEA	102	56
NDPA	130	70
NMOR	116	56
NPYR	100	68
NPIP	114	55
NDBA	116	158, 84
<i>N</i> -methyl-npz	99	100, 56
NDMA-d6	80	50

### 3.4. Sample Solution Preparation

A total of 150 mg of API was weighed into 2 mL centrifuge tubes with filtration equipment and 250  $\mu$ L of the internal standard solution (70 ng/mL NDMA-d6) was added. Then, the solution was centrifuged at 15,000 rpm for 10 min. A clear organic phase was transferred to a sample vial because the selected APIs were practically insoluble in methanol.

### 3.5. Blank Preparation

The blank solution was prepared as described in Section 3.4., but without addition of the sample.

### 3.6. Validation of the GC-MS Method

The GC-MS method for the nine nitrosamines in API was validated according to ICH guidelines. The specificity of the method was examined using a blank solution, standard solution, sample solution, and reference solution (sample solution spiked with standard solution at a 100% concentration level of the specified limit).

The calibration curve was prepared in methanol at concentrations from LOQ to 120% (0.45–21.6 ng/mL) of the specified limit for nitrosamines with respect to sample preparation, the internal standard solution was added in a concentration of 70 ng/mL NDMA-d6. The linearity for every nitrosamine was evaluated by linear regression and statistical evaluation. Linearity experiments included the determination coefficient, *y*-intercept, slope of regression line, and residual sum of squares (*Sr*). The determination coefficients (*R*<sup>2</sup>) of the linear regression for each calibration curve should be  $\geq 0.995$ .



The accuracy of the method was established on three different concentrations: a sample spiked with nitrosamines at 80%, 100%, and 120% of the specification limit; the acceptance criteria were 60–115% and the precision at each concentration should not exceed  $RSD < 21\%$ . In 2020, the European Pharmacopoeia Commission adopted a new general chapter for trace-level analysis of nitrosamine impurities where the recommended maximum allowed limit for nitrosamines in active substances was set at 0.03 ppm (30 ppb) [32]; therefore, in the present study, the same limit was implemented.

Limit of detection (LOD) and limit of quantitation (LOQ) were calculated in two different ways: using peak height (S) and noise height (N), where the LOQ was stated as  $S/N \geq 10$  and LOD as  $S/N \geq 3$ , and using linearity parameters such as the standard deviation of the response ( $\sigma$ ) and the slope of the curve (S); then, LOD and LOQ were determined from the following equations:  $LOD = (3.3 \times \sigma/S)$  and  $LOQ = (10 \times \sigma/S)$ .

The precision of the method was examined as repeatability, intermediate precision, and system precision. Repeatability for all nitrosamines was obtained from analysis of the six sample solutions. Intermediate precision was assessed as a result of the intra-laboratory variations, i.e., a different day or analyst (the RSD should not exceed 21%), while the difference between the results were tested by statistical Horwitz's test. System precision was evaluated by six replicate injections of the standard solution (calculation—Supplementary Note S1).

During validation, the stability of the solutions and robustness of the method also were evaluated. The stability of the solutions was established by measuring the peak area quotient of the nitrosamines and NDMA-d6 for the reference solution, with nitrosamines at a 100% concentration level of the specified limit and for a sample solution of API (cilostazol, sunitinib, and olmesartan medoxomil) spiked with nitrosamines at 100% of the specification level (18.0 ng/mL) after 24 h under autosampler conditions ( $20\text{ }^{\circ}\text{C} \pm 5\text{ }^{\circ}\text{C}$ ). Robustness of the developed method was investigated by injecting the standard solution in different chromatographic conditions: column temperature,  $\pm 5\text{ }^{\circ}\text{C}$ ; rate,  $\pm 1\text{ }^{\circ}\text{C}/\text{min}$ ; and carrier gas pressure,  $\pm 10\%$ .

#### 4. Conclusions

A new, sensitive, simple, environmentally friendly method for extracting nine nitrosamines from active pharmaceutical ingredients, using the high-pressure direct injection GC-MS method, was developed in this study. Based on the performed risk evaluation in the manufacturing process, three APIs were examined. This is the first described method for the determination of nitrosamines in cilostazol and sunitinib malate. The implemented GC-MS method for cilostazol and sunitinib malate showed no influence of the manufacturing process on the development of nitrosamine impurities. The quantitative method for each API was validated according to the requirements of the ICH guideline. The validation of the GC-MS method proved that the method was sensitive, accurate, precise, and suitable for nitrosamine determination in multiple APIs. The results obtained during the validation showed that the LOD and LOQ were in the range from 0.15 to 1.00 ng/mL. In conclusion, the presence of nine nitrosamines was determined as below the acceptable intake limits in selected APIs; however, it is necessary to test other APIs containing specific amine substructures in their structures to ensure the safety of the drugs for patients.

**Supplementary Materials:** The following supporting information can be downloaded at: <https://www.mdpi.com/article/10.3390/ijms232012125/s1>.

**Author Contributions:** Conceptualization, supervision, E.U.S.; methodology, validation, formal analysis, investigation, visualization, A.B.W.; writing—original draft, A.B.W.; writing—review and editing, A.B.W., E.U.S., J.G., M.D.; funding acquisition A.B.W. All authors have read and agreed to the published version of the manuscript.

**Funding:** The project was financially supported by the Ministry of Science and Higher Education: contract number DWD/5/0401/2021.

**Institutional Review Board Statement:** Not applicable.

**Informed Consent Statement:** Not applicable.

**Data Availability Statement:** Not applicable.

**Acknowledgments:** The study was supported by the Polish Ministry of Science and Higher Education, statutory grant no. 841343A. The authors would like to thank Marta Zezula for her support in this study.

**Conflicts of Interest:** The authors declare no conflict of interest.

## References

1. US Food and Drug Administration. FDA Announces Voluntary Recall of Several Medicines Containing Valsartan Following Detection of an Impurity. Available online: <https://www.fda.gov/news-events/press-announcements/fda-announces-voluntary-recall-several-medicines-containing-valsartan-following-detection-impurity> (accessed on 30 August 2022).
2. US Food and Drug Administration. Recalls, Market Withdrawals, & Safety Alerts. Available online: <https://www.fda.gov/safety/recalls-market-withdrawals-safety-alerts> (accessed on 30 August 2022).
3. Tao, X.; Tian, Y.; Liu, W.-H.; Yao, S.; Yin, L. Trace Level Quantification of 4-Methyl-1-Nitrosopiperazin in Rifampicin Capsules by LC-MS/MS. *Front Chem* **2022**, *10*, 834124. [[CrossRef](#)] [[PubMed](#)]
4. European Medicines Agency; Committee for Medicinal Products for Human Use. *ICH Guideline M7 on Assessment and Control of DNA Reactive (Mutagenic) Impurities in Pharmaceuticals to Limit Potential Carcinogenic Risk-Addendum*; European Medicines Agency: Amsterdam, The Netherlands, 2021; Volume 31.
5. International Agency for Research on Cancer. *IARC Publications Website-Overall Evaluations of Carcinogenicity: An Updating of IARC Monographs*; International Agency for Research on Cancer: Lyon, France, 1987; Volume 1–42.
6. Wu, Y.; Levons, J.; Narang, A.S.; Raghavan, K.; Rao, V.M. Reactive Impurities in Excipients: Profiling, Identification and Mitigation of Drug–Excipient Incompatibility. *AAPS PharmSciTech* **2011**, *12*, 1248–1263. [[CrossRef](#)]
7. Yokoo, H.; Yamamoto, E.; Masada, S.; Uchiyama, N.; Tsuji, G.; Hakamatsuka, T.; Demizu, Y.; Izutsu, K.I.; Goda, Y. N-Nitrosodimethylamine (Ndma) Formation from Ranitidine Impurities: Possible Root Causes of the Presence of Ndma in Ranitidine Hydrochloride. *Chem. Pharm. Bull.* **2021**, *69*, 872–876. [[CrossRef](#)]
8. Sedlo, I.; Kolonić, T.; Tomić, S. Presence of Nitrosamine Impurities in Medicinal Products. *Arh. Hig. Rada Toksikol.* **2021**, *72*, 1–5. [[CrossRef](#)] [[PubMed](#)]
9. Bharate, S.S. Critical Analysis of Drug Product Recalls Due to Nitrosamine Impurities. *J. Med. Chem.* **2021**, *64*, 2923–2936. [[CrossRef](#)] [[PubMed](#)]
10. Parr, M.K.; Joseph, J.F. NDMA Impurity in Valsartan and Other Pharmaceutical Products: Analytical Methods for the Determination of N-Nitrosamines. *J. Pharm. Biomed. Anal.* **2019**, *164*, 536–549. [[CrossRef](#)]
11. Liu, J.; Xie, B.; Mai, B.; Cai, Q.; He, R.; Guo, D.; Zhang, Z.; Fan, J.; Zhang, W. Development of a Sensitive and Stable GC-MS/MS Method for Simultaneous Determination of Four N-Nitrosamine Genotoxic Impurities in Sartan Substances. *J. Anal. Sci. Technol.* **2021**, *12*, 3. [[CrossRef](#)]
12. US Food and Drug Administration. *Combined Headspace N-Nitrosodimethylamine (NDMA), N-Nitrosodiethylamine (NDEA), N-Nitrosoethylisopropylamine (NEIPA), and N-Nitrosodiisopropylamine (NDIPA) Impurity Assay by GC-MS/MS*; US Food and Drug Administration: Silver Spring, MD, USA, 2019; pp. 1–7.
13. Tsutsumi, T.; Akiyama, H.; Demizu, Y.; Uchiyama, N.; Masada, S.; Tsuji, G.; Arai, R.; Abe, Y.; Hakamatsuka, T.; Izutsu, K.; et al. Analysis of an Impurity, N-Nitrosodimethylamine, in Valsartan Drug Substances and Associated Products Using GC-MS. *Biol. Pharm. Bull.* **2019**, *42*, 547–551. [[CrossRef](#)]
14. New OMCL Method for Simultaneous Determination of NDMA and NDEA in Sartan. Available online: <https://www.edqm.eu/documents/52006/71923/Ad-hoc-projects-OMCL-Network-ranitidine.pdf/f775fbf3-705e-ce82-d1e1-8eb3c06002d4?t=1628667875462> (accessed on 30 August 2022).
15. Wichitnithad, W.; Sudtanon, O.; Srisunak, P.; Cheewatanakornkool, K.; Nantaphol, S.; Rojsitthisak, P. Development of a Sensitive Headspace Gas Chromatography–Mass Spectrometry Method for the Simultaneous Determination of Nitrosamines in Losartan Active Pharmaceutical Ingredients. *ACS Omega* **2021**, *6*, 11048–11058. [[CrossRef](#)]
16. Lee, D.H.; Hwang, S.H.; Park, S.; Lee, J.; Oh, H.B.; Han, S.B.; Liu, K.-H.; Lee, Y.-M.; Pyo, H.S.; Hong, J. A Solvent-Free Headspace GC/MS Method for Sensitive Screening of N-Nitrosodimethylamine in Drug Products. *Anal. Methods* **2021**, *13*, 3402–3409. [[CrossRef](#)] [[PubMed](#)]
17. Chang, S.-H.; Chang, C.-C.; Wang, L.-J.; Chen, W.-C.; Fan, S.-Y.; Zang, C.-Z.; Hsu, Y.-H.; Lin, M.-C.; Tseng, S.-H.; Wang, D.-Y. A Multi-Analyte LC-MS/MS Method for Screening and Quantification of Nitrosamines in Sartans. *J. Food Drug Anal.* **2020**, *28*, 292–301. [[CrossRef](#)]
18. Yang, J.; Marzan, T.A.; Ye, W.; Sommers, C.D.; Rodriguez, J.D.; Keire, D.A. A Cautionary Tale: Quantitative LC-HRMS Analytical Procedures for the Analysis of N-Nitrosodimethylamine in Metformin. *AAPS J.* **2020**, *22*, 89. [[CrossRef](#)]
19. Schmalder-ripcke, J.; Hannes, C. *Test Method for the Determination of NDMA by LC-MS/MS in Ranitidine Drug Substance and Film Coated Tablets*; EDQM: Strasbourg, France, 2019; pp. 1–7.

20. European Directorate for the Quality of Medicines and HealthCare. LC-MS/MS Method for the Determination of NDEA and NMDA in Valsartan, Irbesartan and Losartan APIs and Finished Dosage Forms. Available online: [https://www.edqm.eu/sites/default/files/medias/fichiers/Sartans/de-by-lc-ms\\_ls.pdf](https://www.edqm.eu/sites/default/files/medias/fichiers/Sartans/de-by-lc-ms_ls.pdf) (accessed on 30 August 2022).
21. Schmidtsdorff, S.; Schmidt, A.H. Simultaneous Detection of Nitrosamines and Other Sartan-Related Impurities in Active Pharmaceutical Ingredients by Supercritical Fluid Chromatography. *J. Pharm. Biomed. Anal.* **2019**, *174*, 151–160. [[CrossRef](#)] [[PubMed](#)]
22. French National Agency for Medicines and Health Products Safety. Determination of NDMA in Valsartan Active Substances and Finished Products by HPLC/UV. Available online: <https://www.edqm.eu/sites/default/files/omcl-method-determination-ndma-valsartan-ansm-september2018.pdf> (accessed on 30 August 2022).
23. Maziarz, M.; Rainville, P. Reliable HPLC/UV Quantification of Nitrosamine Impurities in Valsartan and Ranitidine Drug Substances [TECHNOLOGY BRIEF]. Available online: <https://www.waters.com/webassets/cms/library/docs/720006775en.pdf> (accessed on 30 August 2022).
24. Giménez-Campillo, C.; Pastor-Belda, M.; Campillo, N.; Hernández-Córdoba, M.; Viñas, P. Development of a New Methodology for the Determination of N-Nitrosamines Impurities in Ranitidine Pharmaceuticals Using Microextraction and Gas Chromatography–Mass Spectrometry. *Talanta* **2021**, *223*, 121659. [[CrossRef](#)] [[PubMed](#)]
25. Alshehri, Y.M.; Alghamdi, T.S.; Aldawsari, F.S. HS-SPME-GC-MS as an Alternative Method for NDMA Analysis in Ranitidine Products. *J. Pharm. Biomed. Anal.* **2020**, *191*, 113582. [[CrossRef](#)]
26. Lim, H.H.; Oh, Y.S.; Shin, H.S. Determination of N-Nitrosodimethylamine and N-Nitrosomethylethylamine in Drug Substances and Products of Sartans, Metformin and Ranitidine by Precipitation and Solid Phase Extraction and Gas Chromatography–Tandem Mass Spectrometry. *J. Pharm. Biomed. Anal.* **2020**, *189*, 113460. [[CrossRef](#)]
27. Kao, Y.-T.; Wang, S.-F.; Wu, M.-H.; Her, S.-H.; Yang, Y.-H.; Lee, C.-H.; Lee, H.-F.; Lee, A.-R.; Chang, L.-C.; Pao, L.-H. A Substructure-Based Screening Approach to Uncover N-Nitrosamines in Drug Substances. *J. Food Drug Anal.* **2022**, *30*, 150–162. [[CrossRef](#)] [[PubMed](#)]
28. Sieira, B.J.; Carpinteiro, I.; Rodil, R.; Quintana, J.B.; Cela, R. Determination of N-Nitrosamines by Gas Chromatography Coupled to Quadrupole–Time-of-Flight Mass Spectrometry in Water Samples. *Separations* **2020**, *7*, 3. [[CrossRef](#)]
29. Jurado-Sánchez, B.; Ballesteros, E.; Gallego, M. Comparison of the Sensitivities of Seven N-Nitrosamines in Pre-Screened Waters Using an Automated Preconcentration System and Gas Chromatography with Different Detectors. *J. Chromatogr. A* **2007**, *1154*, 66–73. [[CrossRef](#)] [[PubMed](#)]
30. ChaoYe, S.; HongQuan, C.; SaiFeng, P.; Yun, Z. Determination of Nine Nitrosamines in Water by Solid Phase Extraction and Gas Chromatography–Mass Spectrometry. *J. Environ. Occup. Med.* **2019**, *36*, 1060–1065. [[CrossRef](#)]
31. Association of Official Analytical Chemists. *Appendix F: Guidelines for Standard Method Performance Requirements*; Association of Official Analytical Chemists: Rockville, MD, USA, 2016.
32. European Pharmacopoeia. 2.5.42 N-Nitrosamines in Active Substances; European Pharmacopoeia: Strasbourg, France, 2020.