



Intraductal Delivery and X-ray Visualization of Ethanol-Based Ablative Solution for Prevention and Local Treatment of Breast Cancer in Mouse Models

Elizabeth Kenyon^{1,2}, Erin K. Zaluzec^{1,3}, Katherine Powell^{1,2}, Maximilian Volk^{1,4}, Shatadru Chakravarty^{2,5}, Jeremy Hix^{2,6}, Ripla Arora^{7,8}, Jennifer J. Westerhuis^{9,10}, Matti Kiupel¹¹, Erik M. Shapiro^{2,6}, Lorenzo F. Sempere^{1,2}

¹Precision Health Program, Michigan State University, East Lansing, MI, USA

²Department of Radiology, College of Human Medicine, Michigan State University, East Lansing, MI, USA

³Department of Pharmacology & Toxicology, College of Veterinary Medicine, Michigan State University, East Lansing, MI, USA

⁴College of Osteopathic Medicine, Michigan State University, East Lansing, MI, USA

⁵Advanced Materials Characterization Laboratory/Materials Research Center, Missouri University of Science and Technology, Rolla, MO, USA

⁶Institute for Quantitative (IQ) Health Science and Engineering Advanced Molecular Imaging Facility, Michigan State University, East Lansing, MI, USA

⁷Department of Obstetrics Gynecology and Reproductive Biology, College of Human Medicine, Michigan State University, East Lansing, MI, USA

⁸Institute for Quantitative (IQ) Health Science and Engineering, Michigan State University, East Lansing, MI, USA

⁹Van Andel Research Institute, Grand Rapids, MI, USA

¹⁰Miltenyi Biotec North America, USA

¹¹Veterinary Diagnostic Laboratory, College of Veterinary Medicine, Michigan State University, Lansing, MI, USA

Abstract

Breast cancer is the most prevalent cancer and the second-leading cause of cancer-related death for women in the USA. For high-risk women, prophylactic mastectomy is the most effective primary prevention strategy. Prophylactic mastectomy is an aggressive surgical procedure that completely removes the mammary epithelial cells from which breast cancer arises along with surrounding

Corresponding author: Lorenzo F. Sempere (semperel@msu.edu).

A complete version of this article that includes the video component is available at <http://dx.doi.org/10.3791/63457>.

DISCLOSURES:

The authors have nothing to disclose.

tissue. We seek to develop a minimally invasive intraductal procedure as an alternative to prophylactic mastectomy to locally ablate the mammary epithelial cells before they can become malignant. We and others have developed an intraductal delivery procedure to reach and treat these epithelial cells in rodent models of breast cancer. While the mouse mammary gland with a single non-anastomosed ductal tree opening at the nipple has a much less complex and tortuous architecture than the human breast, chemically induced and genetically engineered mouse models of breast cancer are valuable to produce proof-of-concept studies of new preventative strategies. Here, we describe a procedure for intraductal delivery of an ethanol-based ablative solution containing micro-CT/X-ray tantalum-based contrast agent within the mouse mammary ductal tree for the therapeutic purpose of primary prevention of breast cancer. Intraductal delivery of aqueous reagents (e.g., cytotoxic compounds, siRNAs, AdCre) has been previously described in mouse models. Thus, we focus our protocol description on methodological modifications and unique experimental considerations for optimizing delivery of ethanol, for minimizing local and systemic side effects of ethanol administration, and for *in vivo* visualization of ductal tree filling via micro-CT/fluoroscopy imaging. Visualization of the ductal tree immediately after injection of a contrast-containing solution allows for confirmation of complete filling or unsuccessful outcomes such as underfilling or overfilling. This procedure can be applied for delivery and imaging of other ablative compounds aimed at either preventing tumor formation or locally treating early-stage tumors accessible via the ductal tree.

SUMMARY:

A method of intraductal injection of reagents for an ethanol-based ablative solution to the mouse mammary ductal tree for *in vivo* imaging and breast cancer prevention is described. Injection directly into the nipple opening allows for targeting mammary epithelial cells with minimal collateral tissue damage.

Keywords

ductal tree; intraductal; mammary gland; ductography; chemical ablation; image-guided procedure

INTRODUCTION:

Breast cancer is a common and potentially lethal disease with few options available for prevention¹. The most effective intervention is prophylactic mastectomy; however, only high-risk individuals opt to undergo this procedure as it is a surgery with major life-changing consequences². The procedure completely removes the mammary epithelial cells from which breast cancer arises along with the surrounding tissue. This can result in physical, psychological, and social stress for the individual, and often dissuades individuals from proceeding with this surgical procedure as their first line of primary intervention.

We have demonstrated that delivery of an ablative solution containing 70% ethanol (EtOH) directly into the ductal tree is effective at killing mammary epithelial cells with limited collateral tissue damage and at preventing breast tumors in mouse models³. EtOH has been long used clinically as an ablative or sclerosing agent for local treatment. Percutaneous EtOH injection is used as an ablative agent for unresectable liver tumors, renal and adrenal

neoplasms, and pancreatic cystic tumors⁴⁻⁶; for celiac plexus neurolysis to reduce pain⁷; and for treating breast pseudoaneurysms⁸. Intravascular EtOH injection is used as a sclerosing agent to eliminate swelling and deformation from arteriovenous malformations (AVM), and for cosmetic treatment of spider veins and varicose veins⁹⁻¹³. Like prophylactic mastectomy, the success of prevention with local delivery of an ablative solution hinges on the ability to completely remove all mammary epithelial cells from which cancer could potentially arise. This requires confirmation that the ablative substance has successfully filled the ductal tree, thus contacting all mammary epithelial cells directly. Clinical means for injecting substances within the mammary glands and visualizing them by image-guided fluoroscopy or ductography are readily available^{14,15}; therefore, it will be possible to both deliver and confirm successful delivery when this procedure may warrant evaluation in clinical trials.

Demonstrating the feasibility of this image-guided approach in laboratory animals is a key step in establishing the efficacy and translational feasibility of intraductal (ID) ablation as a preventative measure for breast cancer. In our laboratory, we have developed a method to successfully inject all mammary glands in mice with an ablative solution containing a contrast agent over a course of weekly injections to ensure the animal does not succumb to an overdose of EtOH (Figure 1, Figure 2, refs^{3,16}). This procedure places a 34-gauge needle inside the nipple opening of an isoflurane-anesthetized mouse to inject the test solution. Some key improvements of the procedure include the use of gastight syringes for liquid and gases, injection of higher volumes per ductal trees¹⁷, and extended anti-inflammatory treatment. The preclinical treatment of 5 mg/kg of carprofen, an NSAID, from 2 d before to 7 d after the ID procedure is in line with that of clinical sclerosing therapy for AVM. Typically, after systemic anesthesia the patients receive anti-inflammatory medications, such as NSAIDs, for two days post-procedure that can be extended to mitigate any local inflammation or pain¹². Alcohol intoxication is significantly mitigated by intraperitoneal injection of a 5% sucrose solution in mice. With administration of this sucrose solution, mice can be safely injected with up to 160 μ L of 70% EtOH (up to 4 ductal trees; about 0.4 g/dL of EtOH content in blood); animals fully recovered within 4 h after ID injections. For injection of more than 4 glands in mice and/or higher EtOH concentrations, we perform sequential sessions to allow enough recovery time. Alcohol intoxication in women would be a lesser concern due to the lower proportion of alcohol amount to body weight. Given the number of ductal trees in human breast^{14,15}, about 16, and estimated volume to fill each tree duct^{18,19}, up to 32 mL of 70% EtOH will be administered. This quantity will be much lower than the 50 mL of 95% EtOH administered in other clinical procedures^{4,9}. Intravenous administration of thiamine and glucose solution could be used to further minimize effects of EtOH intoxication, especially in cases where a larger total volume of EtOH may need to be injected and/or for women who have a lower tolerance to alcohol consumption (e.g., allelic variants in alcohol or aldehyde dehydrogenases).

Imaging via micro-CT/fluoroscopy allows us to confirm successful ductal filling of each gland (Figure 1, Figure 2, Figure 3). This can be recorded for future analysis, or assessed in the moment via real-time fluoroscopy imaging, as would be done in clinical application, to limit overall radiation burden imposed on the animal. To further improve specific features of this ablative solution for real time image-guided delivery *in vivo*, we previously compared FDA-approved iodine-containing contrast to a tantalum oxide (TaO_x)-containing

nanoparticle synthesized by the Shapiro lab^{3,16}. TaO_x showed a superior performance as a micro-CT contrast agent for visualizing the initial filling of the ductal tree (Figure 2, Figure 3). TaO_x can be used as a reference contrast to perform a more systematic and longitudinal assessment of other nanoparticle-based blood pool contrast agents (e.g., iodine-, bismuth-, or gold-containing) and compatibility of TaO_x with different concentrations of ethyl cellulose as gelling agent^{20,21}.

PROTOCOL:

All described experiments were conducted under protocols approved by Institutional Animal Care and Use Committee at Michigan State University.

1. Extended anti-inflammatory treatment

- 1.1. Ensure that mice receiving injection solutions containing EtOH or other potential irritants are provided with anti-inflammatory treatment from 2 days before injection to 7 days after injection. The preferred method of dosing is oral delivery through a sucralose gel cup containing carprofen at an appropriate concentration.
- 1.2. Prepare carprofen at required concentration. For this experiment, a working solution of 2 mg/mL (Stock solution is 50 mg/mL) was prepared by diluting in sterile PBS to inject 0.5 mL to achieve a final dose of 1 mg per cup. Add 1% v/v sterile blue food dye (BFD) in place of a fraction of the PBS of the total volume needed for dilution to better visualize complete mixing of the drug into the sucralose gel.
- 1.3. Prepare a cup for carprofen addition as per the manufacturer's recommendation. Unless recommended otherwise, place the cup in 60 °C water bath for 15 min. Remove the cups and dry them off to minimize possibility of contamination.
 - 1.3.1. Use 70% EtOH or an EtOH wipe to clean the surface of the cup lid and allow to dry. Use a syringe to inject the necessary amount of carprofen solution through the lid. For this experiment 500 µL is injected. Cover the injection site with a sticker and shake vigorously for 15 s.
 - 1.3.2. Vortex the cup for 15 s and then store for later use after visually verifying complete mixing. Check the cups for homogenous mixing by looking for the dark blue clumps. Allow the cups to come to room temperature before moving to the refrigerator for storage up to one month.

NOTE: These cups can also be stored at room temperature once dosed if needed but take care to pay attention to drug efficacy guidance from the manufacturer. Dating the sticker helps with keeping track of the injection date without risking a pen or sharp marker puncturing the lid.

- 1.4. When ready to use, wipe down the exterior of the cup with 70% EtOH. Peel off the lid and place the cup into the cage with the mice. Replace the cups every other day or when empty. One cup should be enough for up to 5 mice for 1–2 days.

2. Preoperative preparation

NOTE: This step will occur 2–3 days prior to injection.

- 2.1. Anesthetize the mouse using an isoflurane vaporizer (2–3% isoflurane, 1.5 L/min oxygen) and apply eye lubricant. Maintain anesthesia at 1–3% isoflurane as needed throughout the procedure with a nose cone while carefully monitoring the respiratory rate of the animal.
- 2.2. Apply eye lubricant to the mouse eyes while at the nose cone and then position the mouse on its back. Use a cotton tipped applicator to apply over-the-counter depilatory cream to the nipple area. Rub the cream into the area for 10–30 s with the applicator to help loosen the fur quickly.

NOTE: Do not leave the cream on the animal any longer than necessary and remove completely to avoid burning skin.

- 2.3. After 10–30 s of the application, use warm water on gauze to completely remove the cream and loosen fur from the animal. Perform at least 2–4 rinses of the area with fresh gauze and dry with a dry gauze after final rinse. Check the area of fur removal to ensure good visibility and access to the nipples.
- 2.4. Place the mouse in a separate clean, dry recovery cage on a heating pad. Observe the mouse until fully recovered from anesthesia and return to the home cage.
- 2.5. Provide mice with one cup of sucralose gel solution containing carprofen (1 mg/cup) for pre-dosing of anti-inflammatory agent after recovery. One cup can provide up to five mice with carprofen for up to two days. Check the cup daily and replace as needed, or every other day if not fully consumed.

3. Intraductal injection

NOTE: This step will occur 2–3 days after preoperative preparation.

- 3.1. Prepare 333.3 mM tantalum oxide (TaO_x) stock solution from procured powder formulation as described¹⁶ using phosphate buffered saline (PBS). Use gentle warmth if powder does not fully dissolve into solution.
- 3.2. Make ablative imaging solution by mixing 3 parts stock TaO_x with 7 parts 100% EtOH for a final 70% EtOH 100mM TaO_x solution. Add 1% v/v blue food dye to final solution for aided visualization during injection. Prepare volume based on need for experiment.

NOTE: Gland pairs 1 and 5 can hold up to 30 μL of the solution while all other pairs can be injected with up to 50 μL in 9 weeks of age or older FVB mice.

- 3.3. Anesthetize the mouse with isoflurane as in preoperative preparation and apply eye lubricant. Secure the mouse beneath the stereoscope using tape.
- 3.4. Prepare nipples for the injection by removing any dead skin that covers the nipple opening with fine pointed forceps.

- 3.5. Load 21–51 μL of the injection solution into a 50 μL syringe with a 34 G needle attached. Load an extra 1 μL of the solution to account for the potential leakage of that volume after the needle is removed from the nipple.

NOTE: Volumes above are for the typical procedures presented here. One may inject any volume desired with the understanding that most glands will be overfilled if going above 30 μL or 50 μL in gland pairs 1 and 5, or 2–4, respectively.

- 3.6. Hold the nipple gently with tweezers and insert the needle into the nipple opening until the bevel is fully covered. It may be necessary to pull up the nipple onto the needle rather than pushing the needle into the nipple (Table 1).
- 3.7. Inject the solution at a steady rate of approximately 40 $\mu\text{L}/\text{min}$. Avoid injecting too quickly to ensure that the ductal tree is not damaged. Keep the needle in the nipple for 30 s after the volume is completely injected before removing the needle. This is done to avoid the solution spilling out of the nipple (Figure 2).
- 3.8. While the animal is still anesthetized, inject intraperitoneally 200–250 μL of 5% sucrose solution (up to 10 mL/kg) to minimize the potential effects of alcohol intoxication.

4. micro-CT Imaging

- 4.1. After injecting all desired glands, move the animal swiftly to the micro-CT system and continue to maintain anesthesia using the incorporated isoflurane vaporizer. The animal may be taped to standardize imaging position. For example, taping each hind leg in an extended position helps to keep the leg bones of the animal further from the lower glands of interest in the resultant image.

NOTE: The animals can be imaged using different scan parameters if care is taken to determine an appropriate acceptable lifetime dose of radiation for the animal and the cumulative doses do not exceed this level. Fluoroscopy stills and videos can be generated without performing scans to further reduce radiation exposure (Figure 2).

- 4.2. Position the mouse in appropriate field of view with the fluoroscopy preview function.
- 4.3. Perform TaO_x imaging of the mouse ductal tree with good resolution and opportunity for repeated standard (2 min) acquisition scans. Use the following scan parameters: 90 kVp/88 μA ; field of view (FOV), 36 mm; number of slices, 512; slice thickness, 72 μm ; voxel resolution, 72 μm^3 .
- 4.3.1. Acquire longer (4 or 14 min) high resolution scans for even better resolution in animals that will not be scanned longitudinally using the same parameters.
- 4.4. Following imaging protocol completion, remove the animal from anesthesia and transfer to a separate clean, dry recovery cage on a heating pad. Observe the

animal until fully recovered and then return to home cage. Keep animals injected with ablative solution on carprofen until 7 days after injection.

- 4.5. Make quick renditions of any scanned images within the micro-CT software to better appreciate any contrast leaks or lack of filling (Figure 2).
- 4.6. Perform further formal image analysis for publication or detailed analysis of scans if desired (Figure 3).

NOTE: The primary difference between these methods will be the ability to threshold only the ductal tree (formal rendition) *versus* needing to threshold the entire image (quick rendition). Other measurements and images may be generated using the software packages to best demonstrate the success of ductal tree filling.

5. Image analysis

- 5.1. Make rendering of the injected ductal tree using a specialized software package.
 - 5.1.1 To segment the fat pad (darker compartment compared to peritoneal cavity, femoral muscles, and skin) within which the ductal tree of interest is contained, select spline trace option from the manual menu. Trace the fat pad outline at every third data slice.
 - 5.1.2 Click propagate objects option from the semi-automatic menu to connect all traced and untraced slices into a single object.
- 5.2. Select threshold volume option from the semi-automatic menu to input desired HU range (300–3000 HU is a good starting point) and click threshold rendition button to create a rendition that only displays the contrast (TaO_x) within the ductal tree and eliminates soft tissue.
- 5.3. Toggle view button to rendition as primary to view only the 3D rendition without surrounding structures.

NOTE: Additional software features allow for measurements to be made of the 3D rendition (i.e. length, volume, etc.).

REPRESENTATIVE RESULTS:

Female mice have five pairs of mammary glands with a single ductal tree that opens at the nipple orifice²². At the tips of the developing ductal tree are the terminal end buds (TEBs), proliferative structures that direct growth and branching. After puberty when elongation phase is completed, TEBs regress and become functionally and anatomically indistinguishable from terminal ducts or alveolar buds²³. Terminal ductal lobular units serve a similar function in humans as TEBs do in mice, and are the sites from which breast cancer predominantly arises^{24,25}. We can inject up to 50 μL of 70% EtOH solution to fill the entire ductal tree of thoracic and abdominal mammary glands of 9-week-old FVB/N, NSG, and other mouse strains (Figures 1–3, ref^{3,16}). In a typical experiment, we can inject up to 8 mammary glands with an ablative solution of 70% EtOH and 100 mM TaO_x in two consecutive ID procedures separated by 7 days to allow for animal recovery (Figure 2).

Animals are imaged by micro-CT immediately after the last ID injection to assess successful delivery of the solution to the entire ductal tree (Figure 2). In our experience, the nipples of inguinal glands are suitable for injection in about 60% of animals, and the nipples of cervical glands in about 40%. When suitable, we can inject up to 30 μL of 70% EtOH solution to fill the entire ductal tree of cervical and inguinal mammary glands (Figure 2). FVB and NSG strains generally present more suitable nipples for injection than C57BL/6J or mixed genetic background strains. Whole mount dual staining protocol or 3D confocal microscopy are good orthogonal methods to confirm to what extent the ductal tree was filled (Figure 3). These tissue correlate analyses are compatible with and can be performed after *in vivo* imaging. The obvious limitation of these orthogonal methods is that they require animal termination for tissue collection and analysis, but in an optimization phase of a new ablative formulation they provide independent validation.

DISCUSSION:

Prophylactic mastectomy is currently the most effective intervention for breast cancer, yet it has some serious negative impacts. Local ablation of mammary epithelial cells with an EtOH-based solution is a promising alternative therapy as we demonstrated in a proof-of-concept study on the aggressive FVB-Tg-C3(1)-TAg mouse model of breast cancer³. ID injection of this ablative solution allows for targeting of the mammary epithelial cells from which breast cancer arises with limited collateral damage. Addition of an X-ray contrast agent to the ablative solution allows for enhanced understanding of the effectiveness of the solution at prevention, as we can see whether each ductal tree is successfully filled following injection (Figure 2B). Viewing injected glands by fluoroscopy immediately post-injection mirrors what will likely be done in clinic to confirm successful filling of the ductal tree. Visual confirmation of solution delivery will best inform whether all parts of the tree have been reached in real time. This could allow for further injection to be performed to complete filling at the time or in a future session. It is of great importance that ablative solution reaches all parts of the ductal tree to ensure that all epithelial cells may be accessed for killing (Figure 3). Leaving behind living epithelial cells within the tree would allow for the possibility that breast cancer could still arise. Utilizing contrast in ID injections to image success of the injection could also be useful for other formulations. Other studies have described ID delivery protocols for viral particles (e.g., AdCre, CRISPR guide RNAs), hormones, cytotoxic compounds, siRNAs and/or targeting agents in mice^{3,16,27–32 33–38}, rats^{24,32,39–41}, and rabbits^{42–47}. Independent clinical studies reported successful cannulation of up to eight ducts per breast for local delivery of chemotherapy^{40,48,49}. Visualizing full filling when delivering other solutions aimed at prevention or geared towards treatment would be worthwhile for similar reasons. The knowledge that the solution has reached all branches and terminal ends of the tree will be informative in assessing successful prevention or treatment.

We are not aware of any other intraductal imaging methods in mice^{33,34} or other animal models⁴⁷ that afford the high resolution of TaO_x nanoparticles. Of note, TaO_x in the murine ductal tree outperforms FDA-approved contrast agents for diagnostic ductography^{3,16}. As we continue to assess the ID ablative procedure for its ability to prevent breast cancer, we will be able to determine more precisely from which glands cancer arises through the

added data given through imaging after ID delivery. For instance, one could determine if a gland that was only partially filled is more likely than a non-injected gland to result in tumor formation, which addresses the safety profile and concern of unsuccessful injections on a high-risk woman. This technique has some limitations. This is a relatively challenging mouse technique that requires dexterity and proficiency of the operator to manipulate and successfully cannulate each duct. Each individual injection is an independent event thus unsuccessful injection on one or more glands may compromise result interpretation. Given the size of the murine mammary gland and fragility of the nipple, fluoroscopy or similar image-guidance technique is not available to inform in real-time when to stop infusion. This real time image-guidance will be a requirement for clinical implementation of local delivery of an ablative solution.

ACKNOWLEDGMENTS:

This work was supported, in part, by National Cancer Institute R21 CA226579 and R01 CA258314 grants to LFS and by National Institute of Biomedical Imaging and Bioengineering R01 EB029418 grant to EMS. We would like to thank MSU Institute for Quantitative (IQ) Health Science and Engineering Imaging Core facility for use of their imaging systems and technical expertise. We would like to thank Dr. Danielle Ferguson for reviewing contents of video and figures for adherence to animal welfare guidelines.

REFERENCES:

1. Britt KL, Cuzick J & Phillips KA Key steps for effective breast cancer prevention. *Nat Rev Cancer*. 20 (8), 417–436, (2020). [PubMed: 32528185]
2. Padamsee TJ, Wills CE, Yee LD & Paskett ED Decision making for breast cancer prevention among women at elevated risk. *Breast Cancer Res*. 19 (1), 34, (2017). [PubMed: 28340626]
3. Kenyon E et al. Ductal tree ablation by local delivery of ethanol prevents tumor formation in an aggressive mouse model of breast cancer. *Breast Cancer Res*. 21 (1), 129, (2019). [PubMed: 31779648]
4. Kuang M et al. Ethanol ablation of hepatocellular carcinoma Up to 5.0 cm by using a multipronged injection needle with high-dose strategy. *Radiology*. 253 (2), 552–561, (2009). [PubMed: 19709992]
5. Ansari D & Andersson R Radiofrequency ablation or percutaneous ethanol injection for the treatment of liver tumors. *World J Gastroenterol*. 18 (10), 1003–1008, (2012). [PubMed: 22416173]
6. Zhang WY, Li ZS & Jin ZD Endoscopic ultrasound-guided ethanol ablation therapy for tumors. *World J Gastroenterol*. 19 (22), 3397–3403, (2013). [PubMed: 23801831]
7. Chin M, Chen CL, Chang K, Lee J & Samarasena J Ethanol Ablation of a Peripheral Nerve Sheath Tumor Presenting as a Small Bowel Obstruction. *ACG Case Rep J*. 3 (1), 31–32, (2015). [PubMed: 26504873]
8. Gueng M-K, Chou Y-H, Tiu C-M, Chiou S-Y & Cheng Y-F Pseudoaneurysm of the Breast Treated with Percutaneous Ethanol Injection. *Journal of Medical Ultrasound*. 22 (2), 114–116.
9. Zhang J et al. Comparison between absolute ethanol and bleomycin for the treatment of venous malformation in children. *Exp Ther Med*. 6 (2), 305–309, (2013). [PubMed: 24137179]
10. Wohlgemuth WA et al. Ethanolgel sclerotherapy of venous malformations improves health-related quality-of-life in adults and children - results of a prospective study. *Eur Radiol*. 27 (6), 2482–2488, (2017). [PubMed: 27699470]
11. Steiner F, FitzJohn T & Tan ST Ethanol sclerotherapy for venous malformation. *ANZ J Surg*. 86 (10), 790–795, (2016). [PubMed: 25182969]
12. Sannier K et al. A new sclerosing agent in the treatment of venous malformations. Study on 23 cases. *Interv Neuroradiol*. 10 (2), 113–127, (2004).
13. Domp Martin A et al. Radio-opaque ethylcellulose-ethanol is a safe and efficient sclerosing agent for venous malformations. *Eur Radiol*. 21 (12), 2647–2656, (2011). [PubMed: 21822948]

14. Slawson SH & Johnson BA Ductography: how to and what if? *Radiographics*. 21 (1), 133–150, (2001). [PubMed: 11158649]
15. Sheiman LS & Levesque PH The In's and Out's of Ductography: A Comprehensive Review. *Curr Probl Diagn Radiol*. 45 (1), 61–70, (2016). [PubMed: 26163736]
16. Chakravarty S et al. Tantalum oxide nanoparticles as versatile contrast agents for X-ray computed tomography. *Nanoscale*. 12 (14), 7720–7734, (2020). [PubMed: 32211669]
17. Krause S, Brock A & Ingber DE Intraductal injection for localized drug delivery to the mouse mammary gland. *J Vis Exp*. 80 (issue) (80), (2013).
18. King BL & Love SM The intraductal approach to the breast: raison d'être. *Breast Cancer Res*. 8 (2), 206, (2006). [PubMed: 16677416]
19. Love SM & Barsky SH Anatomy of the nipple and breast ducts revisited. *Cancer*. 101 (9), 1947–1957, (2004). [PubMed: 15382093]
20. Morhard R et al. Development of enhanced ethanol ablation as an alternative to surgery in treatment of superficial solid tumors. *Sci Rep*. 7 (1), 8750, (2017). [PubMed: 28821832]
21. Morhard R et al. Understanding Factors Governing Distribution Volume of Ethyl Cellulose-Ethanol to Optimize Ablative Therapy in the Liver. *IEEE Trans Biomed Eng*. 67 (8), 2337–2348, (2020). [PubMed: 31841399]
22. Hinck L & Silberstein GB Key stages in mammary gland development: the mammary end bud as a motile organ. *Breast Cancer Res*. 7 (6), 245–251, (2005). [PubMed: 16280048]
23. Paine IS & Lewis MT The Terminal End Bud: the Little Engine that Could. *J Mammary Gland Biol Neoplasia*. 22 (2), 93–108, (2017). [PubMed: 28168376]
24. Sivaraman L et al. Effect of selective ablation of proliferating mammary epithelial cells on MNU induced rat mammary tumorigenesis. *Breast Cancer Res Treat*. 73 (1), 75–83, (2002). [PubMed: 12083633]
25. Cardiff RD & Wellings SR The comparative pathology of human and mouse mammary glands. *J Mammary Gland Biol Neoplasia*. 4 (1), 105–122, (1999). [PubMed: 10219910]
26. Arora R et al. Insights from imaging the implanting embryo and the uterine environment in three dimensions. *Development*. 143 (24), 4749–4754, (2016). [PubMed: 27836961]
27. Brock A et al. Silencing HoxA1 by intraductal injection of siRNA lipidoid nanoparticles prevents mammary tumor progression in mice. *Sci Transl Med*. 6 (217), 217ra212, (2014).
28. de Groot JS et al. Intraductal cisplatin treatment in a BRCA-associated breast cancer mouse model attenuates tumor development but leads to systemic tumors in aged female mice. *Oncotarget*. 8 (37), 60750–60763, (2017). [PubMed: 28977823]
29. Wang G et al. Intraductal fulvestrant for Therapy of ERalpha-positive Ductal Carcinoma in Situ (DCIS) of the breast- a Preclinical Study. *Carcinogenesis*. 10.1093/carcin/bgz084, (2019).
30. Yoshida T et al. Effective treatment of ductal carcinoma in situ with a HER-2- targeted alpha-particle emitting radionuclide in a preclinical model of human breast cancer. *Oncotarget*. 7 (22), 33306–33315, (2016). [PubMed: 27119227]
31. Chun YS et al. Intraductally administered pegylated liposomal doxorubicin reduces mammary stem cell function in the mammary gland but in the long term, induces malignant tumors. *Breast Cancer Res Treat*. 135 (1), 201–208, (2012). [PubMed: 22752247]
32. Murata S et al. Ductal access for prevention and therapy of mammary tumors. *Cancer Res*. 66 (2), 638–645, (2006). [PubMed: 16423990]
33. Markiewicz E et al. High resolution 3D MRI of mouse mammary glands with intraductal injection of contrast media. *Magn Reson Imaging*. 33 (1), 161–165, (2015). [PubMed: 25179139]
34. Markiewicz E et al. MRI ductography of contrast agent distribution and leakage in normal mouse mammary ducts and ducts with in situ cancer. *Magn Reson Imaging*. 40 48–52, (2017). [PubMed: 28366759]
35. Annunziato S et al. Comparative oncogenomics identifies combinations of driver genes and drug targets in BRCA1-mutated breast cancer. *Nat Commun*. 10 (1), 397, (2019). [PubMed: 30674894]
36. Rutkowski MR et al. Initiation of metastatic breast carcinoma by targeting of the ductal epithelium with adenovirus-cre: a novel transgenic mouse model of breast cancer. *J Vis Exp*. 10.3791/51171 (85), (2014).

37. Xiang D, Tao L & Li Z Modeling Breast Cancer via an Intraductal Injection of Cre-expressing Adenovirus into the Mouse Mammary Gland. *J Vis Exp*. 10.3791/59502 (148), (2019).
38. Barham W, Sherrill T, Connelly L, Blackwell TS & Yull FE Intraductal injection of LPS as a mouse model of mastitis: signaling visualized via an NF-kappaB reporter transgenic. *J Vis Exp*. 10.3791/4030 (67), e4030, (2012). [PubMed: 22971993]
39. Chun YS et al. Intraductal administration of a polymeric nanoparticle formulation of curcumin (NanoCurc) significantly attenuates incidence of mammary tumors in a rodent chemical carcinogenesis model: Implications for breast cancer chemoprevention in at-risk populations. *Carcinogenesis*. 33 (11), 2242–2249, (2012). [PubMed: 22831956]
40. Stearns V et al. Preclinical and clinical evaluation of intraductally administered agents in early breast cancer. *Sci Transl Med*. 3 (106), 106ra108, (2011).
41. Okugawa H et al. Effect of periductal paclitaxel exposure on the development of MNU-induced mammary carcinoma in female S-D rats. *Breast Cancer Res Treat*. 91 (1), 29–34, (2005). [PubMed: 15868429]
42. Falconer IR The distribution of 131 I- or 125 I-labelled prolactin in rabbit mammary tissue after intravenous or intraductal injection. *J Endocrinol*. 53 (3), 58–59, (1972). [PubMed: 5039234]
43. Fiddler TJ, Birkinshaw M & Falconer IR Effects of intraductal prolactin on some aspects of the ultrastructure and biochemistry of mammary tissue in the pseudopregnant rabbit. *J Endocrinol*. 49 (3), 459–469, (1971). [PubMed: 5104551]
44. Fiddler TJ & Falconer IR The effect of intraductal prolactin on protein and nucleic acid biosynthesis in the rabbit mammary gland. *Biochem J*. 115 (5), 58P–59P, (1969).
45. Bourne RA, Bryant JA & Falconer IR Stimulation of DNA synthesis by prolactin in rabbit mammary tissue. *J Cell Sci*. 14 (1), 105–111, (1974). [PubMed: 4856240]
46. Chadwick A Detection and Assay of Prolactin by the Local Lactogenic Response in the Rabbit. *J Endocrinol*. 27 253–263, (1963). [PubMed: 14079513]
47. Clark A, Bird NK & Brock A Intraductal Delivery to the Rabbit Mammary Gland. *J Vis Exp*. 10.3791/55209 (121), (2017).
48. Mahoney ME et al. Intraductal therapy of ductal carcinoma in situ: a presurgery study. *Clin Breast Cancer*. 13 (4), 280–286, (2013). [PubMed: 23664819]
49. Love SM et al. A feasibility study of the intraductal administration of chemotherapy. *Cancer Prev Res (Phila)*. 6 (1), 51–58, (2013). [PubMed: 23169924]




Figure 1. Workflow of intraductal procedure and image analysis.
Key steps of the ID procedure are highlighted. Please see video for more details.




Figure 2. Successful cannulation and delivery of ablative solution to multiple mammary glands. **A)** Representative nipple variegation in FVB and NSG mouse strains. Long nipples are easier to cannulate than short nipples, whereas too short or vestigial nipples cannot be cannulated. Once cannulated, the size of the nipple does not affect successful intraductal delivery. **B)** Gross anatomical analysis of blue dye in ablative solution provides *ex vivo* evidence of ductal tree filling and delivery success. **C and D)** Real-time fluoroscopy and post-image acquisition 3D micro-CT rendition provide *in vivo* evidence of delivery success. **D)** Successful injection of both abdominal and inguinal glands, and three out of four thoracic glands (fat pad injection in gland #2). Scale bars correspond to 1 mm in images at different magnification.




Figure 3. *In vivo* and *ex vivo* demonstration of ductal tree filling.

70% EtOH/100 mM TaO_x nanoparticles/Evans Blue solution was intraductally injected into the mouse abdominal mammary gland and immediately imaged by micro-CT and processed for dual whole mount staining. Ductal tree was reconstructed using an image analysis software package. The solution entirely fills the carmine alum-stained ductal tree. A separate gland was immunostained for ECadherin (*Cdh1*), cleared using Benzyl Alcohol:Benzyol Benzoate and imaged by confocal microscopy as described²⁶. Ductal tree was reconstructed using image analysis software. Pseudo-coloring rendering of confocal image (i.e., black background to white, green marker signal to magenta) was obtained with image invert function of image editing software. Scale bars correspond to 1 mm in images at different magnification. This figure has been modified from reference³.

Table 1.

Troubleshooting and helpful tips.

Issue	Appearance	Solution
Short nipple (Fig. 2)	Nipple has low profile – hard to grab	It is sometimes easier to hold the skin near the nipple and target the center of the nipple with the needle. The needle will likely dive under the skin. Pulling up slowly may reveal the nipple to be slightly over the tip of the needle and give room to grab and pull it the rest of the way onto the needle. Be very careful when diving below the skin about the angle of the needle. It is easy to inadvertently get a fat pad injection by stabbing at the wrong angle.
Fat nipple	Much larger than other nipples with little peelable dead skin – easily visible without scope	Very easy to get a fat pad injection on these nipples. Be very cautious about angle of the needle when inserting into nipple.
Fat pad injection (Fig. 2)	Swollen around nipple and possibly in nipple itself – easiest to see if color is added to injection solution	If nipple is swelling with first few ul injected, remove needle, and attempt to insert again with more care taken of angle. Begin injection again and watch for further swelling. If swelling continues, abandon attempt. It is very rare to successfully inject a nipple that has started out as a fat pad injection.
Wounds/ scabbing	Open wound or scabbing near injection site of EtOH solution	Apply triple antibiotic ointment to open wounds but leave scabbed wounds alone. Applying ointment to scabs can increase likelihood animal will bother the scab and remove it. Check every 1–2 days until healed depending on severity of wound. Carprofen should be given until healed even if beyond normal window.

Author Manuscript

Author Manuscript

Author Manuscript

Author Manuscript

Materials

Name	Company	Catalog Number	Comments
AnalyzeDirect v12.0	Caliper	n/a	For micro-CT image processing
Carprieve, Carprofen 50 mg/mL	Allivet	50647	For anti-inflammatory treatment
Evans blue	Sigma	E2129-50G	For injection visualization
Hot water bath	Toolots	Yidu_HH-S2	For preparing carprofen cups
Imaris	Bitplane	n/a	For confocal image processing
MediGel Sucralose Cups	ClearH2O	74-02-5022	For delivery of carprofen
Model 1705 RN Syringe, 50 μ L	Hamilton	7655-01	For intraductal injection
Photoshop 2021	Adobe	n/a	For image processing
Quantum GX2 microCT Imaging System	Perkin Elmer	CLS149276	For micro-CT image acquisition
Small Hub RN Needle, 34 gauge, custom (12° bevel angle, 0.375 in, point style 4)	Hamilton	207434	For intraductal injection
Stereo Microscope SZM Series	AmScope	SM-4TPZ-144	For intraductal injection
Sterile blue food dye	McCormick	930641	For injection visualization
Sterile phosphate buffered saline (PBS)	ThermoFisher	14190250	For solution preparation
Stickers	DOT Scientific	DOTSCI-C50	For preparing carprofen cups
Sucrose	Calbiochem	8550-5KG	For intraductal injection
Syringes	Fisher	14-826-79	For preparing carprofen cups
Vortex	VWR	10153-834	For preparing carprofen cups
Warming pump/pad(s)	Braintree Scientific	HTP-1500 120V; AP-R 26E	For intraductal injection/preoperative preparation