scientific reports

Check for updates

Instant quantifcation of sugars OPEN in milk tablets using near‑infrared spectroscopy and chemometric tools

Chanat Thanavanich¹, Nutthatida Phuangsaijai¹, Chanidapha Thiraphatchotiphum¹, **ParichatTheanjumpol2,3 & Sila Kittiwachana1***

Milk tablets are a popular dairy product in many Asian countries. This research aimed to develop an instant and rapid method for determining sucrose and lactose contents in milk tablets using nearinfrared (NIR) spectroscopy. For the quantitative analysis, a training set composed of laboratory-scale milk samples was generated based on a central composite design (CCD) and used to establish partial least squares (PLS) regression for the predictions of sucrose and lactose contents resulting in R² values **of 0.9749 and 0.9987 with the corresponding root mean square error of calibration (RMSEC) values of 1.69 and 0.35. However, the physical diference between the laboratory-scale powder and the fnal product milk tablet samples resulted in spectral deviations that dramatically afected the predictive performance of the PLS models. Therefore, calibration transfer methods called direct standardization (DS) and piecewise direct standardization (PDS) were used to adjust the NIR spectra from the real milk tablet samples before the quantitative prediction. Using high-performance liquid chromatography (HPLC) as a reference method, the developed NIR-chemometric model could be used to instantly predict the sugar contents in real milk tablets by producing root mean square error of prediction** (RMSEP) values for sucrose and lactose of 5.04 and 4.22 with Q^2 values of 0.7973 and 0.9411, **respectively, after the PDS transformation.**

Milk tablets are considered an alternative nutritious snack. This snack not only delivers a pleasant taste of milk but also is a good source of high-quality nutrients such as protein, carbohydrates, and calcium. Along with milk, sugars are among the major ingredients of milk tablets. Often, sucrose is added to sweeten the milk tablets¹. Lactose is the main carbohydrate found in raw milk material². Determination of sugar contents is important for quality control evaluation. In addition, variations in the amount of sugars can infuence the nutritional quality of the finished product contributing to the eating texture of the milk snack^{3,[4](#page-7-3)}.

Traditionally, the quantification of sugar in milk can be achieved using chemical titration methods^{[5](#page-7-4)}. These wet laboratory methods can be simple; however, the experimental procedures are time-consuming and ofen result in large amounts of chemical waste. Several analytical techniques, such as high-performance liquid chromatography $(HPLC)^{6,7}$ $(HPLC)^{6,7}$ $(HPLC)^{6,7}$, electrochemical analysis 8,9 8,9 8,9 8,9 , and nuclear magnetic resonance $(NMR)^{10}$, were also used to accomplish the quantifcation task. Nevertheless, a signifcant drawback was that complicated sample preparation methods were required before the sample detection process. Therefore, these tests are not suitable for manufacturing plants where many samples are obtained from a continuous online production process.

Near-infrared (NIR) spectroscopy investigates the vibrational interaction between samples and electromagnetic radiation in a region of 800–2500 nm. NIR detection has several advantages over other spectroscopic detections in that a large number of samples can be non-destructively measured within a short period without extensive sample preparation. The NIR spectra can be related to interesting chemical properties in the samples using a calibration model called partial least squares (PLS) regression^{11-[13](#page-7-11)}, where the relationship information between the spectral data and the chemical property is mathematically investigated. Tis correlation information can then be used to estimate the chemical properties of unknown samples. For example, the NIR spectrometers were recently applied for compositional analysis of row^{14} and human¹⁵ milk.

1 Department of Chemistry, Faculty of Science, Chiang Mai University, Chiang Mai 50200, Thailand. ²Postharvest Technology Research Center, Faculty of Agriculture, Chiang Mai University, Chiang Mai 50200, Thailand. ³Postharvest Technology Innovation Center, Ministry of Higher Education, Science, Research and Innovation, Bangkok 10400, Thailand. [⊠]email: silacmu@gmail.com

	Concentration of sugars $(\%w/w)$		
Sample names	Sucrose	Lactose	
1. T1 (training)	21.31	36.24	
2. T ₂ (training)	27.11	28.66	
3. T3 (training)	14.20	26.46	
4. I1 (validation)	27.76	25.71	
5. I2 (validation)	27.00	33.40	
6. I3 (validation)	22.31	28.48	
7. E1 (test)	21.78	38.99	
8. E2 (test)	39.04	0.000	
9. E3 (test)	43.69	0.000	
10. E4 (test)	28.04	28.47	
11. E5 (test)	20.06	30.68	
12. E6 (test)	16.72	29.32	
13. E7 (test)	25.48	28.94	

Table 1. Concentrations of sucrose and lactose in milk tablet samples.

Generally, determinations achieved by calibration models are based on the previous knowledge provided by a set of training samples. Accordingly, variations of the test samples could be estimated by those found in the training set to get the optimal predictive results. Hence, the training samples can signifcantly characterize the predictive performance of the calibration models. However, in the production process where variations of the fnal process samples are not expected and should be limited, the fnal products from the manufacturing process may not contain enough variability to establish a training set for developing accurate and robust calibration models. For example, the training samples in powder form that were prepared in a laboratory room had diferent physical structures when compared to the fnal process samples that were eventually pressed into tablets. Sarraguça and Lopes¹⁶ reported that using laboratory-scale samples in the powder form provided more accurate predictive results than the tablet samples produced by the production process. Additionally, Peerapattana et al.^{[17](#page-8-1)} observed that powder samples of mangosteen pericarp retained in glass vials resulted in better predictive results for the prediction of alpha-mangostin content compared to the samples retained in transparent capsules.

Calibration transfer (CT) involves a group of chemometric methods that can be used to minimize inconsistencies from different instrumental measurements¹⁸. After establishing the standardization using the CT methods, a system involving the calibrated model can be developed in which one instrument (primary) can be substituted for another instrument (secondary) and vice versa. For instance, a calibration transfer method called piecewise direct standardization (PDS) was used to minimize the variations between a top-bench spectrometer and a portable NIR detector^{[19,](#page-8-3)20}. As a result, the estimation of the portable NIR measurement can be based on the data-based NIR spectra collected from the top-bench NIR spectrometer. In addition, by following the same methodology, diferent variations among agricultural product samples could be identifed and systematically discarded from unknown samples for adulteration detection^{[21](#page-8-5)}.

Tis research developed an instant method for the detection of the sugar contents in milk tablets based on the NIR detection. In addition, calibration transfer methods were adopted to adjust the diferent variations in the NIR spectra between the laboratory-scale milk samples and the milk tablets obtained from the fnal process. The quantitative analysis was based on multivariate predictions of the PLS models aiming to instantly quantify the concentrations of sucrose and lactose in the milk tablet samples.

Materials and methods

Details of milk samples. A total of 13 diferent milk tablet brands were obtained from local grocery stores in Chiang Mai, Thailand. The relevant details of the milk tablets are summarized in Table [1](#page-1-0). The samples were divided into three groups, namely training (T1–T3), internal validation (I1–I3), and external validation (E1–E7) samples. Each milk tablet was ground into a fne and homogeneous powder using a ceramic mortar and pestle. To generate systematic variations representing the sugar contents in the milk samples, a central composite design (CCD) structure was used^{[22](#page-8-6)}, comprising nine experiments for each sample. For example, for samples T1, T2, and T3, amounts of sucrose (analytical grade,>99% purity, RCI Labscan, Bangkok, Tailand) and lactose (analytical grade,>99% purity, KEMAUS, NSW, Australia) were added to the milk powder according to the coded values of the CCD structure presented in Table [2.](#page-2-0) Then, a combination of the three CCD model samples was used to construct the training set, which resulted in a total of 27 milk samples. The use of the CCD structure was to ensure that the variation in the recorded NIR spectra was related to the concentrations of the sugars in the milk samples and that the number of training samples was sufficient for establishing the prediction models²³. Samples I1, I2, and I3 were used to establish internal validation samples, while the variations in the sugar contents were also generated based on the CCD model. Therefore, 27 additional milk powder samples were used to construct the internal validation set. Samples E1–E7 were used as the external samples to represent the independent test set. These were utilized to evaluate the performance of the calibration models when real samples were introduced.

2

	Code values		Amount of sugar added (%w/w)	
No. of samples	Sucrose	Lactose	Sucrose	Lactose
1	-1	-1	0.000	0.000
\overline{c}	-1	Ω	0.000	15.00
3	-1	1	0.000	30.00
$\overline{4}$	0	-1	15.00	0.000
5	θ	θ	15.00	15.00
6	θ	1	15.00	30.00
7	1	-1	30.00	0.000
8	1	Ω	30.00	15.00
9	1		30.00	30.00

Table 2. The CCD structure for generating the training samples.

It should be noted that two main types of milk tablets were used in this research. Samples E2 and E3 were non-milk-containing tablets or "cheap milk tablets" where artifcial milk favor was added to achieve product satisfaction. On the other hand, the rest of the milk samples were produced from cow's milk as raw material and were referred to as "premium milk tablets".

NIR spectral detection. The NIR spectra of the milk powder (9.00 g) were acquired using a NIR transportation module (width×length×depth: 5.7×29.4×2.0 cm) equipped with the NIRSystem 6500 (Multi-Mode™ Analyzer, Foss, USA) in the range of 400–2500 nm at a 2 nm sampling interval, yielding 1050 data points per spectrum. An average of 64 scans was used for each sample. The samples of the milk tablets were placed inside the NIR transportation module. Layers of the milk tablets were directly attached to the containing glass following the measurement conditions of the powder samples. The milk samples were maintained at a room-controlled temperature of 25 °C for at least 6 h before the NIR detection. Prior to the analysis, the NIR spectra were pretreated by standard normal variate (SNV) to eliminate errors caused by the light scattering during the NIR measurement. Then, they were mean-centered so that the analysis focused on the variance from the data mean rather than the absolute values.

HPLC analysis of sugar determination. Sugar contents in the milk tablet samples were measured using high-performance liquid chromatography (HPLC). For sample preparation, 1.00 g of each ground milk tablet was dissolved in 10 mL ultra-pure water and kept in a water bath (Julabo Labortechnik GMbH, Seelbach, Germany) at 55 °C for 5 min. Then, HPLC-grade acetonitrile was added for protein precipitation^{24,25}. After the denaturation, the sample solution was centrifuged at 10,000 rpm for 5 min. The clear solution was then filtered through a 0.45 µm nylon syringe flter (Agilent Technologies, CA, USA).

The chromatographic analysis of sugar contents in the milk tablets was carried out with a high-performance liquid chromatograph (Agilent 1100 HPLC system, CA, USA) with an Agilent ZORBAX NH₂ column (5 µm, 4.6 mm inner diameter, 150 mm length) operating at 25 °C. Te samples were auto-injected into the HPLC system with an injection volume of 10 µL. A mixture of HPLC-grade acetonitrile and ultra-pure water (75/25%v/v) was used as a mobile phase with a fow rate of 1.00 mL/min. A refractive index detector (RID) was operated at 25 °C. The sugar contents were determined using the external standard calibration curve of the sucrose and lactose standards resulting in R^2 values of 0.9907 and 0.9896, respectively. The concentration values of the sugars in the studied milk tablet samples are summarized in Table [1.](#page-1-0)

Chemometric analysis. *Standardization of NIR spectra using DS and PDS calibration transfers.* Although both forms of the milk samples (tablet and powder) were considered solid, there were diferences, for example, in particle size and tablet compaction pressure. Tese physical variations resulted in signifcant deviations in the recorded NIR spectra[26](#page-8-10). Calibration transfers are multivariate correction methods that can be applied to stabilize variations that may have occurred due to difering instrumental and measurement conditions. In this research, they were used to account for any signal discrepancies between the spectra obtained from the tablet samples and the powder samples. Piecewise direct standardization (PDS) is an extension algorithm of a conventional method called direct standardization $DS^{27,28}$ $DS^{27,28}$ $DS^{27,28}$ $DS^{27,28}$ $DS^{27,28}$. The DS method describes the correlation between the two data matrices (*Xm* and *Xs* referring to master and slave data) by calculating a transformation matrix (*F*) using multiple linear regression models such as MLR, PCR, and PLS:

$$
X_m=X_s\times F
$$

The extension in the PDS algorithm is that each spectral point of the master data $(X_{m,j})$ is specifically related to a spectral subset of the slave data $(X_{s,j})$. The PDS algorithm involves the following steps:

Step 1: Select the spectral points of the master data $(X_{m,j})$ at wavelength *j*. Step 2: Defne the subset spectra of the slave data (*Xs,j*) nearby wavelength *j* to form index *j*−*k* to *j*+*k*

$$
\mathbf{X}_{s,j} = [\mathbf{x}_{s,j-k} \cdot \mathbf{x}_{s,j-k+1}, \ldots, \mathbf{x}_{s,j+k-1} \cdot \mathbf{x}_{s,j+k}]
$$

where *k* is the window size controlling the amount of the spectral data which will be used in the calculation. Step 3: Establish the regression coefficient

$$
X_{m,j}=X_{s,j}\times\boldsymbol{b}_j
$$

where $\bm{b}_{\rm i}$ is a vector containing regression coefficients. Step 4: Generate the transformation matrix (*F*) by organizing the *b^j* into a diagonal matrix

$$
\mathbf{F} = \text{diag}(\mathbf{b}_1^T, \mathbf{b}_2^T, \dots \mathbf{b}_j^T, \dots \mathbf{b}_n^T)
$$

where *n* is the number of spectral channels included. Step 5: Standardize the spectra of unknown samples (X_{sum}) using **F** to obtain the modified spectrum (X_{appS})

 $X_{s,PDS} = F \times X_{s,un}$

In this research, the DS and PDS transformations were used to account for the inconsistencies between the spectra obtained from the powder samples and the tablet samples. Tese transformation methods investigated the correlation between the two datasets. Afer that, the resulting correlation information was applied to adjust the NIR spectra of the milk tablet samples. Consequently, the adjusted data could be compatible when making the prediction using the calibration model established from the NIR spectra of the powder samples without the need to recalibrate the model.

The model optimization was based on a previously published report²¹. The correlation matrices in both DS and PLS were determined using PLS regression which was calculated using the training samples and optimized based on the internal validation samples.

PLS for quantitative analyses. Partial least squares (PLS) regression is among the most potent analysis meth-ods of multivariate calibration models^{[29](#page-8-13)}. The significant advantage of the PLS algorithm is that the variations obtained from both the predictive and response parameters are simultaneously extracted and then used to construct the prediction model. With the use of the PLS model, the correlation between these information blocks could be maximized. In most cases, PLS could successfully offer the optimal predictive performance for the prediction of the NIR spectral data $11,30$ $11,30$ $11,30$.

In this research, the NIR spectra and the sugar contents were, respectively, used as predictive and response parameters for the PLS models. The PLS calculation was done following the procedure described in the previously published literature²⁹. The leave-one-out cross-validation method was applied to identify the optimal number of PLS latent variables³¹. According to Table [1,](#page-1-0) PLS models were developed using training $(T1-T3)$ samples as calibration data. To validate the models, internal validation (I1–I3) and external validation (E1–E7) samples were used for validation and prediction, respectively.

The predictive performance of the PLS models in terms of prediction accuracy was reported by root mean square error of calibration (RMSEC) and root mean square error of prediction (RMSEP). The coefficients of determination for the calibration (R^2) and prediction (Q^2) values were calculated to determine the robustness of the models. In addition, the standard error of the cross-validation (SECV) and the ratio of prediction to devia-tion (RPD) were employed to compare the different predictive performances of the calibration models^{[32](#page-8-16)}. The calculations of the PLS model, PDS calibration transfer, and statistical analyses were implemented using in-house MATLAB scripts (MATLAB, The Math Works Inc., Natick).

Results and discussion

NIR detection. *Exploratory data analysis of NIR spectra data.* Figure [1A](#page-4-0) illustrates the NIR spectra recorded from the milk samples. The corresponding PCA score plot of the NIR spectra was generated to demonstrate the characteristic diferences among the milk samples presented in Fig. [1b](#page-4-0). In this research, the training and validation samples were generated based on the CCD experiments to induce systematic variations due to the sugar concentrations. In the PCA score plot, the samples were scattered across the PCA space, wherein a greater PC1 value represented the sample with high sucrose content. On the other hand, higher PC2 values were associated with the samples with high lactose content.

Consequently, it was clear that the sample organizations were achieved due to the variations in the sugar contents. This demonstrated that the NIR detection could analyze the milk samples with different sugar contents. Most of the test samples were placed within the region of the training samples implying the suitability of the generated training samples. Test sample no. 10 (E10) was located slightly away from the main cluster; however, this could be due to the extracted malt, which was utilized as a favoring agent and resulted in yellow–brown double-layer tablets.

NIR spectral transformation of milk tablet samples. Figure [2](#page-4-1) presents the diference in the NIR spectra between the powder and tablet milk samples. The NIR spectra of milk tablets (E1-E7) were recorded before and after grinding into powder. In Fig. [2](#page-4-1)a, the shapes of the NIR spectra obtained from both physical states were relatively similar. However, the powered milk samples resulted in relatively lower levels of absorbance, which could be related to the particle size of the samples. Tis outcome corresponded to the fnding of a previous report noting that tablet compaction pressure resulted in more intense penetration with higher absorbance intensity of the

4

Figure 1. (**a**) NIR spectra of the studied milk samples (powder) and (**b**) PCA score plot of the corresponding NIR spectral data.

Figure 2. (**a**) NIR spectra of the test milk samples (powder and tablet) and (**b**) PCA score plot of the milk samples after the DS and PDS transformations.

NIR spectra²⁶. The variation that occurred because of the sample's physical conditions could be confirmed in the PCA score plot shown in Fig. [2](#page-4-1)b. From the score plot, the milk samples without any transformation process were separately clustered into two main groups, wherein the tablet and powder samples were placed diferently on the PCA space. The clear separation between the sample clusters implied that the physical differences caused more variation than chemical compositional diferences between milk samples. In other words, in this experiment with the PCA visualization, the physical variation identifed between the powder and the tablet samples could be systematically captured by the frst two PCs.

Afer the spectral data recorded from the tablet samples were transformed using the PDS transformation process, the shape of the spectra was signifcantly changed and became very similar to that of the powder sam-ples, as shown in Fig. [2](#page-4-1)a. The success of the spectral transformation can be confirmed by the PCA score plots shown in Fig. [2](#page-4-1)b, in which the transformed samples were placed nearly exactly in the identical PCA locations as the powder samples. In this case, PDS transformation could efectively remove the variation that occurred due to the physical state of the samples. Commonly, PDS transformation is used to adjust for any diferences in the experimental conditions; for example, the diference observed in the sample detection methods using two diferent instruments, namely primary and secondary. The correlation between the instruments was investigated based on the PLS-NIPAL algorithm during the modeling process³³. As a result, systematic variations in the spectral data obtained from both instruments were extracted concurrently, while also maximizing their covariance values. This correlation information could be used to transform the detection signal from the secondary instrument with respect to the covariation structure of both instruments. Consequently, the transformed signals could be well-suited to the prediction models established from the primary instrument without requiring the prediction model to be regenerated. In this experiment, the reconstruction of the transformed NIR spectra was based on

Model statistics	Sucrose	Lactose
R ²	0.9749	0.9987
Q^2	0.9373	0.9943
RMSEC	1.69	0.35
RMSEP	2.67	0.79
SECV	3.28	5.02
RPD	4.09	12.48

Table 3. Statistical values for the PLS prediction using the training and internal validation milk samples.

the systematic structure, which could be efectively captured by the PCA modeling, as previously discussed. Therefore, PDS could be successfully adapted to the translation of the NIR spectral data of the tablet samples. These results indicated the effective utilization of the PDS transformation model for stabilizing the NIR spectra recorded from samples with diferent physical appearances.

In this research, the DS adjustment which was based on a single calculation of the correlation matrix produced unsatisfactory estimates of the NIR spectra. However, the visualization of this spectral transformation was based on the use of raw NIR spectral data without data pretreatment. Tis global adjustment using DS could be sensitive to the drif in the spectral baseline as visualized in Fig. [2a](#page-4-1) and this demonstrated that the calculation based on the spectral subsection could improve the performance of the calibration transfer process. The effect of this data adjustment on the prediction performance will be discussed in the prediction of the sugar content using the PLS model.

Detection of sugar contents using PLS. *Prediction of internal validation samples.* The NIR spectra of the training samples were used to establish PLS calibration models to estimate the sugar contents in the milk samples. Using the PLS1 algorithm, the calibration model was constructed independently to predict each type of sugar. The validation samples were used to investigate the predictive performance of the developed PLS models, and the predictive results were summarized in Table [3.](#page-5-0)

The PLS model for the prediction of the sucrose content resulted in calibration performance values for RMSEC, RMSEP, R^2 , and Q^2 of 1.69, 2.67, 0.9749, and 0.9373, respectively. At the same time, the prediction of the lactose content generated predictive results for RMSEC, RMSEP, R^2 , and Q^2 values of 0.35, 0.79, 0.9987, and 0.9943, respectively. A relatively small value of RMSECs implied that the PLS successfully ftted the data and that the calibration models adequately modeled the response parameter variations. The PLS models could generate high R^2 and Q^2 values implying that the developed models were stable and could successfully be used to estimate the sugar contents in the prepared powder samples. These results corresponded to the correlation graphs presented in Fig. [3a](#page-6-0),b, where the predicted samples were placed approximately close to the diagonal lines of the graphs suggesting that most of the samples had slight diferences between the reference and the predicted sugar concentrations. The values of Q^2 were slightly lower than the R² values implying that the models were not prone to the overftting problem.

In general, RPD values allow for a comparison of the prediction accuracy of the diferent models. Higher values of RPD would indicate a better predictive ability corresponding to lower bias and SECV values along with higher Q² values. The established model for lactose prediction had a higher RPD value than the prediction of sucrose, indicating that the PLS model had a better predictive performance. These findings suggest that lactose, a naturally occurring sugar in milk, could be used to obtain additional related information from other components, such as the fats and proteins presented in the raw-milk material. In comparison to the sugar naturally found in food products, sucrose, which was added during the preparation process, contained less micronutrient information 34 .

Prediction of sugar content in real milk tablet samples. Table [4](#page-6-1) shows the predictive results of the milk tablets (E1–E7) using the developed PLS models. When the developed PLS models were used to estimate the sugar contents from the NIR spectra directly detected from the tablet samples (P to predict T), higher predictive errors were obtained, resulting in significantly lower Q^2 values (0.7024 and 0.7030 for sucrose and lactose, respectively) and higher RMSEP values (53.44 and 8.14 for sucrose and lactose). The poor predictive results correspond with the PCA score plot shown in Fig. [2b](#page-4-1), where the deviations between the powder and tablet samples were highlighted, noting that the physical structures of the samples could signifcantly afect the NIR measurements.

Signifcant improvements in the predictive results could be obtained afer the NIR spectral data were preprocessed by PDS transformation (P to predict T/P). Figure [4](#page-7-14) illustrates the comparison between the predictive results before and afer the PDS transformation demonstrating the signifcant reduction in the relative errors associated with non-destructive detection. In Table [4,](#page-6-1) the Q^2 values for both sucrose and lactose prediction models increased to 0.7973 and 0.9411, respectively. Te corresponding correlation graphs of the reference and the predicted sugar contents for the PLS models are shown in Fig. [3c](#page-6-0). The DS transformation resulted in the improvement in the prediction of the sucrose content having a reduction of the RMSEP value to 19.27, but it failed to improve the results for the prediction of the lactose content where the RMSEP was increased to 23.73. Tis confrmed that the calculation based on the spectral subset in the PDS method could improve the predictive accuracy of the PLS calibration models.

Figure 3. Correlation graphs for the training and internal validation samples of (**a**) sucrose, (**b**) lactose contents, and (**c**) the real milk tablet samples afer PDS transformation for the predictions of the sucrose and lactose contents.

Table 4. Sugar determinations in real milk tablet samples. *P = powder samples, **T/P = tablet samples adjusted to powder samples, ***T=tablet samples.

The lactose concentrations in samples 8 and 9 were null because they were milk tablets that contained only artificial-milk flavor without the actual cow's milk component. These results indicated that the laboratorydeveloped model could be used to predict real milk tablet samples or the fnal process samples by enabling uncomplicated measurements and real-time analysis of a large number of samples collected during the manufacturing process without the need for sample preparation.

7

Figure 4. Comparison between the prediction results before and afer the DS and PDS transformations of the real milk tablets.

Conclusion

In the manufacturing process, it was not practical to produce modeling samples that had enough variations to generate robust prediction models. In addition, the spectral deviations that occurred while completing the NIR measurements could have affected the prediction accuracy of the developed model. This research demonstrated that using the calibration transfer method widened the utilization ability of the developed calibration models. NIR spectroscopy combined with chemometric analyses can be applied to detect the sugar contents in milk tablets. PDS resulted in an improved level of the predictive performance of the tablet samples. The development process ofered non-destructive, accurate, and rapid techniques for determining sugar contents in real milk tablet samples.

Data availability

The NIR spectroscopic datasets analyzed during the current study are available from the corresponding author upon reasonable request.

Received: 30 August 2022; Accepted: 1 November 2022 Published online: 05 November 2022

References

- 1. RodriguezFurlán, L. T. & Campderrós, M. E. Te combined efects of *Stevia* and sucralose as sugar substitute and inulin as fat mimetic on the physicochemical properties of sugar-free reduced-fat dairy dessert. *Int. J. Gastron. Food Sci.* **10**, 16–23 (2017).
- 2. Gambelli, L. Milk and its sugar-lactose: A picture of evaluation methodologies. *Beverages* **3**, 35. [https://doi.org/10.3390/beverages3](https://doi.org/10.3390/beverages3030035) [030035](https://doi.org/10.3390/beverages3030035) (2017).
- 3. Alcaire, F., Antúnez, L., Vidal, L., Giménez, A. & Ares, G. Aroma-related cross-modal interactions for sugar reduction in milk desserts: Infuence on consumer perception. *Food Res. Int.* **97**, 45–50 (2017).
- 4. Mahato, D. K. *et al.* Sugar reduction in dairy food: An overview with favoured milk as an example. *Foods* **9**, 1400. [https://doi.org/](https://doi.org/10.3390/foods9101400) [10.3390/foods9101400](https://doi.org/10.3390/foods9101400) (2020).
- 5. Shafer, P. A. & Hartmann, A. F. Te iodometric determination of copper and its use in sugar analysis. *J. Biol. Chem.* **45**, 365–390 (1921).
- 6. Mannino, S., Cosio, M. S. & Zimei, P. Microdialysis sampling and high performance liquid chromatography with amperometric detection for sugar analysis in milk products. *Electroanalysis* **8**, 353–355 (1996).
- 7. Rentschler, E. *et al.* Quantifcation of lactulose and epilactose in the presence of lactose in milk using a dual HPLC analysis. *Food Anal. Methods* **9**, 2210–2222 (2016).
- 8. Conzuelo, F. *et al.* An integrated amperometric biosensor for the determination of lactose in milk and dairy products. *J. Agric. Food Chem.* **58**, 7141–7148 (2010).
- 9. Ammam, M. & Fransaer, J. Two-enzyme lactose biosensor based on β-galactosidase and glucose oxidase deposited by AC-electrophoresis: Characteristics and performance for lactose determination in milk. *Sens. Actuators B Chem.* **148**, 583–589 (2010).
- 10. Monakhova, Y. B., Kuballa, T., Leitz, J., Andlauer, C. & Lachenmeier, D. W. NMR spectroscopy as a screening tool to validate nutrition labeling of milk, lactose-free milk, and milk substitutes based on soy and grains. *Dairy Sci. Technol.* **92**, 109–120 (2012).
- 11. Fernández-Espinosa, A. J. Combining PLS regression with portable NIR spectroscopy to on-line monitor quality parameters in intact olives for determining optimal harvesting time. *Talanta* **148**, 216–228 (2016).
- 12. Oliveira, M. M., Cruz-Tirado, J. P., Roque, J. V., Teóflo, R. F. & Barbin, D. F. Portable near-infrared spectroscopy for rapid authentication of adulterated paprika powder. *J. Food Compos. Anal.* **87**, 103403 (2020).
- 13. Rahman, A., Wang, S., Yan, J. & Xu, H. Intact macadamia nut quality assessment using near-infrared spectroscopy and multivariate analysis. *J. Food Compos. Anal.* **102**, 104033 (2021).
- 14. Riu, J., Gorla, G., Chakif, D., Boqué, R. & Giussani, B. Rapid analysis of milk using low-cost pocket-size NIR spectrometers and multivariate analysis. *Foods* **9**, 1090 (2020).
- 15. Jorge dos Santos, V. *et al.* MicroNIR spectroscopy and multivariate calibration in the proximal composition determination of human milk. *LWT* **147**, 111645 (2021).
- 16. Sarraguça, M. C. & Lopes, J. A. Quality control of pharmaceuticals with NIR: From lab to process line. *Vib. Spectrosc.* **49**, 204–210 (2009)
- 17. Peerapattana, J., Otsuka, K. & Otsuka, M. Application of NIR spectroscopy for the quality control of mangosteen pericarp powder: Quantitative analysis of alpha-mangostin in mangosteen pericarp powder and capsule. *J. Nat. Med.* **67**, 452–459 (2013).
- 18. Liang, C., Yuan, H. F., Zhao, Z., Song, C. F. & Wang, J. J. A new multivariate calibration model transfer method of near-infrared spectral analysis. *Chemom. Intell. Lab. Syst.* **153**, 51–57 (2016).
- 19. Eliaerts, J. *et al.* Evaluation of a calibration transfer between a bench top and portable Mid-InfraRed spectrometer for cocaine classifcation and quantifcation. *Talanta* **209**, 120481 (2020).
- 20. McVey, C., McGrath, T. F., Haughey, S. A. & Elliott, C. T. A rapid food chain approach for authenticity screening: The development, validation and transferability of a chemometric model using two handheld near infrared spectroscopy (NIRS) devices. *Talanta* **222**, 121533 (2021).
- 21. Wongsaipun, S., Teanjumpol, P. & Kittiwachana, S. Development of a universal calibration model for quantifcation of adulteration in Tai jasmine rice using near-infrared spectroscopy. *Food Anal. Methods* **14**, 997–1010 (2021).
- 22. Srithongkul, C., Krongchai, C., Santasup, C. & Kittiwachana, S. An investigation of the efect of operational conditions on a sequential extraction procedure for arsenic in soil in Tailand. *Chemosphere* **242**, 125230 (2020).
- 23. Bezerra, M. A., Santelli, R. E., Oliveira, E. P., Villar, L. S. & Escaleira, L. A. Response surface methodology (RSM) as a tool for optimization in analytical chemistry. *Talanta* **76**, 965–977.<https://doi.org/10.1016/j.talanta.2008.05.019> (2008).
- 24. Dai, X. *et al.* Accurate analysis of urea in milk and milk powder by isotope dilution gas chromatography-mass spectrometry. *J. Chromatogr. B Analyt. Technol. Biomed. Life Sci.* **878**, 1634–1638 (2010).
- 25. Chiavelli, L. U. R. *et al.* Optimization of milk sample cleanup using response surface methodology. *Food Anal. Methods* **13**, 166–175 (2020)
- 26. Duy, P. K., Chun, S., Lee, Y. & Chung, H. Investigation of the particle size-dependent near-infrared spectral features of binary mixture samples in conjunction with Monte Carlo simulation and the infuence of particle size on the accuracy of quantitative analysis. *Analyst* **143**, 4306–4315 (2018).
- 27. Wang, Y., Veltkamp, D. J. & Kowalski, B. R. Multivariate instrument standardization. *Anal. Chem.* **63**, 2750–2756 (1991).
- 28. Wulfert, F., Kok, W. T., de Noord, O. E. & Smilde, A. K. Correction of temperature-induced spectral variation by Continuous Piecewise Direct Standardization. *Anal. Chem.* **72**, 1639–1644 (2000).
- 29. Brereton, R. G. *Chemometrics: Data Analysis for the Laboratory and Chemical Plant* (Wiley, 2003).
- 30. Mustorgi, E. *et al.* A chemometric strategy to evaluate the comparability of PLS models obtained from quartz cuvettes and disposable glass vials in the determination of extra virgin olive oil quality parameters by NIR spectroscopy. *Chemom. Intell. Lab. Syst.* **199**, 103974 (2020).
- 31. Kittiwachana, S. *et al.* One class classifers for process monitoring illustrated by the application to online HPLC of a continuous process. *J. Chemom.* **24**, 96–110 (2010).
- 32. Wongsaipun, S., Krongchai, C., Jakmunee, J. & Kittiwachana, S. Rice grain freshness measurement using rapid visco analyzer and chemometrics. *Food Anal. Methods* **11**, 613–623 (2018).
- 33. Bouveresse, E. & Massart, D. L. Standardisation of near-infrared spectrometric instruments: A review. *Vib. Spectrosc.* **11**, 3–15. [https://doi.org/10.1016/0924-2031\(95\)00055-0](https://doi.org/10.1016/0924-2031(95)00055-0) (1996).
- 34. Erickson, J. & Slavin, J. Total, added, and free sugars: Are restrictive guidelines science—Based or achievable?. *Nutrients* **7**, 2866– 2878 (2015).

Acknowledgements

This project is funded by National Research Council of Thailand (NRCT) and Chiang Mai University: N42A650304. Postharvest Technology Innovation Centre, Office of the Higher Education Commission, Bangkok, Tailand is also gratefully acknowledged.

Author contributions

C.T. carried out experiments, performed chemometric calculations, and wrote the original manuscript; N.P., C.T. and P.T. carried out experiments; S.K. conceptualized, supervised the funding of this work, reviewed, and edited the manuscript. All authors reviewed the manuscript and contributed to the fnal version.

Competing interests

The authors declare no competing interests.

Additional information

Correspondence and requests for materials should be addressed to S.K.

Reprints and permissions information is available at [www.nature.com/reprints.](www.nature.com/reprints)

Publisher's note Springer Nature remains neutral with regard to jurisdictional claims in published maps and institutional afliations.

Open Access Tis article is licensed under a Creative Commons Attribution 4.0 International \odot License, which permits use, sharing, adaptation, distribution and reproduction in any medium or format, as long as you give appropriate credit to the original author(s) and the source, provide a link to the Creative Commons licence, and indicate if changes were made. The images or other third party material in this article are included in the article's Creative Commons licence, unless indicated otherwise in a credit line to the material. If material is not included in the article's Creative Commons licence and your intended use is not permitted by statutory regulation or exceeds the permitted use, you will need to obtain permission directly from the copyright holder. To view a copy of this licence, visit<http://creativecommons.org/licenses/by/4.0/>.

 \circ The Author(s) 2022