



Since January 2020 Elsevier has created a COVID-19 resource centre with free information in English and Mandarin on the novel coronavirus COVID-19. The COVID-19 resource centre is hosted on Elsevier Connect, the company's public news and information website.

Elsevier hereby grants permission to make all its COVID-19-related research that is available on the COVID-19 resource centre - including this research content - immediately available in PubMed Central and other publicly funded repositories, such as the WHO COVID database with rights for unrestricted research re-use and analyses in any form or by any means with acknowledgement of the original source. These permissions are granted for free by Elsevier for as long as the COVID-19 resource centre remains active.



# Clinical performance evaluation of BD SARS-CoV-2 reagents for BD MAX<sup>TM</sup> System in asymptomatic individuals

Karen Yanson<sup>\*\*</sup>, William Laviers, Faten Suhaidi, Zachary Greeley, Courtney Merryman, Reda Proctor, Dominique Hall, Lori Neely<sup>\*</sup>

Becton, Dickinson and Company, BD Life Sciences – Integrated Diagnostic Solutions, Sparks, MD, USA

## ARTICLE INFO

### Article history:

Received 5 August 2022

Revised in revised form 14 October 2022

Accepted 11 November 2022

Available online 17 November 2022

### Keywords:

COVID-19  
SARS-CoV-2  
Asymptomatic  
Reagents  
BD MAX

## ABSTRACT

Transmission by asymptomatic individuals is a persistent hurdle in the effort to control the spread of SARS-CoV-2. Therefore, it is essential to continue developing assays and evaluate their performance for detection of SARS-CoV-2 in individuals without COVID-19 symptoms. In this study, 223 nasopharyngeal swab specimens collected from COVID-19 asymptomatic individuals were tested using the BD SARS-CoV-2 (RT-PCR-based) reagents for the BD MAX<sup>TM</sup> System and compared with results obtained with the Biomerieux BioFire<sup>®</sup> Respiratory RT-PCR Panel. Positive and negative percent agreements of 100% (95% CI, 84.5%–100%) and 99.0% (95% CI, 96.5%–99.7%), respectively, were observed for the BD SARS-CoV-2 assay. These results demonstrate the effectiveness of the BD SARS-CoV-2 assay for detecting SARS-CoV-2 in asymptomatic individuals and suggest that this assay can facilitate optimized case surveillance and infection control efforts. Investigations using larger sample sizes of asymptomatic individuals would be beneficial to support the findings in this study.

© 2022 The Author(s). Published by Elsevier Inc. This is an open access article under the CC BY-NC-ND license (<http://creativecommons.org/licenses/by-nc-nd/4.0/>)

## 1. Introduction

As the reactive strategies implemented during the SARS-CoV-2 pandemic to curb the spread of the virus are delivering their expected outcomes and COVID-19 becomes an endemic disease in many parts of the world, testing will continue to play a major role in surveillance and infection control strategies [1]. Since those activities require the availability of rapid and accurate diagnostic methods [2], it is expected that molecular testing for SARS-CoV-2, in tandem with other mitigation measures such as the isolation of infected individuals and continued vaccination efforts, will continue to play a crucial role in preventing COVID-19 outbreaks in communities and workplaces in the foreseeable future [3]. Given that asymptomatic transmission frequently occurs [4–9], and that as many as 40% of confirmed cases (those receiving a screening or confirmatory polymerase chain reaction [PCR] test) involve asymptomatic individuals [10,11], it is crucial to continue developing and validating highly sensitive SARS-CoV-2 assays to ensure that asymptomatic cases are identified and properly managed.

Since the beginning of the pandemic in 2019, the World Health Organization (WHO) and the U.S. Food and Drug Administration

(FDA) have issued emergency use (Emergency Use Listing [EUL] and Emergency Use Authorization [EUA], respectively) for the development of *in vitro* diagnostic assays to detect the SARS-CoV-2 virus, including molecular and antigen testing systems [12–14]. At the onset of the pandemic, however, limited resources (in terms of health care providers' (HCP) effort, laboratory supplies, and assay reagents) and overwhelming demand for SARS-CoV-2 testing resulted in tests being primarily performed with symptomatic individuals. As the understanding of SARS-CoV-2 transmission grew and the role of asymptomatic individuals in viral transmission became established, however, the validation and availability of molecular testing for asymptomatic individuals gained importance. In June 2020, a guidance document was released by the FDA to inform tests developers on the validation of COVID-19 assays intended for the screening of individuals without symptoms or other reasons to suspect infection with the virus, and for pooled sample testing [15]. In March 2021, the FDA also released supplemental guidance for the development of assays involving a serial testing indication for individuals without symptoms or other reasons to suspect COVID-19 [16].

The BD SARS-CoV-2 Reagents for BD MAX<sup>TM</sup> System assay (“BD SARS-CoV-2 assay;” Becton, Dickinson and Company; BD Life Sciences – Integrated Diagnostic Solutions, Sparks, MD) is a molecular testing system using reverse transcriptase-polymerase chain reaction (RT-PCR) that was initially authorized under EUA for the detection of SARS-CoV-2 RNA from individuals suspected of COVID-19 [17]. This study evaluated the performance of this assay in a population

\* Corresponding author. Tel: 410-316-3616; fax: 410-316-3690.

\*\* Corresponding author. Tel: 410-316-4793; fax: 410-316-4041.

E-mail addresses: [Karen\\_Yanson@bd.com](mailto:Karen_Yanson@bd.com) (K. Yanson), [Lori.Neely@BD.com](mailto:Lori.Neely@BD.com) (L. Neely).

without symptoms or other reasons to suspect COVID-19 infection by comparing its results to those obtained with the Biomerieux BioFire® Respiratory Panel 2.1 (“BioFire SARS-CoV-2 assay;” Biomerieux, BioFire Diagnostics, Salt Lake City, UT).

## 2. Methods and materials

### 2.1. Specimen collection and storage

Clinical nasopharyngeal (NP) specimens were consecutively collected from 3 clinical sites and one specimen vendor in the US from October 2020 to October 2021 to meet the required sample size, based on the FDA EUA guidance documentation for developers of molecular diagnostic tests [18] (July 2020 & October 2021 versions). These specimens were obtained from asymptomatic male and female individuals 1 year old and over who elected to undergo SARS-CoV-2 testing for any reasons other than (1) confirmed exposure to the virus or (2) suspected active infection. The NP specimens were stored in universal viral transport media (BD Universal Viral Transport, Becton, Dickinson and Company; BD Life Sciences – Integrated Diagnostic Solutions, Sparks, MD) or Viral Transport Media (VTM; COPAN Diagnostics Inc., Murrieta, CA) and frozen at -70°C.

### 2.2. Participant consent statement

For those specimens collected prospectively, no study procedures were performed without an informed consent process and signature of a consent form. Institutional review board approval of the protocol was received from Advarra Institutional Review Board prior to study initiation and de-identified specimens were used for testing. De-identified remnant specimens were obtained with Ethical & Independent Review Services or Western Institutional Review Board approval with waiver of informed consent. The study was also conducted according to the principles set forth by the Declaration of Helsinki and Good Clinical Practices.

### 2.3. Testing procedure

Deidentified clinical specimens were allowed to thaw at room temperature before aliquoting. For each clinical specimen, 3 aliquots were prepared according to manufacturers’ instructions. One aliquot (750uL of pipetted sample) was prepared for testing with the BD MAX SARS-CoV-2 assay, which utilizes multiplexed primers and probes that are designed to amplify 2 unique regions of the SARS-CoV-2 nucleocapsid (N) gene, N1 and N2, and the human ribonuclease P (RNase P) gene. Another aliquot (300uL of pipetted sample) was prepared for testing with an FDA De Novo granted RT-PCR assay, the Biomerieux BioFire® Respiratory Panel 2.1 (“BioFire SARS-CoV-2 assay;” Biomerieux, BioFire Diagnostics, Salt Lake City, UT), which utilizes multiplexed primers and probes that are designed to amplify unique regions of the SARS-CoV-2 spike and membrane proteins. The third aliquot (400uL of pipetted sample) was stored frozen for future discordant testing with the EUA Cepheid Xpert® Xpress SARS-CoV-2/

Flu/RSV assay (“Xpert SARS-CoV-2 assay,” Cepheid, Sunnyvale, CA), if necessary. All initial testing was performed internally at BD Life Sciences – Integrated Diagnostic Solutions, however, discordant samples were shipped on dry ice to an external laboratory for testing with the Xpert SARS-CoV-2 assay. All testing, inclusive of input sample volume, was conducted as described in each test’s respective *Instructions for Use* document. The study team was blind to the original reported standard-of-care result when performing the tests. The standard-of-care result was only utilized to ensure the sample cohort would have sufficient numbers of positive and negative samples and was not used in the assessment of agreement between the test under evaluation and the comparator methods. All environmental monitoring swabs generated negative results for all SARS-CoV-2 targets before any testing occurred and external processing controls yielded the expected results on each day of validation testing.

### 2.4. Data analysis

Test results from the BioFire SARS-CoV-2 assay were used as the comparator. Specimens are considered positive for SARS-CoV-2 by the BD MAX SARS-CoV-2 assay (index test) if either the N1 or N2 gene target are detected. The positive percent agreement (PPA) and the negative percent agreement (NPA) between the BD MAX SARS-CoV-2 assay results and comparator were calculated as [(index and comparator positive) / total comparator positive] and [(index and comparator negative) / total comparator negative], respectively [19]. The 95% confidence intervals (CI) were computed based on the Wilson Score method [20]. Acceptance criteria for the BD MAX SARS-CoV-2 assay for FDA’s Emergency Use Authorization for asymptomatic SARS-CoV-2 testing were ≥95% for PPA (lower bound of the two-sided 95% confidence interval >76%) and ≥98% NPA (lower bound of the two-sided 95% confidence interval >95%) [12]. This manuscript was prepared according to STARD guidelines for diagnostic accuracy studies reporting [21].

## 3. Results

A total of 224 specimens meeting the eligibility criteria were collected; one specimen was excluded due to insufficient volume to test, bringing the number of evaluable specimens to 223. Of those evaluable specimens, 36.8% (82/223) were from male and 63.2% (141/223) from female participants, and the following age distribution was observed: 3.1% (7/223) from participants <18 years old, 39.0% (87/223) for participants 18 to 29 years old, 26.9% (60/223) for those 30 to 39, 27.8% (62/223) for participants 40 to 65 years of age, and 3.1% (7/223) for those over the age of 65. Twenty-one (21) specimens were positive and 200 were negative by both BD MAX SARS-CoV-2 and the comparator method (positivity rate of 9.4% for both assays). There were 2 BD MAX SARS-CoV-2-positive but comparator negative specimens. The PPA was 100% (95% CI, 84.5%–100%) and the NPA was 99.0% (95% CI, 96.5%–99.7%) (Table 1). Upon testing of the 2 discordant specimens by a third method, the Xpert SARS-CoV-2 assay, it was found that one BD MAX SARS-CoV-2 positive specimen was also

**Table 1**  
Performance of the BD MAX SARS-CoV-2 assay compared to BioFire SARS-CoV-2 assay.

		BioFire		Total
		Positive	Negative	
MAX	Positive	21	2	23
	Negative	0	200	200
	Total	21	202	223
PPA 100% (CI, 84.5%, 100%)				
NPA 99.0% (CI, 96.5%, 99.7%)				

Abbreviations: PPA = positive percent agreement; NPA = negative percent agreement; CI = confidence interval.

**Table 2**

Discordant specimens tested by Xpert SARS-CoV-2 assay.

Specimen ID	MAX N1 Ct	MAX N2 Ct	BioFire	Xpert (Ct)
1	Pos (35.3)	Neg (-1)	Neg	Pos (42.4)
2	Neg (-1)	Pos (34.2)	Neg	Neg (0.0)

Abbreviations: Ct = threshold cycle.

positive by Xpert SARS-CoV-2 and 1 specimen positive by BD MAX SARS-CoV-2 was negative by Xpert SARS-CoV-2 (Table 2). Three (3) invalid test results were observed when executing the BioFire RP 2.1 test. Per BioFire IFU, these invalid samples were subjected to a repeat run and all were resolved, that is, upon repeat testing, 2 results were negative for SARS-CoV-2 and one was positive. There were no invalid results observed on the BD MAX assay.

#### 4. Discussion

The BD MAX SARS-CoV-2 assay showed a 100% (21/21) (95% confidence interval (84.5%–100%) PPA and 99.0% (200/202) (95% confidence interval (96.5%–99.7%) NPA, compared to the BioFire RP2.1 SARS-CoV-2 assay. Of the 2 BD MAX SARS-CoV-2 assay-positive/ BioFire RP2.1 SARS-CoV-2 assay-negative specimens, one was positive and one was negative by the discordant test.

Although the precise percentage of the COVID-19 population that remains asymptomatic through the course of the infection has not been determined [6-8,22,23], it is clear that this segment of the population contributes to the spread of SARS-CoV-2 and may exhibit viral loads comparable to those found in individuals with symptoms [6,23-25]. The breakthrough infections that occur among vaccinated individuals are also, in part, likely due to transmission from the asymptomatic population [26]. In order to effectively mitigate SARS-CoV-2 transmission, it may be beneficial to test asymptomatic individuals – especially those with contact or within close proximity to positive individuals, or those engaging in social activities (travel, sport venues, school, etc.) with large crowds or where physical proximity to others is required.

In response to the COVID-19 pandemic, substantial effort has been dedicated to developing sensitive diagnostic assays to detect the SARS-CoV-2 virus. With emerging SARS-CoV-2 variants, molecular testing platforms, particularly those detecting multiple genome targets, exhibit greater sensitivity and offer more robustness than antibody tests [27]. While it is already known that molecular testing is a highly sensitive SARS-CoV-2 detection method when used with symptomatic individuals, our study confirms that such molecular testing also performs well to detect the virus in asymptomatic individuals. However, the results obtained should be considered in relation to the limited number of positive samples evaluated, and as such, conducting further studies involving larger sample sizes would be useful.

#### 5. Conclusion

The BD MAX SARS-CoV-2 assay met FDA's acceptance criteria for EUA SARS-CoV-2 detection method in asymptomatic individuals and displayed a high clinical sensitivity and specificity for the detection of SARS-CoV-2, even in individuals without symptoms or other reasons to suspect a COVID-19 infection. Investigations using larger sample sizes of asymptomatic individuals would be beneficial to support the findings in this study.

#### Authors' Contributions

*Karen Yanson*: Conceptualization, resources, writing review and editing, project administration; *William Laviers*: Methodology, investigation, validation, formal analysis, data curation; *Faten Suhaidi*:

Investigation; *Zachary Greeley*: Investigation; *Courtney Merryman*: Investigation; *Reda Proctor*: Investigation; *Dominique Hall*: Investigation; and *Lori Neely*: Supervision, methodology, writing review and editing, resources.

#### Declaration of competing interests

KY, WL, FS, ZG, CM, RP, DH, and LN are employees of the study sponsor, Becton, Dickinson and Company; these authors have no other potential conflicts of interest to declare.

#### Acknowledgments

We thank H el ene Tanguay, Ph.D. MWC (Becton, Dickinson and Company, BD Life Sciences – Integrated Diagnostic Solutions) for input on the content of this manuscript and editorial assistance. We are also grateful to the study participants who allowed this work to be performed.

#### Funding

This study was funded by Becton, Dickinson and Company, BD Life Sciences – Integrated Diagnostic Solutions.

#### References

- [1] The Lancet Infectious, D. Transitioning to endemicity with COVID-19 research. *Lancet Infect Dis* 2022;22:297.
- [2] Peeling RW, Heymann DL, Teo YY, Garcia PJ. Diagnostics for COVID-19: moving from pandemic response to control. *Lancet* 2022;399:757–68.
- [3] Miller MB, Ooi EE, Rhoads DD, Kulldorff M, Anderson DE, Lee H, et al. As omicron takes hold and other new variants arise, COVID-19 testing remains the universally agreed tool to effect transition from pandemic to endemic state. *Front Public Health* 2022;10:883066.
- [4] Furukawa NW, Brooks JT, Sobel J. Evidence supporting transmission of severe acute respiratory syndrome coronavirus 2 while presymptomatic or asymptomatic. *Emerg Infect Dis* 2020;26:e1–6.
- [5] Kronbichler A, Kresse D, Yoon S, Lee KH, Effenberger M, Shin JI. Asymptomatic patients as a source of COVID-19 infections: a systematic review and meta-analysis. *Int J Infect Dis* 2020;98:180–6.
- [6] Nikolai LA, Meyer CG, Kremsner PG, Velavan TP. Asymptomatic SARS Coronavirus 2 infection: invisible yet invincible. *Int J Infect Dis* 2020;100:112–6.
- [7] Al-Sadeq DW, Nasrallah GK. The incidence of the novel coronavirus SARS-CoV-2 among asymptomatic patients: a systematic review. *Int J Infect Dis* 2020;98:372–80.
- [8] Han D, Li R, Han Y, Zhang R, Li J. COVID-19: Insight into the asymptomatic SARS-CoV-2 infection and transmission. *Int J Biol Sci* 2020;16:2803–11.
- [9] Sayampanathan AA, Heng CS, Pin PH, Pang J, Leong TY, Lee JY. Infectivity of asymptomatic vs symptomatic COVID-19. *Lancet* 2021;397:93–4.
- [10] Ma Q, Liu J, Liu Q, Kang L, Liu R, Jing W, et al. Global percentage of asymptomatic SARS-CoV-2 infections among the tested population and individuals with confirmed COVID-19 diagnosis: a systematic review and meta-analysis. *JAMA Netw Open* 2021;4:e2137257.
- [11] Parikh A, Cooper L, Frogel D, Le Bengier K, Cooper CK, Parvu V. Large-scale SARS-CoV-2 antigen testing with real-world specimens. *Front Public Health* 2022;10:836328.
- [12] U.S. Food and Drug Administration. Coronavirus Disease 2019 (COVID-19) emergency use authorization for medical devices—in vitro diagnostics EUAs. Available from: <https://www.fda.gov/medical-devices/coronavirus-disease-2019-covid-19-emergency-use-authorizations-medical-devices/in-vitro-diagnostics-euas#individual-molecular>. Accessed October 5, 2022.
- [13] WHO. Global research collaboration for infectious disease preparedness. COVID 19: Public Health Emergency of International Concern (PHEIC). global research and innovation forum: towards a research roadmap. Feb. 11–12, 2020. Available from: [https://www.who.int/blueprint/priority-diseases/key-action/Global\\_Research\\_Forum\\_FINAL\\_VERSION\\_for\\_web\\_14\\_feb\\_2020.pdf?ua=1](https://www.who.int/blueprint/priority-diseases/key-action/Global_Research_Forum_FINAL_VERSION_for_web_14_feb_2020.pdf?ua=1). Accessed October 5, 2022.

- [14] U.S. Food and Drug Administration. *Policy for diagnostic tests for coronavirus disease-2019 during the public health emergency*. Available from: [www.fda.gov/regulatory-information/search-fda-guidance-documents/policy-diagnostic-tests-coronavirus-disease-2019-during-public-health-emergency](https://www.fda.gov/regulatory-information/search-fda-guidance-documents/policy-diagnostic-tests-coronavirus-disease-2019-during-public-health-emergency). Accessed February 23, 2021.
- [15] U.S. Food and Drug Administration. *Coronavirus (COVID-19) update: facilitating diagnostic test availability for asymptomatic testing and sample pooling*. FDA statement, 2020. Available from: <https://www.fda.gov/news-events/press-announcements/coronavirus-covid-19-update-facilitating-diagnostic-test-availability-asymptomatic-testing-and>. Accessed October 5, 2022.
- [16] U.S. Food and Drug Administration. *Coronavirus (COVID-19) Update: FDA takes steps to streamline path for COVID-19 screening tools, provides information to help groups establishing testing programs*. FDA statement, 2021. Available from: <https://www.fda.gov/news-events/press-announcements/coronavirus-covid-19-update-fda-takes-steps-streamline-path-covid-19-screening-tools-provides>. Accessed October 5, 2022.
- [17] BD SARS-CoV-2 Reagents for BD MAX™ System [package insert]. Sparks-Glencoe, MD, USA: Becton, Dickinson and Company; 2020.
- [18] U.S. Food and Drug Administration. *Template for developers of molecular diagnostic tests*. Oct. 6, 2021. Available from: <https://www.fda.gov/medical-devices/coronavirus-disease-2019-covid-19-emergency-use-authorizations-medical-devices/in-vitro-diagnostics-euas# covid19ivdTemplates>. Accessed May 13, 2022.
- [19] U.S. Food and Drug Administration. *Statistical guidance on reporting results from studies evaluating diagnostic tests - guidance for industry and FDA staff*. March 2007. Available from: <https://www.fda.gov/regulatory-information/search-fda-guidance-documents/statistical-guidance-reporting-results-studies-evaluating-diagnostic-tests-guidance-industry-and-fda>. Accessed September 26, 2022.
- [20] Wilson EB. Probable inference, the law of succession, and statistical inference. *J Am Statist Assoc* 1927;22:209–12.
- [21] Bossuyt PM, Reitsma JB, Bruns DE, Gatsonis CA, Glasziou PP, Irwig L, et al. STARD 2015: an updated list of essential items for reporting diagnostic accuracy studies. *Radiology* 2015;277:826–32.
- [22] Byambasuren O, Cardona M, Bell K, Clark J, McLaws ML, Glasziou P. Estimating the extent of asymptomatic COVID-19 and its potential for community transmission: systematic review and meta-analysis. *Official J Assoc Med Microbiol Infect Dis Canada* 2020;5(4) Vol. 5, No. 4. doi: 10.3138/jammi-2020-0030.
- [23] Gao Z, Xu Y, Sun C, Wang X, Guo Y, Qiu S, et al. A systematic review of asymptomatic infections with COVID-19. *J Microbiol Immunol Infect* 2021;54:12–6.
- [24] Lee S, Kim T, Lee E, Lee C, Kim H, Rhee H, et al. Clinical course and molecular viral shedding among asymptomatic and symptomatic patients with SARS-CoV-2 infection in a community treatment center in the Republic of Korea. *JAMA Intern Med* 2020;180:1447–52.
- [25] Van Vinh Chau N, Lam VT, Dung NT, Yen LM, Minh NNQ, Hung LM, et al. The natural history and transmission potential of asymptomatic severe acute respiratory syndrome coronavirus 2 infection. *Clin Infect Dis* 2020;71:2679–87.
- [26] CDC COVID-19 Vaccine Breakthrough Case Investigation Team. COVID-19 vaccine breakthrough infections reported to CDC - United States, January 1–April 30, 2021. *MMWR Morb Mortal Wkly Rep* 2021;70:792–3.
- [27] Yanson K, Laviers W, Neely L, Lockamy E, Castillo-Hernandez LC, Oldfield C, et al. Performance evaluation of the BD SARS-CoV-2 reagents for the BD MAX system. *J Clin Microbiol* 2021;59:e0101921.