



Published in final edited form as:

Nat Protoc. 2022 March ; 17(3): 748–780. doi:10.1038/s41596-021-00665-4.

## Synthesis of siRNA nanoparticles to silence plaque-destabilizing gene in atherosclerotic lesional macrophages

Xiangang Huang<sup>1,†</sup>, Chuang Liu<sup>1,†</sup>, Na Kong<sup>1</sup>, Yufen Xiao<sup>1</sup>, Arif Yurdagul Jr.<sup>2</sup>, Ira Tabas<sup>3,4,5</sup>, Wei Tao<sup>1,\*</sup>

<sup>1</sup>Center for Nanomedicine and Department of Anesthesiology, Brigham and Women's Hospital, Harvard Medical School, Boston, MA 02115, USA.

<sup>2</sup>Department of Molecular and Cellular Physiology, LSU Health Shreveport, Shreveport, LA, 71130, USA.

<sup>3</sup>Department of Medicine, Columbia University Irving Medical Center, New York, NY 10032, USA.

<sup>4</sup>Department of Physiology and Cellular Biophysics, Columbia University Irving Medical Center, New York, NY 10032, USA.

<sup>5</sup>Department of Pathology and Cell Biology, Columbia University Irving Medical Center, New York, NY 10032, USA.

### Abstract

Macrophages in atherosclerotic lesions promote plaque progression and are an attractive therapeutic target in cardiovascular research. Here we present a protocol for synthesis of small interfering RNA (siRNA) nanoparticles (NP) that target lesional macrophages as a potential treatment for atherosclerosis. Ca<sup>2+</sup>/calmodulin-dependent protein kinase  $\gamma$  (CaMKII $\gamma$ ) activity in macrophages of advanced human and mouse atherosclerotic plaques drives necrosis via down-regulating the expression of the efferocytosis receptor MerTK. Therefore, selective inhibition of CaMKII $\gamma$  in lesional macrophages holds great promise for the treatment of advanced atherosclerosis. We have recently developed a siRNA NP platform that can selectively silence CaMKII $\gamma$  in macrophages, resulting in increased plaque stability. We provide a detailed protocol for the synthesis of NP components, the preparation and characterization (physicochemical and *in vitro*) of siRNA NPs, and the evaluation of *in vivo* therapeutic effects of siRNA NPs and their biocompatibility in atherosclerotic mice. Our siRNA-loaded polymer-lipid hybrid NPs are constructed via a robust self-assembly method, exhibiting excellent *in vivo* features for systemic siRNA delivery. Following this protocol, it takes 3-5 d to prepare the siRNA NPs, 8-10 d to characterize the NPs, and 4-5 weeks to evaluate their therapeutic effects in established

\*Corresponding author. wtao@bwh.harvard.edu.

†These authors contributed equally to this work.

#### Author Contributions

Conceptualization, W.T., I.T., A.Y.J.; writing and editing, W.T., I.T., A.Y.J., X.H., C.L., N.K., Y.X.; all the authors read and approved the final manuscript.

#### Competing interests

The authors declare no competing interests.

#### Supplementary information

Supplementary information is available for this paper.

atherosclerotic mice. By changing the RNA molecules loaded in the NPs, lesional macrophages can be targeted for the exploration and validation of new targets/pathways in atherosclerosis.

---

## Introduction

Macrophage activation has been implicated in the pathogenesis of many human diseases, including atherosclerosis, obesity, diabetes, cancer, and skin diseases.<sup>1</sup> Among them, atherosclerosis, an inflammatory disease involving arterial plaque deposition which increases the risk of cardiovascular disease and stroke, is highly correlated with macrophage activation.<sup>2,3</sup> Therefore, inhibiting the activation of pro-atherogenic pathways in macrophages holds great promise for the treatment of atherosclerotic vascular diseases.<sup>4</sup> However, to date, only suboptimal therapeutic effects have been achieved due to the lack of therapeutics or treatments with high specificity and efficacy.

RNA interference (RNAi) is a powerful method capable of selectively silencing any gene of interest, showing great potential to treat various diseases.<sup>5-10</sup> The success of RNAi therapy largely depends on the efficient *in vivo* delivery of therapeutic small interfering RNA (siRNA) to the target cells. To reach the RNAi machinery in the cytoplasm, siRNA molecules need to overcome several physiological barriers, including being degraded by nucleases, being rapidly cleared by the mononuclear phagocyte system (MPS) or renal filtration,<sup>11</sup> poor cellular uptake, and inefficient endosomal escape.<sup>8</sup> Although the rapid development of siRNA technology has yielded relatively stable siRNA with reduced side effects through chemical modification and sequence optimization,<sup>12-14</sup> the efficient delivery of siRNA into target cells is still challenging. Therefore, it is highly desirable to develop a safe and efficient delivery vehicle that can protect the siRNA from enzymatic degradation, deliver it to target cells and facilitate its cellular uptake and endosomal escape.

To this end, various delivery systems have been developed for more effective delivery of siRNA.<sup>15</sup> Many of them have been investigated in clinical trials for the treatment of different diseases, including cancer, genetic disorders, fibrosis, hypercholesterolemia, and viral infections.<sup>16</sup> The research interests of the scientific community in siRNA-based therapy have been fueled by the first U.S. Food and Drug Administration (FDA)-approved nano-formulated siRNA drug (Patisiran, Alnylam Pharmaceuticals), treating a rare genetic disorder.<sup>17</sup> Although several polymeric or lipid nanoparticle (NP)-based siRNA delivery systems have been developed for the treatment of vascular diseases,<sup>18-24</sup> siRNA-based therapy via targeting macrophages in atherosclerotic lesions for atherosclerosis treatment is still in its infancy. It is worth mentioning that Leuschner et al. pioneered the siRNA-based therapy using cationic lipid NPs for *in vivo* treatment of atherosclerosis via targeting macrophages.<sup>25</sup> While these polymeric or lipid NP-based systems have shown satisfying therapeutic effects, they might not be applicable in the treatments where higher gene therapy efficacy are required. Therefore, novel strategies which can further enhance the protection of siRNA and prolong its *in vivo* circulation time are needed.

## Development of the Protocol

To address outstanding limitations in nucleic acid delivery to macrophages, we developed a macrophage-targeting siRNA NP platform for the treatment of atherosclerosis via selectively silencing a functional gene in macrophages of atherosclerotic lesions.<sup>26</sup> The macrophage gene *Camk2g*, encoding a calcium-activated kinase named CaMKII $\gamma$  (Ca<sup>2+</sup>/calmodulin-dependent protein kinase  $\gamma$ ), was chosen as the target based on our previous work. We demonstrated that the activation of lesional macrophage CaMKII $\gamma$  in advanced human and mouse atherosclerotic plaques was a major stimulus for the formation of necrotic lesions covered with thin fibrous caps<sup>27</sup>. Efferocytosis is a specialized phagocytic process for the removal of dead cells.<sup>28</sup> Mechanistic studies revealed that the CaMKII $\gamma$  suppresses the activation of transcription factor 6 (ATF6)/liver X receptor alpha (LXR $\alpha$ )/c-Mer proto-oncogene tyrosine kinase (MerTK) pathway that usually promotes the efferocytosis. Consequently, the inefficient efferocytosis causes uncleared cells secondary necrosis, which eventually results in plaque necrosis and fibrous cap thinning (Fig. 1).<sup>27</sup> Thin-capped, necrotic lesions destabilize the plaques in humans, leading to harmful consequences such as plaque rupture, acute luminal thrombosis, and various acute atherothrombotic vascular events.<sup>29</sup>

Our siRNA NPs are rationally designed and constructed. First, a cationic lipid is used to complex the siRNA via electrostatic interactions. Then, a hydrophobic poly(lactic-co-glycolic) acid (PLGA) polymer further interacts with the siRNA/cationic lipid complex via hydrophobic interactions, followed by PEGylation of the complex using lipid-polyethylene glycol (lipid-PEG), yielding the siRNA NPs. To enable the targeted delivery of siRNA to macrophages, the macrophage-targeting peptide ligand S2P is conjugated to the siRNA NPs through the lipid-PEG. In a fat-fed *Ldlr*<sup>-/-</sup> mouse model with established atherosclerosis, the administration of the S2P-siCamk2g NPs can silence the CaMKII $\gamma$  of the lesional macrophage and subsequently enhance efferocytosis, resulting in thicker fibrous cap and less atherosclerotic plaque necrosis. Through this formulation, we developed a promising siRNA NP platform for the treatment of atherosclerosis via genetic manipulation of a specific molecule in lesional macrophages, providing a new therapeutic strategy for atherosclerosis and other-related vascular diseases.

Here we outline the procedures for the preparation and characterization of siRNA NPs, and the evaluation of *in vivo* therapeutic effects of siRNA NPs in atherosclerotic mice. The NP components cationic lipid (G0-C14) and peptide-modified lipid-PEG (DSPE-PEG-S2P) can be easily synthesized by a one-step reaction. The siRNA-loaded polymer-lipid hybrid NPs are conveniently achieved by dropwise addition of the siRNA complex into the lipid-PEG aqueous solution at room temperature (25 °C). After quick purification by ultrafiltration, the resulting siRNA NPs can be readily used for *in vitro* and *in vivo* evaluation of gene silencing efficiency and therapeutic effects for atherosclerosis. We anticipate that these procedures will be instrumental for researchers seeking to build siRNA delivery systems to treat atherosclerosis and other diseases.

## Overview of the procedure

The activation of lesional macrophages is one of the causes of the instability of the atherosclerotic plaques. We present here a robust strategy for the treatment of atherosclerosis in mice using lesional macrophage-targeting siRNA NPs. This protocol contains five main stages (Fig. 2). We first describe the synthesis and purification of NP components, including cationic lipid (G0-C14) and targeted lipid-PEG (DSPE-PEG-S2P) (Fig. 2a). Targeting is achieved by modifying the lipid with a peptide sequence that binds to/interacts with the macrophage receptor stabilin-2. Nuclear magnetic resonance (NMR) is used to confirm the structure of the compounds. We then provide a step-by-step nano-precipitation approach for the preparation of targeted siRNA NPs (Fig. 2b). The preparation and optimization of the macrophage-targeting siRNA NPs are key features of the protocol. Five components, including the synthesized (G0-C14, DSPE-PEG-S2P) and commercially available (siRNA, PLGA, DSPE-PEG) compounds, are used to construct the siRNA NPs. The ratio of components, the formulation, the stirring speed, and the incubation time significantly impact the properties of the resulting NPs. For example, if the stirring speed is much lower than the recommended speed (1000 rpm), large NPs with low siRNA encapsulation efficiency will be formed. All the optimized parameters are discussed in this protocol. To confirm the successful formation of siRNA NPs, their physicochemical properties are fully characterized (Fig. 2c). Dynamic light scattering (DLS) is used to determine the size and colloidal stability of siRNA NPs. Transmission electron microscopy (TEM) is employed to characterize the morphology of siRNA NPs. Gel electrophoresis is used to investigate the stability of siRNA encapsulated by NPs. The resulting siRNA NPs are further characterized *in vitro* (Fig. 2d). Confocal microscopy is used to evaluate the endosomal escape ability of siRNA NPs by measuring the fluorescent signals of dye-labeled siRNA. A plate reader is implemented to quantify the gene silencing efficiency of the siRNA targeting luciferase (siLuc) NPs by measuring the bioluminescence signals. A plate reader is also used to assess the *in vitro* toxicity of siRNA NPs (the alamarBlue cell viability reagent). Finally, we evaluate the *in vivo* therapeutic effect of siCamk2g NPs in a WD-fed *Ldlr*<sup>-/-</sup> mouse model with established atherosclerosis (Fig. 2e). To confirm that the siCamk2g NPs can improve the plaque stability, microscopy is employed to determine the fibrous cap thickness and necrotic core area (picosirius red staining). To further demonstrate that the enhanced plaque stability is ascribed to the promoted efferocytosis, fluorescence microscopy is implemented to characterize the macrophage-associated apoptotic cells (TUNEL assay). In addition, we evaluate the *in vivo* toxicity of the siRNA NPs using H&E staining and hematological analysis. Throughout the procedure (Fig. 2), we discuss the parameters that need to be optimized and how they influence the properties of the siRNA NPs, and how to troubleshoot any problems encountered when following the procedures.

## Applications of the method

The key features of this protocol are the preparation, optimization and characterization of the siRNA-loaded polymer-lipid hybrid NPs, which have various applications in the biomedical field. The first application of this siRNA NP platform is the treatment of atherosclerosis which is described in this protocol using a mouse model. It has been previously established that the activation of lesional macrophage CaMKII $\gamma$  in advanced human and mouse atherosclerotic plaques is a major stimulus for the formation of necrotic

lesions covered with thin fibrous caps.<sup>27</sup> Therefore, the macrophage gene *Camk2g*, encoding CaMKII $\gamma$ , was chosen as the target for RNAi therapy. After administration to the mouse model with atherosclerosis, the targeted S2P-siCamk2g NPs were able to selectively silence the CaMKII $\gamma$  of the lesional macrophage and subsequently enhance efferocytosis, resulting in thicker fibrous cap and less atherosclerotic plaque necrosis.<sup>26</sup> These siRNA NPs hold great promise for the treatment of atherosclerosis as demonstrated in a mouse model. For the treatment of other diseases, potential users can obtain tailored siRNA NP platforms conveniently by replacing the siRNA used in this protocol with another functional siRNA sequence. For example, this siRNA NP platform has also been used successfully for the treatment of non-small cell lung cancer (NSCLC) in a mouse model as demonstrated by another study in our group, in which, Prohibitin 1 (PHB1) was chosen as the target to silence. The siPHB1-loaded NPs were <100 nm in diameter, had long blood circulation time ( $t_{1/2}$  of ~8.1 h), had high tumor accumulation (~5.2 mg/kg), and could effectively inhibit tumor growth in mice bearing NCI-H460 and A549 tumors.<sup>6</sup> In addition to siRNA delivery, the present NP platform can also be used for the delivery of large nucleic acid molecules such as mRNA with proper optimization.<sup>30,31</sup>

The wide deployment of two highly effective NP-based mRNA COVID-19 vaccines (mRNA-1273, Moderna<sup>32</sup>; and BNT162b2, Pfizer<sup>33</sup>) have fueled the research interests in NP-based mRNA delivery.<sup>34-37</sup> It is worth noting that the NPs used in these two vaccines are chemically dissimilar to our NPs. For example, they contain ionizable lipids rather than cationic lipids, and they do not contain PLGA or targeting ligands. Despite that, the application of our polymer-lipid hybrid NP platform for mRNA delivery for cancer therapy has been demonstrated in mouse models in two other studies by our group.<sup>30,31</sup> Phosphatase and tensin homologue deleted on chromosome 10 (*PTEN*) is a well-defined tumor suppressor gene that is lost or mutated in various human cancers.<sup>38</sup> The restoration of functional *PTEN* in *PTEN*-null tumor cells has been shown to suppress, at least temporarily, the growth of tumors. Indeed, the *PTEN*-mRNA can be efficiently incorporated into our polymer-lipid hybrid NPs with slight optimization to the protocol presented here. The resulting *PTEN*-mRNA-loaded NPs exhibit good stability in physiological environments and achieve high *PTEN*-mRNA transfection efficiency in prostate cancer cells without eliciting significant cytotoxicity. Furthermore, in several mouse models of prostate cancer, the systemic administration of the *PTEN*-mRNA-loaded NPs effectively inhibited the tumor growth in the mice.<sup>30</sup> Similarly, another tumor suppressor gene named *p53* can also be employed as the therapeutic target for cancer therapy. We demonstrated that *p53*-mRNA-loaded polymer-lipid hybrid NPs can inhibit the tumor growth of mice bearing different tumor models, and also sensitize the tumors to mTOR inhibition.<sup>31</sup> Several recently published reviews have discussed the potential of our polymer-lipid hybrid NPs as an efficient siRNA or mRNA delivery platform for the treatment of different diseases.<sup>39-44</sup>

Overall, the polymer-lipid hybrid NP platforms generated by following our protocol are able to efficiently deliver siRNA and mRNA to target sites *in vivo* for the treatments of atherosclerosis and cancer, as validated by the studies of our group.<sup>6,26,30,31</sup> The NPs we prepared here have been modified with the targeting peptide S2P, a peptide that recognizes the macrophage receptor stabilin-2, allowing for the specific targeting of lesional macrophages. Besides the application for the treatment of atherosclerosis, our polymer-

lipid hybrid NPs without any targeting ligand have been reported to show significant accumulations in tumors, holding great promise for cancer therapy.<sup>6</sup> Proper modifications are required before these NPs can be used for other applications. We anticipate that this protocol will be of interest to researchers who are developing RNA NPs to treat atherosclerosis, cancer, and beyond. The techniques involved in this protocol can also be used by specialists in the fields of chemistry, nanotechnology, cardiovascular diseases, biology, RNA therapeutics, and pharmacy.

### Comparison with other methods

The present protocol describes the construction of a siRNA-loaded polymer-lipid hybrid NP platform for the treatment of atherosclerosis in mice via the regulation of lesional macrophage activation. Although our siRNA NPs-based therapeutic strategy for atherosclerosis has been rarely reported, a variety of other siRNA NP platforms have been developed for the treatment of a wide range of diseases.<sup>15</sup> One of the widely investigated siRNA NP platforms is the polymeric NP. For example, a PEGylated poly guanidinium (PEG-pGu)-based siRNA NP platform has been developed for Alzheimer's disease (AD) therapy.<sup>45</sup> To increase the blood-brain barrier (BBB) penetration, a galactose (Gal) ligand was introduced to the siRNA NPs. In a transgenic AD mouse model, the administrated Gal-siRNA NPs were able to cross the BBB via glycemia-controlled glucose transporter-1 (Glut1)-mediated transport, improving the cognitive behavior of the mice. One concern of this cationic polymer system is the potential toxicity. To address this issue, a degradable influenza virus-inspired cationic polymer system was prepared for siRNA delivery.<sup>46</sup> Upon internalization by the cells, the resulting siRNA NPs could quickly escape from the endosome through an influenza virus mimicking endosomal escape mechanism, releasing the siRNA into the cytosol via the degradation of the polymer. Despite the fact that the toxicity issues of cationic polymers can be alleviated by the incorporation of degradable bonds, the instability of the polymeric NPs can be problematic when administrated *in vivo*. The only driving force for the siRNA/cationic polymer complex is electrostatic interactions, which can be largely weakened in the ionic strength condition of the body.

Another siRNA NP delivery platform is the DNA nanostructure.<sup>47,48</sup> As one of the well-characterized DNA nanostructures, DNA tetrahedron has been widely used for the delivery of therapeutics including siRNA. For instance, a self-assembled DNA tetrahedron was prepared for *in vivo* siRNA delivery.<sup>49</sup> The siRNA cargos were conveniently incorporated onto one of the six edges of the DNA tetrahedrons through one-step self-assembly. Importantly, different cancer-targeting ligands could be easily conjugated to surfaces of the DNA tetrahedrons via DNA hybridization. The high programmability of the DNA strands allowed the precise control of the size of the tetrahedrons, spatial orientation, and density of the ligands. The siRNA-loaded DNA tetrahedral NPs could deliver the therapeutic siRNA to the tumor and silence the target genes. A similar siRNA-loaded DNA tetrahedral NP platform has also been developed for targeted delivery of therapeutic siRNA to the kidney of the mice, treating acute kidney injury (AKI).<sup>50</sup> In addition, a DNA nanosuitcase platform<sup>51</sup> and a DNA nanotube platform<sup>52</sup> have also been constructed for siRNA delivery. A major advantage of all these DNA nanostructure-based NP platforms is their negligible toxicity owing to the inherent biocompatibility of DNA materials. However, these DNA

nanostructures do not provide enough protection for the siRNA cargos, making them vulnerable to enzymatic degradation.

Among all the siRNA NP platforms, the lipid NP is the most clinically advanced platform. Indeed, the first FDA-approved siRNA therapeutic (Patisiran) is based on lipid NP formulations.<sup>17</sup> This type of lipid NP is typically composed of: the ionizable lipid, phospholipid, cholesterol, lipid-PEG, and siRNA.<sup>53</sup> Homogenous lipid NPs with a diameter ~50-100 nm can be formed by rapidly mixing siRNAs with other components using a microfluidic mixing technique.<sup>54</sup> Notably, these lipid NPs are neutral at pH 7.4, but protonated in the acidic environment of endosomes, improving the endosomal escape of siRNA via the fusion of lipids with the membrane of endosomes. Despite the successful applications of siRNA lipid NPs,<sup>55,56</sup> specific designs might be required to improve their delivery efficiency for different diseases.<sup>57</sup> For example, for specific cell types-related diseases such as macrophage-related atherosclerosis, a macrophage-targeting ligand can be introduced to the NPs to improve their therapeutic efficiency. Recently, a well-designed selective organ targeting (SORT) NP platform shows its capability to specifically deliver the therapeutics (Cas9 mRNA, single guide RNA) to lung, spleen and liver to enable efficient gene editing, holding great promise for the treatment of diverse organs-related diseases.<sup>58,59</sup>

Here, we provide a detailed protocol for the development of a polymer-lipid hybrid NP platform for efficient *in vivo* delivery of siRNA.<sup>6,26</sup> To date, siRNA-loaded polymer-lipid hybrid NPs have been typically synthesized by solvent evaporation or double emulsion methods.<sup>60-63</sup> However, our siRNA-loaded polymer-lipid hybrid NP is prepared through a simple nano-precipitation method, and has a unique structure consisting of a siRNA/cationic lipid/PLGA polymer core and a lipid-PEG shell. Compared to the ionizable lipid NP platform, we have introduced an extra hydrophobic polymer core to our platform, aiming to further protect the unstable nucleic acid cargo and prolong the blood circulation time. We have previously prepared different polymer-lipid hybrid NPs with several different cationic lipids, and the most efficient NP was prepared with G0-C14 lipid which contains five tertiary amines and seven hydrocarbon tails.<sup>6</sup> We use the lowest concentration of G0-C14 (siRNA/G0-C14 weight ratio, 1:15) that is able to condense the siRNA to prepare our siRNA NPs, while no significant *in vitro* and *in vivo* toxicities of the NPs are observed. The tertiary amines can provide sufficient positive charges to condense the negatively charged nucleic acids, and the hydrocarbon tails can bind to the hydrophobic polymer via hydrophobic interaction. Therefore, the amines and hydrocarbon tails are pivotal for the construction of our platform, whereas ionizable lipids do not have sufficient cationic charge to be suitable for our platform.

Our polymer-lipid hybrid NP platform has several advantages over the other *in vivo* delivery platform discussed above. First, the NPs have a long blood circulation time ( $t_{1/2}$  of ~8.1 h compared to  $t_{1/2}$  of ~4 h for a MC3-based lipid NP<sup>64</sup>). The formation of the NPs is driven by the electrostatic interactions between the anionic nucleic acids and cationic lipids, and the hydrophobic interaction between the hydrophobic PLGA and lipid tails. Utilizing both electrostatic and hydrophobic forces produces NPs with high *in vivo* stability.<sup>6</sup> Second, targeting ligands can be easily introduced to the NPs. Given the multicomponent structure of the NPs, different targeting ligands can be conveniently introduced to the lipid-PEG shell

without changing the main structure of the NPs. The ligand modification is described in this protocol.<sup>26</sup> In addition, to realize the on-demand release of the therapeutics, the hydrophobic PLGA polymer of the NPs can be replaced by a stimuli-responsive hydrophobic polymer. For example, we developed a redox responsive NP platform by replacing the PLGA polymer with a reduction-responsive polymer.<sup>31</sup> Taken together, our polymer-lipid hybrid NP platform is highly flexible (to be engineered with different functions) and easy-to-prepare, holding great promise for efficient delivery of siRNA and other types of nucleic acids.

## Experimental design

Here we describe the design and rationale for each procedure in this protocol.

**Synthesis of cationic lipid and ligand (Steps 1-24)**—The cationic lipid termed G0-C14 is synthesized by a unique synthetic strategy involving effective ring-opening of epoxide (1,2-epoxytetradecane) by amine compound (poly (amido amide) (PAMAM) dendrimer G0) (Fig. 3a).<sup>60</sup> This synthetic strategy requires no solvent and protection/deprotection steps, making it highly suitable for high-throughput screening of cationic lipids for superior nucleic acid delivery.<sup>65</sup> The generation 0 of PAMAM dendrimer is employed for the synthesis of cationic lipids based on its capability to condense the siRNA and its lower cytotoxicity than higher generations. In addition, we have previously synthesized cationic lipids using different cationic amino molecules, including G0, branched polyethylenimine (Mw ~ 800) (PEI), polypropylenimine tetramine dendrimer, generation 1 (DAB) and diethylene triamine. Among them, the G0-lipid shows the highest siRNA delivery efficiency.<sup>6</sup> We chose 1,2-epoxytetradecane for the synthesis of cationic lipid because lipid-like compounds with 14 carbons have been shown to be optimal for siRNA delivery.<sup>65</sup> The maximal amounts of lipid tails that can be conjugated to G0 are eight according to the four free amines of G0. To obtain G0-C14 with one less tail than the fully conjugated G0-C14, seven equivalents (eq.) of lipid tails (1,2-epoxytetradecane) are reacted with one eq. of G0. The obtained crude product can be purified by silica gel chromatography with gradient elution ( $\text{CH}_2\text{Cl}_2 \rightarrow 75:22:3 \text{ CH}_2\text{Cl}_2/\text{MeOH}/\text{NH}_4\text{OH}$ ). The chemical structure of the G0-C14 can be determined by  $^1\text{H}$  NMR (The purified G0-C14 contains about 7 lipid tails as calculated from the  $^1\text{H}$  NMR spectrum<sup>60</sup>).

The ligand-modified lipid-PEG (DSPE-PEG-S2P) is prepared by an efficient and widely used click reaction between the thiol group of the targeting peptide sequence: CRTLTVRKC (S2P) and the maleimide (Mal) group of Mal-modified DSPE-PEG (DSPE-PEG-Mal) (Fig. 3b). We use thiol-maleimide click chemistry for two reasons. First, the thiol group can be easily and conveniently introduced during the custom synthesis of the peptide. Hence, the peptide ligand can be conjugated to the lipid-PEG in one step. Second, the thiol and maleimide groups are stable in both organic and aqueous solutions, leading to high reaction efficiency. The resulting crude DSPE-PEG-S2P can be purified through washes with cold methanol. The successful conjugation of S2P to DSPE-PEG can be validated by  $^1\text{H}$  NMR. If a user wants to target another cell type or organ, we recommend using an alternative cysteine-terminated peptide or other targeting ligand in place of S2P.



**Preparation of nucleic acid-loaded NPs (Steps 25-52)**—The siRNA NPs are prepared by a robust nano-precipitation method that is well-established by our group. Due to the vulnerability of free siRNA to nucleases, the siRNA needs to be fully protected before it reaches its intended destination. To this end, our NPs are designed to have a core-shell structure, where the siRNA is loaded in the core and protected from enzymatic degradation. To fabricate the siRNA-containing core, we first condense the siRNA with cationic lipid G0-C14 through electrostatic interaction between the negatively charged siRNA and the positively charged G0-C14. However, the siRNA/G0-C14 complex might not be stable because the electrostatic interaction may be easily weakened by the high ionic strength condition of the body. Therefore, a hydrophobic polymer PLGA is used to cover the siRNA/G0-C14 complex via the hydrophobic interaction between the hydrophobic lipid tails of G0-C14 and the hydrophobic PLGA. It is worth noting that the cationic lipid G0-C14 plays a pivotal role in the fabrication of siRNA NPs. The cationic part of the G0-C14 should have enough positive charges to efficiently condense the negatively charged siRNA and the lipid part of the G0-C14 should have enough hydrophobicity to interact with the hydrophobic PLGA. To reduce the non-specific cellular uptake of the NPs and therefore prolong their blood circulation time *in vivo*, many neutral hydrophilic polymers have been used to modify the surfaces of the NPs.<sup>66</sup> Among them, PEG is the most popular polymer due to its low toxicity and simple chemistry. Based on this, we introduced a PEG shell to cover the siRNA/G0-C14/PLGA core to increase the stability and blood circulation time of the siRNA NPs. The PEGylation process can be easily conducted by dropwise addition of the organic solution containing the siRNA/G0-C14/PLGA core to the aqueous solution containing the DSPE-PEG shell. For the introduction of targeting ligand such as S2P to the NPs, researchers just need to add the DSPE-PEG-S2P (or other targeting ligand-modified lipid) to the aqueous solution during the preparation of NPs. The targeting capability of the siRNA NPs is enabled by the S2P ligand, therefore, a high DSPE-PEG-S2P/DSPE-PEG ratio of 1:1 is recommended to maximize the targeting effect of the siRNA NPs to macrophages. Here, the DSPE-PEG-S2P/DSPE-PEG ratio of 1:1 was used for all *in vitro* and *in vivo* experiments. An optimal G0-C14/siRNA weight ratio was used to ensure the condensation of the siRNA by G0-C14; an optimal PLGA/siRNA weight ratio was used to ensure the protection of siRNA by PLGA; an optimal lipid-PEG/siRNA weight ratio was used to ensure the formation of stable siRNA NPs. In addition, an optimal stirring speed was used to achieve NPs with high siRNA encapsulation efficiency and proper size; an optimal incubation time was used to make sure the formation of siRNA NPs and avoid the degradation of the siRNA.

To remove the small amount of organic solvent in the NP solutions, the NPs are further washed several times with water using an ultrafiltration purification method. We chose the ultrafiltration purification method based on the following considerations. First, it saves a lot of time since the purification can normally be done in one or two hours. Short operation time lowers the risk of degradation of siRNA during the purification process, and is especially important if the NP platform is used to deliver very unstable RNA therapeutics such as mRNA. Second, the concentration of the NPs can be easily adjusted. The ultrafiltration can result in a highly concentrated solution of NPs, which can be easily redispersed in any buffers to obtain NPs with any desired concentration according to different needs. This

is quite useful when the NPs are applied *in vivo*, where high concentrations of NPs are normally required. In addition, our siRNA NP solution can be frozen and kept in the freezer for a few weeks, and then thawed for use. The thawed siRNA NPs show no significant changes in activity, size and polydispersity index (PDI). We used a 20-mL glass vial to produce maximal 400  $\mu$ L of NPs (siRNA 5 $\mu$ M). If larger volume of NPs (e.g. 2 mL) is required, the researchers can easily produce 2 mL of NPs by using 5 glass vials (20-mL).

Our polymer-lipid hybrid NPs can also be used as efficient delivery platforms for mRNA.<sup>30,31</sup> However, slight modifications are necessary for the preparation of mRNA NPs compared to the siRNA NP formulation presented here. The mRNA molecules are much larger than the siRNA molecules, and therefore the mRNA/G0-C14 weight ratio needs to be optimized to ensure the sufficient condensation of the mRNA by the cationic lipids. Compared to the preparation of siRNA NPs, where the siRNA/G0-C14 ratio of 1:10 is required, the preparation of mRNA NPs needs a larger amount of cationic G0-C14 lipids (mRNA/G0-C14 weight ratio, 1:15). In addition, the stability of large and single-strand mRNA is much lower than the small and double-strand siRNA. Therefore, the working area and all the tools need to be cleaned carefully with RNaseZAP to remove all the RNase contamination before starting to work with mRNA. Meanwhile, it is important to maintain the mRNA and mRNA NPs at a low temperature of 4 °C or below since mRNA is prone to degradation at room temperature.

#### **Dynamic light scattering (DLS) characterization of siRNA NPs (Steps 53-58)—**

It has been previously reported that particles with diameters greater than ~200 nm cannot pass the spleen's fenestrations which are typically below 200-500 nm in width.<sup>67,68</sup> Thus, NPs with diameters in the range of 30-200 nm are ideal for systemic *in vivo* delivery because they can remain in circulation and avoid rapid renal clearance. To form siRNA NPs with ideal size, dynamic light scattering (DLS) is employed to monitor the size of each NP formulation. Most of the DLS instruments are able to detect NPs with sizes from 10 nm to 1 $\mu$ m within a few minutes. Normally, 1.5 mL of sample is enough for detection, and it can be recycled for other tests. DLS can also be used to measure the zeta potential of the NPs. Any surface modifications of the NPs can lead to changes in zeta potential. Thus, by monitoring the changes in zeta potentials of the NPs, we can determine whether the surface modifications and alterations in surface charge are successful. If different users need to reproduce the NPs according to this protocol, DLS can be used to quickly confirm the successful production of the NPs.

#### **Transmission electron microscope (TEM) characterization of siRNA NPs (Steps 59-68)—**

Based on our design for the siRNA NPs, a core-shell structure is important with regard to protecting the siRNA from enzymatic degradation and prolonging the blood circulation time. The transmission electron microscope (TEM) can clearly characterize the morphology and structure of the NPs. We thus use the TEM to confirm the core-shell structure of the siRNA NPs. In order to improve the contrast of the sample in TEM images, 1 % (wt/vol) of uranyl acetate is used to provide negative stain contrast in TEM images. The TEM images of the siRNA NPs in this protocol are produced using a JEOL 1200EX electron microscope at 80 kV.

**Gel electrophoresis of siRNA NPs (Steps 69-91)**—Our polymer-lipid hybrid NPs are designed to protect the loaded siRNA from enzymatic degradation. To confirm the protective effect of the NPs, the serum stability of siRNA NPs is assessed using agarose gel electrophoresis by monitoring the siRNA band. We incubate the siRNA-loaded NPs in a solution containing 50% (vol/vol) serum at 37 °C, a condition mimicking the *in vivo* environments, for different time periods. Theoretically, the enzymes in the serum can quickly degrade the siRNA at 37 °C upon direct exposure, resulting in the fading or disappearance of the siRNA band in agarose gel electrophoresis. In order to directly compare the siRNA that is loaded in the NPs with the free siRNA, the siRNA in the NPs is isolated before loading onto the agarose gel. To isolate the siRNA from the NPs, siRNA NPs are first pelleted using Amicon Ultra-15 centrifugal filter units (MWCO, 100 kDa). Then, the pellet is dissolved in chloroform to destroy the core-shell structure of the NPs, releasing the incorporated siRNA which is extracted from chloroform using 0.5 M NaCl containing 0.1% (wt/vol) SDS.<sup>6</sup> The siRNA protection efficacy of siRNA by the polymer-lipid hybrid NPs is determined by gel electrophoresis using ethidium bromide-infused E-Gel 4% (wt/vol) agarose gels which are imaged under ultraviolet (UV) light.

**Determination of siRNA encapsulation efficiency (EE) (Steps 92-101)**—Fluorescence detection is a highly sensitive technique widely used for the detection of compounds in solution even at very low concentrations. To determine the siRNA encapsulation efficiency (EE) of the NPs, we label the siRNA with a fluorescent dye Dy647 (Dy647-siRNA). Alternatively, the siRNA EE of the NPs can be directly determined using ribogreen, one of the most sensitive fluorescent dyes for the detection of RNA in solution. Instead of ribogreen, we used dye-labeled siRNA (Dy647-siRNA) as this was also required for other experiments. The Dy647-siRNA is loaded into the NPs following the approaches we described above. The unloaded Dy647-siRNA can be easily removed in the ultrafiltration purification step. The encapsulation of fluorescent molecules into the cores of the NPs can result in changes in their fluorescence signals. Therefore, the Dy647-siRNA needs to be released from the NPs before it can be precisely determined by fluorescence detection. To this end, Dy647-siRNA-loaded NPs are dissolved in organic solvent dimethyl sulfoxide (DMSO), a solvent that can destroy the NPs by dissolving all the NP components, to release the encapsulated Dy647-siRNA. The fluorescence of the free siRNA can then be measured using a fluorescent plate reader, such as the Infinite M200 PRO spectrophotometer. The initial amount of Dy647-siRNA loaded into the NPs is also dissolved in DMSO and its fluorescence is measured. The fluorescence of only DMSO is also measured and used as a blank control. All the measurements are conducted in triplicate. The siRNA EE of NPs can be calculated according to the following equation:

$$\% EE = 100 \times (\text{sample} - \text{blank}) / (\text{load} - \text{blank}) .$$

**Evaluation of the stability of siRNA NPs (Steps 102-108)**—The stability of the siRNA NPs in biological media is critical to their long blood circulation, which allows the accumulation of siRNA NPs at lesional sites. To determine the stability of NPs *in vitro*, NPs are incubated with PBS containing 10% (vol/vol) serum (mimicking the *in vivo* environments) and their sizes are measured by DLS after different incubation time periods

(up to 3 days). All the measurements are conducted in triplicates. The siRNA NPs after incubation for different time periods are considered to be stable at these conditions if their sizes show no significant changes compared to that of the freshly prepared NPs.

**Characterization of endosomal escape of siRNA NPs (Steps 109-121)**—The inefficient endosomal escape of siRNA NPs is reported to be a limiting factor for efficient siRNA delivery.<sup>57,69</sup> Therefore, the endosomal escape ability of NPs can increase siRNA delivery efficiency. To investigate the endosomal escape ability of our siRNA NPs, we first label NPs via labeling the encapsulated siRNA with a fluorescent dye Dy647 (red). So, we can track the locations of siRNA NPs in the cells (we used HeLa cells) by monitoring the red fluorescence. Then, we use LysoTracker Green (green) to stain the late endosome whose locations in cells can be tracked via monitoring the green fluorescence. Upon incubation with Dy647-siRNA NPs for different time periods (1 or 4 h), the HeLa cells are first stained with LysoTracker Green for late endosome staining and then stained with Hoechst 33342 (blue) for nuclei staining. The cells are observed and imaged using an FV1000 laser scanning confocal microscope. The overlap of the red and green fluorescence signals of cells suggests that siRNA NPs are trapped inside the endosomes. The successful endosomal escape of siRNA NPs is defined by the separation of the red and green fluorescence signals of the cells.

***In vitro* evaluation of gene silencing by siRNA NPs (Steps 122-131)**—Before the evaluation of siRNA NPs *in vivo*, we assess the gene silencing ability of siRNA NPs *in vitro*. To this end, one of the most widely used reporter gene systems based on luciferase-expressing HeLa (HeLa-Luc) cells is used for gene silencing assessment. Normally, the gene silencing efficiency of siRNA NPs can be determined by measuring the protein expression levels using western blot analysis, which is relatively costly and time-consuming. The advantage of the HeLa-Luc system is that the gene silencing efficiency of the siRNA NPs can be quickly determined by measuring the bioluminescence using a plate reader or an In-Vivo Xtreme imaging system with a charge-coupled device (CCD) camera. A luciferase substrate must be added to the cells and incubated for 5 min before the measurements. Both quantitative data and images can be generated and used for gene silencing efficiency analysis. The gene silencing of our siRNA NPs is evaluated by incubating the siLuc NPs with HeLa-Luc cells for 24 h in 12-well plates using RPMI 1640 medium containing 10% (vol/vol) FBS as the culture medium. Upon 24 h incubation, the bioluminescence of the cells is measured using an In-Vivo Xtreme imaging system with a CCD camera. The bioluminescence signal of cells treated with PBS is measured and used as a PBS control. The treatment of PBS should not lead to any silencing effect of the *Luc* gene. All the measurements are conducted in triplicate. To quantify the gene silencing efficiency, the average radiance (photons per second per cm<sup>2</sup> per steradian) within regions of interest is quantified using Bruker MI SE software.

***In vitro* viability, apoptosis, and proliferation assessment of siRNA NPs (Steps 132-164)**—Biocompatibility is a prerequisite in the development of NPs as safe *in vivo* siRNA delivery systems. In addition to the *in vivo* safety issues, the toxicity of the delivery systems can also interfere with the *in vitro* evaluations. For example, a false gene silencing

effect can be achieved after incubating toxic siLuc NPs with HeLa-Luc cells. This can be ascribed to the fact that the toxicity of siLuc NPs can also decrease Luc expression by suppressing the growth of HeLa-Luc cells. To evaluate the safety of our siRNA NPs, we first investigate their *in vitro* toxicity. Different cells, including HeLa-Luc, RAW 264.7, and HEK-293 cells, are incubated with various amounts of siLuc NPs in 96-well plates for 24 h. After which, the medium of the cells is replaced with a fresh medium and incubated for another 48 h. After the addition of the alamarBlue cell viability and proliferation reagent, the viability of the cells can be determined by measuring the absorbance at 570 nm using a plate reader. The absorbance of alamarBlue-treated medium is also measured and used as a blank control. The absorbance of cells treated with PBS is measured and used as a PBS control. The treatment of PBS should not lead to any cytotoxicity in the cells. All the measurements are conducted in 5 replicates. The viability of NP-treated cells can be calculated according to the following equation:

$$\% \text{ Cell viability} = 100 \times (\text{sample} - \text{blank}) / (\text{PBS} - \text{blank}) .$$

Some reagents can induce cell apoptosis, especially early apoptosis, which may not be detected by a cell viability reagent. We therefore also assess the apoptosis of cells treated with our siRNA NPs. RAW 264.7 cells are incubated with siLuc NPs containing 50 nM siLuc in 24-well plates for 24 h. After which, the medium of the cells is replaced with a fresh medium and incubated for another 48 h. After the double staining with Annexin V-FITC and propidium iodide (PI) (Annexin V-FITC Apoptosis Detection kit), the apoptosis of the cells can be determined by flow cytometry. The apoptosis of cells treated with PBS is measured and used as a PBS control. The treatment of PBS should not lead to any significant apoptosis to the cells. All the measurements are conducted in triplicates. The early apoptosis of cells is defined by Annexin V-FITC positive and PI negative; and the late apoptosis or dead cells are defined by both FITC Annexin V and PI positive.

To investigate whether the siRNA NPs have an adverse effect on the growth of cells that is not detectable in the short term, we evaluate the proliferation of siRNA NPs-treated cells over 7 days. RAW 264.7 cells are incubated with siLuc NPs containing 50 nM siLuc in 96-well plates for 24 h. After which, the medium of the cells is replaced with a fresh medium and incubated for another 7 days. After the addition of the alamarBlue cell viability and proliferation reagent, the proliferation of cells can be determined by measuring the absorbance at 570 nm using a plate reader. Notably, alamarBlue is a non-toxic agent that can be applied repeatedly on the same samples, allowing continuous monitoring of the samples.<sup>70</sup> Thus, the proliferation of cells is monitored every day and the medium containing alamarBlue reagent is replaced by a fresh medium after each measurement. The absorbance of alamarBlue-treated medium is also measured and used as a blank control. The absorbance of cells treated with PBS is measured and used as a PBS control. The treatment of PBS should not lead to any inhibition to the growth of cells. All the measurements are conducted in 5 replicates. The proliferation of the NP-treated cells can be calculated according to the following equation:

$$\% \text{ Relative cell proliferation} = 100 \times (\text{sample}_{\text{day } n} - \text{sample}_{\text{day } 1}) / (\text{sample}_{\text{day } 1} - \text{blank}), n = 1-7 .$$

***In vivo* evaluation of the therapeutic effect of siRNA NPs (Steps 165-189)**—To assess the *in vivo* therapeutic effect of siRNA NPs for the treatment of atherosclerosis, we first establish a mouse model of atherosclerosis. It is well-documented that the *Ldlr* gene-deficient mice can easily develop atherosclerotic plaques without any surgeries or injections by only feeding the mice with a high-fat Western diet (WD) for 6 to 8 weeks. Therefore, eight- to ten-week-old male mice with deficient *Ldlr* gene (*Ldlr*<sup>-/-</sup>) are obtained from Jackson laboratory, and fed with WD for 8 weeks to develop the atherosclerotic plaques. It is worth noting that the influence of sex may need to be considered in many mouse models for human diseases.<sup>71</sup> We have used eight- to ten-week-old male mice with deficient *Ldlr* gene (*Ldlr*<sup>-/-</sup>) for the evaluation of the therapeutic effect of siRNA NPs. However, we speculate that this procedure would also be applicable to female mice since all the atherosclerosis-related phenotypes, including cholesterol, high-density lipoprotein (HDL), low-density lipoprotein (LDL), triglycerides, and Apolipoprotein B (ApoB), elevate similarly both in male and female *Ldlr* knockout mice as demonstrated by Jackson laboratory (supplier of lab mice). Therefore, sex should not have a significant influence on the results generated by this protocol. The mice are randomly divided into 3 groups (with 7-9 mice per group) and injected with two doses of PBS, control siRNA NPs or siCamk2g NPs per week, respectively. After 4 weeks of treatment, the mice are euthanized with CO<sub>2</sub> and their hearts with attached aortic roots are either fixed in 4% (wt/vol) paraformaldehyde followed by paraffin embedding or treated with optical cutting temperature (OCT) compound and frozen. Sections of the aortic root with different sizes are isolated from the paraffin-embedded or frozen samples for the analysis of atherosclerotic lesions.

According to the proposed therapeutic mechanism of our siCamk2g NPs (Fig. 1), the treatment of siCamk2g NPs can knock down the expression of CaMKII $\gamma$ , resulting in the activation of the MerTK pathway in macrophages of the atherosclerotic lesions. To validate this mechanism, the expressions of CaMKII $\gamma$  and MerTK in lesional macrophages are characterized by immunofluorescence imaging. After blocking with bovine serum albumin (BSA), the obtained aortic root sections are first stained with primary anti-Mac2 and anti-MerTK and subsequently stained with fluorescently labeled secondary antibodies after PBS rinses. The nuclei of the macrophages are further stained with 4',6-diamidino-2-phenylindole phenylindole (DAPI). The stained sections are imaged by a Zeiss epifluorescence microscope and the images are analyzed by FIJI software. For statistical analysis, n = 7 to 9 mice per group.

Then, to further validate that the activation of the MerTK pathway can promote the efferocytosis of the lesional macrophages, the apoptotic cells in aortic root sections are detected by a TUNEL (terminal deoxynucleotidyl transferase-mediated deoxyuridine triphosphate nick-end labeling) assay, a method for detecting apoptotic cells, using immunofluorescence imaging. The aortic root sections are first stained with TUNEL staining reagents and washed with PBS. Then, the sections are blocked with BSA, stained with primary anti-Mac2 antibody, incubated with fluorescently labeled secondary antibodies and stained with DAPI. The stained sections are imaged using a Zeiss epifluorescence microscope and the images are analyzed using FIJI software. The apoptotic cells (TUNEL<sup>+</sup>) should be surrounded by Mac2<sup>+</sup> macrophages when efferocytosis occurs. In contrast,

isolated apoptotic cells have condensed nucleus and low antibody reactivity and are not associated with neighboring macrophages. Thus, the lesional efferocytosis can be quantified as the ratio of TUNEL<sup>+</sup> cells (red) in contact with Mac2<sup>+</sup> macrophages (green) to isolated TUNEL<sup>+</sup> cells. For statistical analysis, n = 7 to 9 mice per group.

Finally, to confirm that the promoted efferocytosis can improve the plaque stability by reducing plaque necrosis and inhibiting fibrous cap thinning, the necrotic core area and fibrous cap thickness of the plaques in aortic root sections are determined by picosirius red staining-based microscopy imaging. The picosirius red staining is a widely used histological method to visualize collagen in paraffin-embedded tissue sections. Thus, the collagen of aortic root sections is stained with picosirius red and imaged with light microscopy. The necrotic core area can be obtained directly from the images; the fibrous (collagen) cap thickness is quantified from one midpoint and two shoulder regions of the plaque and normalized to the lesion size.<sup>72</sup> For statistical analysis, n = 7 to 9 mice per group.

***In vivo* evaluation of the toxicity of the siRNA NPs (Steps 190-198)**—The biosafety of the siRNA NPs is essential for their use in *in vivo*. We therefore investigate the *in vivo* toxicity of the siRNA NPs under the same conditions (e.g. dose, frequency) as the therapeutic effect assay. Six-week-old female BALB/c mice obtained from Jackson laboratory are injected intravenously with two doses of PBS, or S2P<sub>50</sub> siCamk2g NPs (1 nmol siRNA per injection) per week for 4 weeks. The mice are euthanized by CO<sub>2</sub> inhalation three days after the last administration. To evaluate whether the siRNA NPs can cause any damage to the organs of the mice, the major organs are fixed in 4% (wt/vol) paraformaldehyde followed by paraffin embedding. Sections with different sizes are isolated from the paraffin-embedded samples and subjected to H&E staining. In addition, the blood from the mice is collected and various hematological parameters are evaluated, including alanine aminotransferase (ALT), aspartate aminotransferase (AST), blood urea nitrogen (BUN), hematocrit (HCT), hemoglobin (Hb), lymphocyte count (LY), mean corpuscular volume (MCV), mean corpuscular hemoglobin (MCH), mean corpuscular hemoglobin concentration (MCHC), neutrophil count (NPC), red blood cell count (RBC), and white blood cell count (WBC). For statistical analysis, n = 3 mice per group.

**Recommended controls in the procedures**—We recommend preparing the following control samples:

- Step 71 Untreated extracted siRNA control: extracted siRNA from siRNA S2P<sub>50</sub> NPs without incubation with serum.
- Step 73 Untreated siRNA control: free siRNA without incubation with serum (incubate with PBS instead).
- Step 99 Blank control: transfer 90  $\mu$ L of 95% (vol/vol) DMSO solution to the 96-well plate.
- Step 107 Untreated siRNA NPs control: the freshly prepared siLuc NPs without incubation with PBS containing 10% (vol/vol) serum.

- Step 124 PBS control (or 0 nM siLuc): for the bioluminescence signal of cells treated with PBS.
- Step 132 Blank control: the absorbance of medium without any cells.
- Step 134 PBS control: the absorbance of cells treated with PBS.
- Step 145 PBS control: the apoptosis of cells treated with PBS.
- Step 151 Blank control: unstained cells; AV-FITC control: cells stained with AV-FITC (no PI); PI control: cells stained with PI (no AV-FITC).
- Step 153 Blank control: the absorbance of medium without cells.
- Step 176 PBS treated control: gene silencing of sections from mice treated by PBS; siCtrl NPs treated control: gene silencing of sections from mice treated by S2P<sub>50</sub> siCtrl NPs.
- Step 185 PBS treated control: efferocytosis of sections from mice treated by PBS; siCtrl NPs treated control: efferocytosis of sections from mice treated by S2P<sub>50</sub> siCtrl NPs.
- Step 189 PBS treated control: plaque stability of sections from mice treated by PBS; siCtrl NPs treated control: plaque stability of sections from mice treated by S2P<sub>50</sub> siCtrl NPs.
- Step 198 PBS treated control: H&E staining of sections from mice treated by PBS; hematological parameters of the blood from mice treated by PBS.

### Expertise needed to implement the protocol

This protocol describes all the procedures for each experiment in detail and is easy to follow. Any researchers trained in nanotechnology, gene delivery, or biochemistry should be able to implement the protocol.

### Limitations

The highly stable lipid-PEG shell of the NPs can reduce their cellular uptake, making them less suitable for *in vitro* or *ex vivo* applications, especially in situations where fast gene transfection is required. One solution to this issue is the introduction of targeting ligands. As demonstrated in this protocol, an S2P ligand is employed to improve the specific uptake of NPs by macrophages.<sup>26</sup> Another solution to this issue is the replacement of lipid-PEG. We have previously established that lipid-PEGs with long and saturated lipophilic tails had a long dissociation half-life ( $t_{1/2}$ ) when presented on NPs.<sup>73</sup> Therefore, the cellular uptake efficiency of the NPs can be tuned by using lipid-PEGs with short or unsaturated lipophilic tails. Despite the encouraging preclinical results that have been achieved using our polymer-lipid hybrid NP platform, the translation of our platform to clinical application needs further optimization. For example, the cationic lipids in our platform are biocompatible with most of the applications. However, in settings that require a high concentration of cationic lipids for NP synthesis, toxicity becomes an issue. Although negligible cytotoxicity was observed in cells treated with our NPs (siRNA 6.25-50 nM), we speculated that cytotoxicity > 20 % could be observed with siRNA doses > 50 nM (Supplementary Fig. 4a). This issue might be



addressed by developing novel stimuli-responsive cationic lipids which can be degraded in the presence of various internal stimuli.<sup>74</sup>

## Materials

### Biological Materials

- Mouse model of interest. ! **CAUTION** All experiments involving live animals must conform to the requirements of your institutional animal care and use committee (IACUC) and be performed strictly following protocols that are approved by IACUC. All our procedures were conducted in accordance with the protocols reviewed and approved by IACUC at Harvard Medical School and Columbia University Irving Medical Center. To produce the results shown here, we used:
  - *Ldlr*<sup>-/-</sup> mice, eight- to ten-week-old male mice (The Jackson Laboratory, B6.129S7-*Ldlr*<sup>tm1Her</sup>/J, Stock No: 002207)
  - BALB/c mice, six-week-old female BALB/c mice (The Jackson Laboratory, BALB/cJ, Stock No: 000651)
- Cells of interest. ! **CAUTION** It is important to make sure that all the cell lines used in your research have no issues of cell line misidentification and cross-contamination. Our cell lines used in all the procedures had no issues of cell line misidentification and cross-contamination according to the latest report provided by International Cell Line Authentication Committee (ICLAC). <http://iclac.org/databases/cross-contaminations/>. ! **CAUTION** All the cell lines used in your research should be tested periodically to ensure they are authentic and not infected with mycoplasma. Our cell lines used in all the procedures were free of mycoplasma as determined by Mycoplasma Detection Kit (Thermo Fisher Scientific, M7006). To produce the results shown here, we used:
  - HeLa-Luc cells (a HeLa cell line stably transfected with a firefly luciferase reporter gene, Sigma-Aldrich, cat. no. 11033106). <https://scicrunch.org/resolver/RRID><sup>75</sup>: CVCL\_2939. (ECACC Cat# 11033106, RRID:CVCL\_2939)
  - RAW 264.7 (ATCC, TIB-71<sup>TM</sup>). <https://scicrunch.org/resolver/RRID>: CVCL\_0493. (ATCC Cat# TIB-71, RRID:CVCL\_0493)
  - HEK-293 (ATCC, CRL-1573<sup>TM</sup>). <https://scicrunch.org/resolver/RRID>: CVCL\_0045. (ATCC Cat# CRL-1573, RRID:CVCL\_0045)

### Reagents

- G0 (amine) (Sigma-Aldrich, cat. no. 412368, CAS. no. 155773-72-1)
- Methanol (MeOH, Sigma-Aldrich, cat. no. 322415, CAS. no. 67-56-1). ! **CAUTION** Methanol is vaporizable and harmful to naked eyes. Work with personal protective equipment when handling methanol (e.g. safety goggles and face shields).

- 1,2-epoxytetradecane (epoxide) (Sigma-Aldrich, cat. no. 260266, CAS. no. 3234-28-4). ! **CAUTION** 1,2-epoxytetradecane can cause severe damage to the respiratory system. Perform all operations under a chemical fume hood.
- Dichloromethane (CH<sub>2</sub>Cl<sub>2</sub>, Sigma-Aldrich, cat. no. 34856, CAS. no. 75-09-2). ! **CAUTION** Transferring dichloromethane that involves cannulas and needles (syringe) should be careful. Piercing the needles containing a residual amount of dichloromethane into tissues can cause severe injuries.
- Ammonium hydroxide solution (NH<sub>4</sub>OH, Sigma-Aldrich, cat. no. 338818). ! **CAUTION** Ammonium hydroxide solution can cause damage to the respiratory system. Handle this reagent under a chemical fume hood.
- Silica gel (high-purity grade, average pore size 60 Å, 70-230 mesh, 63-200 µm, cat. no. 60741-1KG)
- Silica-gel-coated plates (Millipore Sigma, cat. no. MC1156840001)
- DSPE-PEG-Mal (Nanocs Inc., cat. no. PG2-DSML-3k, 3400 Da)
- DSPE-PEG (Avanti Polar Lipids, cat. no. 880320, 3400 Da)
- Targeting peptide CRTLTVRKC (S2P) (GLS Biochem Systems Inc., cat. no. 65235)
- Nuclease-free water (HyPure water) (Fisher Scientific, cat. no. SH3053801)
- RNaseZAP (Fisher Scientific, cat. no. AM9785)
- Acetone (Sigma-Aldrich, cat. no. 650501, CAS. no. 67-64-1). ! **CAUTION** Flammable and vaporizable acetone is harmful to the central nervous system. Work with acetone under a fume hood.
- *N,N*-Dimethylformamide (DMF, Sigma-Aldrich, cat. no. 227056, CAS. no. 68-12-2) ! **CAUTION** *N,N*-Dimethylformamide may cause abdominal pain, headache, skin problems, constipation, nausea and vomiting. Wear personal protective equipment (e.g. eyeshields, face shields, lab coats and gloves) during operation.
- 50:50 Poly(DL-lactide-*co*-glycolide) (PLGA, LACTEL Absorbable Polymers, cat. no. B6013-2, Inherent Viscosity dL/g: 0.55 - 0.75)
- 1X PBS (Fisher Scientific, cat. no. SH3025601, pH 7.4)
- NaOH solution (1M, Sigma-Aldrich, cat. no. 1091371000 )
- HCl solution (1M, Sigma-Aldrich, cat. no. 1506961000)
- UAR-EMS Uranyl acetate (Electron Microscopy Sciences, cat. no. 22405)
- Fetal bovine serum (FBS, Thermo Fisher Scientific, cat. no. 26140095)
- Chloroform (CHCl<sub>3</sub>, Sigma-Aldrich, cat. no. C2432, CAS. no. 67-66-3). ! **CAUTION** Chloroform is toxic and can cause severe health problems, including

central nervous system damage, irritation of skin and eyes, etc. Perform all operations under a chemical fume hood.

- Ethidium bromide-infused E-Gel 4% (wt/vol) agarose gel (Thermo Fisher Scientific, cat. no. G501804)
- 10X BlueJuice Gel Loading Buffer (Thermo Fisher Scientific, cat. no. 10816015)
- Dimethyl sulfoxide (DMSO, Sigma-Aldrich, cat. no. 276855, CAS. no. 67-68-5).
- RPMI 1640 medium (Thermo Fisher Scientific, cat. no. 11875119)
- Penicillin-Streptomycin (10,000 U/mL Penicillin, 10,000 µg/mL Streptomycin, Thermo Fisher Scientific, cat. no. 15140122)
- Hoechst 33342 (Thermo Fisher Scientific, cat. no. H1399)
- LysoTracker™ Green (Thermo Fisher Scientific, cat. no. L7526)
- Steady-Glo Luciferase reagent (Promega Corporation, cat. no. E2520)
- Dulbecco's Modified Eagle Medium (DMEM, Thermo Fisher Scientific, cat. no. 11965092)
- AlamarBlue cell viability reagent (Thermo Fisher Scientific, cat. no. DAL1100)
- Annexin V-FITC Apoptosis Detection kit (Thermo Fisher Scientific, cat. no. V13245)
- “Western” purified atherogenic diet (Envigo, TD.88137 Adjusted calories diet, 42% kcal from fat)
- Paraformaldehyde Solution, 4% (wt/vol) in PBS (Fisher Scientific, cat. no. AAJ19943K2)
- Bovine serum albumin (BSA, Sigma-Aldrich, cat. no. A2153, CAS. no. 9048-46-8)
- 10X Tris Buffered Saline (10X TPS, Bio-Rad Laboratories, Inc. cat. no. 1706435)
- Tween 20 (Sigma-Aldrich, cat. no. P2287)
- Goat serum (Sigma-Aldrich, cat. no. G9023)
- Triton™ X-100 (Sigma-Aldrich, cat. no. X100)
- 4',6-diamidino-2-phenylindole phenylindole (DAPI, Sigma-Aldrich, cat. no. D8417, CAS. no. 28718-90-3)
- Picrosirius red (Polysciences, cat. no. 24901A)
- Anti-luciferase RNA (siLuc), Dy647-labelled siLuc (Dy647-siLuc), and anti-CaMKIIγ RNA (siCamk2g) (Dharmacon, custom siRNA, siRNA sequence can be found in Table 1).

## Equipment

- Heating plate (Isotemp™ Hot Plate, Fisher Scientific, cat. no. HP88857296). ! **CAUTION** Avoid contacting heated part of the plate during heating operation, which can cause severe burns.
- 25-mL round bottom flask (Chemglass, cat. no. CG-618-01)
- 500-mL round-bottom flask (Chemglass, cat. no. CG-618-09)
- Blood collection tube (Greiner bio-one, cat. no. 454012)
- Amicon Ultra-15 centrifugal filter units (Sigma-Aldrich, cat. no. UFC9003, Mw cutoff, 3 kDa)
- 10-mL glass vial (Chemglass, cat. no. CG-4908-05).
- 20-mL glass vial (Chemglass, cat. no. CG-4908-01)
- Gel loading tips (Fisher Scientific, cat. no. 02-707-139)
- Carbon support copper grids (Electron Microscopy Sciences, cat. no. FCF200-CU)
- 96-well plate (Fisher Scientific, cat. no. CLS3595)
- Nunc glass-bottom dish (Thermo Fisher Scientific, cat. no. 150680)
- 12-well plate (Fisher Scientific, cat. no. CLS3512)
- 6-well plate (Fisher Scientific, cat. no. CLS3516)
- Syringes (BD, 1/2 cc insulin syringe U-100)
- Oil bath
- Rotary evaporator (Buchi Labortechnik AG)
- Centrifuge (Thermo Scientific Sorvall Legend X1R Centrifuge cat. no. 75004263)
- Orion Star™ A221 Portable pH Meter (Sigma-Aldrich, cat. no. STARA221)
- Nanodrop spectrophotometer (Thermo Fisher Scientific, NanoDrop™ 1000 Spectrophotometer)
- Sonication bath (Fisher Scientific, cat. no. FB11207)
- Digital Vortexer (Fisher Scientific, cat. no. 0215370)
- DLS measurement (Brookhaven Instruments Corporation)
- Transmission electron microscopes (JEOL 1200EX)
- E-Gel Power Snap Electrophoresis Device (Thermo Fisher Scientific, cat. no. G8100)
- Confocal laser scanning microscope (Olympus FLUOVIEW FV1000)
- Flow cytometer (BD LSR Fortessa)

- Bioluminescence imaging system (Bruker In-Vivo Extreme II Optical/X-ray system)
- Plate reader (Tecan, Infinite<sup>®</sup> M200 PRO)
- GraphPad (<https://www.graphpad.com>)
- FIJI software (<https://imagej.net/Fiji>).

### Reagent setup

- HEPES solution (Thermo Fisher Scientific, 1 M, pH 7.0-7.6, cat. no. 15630080):  
To adjust the pH of HEPES solution to pH 6.5, introduce HCl while monitoring the pH with a pH meter. If the pH goes too low, raise it by carefully adding NaOH. ! **CAUTION** Users should exercise extreme caution when working with NaOH and HCl. Wear gloves, eye protection, face shields, and lab coats during operation. The buffer can be stored at 4 °C for up to 6 months.
- RPMI 1640 cell culture medium:  
To prepare the cell culture medium, add 50 mL FBS and 5 mL Penicillin-Streptomycin into 450 mL RPMI 1640 medium. The as-prepared cell culture medium contains 10% (vol/vol) FBS, 100 U/mL penicillin and 100 µg/mL streptomycin. The medium can be stored at 4 °C for up to 6 weeks.
- 1X Tris Buffered Saline with 0.1% (vol/vol) Tween 20:  
Dilute 100 mL 10X Tris Buffered Saline with 900 mL deionized water, and add 1 mL Tween 20. A brief sonification will help to homogeneously mix the as-prepared 1X Tris Buffered Saline with 0.1% (vol/vol) Tween 20. The buffer can be stored at 4 °C for up to 3 months.
- Blocking buffer (1× PBS containing 5% (vol/vol) normal serum and 0.3% (vol/vol) Triton<sup>™</sup> X-100):  
To prepare a 10 mL blocking buffer, add 0.5 mL normal goat serum into 9.5 mL 1× PBS. While stirring, add 30 µl Triton<sup>™</sup> X-100 into the above solution. The buffer should be freshly prepared before use.
- Antibody Dilution Buffer (1× PBS containing 1% (wt/vol) BSA and 0.3% (vol/vol) Triton X-100):  
To prepare a 10 mL antibody dilution buffer, add 30 µl Triton<sup>™</sup> X-100 into 10 mL 1× PBS. Then add 0.1 g BSA to the solution. Make sure the solution is homogeneously mixed. The buffer should be freshly prepared before use.
- Anti-CaMKIIγ (1:400, Novus Biologicals, cat. no. NBP2-15685). [https://scicrunch.org/resolver/RRID: AB\\_2892990](https://scicrunch.org/resolver/RRID:AB_2892990). (Novus, Cat# NBP2-15685, RRID:AB\_2892990)  
  
Dilute the primary antibody with antibody dilution buffer. The buffer should be freshly prepared before use.

- Anti-Mac2 antibody (1:10,000, Cedarlane Labs, cat. no. CL8942AP). [https://scicrunch.org/resolver/RRID: AB\\_2814900](https://scicrunch.org/resolver/RRID:AB_2814900). (Fluidigm Cat# 3153026B, RRID:AB\_2814900)  
Dilute the primary antibody with antibody dilution buffer. The buffer should be freshly prepared before use.
- Anti-MerTK (1:500, R&D Systems, cat. no. AF591). [https://scicrunch.org/resolver/RRID: AB\\_2098565](https://scicrunch.org/resolver/RRID:AB_2098565). (R and D Systems Cat# AF591, RRID:AB\_2098565).  
Dilute the primary antibody with antibody dilution buffer. The buffer should be freshly prepared before use.
- Alexa Fluor 647 Goat-anti Mouse IgG (Life Technologies, A-28181). [https://scicrunch.org/resolver/RRID: AB\\_2536165](https://scicrunch.org/resolver/RRID:AB_2536165). (Thermo Fisher Scientific Cat# A28181, RRID: AB\_2536165).  
Dilute the second antibody with antibody dilution buffer. The buffer should be freshly prepared before use.
- Mouse model of atherosclerosis:  
To establish a mouse model of atherosclerosis, purchase eight- to ten-week-old male mice with deficient Ldlr gene (*Ldlr*<sup>-/-</sup>) from Jackson laboratory. Then, feed the mice with a high-fat Western diet (WD) for 8 weeks to develop atherosclerotic plaques.<sup>26</sup>

## Procedure

### Synthesis of cationic lipid (G0-C14) ● Timing ~ 3 d

**CRITICAL** This procedure should result in about 1.9 g of purified cationic lipid.

1. Weigh 1.0 g (1.936 mmol, 1 eq.) of G0 (amine) and add it to a 25-mL round bottom flask.  
**▲ CRITICAL STEP** The G0 is 20 % (wt/vol) in methanol, so 1.0 g of G0 is equal to  $1.0/0.2 = 5.0$  g G0 solution.
2. Weigh 2.878 g (13.552 mmol, 7 eq.) of 1,2-epoxytetradecane (epoxide) and add it to the round bottom flask containing G0.  
**▲ CRITICAL STEP** The 2-epoxytetradecane has a purity of 85% (wt/wt), so 2.878 g of 1,2-epoxytetradecane is equal to  $2.878/0.85 = 3.386$  g crude 1,2-epoxytetradecane. The ratio of the amines (G0) and hydrocarbon tails (C14) are critical for the construction of our platform. G0 has four primary amines that can maximally react with 8 lipid tails (1,2-epoxytetradecane). The feeding ratio of epoxide and G0 is 7:1 here because we want to form products with 7 lipid tails as these materials display high siRNA delivery efficiency.
3. Add a magnetic bar to the round bottom flask and equip the flask with a drying tube.

4. Place the flask in an oil bath supported by a hot-plate magnetic stirrer and start stirring the reaction solution.
5. Set the heating temperature to 90 °C and start heating the solution.
6. After 30 min, when all of the solvent (methanol) evaporates, replace the drying tube with an argon balloon.

▲ **CRITICAL STEP** Argon protection is mandatory due to the involvement of amine in the reaction.

7. Stir the reaction solution at 90 °C for 48 h under argon atmosphere.

▲ **CRITICAL STEP** The viscosity of the reaction mixture increases with the evaporation of the methanol. Ensure your magnetic stir bar is large enough so that it continues stirring in the increased viscosity solution.

8. After the completion of the reaction, remove the balloon and the magnetic stir bar.

**CRITICAL STEP** This is a solvent-free reaction, no need to remove solvent after the reaction.

9. Add 2.0 mL of CH<sub>2</sub>Cl<sub>2</sub> to dissolve the reaction mixture.
10. Purify the product by flash silica gel column chromatography (silica gel, high-purity grade). Load the crude product solution to the silica column (diameter 40 mm, height 400 mm), remove the unreacted starting materials using CH<sub>2</sub>Cl<sub>2</sub>, and then elute the product using 75:22:3 CH<sub>2</sub>Cl<sub>2</sub>/MeOH/NH<sub>4</sub>OH. The purified G0-C14 contains 7 lipid tails as calculated from the <sup>1</sup>H NMR spectrum<sup>60</sup>

? TROUBLESHOOTING

11. Collect the eluents in, 20-mL fractions using 30-mL tubes, around 25 fractions should be collected. Monitor the fractions by thin-layer chromatography (TLC) using silica-gel-coated plates and 75:22:3 CH<sub>2</sub>Cl<sub>2</sub>/MeOH/NH<sub>4</sub>OH as eluent (R<sub>f</sub> is about 0.3 for product). Visualize the eluted product on the TLC plate by I<sub>2</sub>.
12. Product is expected to appear between fraction No.5 and No.10. Combine and transfer all the fractions containing pure product to a clean 500-mL round-bottom flask. Remove all the solvents using a rotary evaporator at 150 mbar and 25 °C in a water bath. The maximum temperature can be used is 45 °C and it normally takes 1-2 h to get the oily product.
13. Place the flask under a high vacuum for 24 h to remove any residual solvent.
14. Weigh the pure G0-C14, calculate the yield and store it at -20 °C.

? TROUBLESHOOTING

■ **PAUSE POINT** The pure G0-C14 can be store at -20 °C for at least 1 year.

### Synthesis of ligand-modified lipid-PEG (DSPE-PEG-S2P) ● Timing ~ 3 d

**CRITICAL** This procedure should result in about 105 mg of purified ligand-modified lipid.

- 15 Weigh 200.0 mg (0.059 mmol, 1.0 eq.) of DSPE-PEG-Mal and add it to a 10-mL round bottom flask.
- 16 Weigh 63.7 mg (0.059 mmol, 1.0 eq.) of targeting peptide CRTLTVRKC (S2P) and add it to the round bottom flask containing DSPE-PEG-Mal.
- 17 Add 2 mL of 1 × HEPES buffer (pH 6.5) to dissolve the starting materials.  
**CRITICAL STEP** Controlling the pH by HEPES buffer is important for the dissolution of S2P peptide and maximal reaction efficiency.
- 18 Add a magnetic bar to the round bottom flask and equip the flask with an argon balloon.
- 19 Place the flask on a hot-plate magnetic stirrer and start stirring the reaction solution.
- 20 Stir the reaction solution at room temperature for 48 h under argon atmosphere.
- 21 After the completion of the reaction, remove the balloon and the magnetic stir bar.
- 22 Pour the reaction solution into 30 mL of cold methanol in a 50-mL falcon tube to precipitate the product. Collect the precipitates by centrifugation (2305 × g, 5 min, 4 °C) and wash the precipitate 2 times with 15 mL of cold methanol.  
**▲ CRITICAL STEP** Alternatively, purify the product with ultrafiltration using Amicon Ultra-15 centrifugal filter units (Mw cutoff, 3 kDa). After 5 rounds of washes, collect the filtrates containing pure product and remove the water by lyophilization.
- 23 Collect and transfer the precipitates to 2-mL Eppendorf tubes. Place the tubes under high vacuum for 24 h to remove any residual solvent.
- 24 Weigh the pure DSPE-PEG-S2P, calculate the yield and store it at −20 °C.

#### ? TROUBLESHOOTING

■ **PAUSE POINT** The pure DSPE-PEG-S2P can be stored at −20 °C for at least 1 year.

#### Preparation of siRNA stock solution ● Timing ~ 1 h

- 25 Order the anti-luciferase RNA (siLuc), Dy647-labelled siLuc (Dy647-siLuc), and anti-CaMKII $\gamma$  RNA (siCamk2g) (Table 1). We recommend ordering 1.0  $\mu$ mol of each siRNA in order to perform the following in vitro and in vivo procedures.
- 26 Upon receiving the siRNA, spin down the lyophilized siRNA product.  
**CRITICAL STEP** Normally, the lyophilized siRNAs are light and have very small volumes, thus they can be easily flushed out from the tubes by airflow when they are near the upper parts of the tubes. So, do not open the tubes before spin-down.



- 27 To make a siRNA stock solution of 100  $\mu\text{M}$ , add 0.75 mL of nuclease-free water (HyPure water) to a tube containing 75 nmol siRNA.
- 28 Place the tube on a vortex for 30 s to dissolve the siRNA and spin down the solution. Repeat this process 3 times.
- ▲ **CRITICAL STEP** It takes time to dissolve the siRNA in water, so, make sure the siRNA is completely dissolved before it is used for further experiment. Otherwise, the concentration of siRNA cannot be precisely determined.
- CRITICAL STEP** Do not use a sonicator to facilitate the dissolution of the siRNA, because it can destroy the siRNA.
- 29 Transfer 1  $\mu\text{L}$  of the obtained siRNA solution to a Nanodrop spectrophotometer to determine the concentration of the siRNA solution.
- 30 Aliquot the siRNA solution into 50- $\mu\text{L}$  fractions and store them at  $-20\text{ }^{\circ}\text{C}$ .
- ▲ **CRITICAL STEP** It's important to aliquot the siRNA solution to minimize repeated freeze-thaw cycles, which can degrade the siRNA. Notably, we have not observed significant degradation of siRNA with 3 freeze-thaw cycles.
- **PAUSE POINT** The siRNA aliquots can be store at  $-20\text{ }^{\circ}\text{C}$  for 1 year.

#### Preparation of siRNA-loaded NPs ● Timing ~ 4 h

**CRITICAL** This procedure should produce about 400  $\mu\text{L}$  of purified siRNA loaded NPs (5 $\mu\text{M}$  siRNA content). This scale is sufficient for most of the downstream characterization and in vitro and in vivo experiments. To perform downstream gel electrophoresis studies (Steps 69-91), this procedure will need to be scaled 7 times of siRNA NPs. For serum stability studies (Steps 102-108), this procedure will need to be scaled 2 times. For in vivo studies (Steps 165-170), this procedure will need to be scaled 5 times. See the Experimental Design section for advice on scaling the reaction.

- 31 Clean the working area and all the tools with RNaseZAP to remove all the RNase contamination.
- ▲ **CRITICAL STEP** This is especially important if the NPs are used to deliver more unstable nucleic acid cargoes, such as mRNA.
- 32 Weigh 10 mg of G0-C14 (from Step 14) and add it to a 10-mL glass vial.
- 33 To make an acetone solution containing 2.5 mg/mL G0-C14, add 4 mL of acetone to a glass vial and place the vial in a sonication bath for 10 min to dissolve the G0-C14.
- ▲ **CRITICAL STEP** The dissolution of G0-C14 in acetone is very slow even under sonication, so be patient to make sure that a clear solution is obtained (this usually takes 10 min). It is very difficult or impossible to dissolve the G0-C14 when its concentration is higher than 2.5 mg/mL. So the maximum G0-C14 concentration is 2.5 mg/mL. In another study of our group, we have shown that acetone can be replaced by dimethylformamide (DMF).<sup>31</sup>

- 34 Weigh 20 mg of PLGA and add it to a 10-mL glass vial.
- 35 To make an acetone solution containing 5 mg/mL PLGA, add 4 mL of acetone to the glass vial and place the vial in a sonication bath for 3 min to dissolve the PLGA.
- 36 Add a magnetic bar to a 20-mL glass vial and rinse the vial with 10 mL of HyPure water, repeat the rinse 3 times.
- ▲ **CRITICAL STEP** This wash step can reduce the RNase contamination of the vials and the magnetic bar.
- 37 Weigh 1 mg of hybrid lipid-PEGs containing different ratios of DSPE-PEG-S2P (from Step 24)/DSPE-PEG and add it to the 20-mL glass vial. Select the composition for your intended application using Table 2. We recommend 1:1 ratio for atherosclerotic lesional macrophages-targeting application. The targeting capability of the siRNA NPs is enabled by the S2P ligand, therefore, a high DSPE-PEG-S2P/DSPE-PEG ratio of 1:1 is recommended to maximize the targeting effect of the siRNA NPs to macrophages. Here, the DSPE-PEG-S2P/DSPE-PEG ratio of 1:1 was used for all in vitro and in vivo experiments.
- 38 To make an aqueous solution containing 0.1 mg/mL lipid-PEGs, add 10 mL of HyPure water to the glass vial and place the vial on a vortex for 30 s to dissolve the lipid-PEGs.
- 39 Place the glass vial on a magnetic stirrer and start stirring the lipid-PEG solution at a stirring speed of 600 rpm.
- 40 Thaw the siRNA aliquots on the ice. One aliquot is required for the preparation of 400  $\mu$ L of siRNA loaded NPs (5 $\mu$ M siRNA content).
- CRITICAL STEP** Thawing siRNA aliquots on ice can minimize the degradation of siRNA. Do not thaw the siRNA aliquots at 37 °C water bath.
- 41 Transfer 20  $\mu$ L of 100  $\mu$ M stock siRNA solution (2 nmol) to a 1.7 mL Eppendorf tube.
- ▲ **CRITICAL STEP** To prepare NPs containing different types of siRNA, just use different siRNA sequences (e.g. siLuc, Dy647-siLuc, siCamk2g) during the preparation of NPs. No further modifications of the procedures are required.
- 42 Add 100  $\mu$ L of 2.5 mg/mL G0-C14 acetone solution (250  $\mu$ g, from Step 33) to the tube to complex the siRNA. Mix by gentle pipetting (10 s) to form siRNA/G0-C14 complexes via the electrostatic interactions.
- ? TROUBLESHOOTING
- 43 Add 500  $\mu$ L of 5 mg/mL PLGA acetone solution (2.5 mg, from Step 35) to the tube. Mix by gentle pipetting (10 s) to achieve a homogeneous siRNA/G0-C14/PLGA solution.

- 44 Immediately add the siRNA/G0-C14/PLGA solution dropwise to the 20-mL vial containing 10 mL of the lipid-PEG solution (from Step 37) under stirring (600 rpm).
- ▲ **CRITICAL STEP** The acetone solution should be added to the aqueous solution dropwise to allow the assembly of the core-shell structure of the polymer-lipid hybrid NPs. Use the gel loading tips to add the solution, because the speed of adding solutions can be more precisely controlled by this type of tips compared to the normal tips. In addition, immerse the tips in the aqueous solution during the addition instead of on the top of it.
- 45 After the completion of the addition, raise the stirring speed to 1000 rpm and stir the solution at room temperature for 1 h.
- ▲ **CRITICAL STEP** The stirring speed is important because vigorous stirring can facilitate the formation of the polymer-lipid hybrid NPs. In another study of our group, we have shown that the stirring time can be shortened to 30 min to minimize the degradation of the nucleic acid cargos at room temperature.<sup>30</sup>
- 46 Rinse the Amicon Ultra-15 centrifugal filter unit (MW cutoff, 100 kDa) with 10 mL of HyPure water, repeat this rinse process 3 times.
- 47 After the completion of the assembly of siRNA NPs, transfer the NP solution (from Step 45) to the Amicon tube.
- 48 Concentrate the NP solutions by centrifugation for 20 min at  $3011 \times g$ , 4 °C.
- 49 Wash the NPs by adding 10 mL of HyPure water to the Amicon tube, redisperse the NPs by pipetting. Concentrate the NP solutions by centrifugation for 20 min at  $3011 \times g$ , 4 °C.
- 50 Repeat the wash step 3 times.
- ▲ **CRITICAL STEP** The centrifugation requires a longer time after the first wash, so increase the time for centrifugation accordingly to make sure that most of the solution is removed by filtration at each wash.
- 51 After the final wash, add 400  $\mu\text{L}$  of PBS (pH 7.4) to redisperse the siRNA NPs. This should result in an siRNA concentration of  $\sim 5 \mu\text{M}$  (RNA encapsulation efficiency is measured in Steps 92-101).
- ▲ **CRITICAL STEP** Control the final concentration of siRNA NP by adding a different volume of PBS according to different needs. For example, 400  $\mu\text{L}$  of siRNA NPs with a high siRNA concentration (20  $\mu\text{M}$  of siRNA) are required for Step 92. To this end, add 100  $\mu\text{L}$  of PBS (pH 7.4) to redisperse the siRNA NPs to obtain 100  $\mu\text{L}$  of siRNA NPs with  $\sim 20 \mu\text{M}$  siRNA, combine 4 batches of this solution to obtain 400  $\mu\text{L}$  of siRNA NPs (20  $\mu\text{M}$  of siRNA content).
- 52 Transfer the siRNA NPs (siRNA concentration, 5  $\mu\text{M}$ , 400  $\mu\text{L}$ ) to a 1.7 mL Eppendorf tube and put it on ice immediately.

■ **PAUSE POINT** The siRNA NPs with siRNA concentration of 5  $\mu\text{M}$  or higher can be used freshly prepared or stored at  $-20\text{ }^{\circ}\text{C}$  for at least 1 month. One freeze-thaw cycle has no significant influence on the effectiveness of the siRNA NPs, however, repeated freeze-thaw cycles should be avoided. If mRNA NPs are prepared, they can be used freshly prepared or stored at  $-80\text{ }^{\circ}\text{C}$  for two weeks.

#### Dynamic light scattering (DLS) characterization of siRNA NPs ● Timing ~ 5 h

- 53 Prepare 400  $\mu\text{L}$  of each siLuc NPs (S2P<sub>0</sub>, and S2P<sub>50</sub> NPs, siLuc 5  $\mu\text{M}$ ) according to the procedures for the preparation of NPs (Steps 31-52 using the formulations specified in Table 2).
- 54 Add 150  $\mu\text{L}$  of S2P<sub>0</sub>, or S2P<sub>50</sub> NPs (siLuc, 5  $\mu\text{M}$ ) to a 1.7 mL Eppendorf tube.
- 55 Add 1.35 mL of PBS (pH 7.4) to the tube to dilute the siLuc NPs.
- 56 Vortex and spin down the solution to achieve 1.5 mL of siLuc NPs (siLuc, 0.5  $\mu\text{M}$ ).  
**▲ CRITICAL STEP** The volume of the siLuc NPs has to be enough for the DLS measurement, otherwise unreliable results could be obtained.
- 57 Transfer the siLuc NPs to a cuvette and measure the size (intensity and distribution) of the siLuc NPs by DLS.  
? TROUBLESHOOTING
- 58 After the size measurement, recover the sample and measure its zeta potential.

#### Transmission electron microscope (TEM) characterization of siRNA NPs ● Timing ~ 6 h

- 59 Prepare 400  $\mu\text{L}$  of each siLuc NPs (S2P<sub>0</sub>, S2P<sub>10</sub>, S2P<sub>30</sub> and S2P<sub>50</sub> NPs; siLuc 5  $\mu\text{M}$ ) according to the procedure for the preparation of NPs (Steps 31-52 using the formulations specified in Table 2).
- 60 Put a filter paper on the bench.
- 61 Mount the carbon-coated copper grids in a tweezer and rest the tweezer on the filter paper.
- 62 Place 20  $\mu\text{L}$  of S2P<sub>0</sub>, S2P<sub>10</sub>, S2P<sub>30</sub>, or S2P<sub>50</sub> NPs on the grid and incubate at room temperature for 10 min.  
**▲ CRITICAL STEP** The volume can be adjusted in the range of 3-20  $\mu\text{L}$  according to the concentration of the NPs.
- 63 After sample incubation, blot the sample with filter paper.
- 64 Wash the sample by placing 3  $\mu\text{L}$  of ultrapure water on the grid.  
**▲ CRITICAL STEP** A second wash can be applied if too much salt is observed in TEM imaging.
- 65 Blot the water with filter paper.

- 66 Immediately place 3  $\mu\text{L}$  of stain solution (UAR-EMS Uranyl acetate replacement stain) on the grid for 10 s and blot the stain solution with filter paper.
- 67 Dry the sample under air for more than 1 h.
- 68 When the sample is completely dry, image the sample using a JEOL 1200EX electron microscope at 80 kV.

? TROUBLESHOOTING

### Gel electrophoresis of siRNA NPs ● Timing ~ 16 h

- 69 Prepare 2.8 mL of siLuc NPs (siLuc S2P<sub>50</sub> NPs; siLuc 5  $\mu\text{M}$ ) according to the procedure for the preparation of NPs (Steps 31-52 using the formulations specified in Table 2. This is a 7-fold larger scale than reported in the above steps).
- 70 Prepare 5 Eppendorf tubes (1.7 mL), add 0.5 mL of fetal bovine serum (FBS) to each tube.
- 71 Add 0.5 mL of siLuc S2P<sub>50</sub> NPs (siLuc 5  $\mu\text{M}$ ) to each of the 5 tubes containing serum. Prepare an untreated extracted siRNA control: extracted siRNA from siRNA S2P<sub>50</sub> NPs without incubation with serum.
- 72 Mix the solution by vortexing and spin down the solution.  
**▲ CRITICAL STEP** Use a solution containing as high as 50% (vol/vol) serum to mimic the *in vivo* environments.
- 73 Similarly, prepare 5 Eppendorf tubes (1.7 mL), add 0.5 mL of fetal bovine serum (FBS) to each tube. Add 0.5 mL of free siLuc PBS solution (siLuc 5  $\mu\text{M}$ ) to each of the 5 tubes containing serum. Mix the solution by vortexing and spin down the solution. Prepare an untreated siRNA control: free siRNA without incubation with serum (incubate with PBS instead).
- 74 Incubate all the 5 siLuc S2P<sub>50</sub> NP samples and 5 free siLuc samples at 37°C in a water bath.
- 75 After incubation for different time periods (0, 3, 6, 9, or 12 h), take out one siLuc S2P<sub>50</sub> NP sample and one siLuc sample at each time point.
- 76 Transfer the NP samples to Amicon Ultra-15 centrifugal filter units (MWCO, 100 kDa), concentrate the samples by centrifugation at 3011  $\times$  g, 4 °C for 15 min.  
**▲ CRITICAL STEP** Adjust the time for centrifugation accordingly if necessary.
- 77 Transfer the concentrated NP samples to 1.7 mL Eppendorf tubes.
- 78 Add 0.5 mL of chloroform to each tube to dissolve the concentrated NP samples.

**▲ CRITICAL STEP** It's critical to keep the volume of water in NP samples as low as possible, so the chloroform can efficiently destroy the core-shell structure of the NPs and release the encapsulated siRNA.

- 79 Extract the siLuc by adding 0.5 M NaCl containing 0.1% (wt/vol) SDS to each tube containing chloroform dissolved NP samples.
- 80 Vortex and spin down the solutions.
- 81 Collect the upper aqueous layers.
- 82 Repeat the extraction process 3 times.
- 83 Combine and transfer all the aqueous layers to 2 mL Eppendorf tubes.
- 84 Transfer the free siLuc samples to 2 mL Eppendorf tubes and dilute them to the same volume as the NP samples.
- 85 Store all the 5 NP samples and 5 free siLuc samples at  $-20^{\circ}\text{C}$ .

**■ PAUSE POINT** The samples can be stored at  $-20^{\circ}\text{C}$  for at least 2 weeks.

- 86 Unpack the ethidium bromide-infused E-Gel 4% (wt/vol) agarose gel and load the gel into a horizontal electrophoresis apparatus.
- 87 Add 20  $\mu\text{L}$  of each of the 10 siRNA samples to 10 Eppendorf tubes (0.6 mL) along with control samples.
- 88 Add 2  $\mu\text{L}$  of loading buffer to each tube, mix them by vortexing and spin down.
- 89 Load 20  $\mu\text{L}$  of each of the 10 siRNA samples to the wells of the agarose gel.
- 90 Run the gel electrophoresis for 30 min using E-Gel Power Snap System (Voltage 48 V).
- 91 After the completion of gel electrophoresis, visualize the agarose gel in an ethidium channel using a Syngene PXi imager.

? TROUBLESHOOTING

#### Determination of siRNA encapsulation efficiency (EE) ● Timing ~ 4 h

- 92 Prepare 400  $\mu\text{L}$  of Dy647-siLuc NPs (S2P<sub>0</sub>, and S2P<sub>50</sub> NPs; Dy647-siLuc 20  $\mu\text{M}$ ) according to the procedure for the preparation of NPs (Steps 31-52 using the formulations specified in Table 2).
 

**▲ CRITICAL STEP** Conduct all the handling with Dy647-siLuc in the dark since the Dy647 is very sensitive to light exposure.
- 93 Add 5  $\mu\text{L}$  of Dy647-siLuc S2P<sub>0</sub> or S2P<sub>50</sub> NPs to a 0.6 mL Eppendorf tube.
- 94 Add 95  $\mu\text{L}$  of DMSO to the tube to dissolve the NPs and release the encapsulated Dy647-siLuc.
- 95 Mix the solution by vortexing and spin down the solution (final feeding Dy647-siLuc, 1  $\mu\text{M}$ ).

**▲ CRITICAL STEP** It is critical to completely dissolve the NPs in DMSO because the incomplete dissolution of the NPs can lead to the decrease of siRNA EE. Repeat the vortexing process if necessary.

- 96 Similarly, prepare a free Dy647-siLuc solution (Dy647-siLuc, 20  $\mu$ M). Add 5  $\mu$ L of free Dy647-siLuc to a 0.6 mL Eppendorf tube containing 95  $\mu$ L of DMSO. Mix the solution by vortexing and spin down the solution (final free Dy647-siLuc, 1  $\mu$ M).
- 97 Prepare all the samples in 6 replicates.
- 98 Transfer 90  $\mu$ L of each sample to a 96-well plate.
- 99 Prepare a Blank control: Transfer 90  $\mu$ L of 95% (vol/vol) DMSO solution to the 96-well plate.
- 100 Measure the fluorescence intensities of all the samples using an Infinite M200 PRO spectrophotometer (ex. 645 nm, em. 680 nm).
- 101 The siRNA EE of NPs can be calculated according to the following equation:

$$\% \text{ EE} = 100 \times (\text{sample} - \text{blank}) / (\text{feeding} - \text{blank}) .$$

#### ? TROUBLESHOOTING

#### Evaluation of the stability of siRNA NPs ● Timing ~ 3 d

- 102 Prepare 800  $\mu$ L of siLuc NPs (S2P<sub>0</sub>, and S2P<sub>50</sub> NPs; siLuc, 5  $\mu$ M) according to the procedure for the preparation of NPs (Steps 31-52 using the formulations specified in Table 2. The procedure will need to be doubled to produce the required amount of NPs for this assay).
- 103 Add 200  $\mu$ L of siLuc S2P<sub>0</sub> or S2P<sub>50</sub> NPs to a 2 mL Eppendorf tube. Prepare all the samples in triplicates.
- 104 Add 1.6 mL of PBS and 200  $\mu$ L of FBS to each of the tubes to dilute the siLuc NPs.
- 105 Mix the solution by vortexing and spin down the solution.
- 106 Incubate the siLuc NP solutions at 37 °C over 3 days. At different time points (0, 12, 24, 36, 48, 60, and 72 h), take out the siLuc NPs and transfer 1.5 mL of them to a cuvette.
- 107 Prepare an untreated siRNA NP control in triplicate: the freshly prepared siLuc NPs without incubation with PBS containing 10% (vol/vol) serum.
- 108 Monitor the sizes of the NPs by DLS, recover the samples, and then continue incubating them at 37 °C.

**▲ CRITICAL STEP** It is normal that the NPs in FBS show an extra 20 nm peak compared to the NPs in PBS, because this peak belongs to FBS.

### Characterization of the endosomal escape of siRNA NPs ● Timing ~ 30 h

- 109** Seed 1 mL of HeLa-Luc cells in Nunc glass-bottom dishes at a density of  $2 \times 10^5$  cells per well (complete RPMI 1640 medium). Incubate the cells at 37 °C in a humidified incubator containing 5% (vol/vol) CO<sub>2</sub> for 24 h.
- 110** Prepare 400 µL of Dy647-siLuc NPs (S2P<sub>50</sub> NPs; Dy647-siLuc 5 µM) according to the procedure for the preparation of NPs (Steps 31-52 using the formulations specified in Table 2)
- 111** Add 10 µL of Dy647-siLuc NPs to a 1.7 mL Eppendorf tube. Prepare the samples in duplicates.
- 112** Add 1 mL of RPMI 1640 medium to each of the tubes to dilute the NPs.
- 113** Mix the solution by pipetting gently.
- 114** Remove the cell medium (from Step 109), add 1 mL of medium containing Dy647-siLuc NPs to each of the two duplicate wells, respectively. The final concentration of Dy647-siLuc in each well is 50 nM.
- 115** Put the cells back in the incubator.
- 116** Remove the NPs after 1 h or 4 h incubation. Wash the cells 3 times with PBS.
- 117** Stain the late endosome of the cells by adding 1 mL of medium containing 50 µM Hoechst 33342 (stock solution 1 mM) and incubate for 30 min.
- ▲ **CRITICAL STEP** The incubation time can be varied from 30 min to 2 h depending on the cell type and density.
- 118** Remove the staining solution and rinse the cells 3 times with PBS.
- 119** Stain the nuclei of the cells by adding 1 mL of medium containing 5 µg/mL Hoechst 33342 (stock solution 10 mg/mL) and incubate for 20 min.
- ▲ **CRITICAL STEP** The nuclei can be easily over-stained, so decrease the incubation time if necessary.
- 120** Remove the staining solution and rinse the cells 3 times with PBS.
- 121** Add 1 mL of fresh medium to the cells, and then observe and image the cells using an FV1000 laser scanning confocal microscope.

? TROUBLESHOOTING

### *In vitro* evaluation of gene silencing by siRNA NPs ● Timing ~ 2 d

- 122** Seed 1 mL of HeLa-Luc cells in a 12-well plate at a density of  $6 \times 10^4$  cells per well (complete RPMI 1640 medium). Incubate the cells at 37 °C in a humidified incubator containing 5% (vol/vol) CO<sub>2</sub> for 24 h (or until 50-70% confluent).
- 123** Prepare 400 µL of siLuc NPs (S2P<sub>50</sub> NPs; Dy647-siLuc 5 µM) according to the procedure for the preparation of NPs (Steps 31-52 using the formulations specified in Table 2.)



- 124** Add 2.5, 5, and 10  $\mu\text{L}$  of Dy647-siLuc NPs to 3 Eppendorf tubes (1.7 mL), respectively. Prepare all the samples in triplicates. Prepare a PBS control (or 0 nM siLuc): for the bioluminescence signal of cells treated with PBS.
- 125** Add 1 mL of RPMI 1640 medium to each of the tubes to dilute the NPs.
- 126** Mix the solution by pipetting gently.
- 127** Remove the cell medium, add 1 mL of the medium containing Dy647-siLuc NPs (or the PBS control) to the cells (in triplicate). Final concentrations of Dy647-siLuc in the wells are 0, 12.5, 25, and 50 nM, respectively.
- ▲ CRITICAL STEP** For optimal gene silencing efficiency, make sure the cells are 50-70% confluent at the time of adding NPs.
- 128** Put the cells back in the incubator.
- 129** Remove the NPs after 24 h incubation. Add 1 mL of medium containing Steady-Glo Luciferase reagent (substrate) to the cells and keep the cells in the incubator for 5 min.
- ▲ CRITICAL STEP** Steady-Glo Luciferase reagents can easily lose their activities through repeated freeze-thaw cycles, divide the reagent into aliquots and always protect from light.
- 130** Determine the luciferase expression of the cells by measuring their bioluminescence signals using an In-Vivo Xtreme imaging system with a charge-coupled device (CCD) camera.
- 131** Determine the gene silencing efficiency by quantifying the average radiance (photons per second per  $\text{cm}^2$  per steradian) within regions of interest using Bruker MI SE software.

? TROUBLESHOOTING

***In vitro* viability assessment of siRNA NPs ● Timing ~ 4 d**

- 132** Seed 100  $\mu\text{L}$  of HeLa-Luc (RPMI 1640 medium), RAW 264.7 (DMEM medium), and HEK-293 (EMEM medium) cells in three 96-well plates at a density of  $5 \times 10^3$  cells per well, respectively. Incubate the cells at 37 °C in a humidified incubator containing 5% (vol/vol)  $\text{CO}_2$  for 24 h. Prepare a Blank control – wells without cells: to measure the absorbance of medium without any cells. Prepare enough wells to run 5 replicates of each sample.
- 133** Prepare 400  $\mu\text{L}$  of siLuc NPs (S2P<sub>50</sub> NPs; siLuc 5  $\mu\text{M}$ ) according to the procedure for the preparation of NPs (Steps 31-52 using the formulations specified in Table 2)
- 134** Add 2.5, 5, 10, and 20  $\mu\text{L}$  of siLuc NPs to 4 Eppendorf tubes (2 mL), respectively. Add 2 mL of RPMI 1640 medium to each of the tubes to dilute the NPs. Prepare a PBS control: to measure the absorbance of cells treated with PBS.

- 135 Mix the solution by pipetting gently.
- 136 Remove the cell medium, add 100  $\mu\text{L}$  of the medium containing siLuc NPs (or PBS control) to the wells of each plate. Prepare 5 replicate wells for each NP sample. The final concentrations of siLuc in the wells are 0, 6.25, 12.5, 25, and 50 nM, respectively.
- 137 Put the cells back in the incubator.
- 138 Remove the NPs after 24 h incubation. Add fresh medium to the cells and incubate for another 48 h.
- ▲ CRITICAL STEP** The evaporation of cell medium in the outer wells of the 96 well plate is significant after 48 h incubation, which can influence the growth of cells. So, we recommend only using the inner wells of the plate.
- 139 Remove the cell medium and add 100  $\mu\text{L}$  of medium containing alamarBlue cell viability and proliferation reagents (add 10  $\mu\text{L}$  alamarBlue reagents to 90  $\mu\text{L}$  of medium) to the cells.
- 140 Put the cells back in the incubator and incubate for 1 h.
- 141 Determine the viability of the cells by measuring the absorbance at 570 nm using a plate reader.
- 142 Calculate the viability of NP-treated cells according to the following equation:

$$\% \text{ Cell viability} = 100 \times (\text{sample} - \text{blank}) / (\text{PBS} - \text{blank}).$$

#### ***In vitro* apoptosis assessment of siRNA NPs ● Timing ~ 4 d**

- 143 Seed 0.5 mL of RAW 264.7 (DMEM medium) cells in a 24-well plate at a density of  $3 \times 10^4$  cells per well. Incubate the cells at 37 °C in a humidified incubator containing 5% (vol/vol) CO<sub>2</sub> for 24 h.
- 144 Prepare 400  $\mu\text{L}$  of siLuc NPs (S2P<sub>0</sub> NPs, S2P<sub>50</sub> NPs; siLuc 5  $\mu\text{M}$ ) according to the procedure for the preparation of NPs (Steps 31-52 using the formulations specified in Table 2)
- 145 Add 20  $\mu\text{L}$  of S2P<sub>0</sub> NPs and S2P<sub>50</sub> NPs to 2 Eppendorf tubes (2 mL), respectively. Add 2 mL of DMEM medium to each of the tubes to dilute the NPs. Prepare a PBS control: to later measure the apoptosis of cells treated with only PBS.
- 146 Mix the solution by pipetting gently.
- 147 Remove the cell medium, add 0.5 mL of medium containing siLuc NPs (or PBS control) to the cells. Prepare 3 replicates for each NP sample. The final concentration of siLuc in each well is 50 nM.
- 148 Put the cells back in the incubator.

- 149 Remove the NPs after 24 h incubation. Add fresh medium to the cells and incubate for another 48 h.
- 150 Remove the cell medium, wash once with PBS, detach cells with trypsin, wash and collect the cells by centrifugation ( $1000 \times g$ ,  $4 \text{ }^{\circ}\text{C}$ , 5 min).
- 151 Double stain the cells with Annexin V (AV)-FITC and propidium iodide (PI) (Annexin V-FITC Apoptosis Detection kit). Prepare a Blank control: Unstained cells. Prepare a AV-FITC control: Cells stained with AV-FITC (no PI). Prepare a PI control: Cells stained with PI (no AV-FITC).

**▲ CRITICAL STEP** Keep the stained cells in dark or cover them with aluminum foil due to the light sensitivity of AV-FITC and PI.

- 152 Determine the degree of apoptosis of the cells by flow cytometry. Analyze AV-FITC staining using the FITC signal detector and PI staining using the phycoerythrin (PE) emission signal detector. For the gating strategy for flow cytometry: first, set a live gate (gate 1) in the FSC (forward scatter)/SSC (side scatter) plot to exclude cell debris. Then, use the cell populations from gate 1 for the downstream AV/PI display. Flow cytometry density plots show AV (X-axis) and PI (Y-axis) staining of cells. The right lower quadrant represents AV positive/PI negative staining indicating early apoptosis. The right upper quadrant represents both high AV and PI staining indicating late apoptosis. The left upper quadrant represents low AV and high PI staining indicating necrosis. The left lower quadrant indicates viable cells.

? TROUBLESHOOTING

### ***In vitro* proliferation assessment of siRNA NPs ● Timing ~ 8 d**

- 153 Seed 2 mL of RAW 264.7 (DMEM medium) cells in a 6-well plate at a density of  $2 \times 10^4$  cells per well. Incubate the cells at  $37 \text{ }^{\circ}\text{C}$  in a humidified incubator containing 5% (vol/vol)  $\text{CO}_2$  for 24 h. Prepare enough cell-seeded wells to run triplicate of each sample. Leave several wells empty as a Blank control: to later measure the absorbance of medium without cells.
- 154 Prepare 400  $\mu\text{L}$  of siLuc NPs (S2P<sub>0</sub> NPs, S2P<sub>50</sub> NPs; siLuc 5  $\mu\text{M}$ ) according to the procedure for the preparation of NPs (Steps 31-52 using the formulations specified in Table 2.)
- 155 Add 80  $\mu\text{L}$  of S2P<sub>0</sub> NPs and S2P<sub>50</sub> NPs to 2 falcon tubes (15 mL), respectively. Add 8 mL of DMEM medium to each of the tubes to dilute the NPs.
- 156 Mix the solution by pipetting gently.
- 157 Remove the cell medium, add 2 mL of medium containing siLuc NPs (or PBS control) to the cells. Prepare 3 replicates for each NP sample. The final concentration of siLuc in each well is 50 nM.
- 158 Put the cells back in the incubator.

- 159** Remove the NPs after 24 h incubation. Add fresh medium to the cells and incubate for another 7 days.
- 160** At each time point (day 1, 2, 3, 4, 5, 6, and 7), remove the cell medium and add 2 mL of medium containing alamarBlue cell viability and proliferation reagents to the cells.
- ▲ CRITICAL STEP** Handle the alamarBlue reagents in dark due to their light sensitivity.
- 161** Put the cells back to the incubator and incubate for 1 h.
- 162** Determine the proliferation of the cells by measuring the absorbance at 570 nm using a plate reader.
- 163** Remove the medium containing alamarBlue reagents after each measurement, and add fresh medium to cells to allow the continued growth of the cells.
- 164** Calculate the proliferation of the cells according to the following equation:

$$\% \text{ Relative cell proliferation} = 100 \times (\text{sample}_{\text{day } n} - \text{sample}_{\text{day } 1}) / (\text{sample}_{\text{day } 1} - \text{blank}), n = 1-7.$$

#### Treatment of WD-fed *Ldlr*<sup>-/-</sup> mice with S2P<sub>50</sub> siCamk2g NPs ● Timing ~ 4 weeks

- 165** Prepare 2 mL of S2P<sub>50</sub> siCamk2g NPs and S2P<sub>50</sub> siCtrl NPs (siRNA 5 μM) according to the procedure for the preparation of NPs (2 mL of NPs is required for single-dose injection for 7~9 mice). The siCtrl NPs is the control NPs prepared using non-therapeutic siLuc. Follow Steps 31-52 using the formulations specified in Table 2. This is 5-fold larger scale than reported in the above steps.
- ? TROUBLESHOOTING
- 166** After establishing the mouse model of atherosclerosis (see Reagent setup section), randomly divide the mice into 3 groups and inject intravenously two doses of PBS, S2P<sub>50</sub> siCtrl NPs (1 nmol siRNA per injection, 200 μL), or S2P<sub>50</sub> siCamk2g NPs (1 nmol siRNA per injection, 200 μL) per week to the mice, respectively. In each week, inject the second dose 3 or 4 days after the injection of the first dose. For statistical analysis, n = 7 to 9 mice per group are required per assay.
- ? TROUBLESHOOTING
- 167** Maintain the mice on the WD during the treatment period.
- 168** After 4 weeks of treatment, euthanize the mice with CO<sub>2</sub>, collect the blood by left ventricular puncture, and perfuse the vasculature with cold PBS.
- 169** Fix the heart with attached aortic roots in 4% (wt/vol) paraformaldehyde followed by paraffin embedding or treat it with optical cutting temperature (OCT) compound and freeze it.

- 170 Isolate sections of the aortic root with sizes of 8- or 6- $\mu$ m from the acetone-fixed/paraffin-embedded or frozen samples for the analysis of atherosclerotic lesions.

### **Silence of CaMKII $\gamma$ in lesional macrophages of WD-fed *Ldlr*<sup>-/-</sup> mice ● Timing ~ 2 d**

- 171 Block the obtained aortic root sections with 1% (wt/vol) bovine serum albumin (BSA) at room temperature for 1 h. For statistical analysis, n = 7 to 9 sections per group are required per assay.
- 172 Incubate the sections with primary anti-CaMKII $\gamma$  (1:400) and anti-Mac2 (1:10,000), or with anti-MerTK (1:500) at 4°C overnight.
- ▲ **CRITICAL STEP** To maintain the optimal activity of the primary antibody, handle the antibodies on ice.
- 173 Rinse the sections with PBS three times.
- 174 Incubate the sections with fluorescently labeled secondary antibodies for 2 h at room temperature.
- 175 After rinsing the sections with PBS three times, incubate the sections with 4',6-diamidino-2-phenylindole phenylindole (DAPI) for 20 min at room temperature.
- ▲ **CRITICAL STEP** Keep the stained sections in the dark or cover them with aluminum foil due to the light sensitivity of the dye-labelled antibodies and DAPI.
- 176 Rinse the sections (n = 7 to 9 sections per group) with PBS three times. Observe and image the sections using a Zeiss epifluorescence microscope, and analyze the images using FIJI software. The PBS treated control is used to investigate gene silencing of sections from mice treated by PBS. The siCtrl NPs treated control is used to investigate gene silencing of sections from mice treated by S2P<sub>50</sub> siCtrl NPs.

? TROUBLESHOOTING

### **Efferocytosis in lesional macrophages of WD-fed *Ldlr*<sup>-/-</sup> mice ● Timing ~ 2 d**

- 177 Incubate the Acetone-fixed aortic root sections (from Step 170) with TUNEL staining reagents at 37°C for 1 h. For statistical analysis, n = 7 to 9 sections per group are required per assay.
- 178 Wash the sections three times with PBS.
- 179 Block the sections with 1% (wt/vol) BSA at room temperature for 1 h.
- 180 Incubate the sections with primary anti-Mac2 (1:10,000) at 4°C overnight.
- 181 Rinse the sections with PBS three times.
- 182 Incubate the sections with fluorescently labeled secondary antibodies for 2 h at room temperature.

- 183** After rinsing the sections with PBS three times, incubate the sections with DAPI for 20 min at room temperature.
- ▲ CRITICAL STEP** Keep the stained sections in the dark or cover them with aluminum foil.
- 184** Rinse the sections with PBS three times. Observe and image the sections using a Zeiss epifluorescence microscope, and analyze the images using FIJI software.
- 185** Quantify the lesional efferocytosis as the ratio of TUNEL<sup>+</sup> cells (red) in contact with Mac2<sup>+</sup> macrophages (green) to isolated TUNEL<sup>+</sup> cells. The PBS treated control is used to investigate efferocytosis of sections from mice treated by PBS. The siCtrl NPs treated control is used to investigate efferocytosis of sections from mice treated by S2P<sub>50</sub> siCtrl NPs.

#### **Plaque stability of WD-fed *Ldlr*<sup>-/-</sup> mice ● Timing ~ 2 d**

- 186** Stain the collagen of aortic root sections (from Step 170) with picosirius red according to the manufacturer's instructions. For statistical analysis, n = 7 to 9 sections per group are required per assay.
- 187** Observe and image the sections with light microscopy.
- 188** Obtain the necrotic core area by measuring the area from the images.
- 189** Quantify the fibrous (collagen) cap thickness from one midpoint and two shoulder regions of the plaque and normalize it to the lesion size (Figure 5h, g, i, j). The PBS treated control is used to investigate plaque stability of sections from mice treated by PBS. The siCtrl NPs treated control is used to investigate plaque stability of sections from mice treated by S2P<sub>50</sub> siCtrl NPs.

#### ***In vivo* evaluation of the toxicity of S2P<sub>50</sub> siCamk2g NPs ● Timing ~ 5 weeks**

- 190** Purchase six-week-old female BALB/c mice from Jackson laboratory.
- 191** Conduct all the procedures under the protocols approved by the Institutional Animal Care and Use Committees at your institution.
- 192** Prepare 800 μL of S2P<sub>50</sub> siCamk2g NPs (siRNA 5 μM) according to the procedure for the preparation of NPs. Follow Steps 31-52 using the formulations specified in Table 2. This is 2-fold larger scale than reported in the above steps.
- 193** Randomly divide the mice into 2 groups and inject intravenously two doses of PBS or S2P<sub>50</sub> siCamk2g NPs (1 nmol siRNA per injection, 200 μL) per week to the mice, respectively. For statistical analysis, n = 3 mice per group.
- 194** After 4 weeks of administration, euthanize the mice with CO<sub>2</sub> three days after the last administration, collect the blood by left ventricular puncture.
- ▲ CRITICAL STEP** After blood collection, immediately turn the blood collection tubes upside down a few times to mix blood with the anticoagulants.

- 195** Collect the major organs and fix them in 4% (wt/vol) paraformaldehyde followed by paraffin embedding.
- 196** Isolate sections from the paraffin-embedded organs, stain them with hematoxylin solution and eosin Y solution (H&E) according to the manufacturer's instructions.
- 197** Observe and image the sections with microscopy.
- 198** Evaluate various hematological parameters of the blood, including alanine aminotransferase (ALT), aspartate aminotransferase (AST), blood urea nitrogen (BUN), hematocrit (HCT), hemoglobin (Hb), lymphocyte count (LY), mean corpuscular volume (MCV), mean corpuscular hemoglobin (MCH), mean corpuscular hemoglobin concentration (MCHC), neutrophil count (NPC), red blood cell count (RBC), and white blood cell count (WBC). PBS treated control is used to investigate H&E staining of sections from mice treated by PBS; hematological parameters of the blood from mice treated by PBS.

## Timing

- Steps 1-14, synthesis of cationic lipid (G0-C14): 3 d
- Steps 15-24, synthesis of ligand-modified lipid-PEG (DSPE-PEG-S2P): 3 d
- Steps 25-30, preparation of siRNA stock solution: 1 h
- Steps 31-52, preparation of siRNA-loaded NPs: 4 h
- Steps 53-58, dynamic light scattering (DLS) characterization of siRNA NPs: 5 h
- Steps 59-68, transmission electron microscope (TEM) characterization of siRNA NPs: 6 h
- Steps 69-91, gel electrophoresis of siRNA NPs: 16 h
- Steps 92-101, determination of the siRNA encapsulation efficiency (EE): 4h
- Steps 102-108, evaluation of the stability of siRNA NPs: 3 d
- Steps 109-121, characterization of the endosomal escape of siRNA NPs: 30 h
- Steps 122-131, *in vitro* evaluation of gene silencing by siRNA NPs: 2 d
- Steps 132-142, *in vitro* viability assessment of siRNA NPs: 4 d
- Steps 143-152, *in vitro* apoptosis assessment of siRNA NPs: 4 d
- Steps 153-164, *in vitro* proliferation assessment of siRNA NPs: 8 d
- Steps 165-170, treatment of WD-fed *Ldlr*<sup>-/-</sup> mice with S2P<sub>50</sub> siCamk2g NPs: 4 weeks

Steps 171-176, the silence of CaMKII $\gamma$  in lesional macrophages of WD-fed *Ldlr*<sup>-/-</sup> mice: 2 d

Steps 177-185, efferocytosis evaluation in lesional macrophages of WD-fed *Ldlr*<sup>-/-</sup> mice: 2 d

Steps 186-189, plaque stability evaluation of WD-fed *Ldlr*<sup>-/-</sup> mice: 2 d

Steps 190-198, *in vivo* evaluation of the toxicity of S2P50 siCamk2g NPs: 5 weeks

## Troubleshooting

Troubleshooting advice can be found in Table 3

## Anticipated results

Following this protocol, researchers will be able to synthesize the NP components, prepare the siRNA-loaded polymer-lipid hybrid NPs and characterize the NPs. We also describe the *in vivo* treatment of atherosclerosis in a mouse model as one of the applications of our polymer-lipid hybrid NPs. The results generated by this protocol can be anticipated based on our previous studies.<sup>26,60</sup> For the synthesis of G0-C14, a light green solid with a yield of about 50% was achieved from 1.0 g of G0 starting material. The conjugation of DSPE-PEG-Mal (200 mg) with S2P (63.7 mg) produced a white powder (DSPE-PEG-S2P) with a yield of about 30%. The structures of G0-C14 and DSPE-PEG-S2P were confirmed by <sup>1</sup>H NMR analysis, which can be found in Supplementary Fig. 1.

The physicochemical properties of the siRNA NPs were characterized by DLS and TEM. The non-targeted siRNA NPs (S2P<sub>0</sub> NPs) had a diameter of 108 ± 2.8 nm (intensity size distribution) and a zeta potential of -7.43 ± 0.94 mV (Fig. 4a, b). As expected, the conjugation of S2P ligand to the NPs resulted in targeted S2P<sub>50</sub> NPs with a slightly increased diameter of 116.2 ± 2.5 nm (intensity size distribution) and zeta potential of -2.66 ± 1.01 mV. TEM images revealed that both the non-targeted (S2P<sub>0</sub> NPs) and targeted (S2P<sub>50</sub> NPs) NPs had spherical morphologies (Fig. 4d). It is worth noting that both S2P<sub>0</sub> and S2P<sub>50</sub> NPs had similar siRNA EE of about 50% (Fig. 4c). The stability of the siRNA NPs in biological media is critical to their long blood circulation, which allows the accumulation of siRNA NPs at lesional sites. We demonstrated that the S2P<sub>0</sub> and S2P<sub>50</sub> NPs were stable in PBS solution containing 10% (vol/vol) serum at 37 °C for at least 3 days as evidenced by the negligible changes in their particle size (Supplementary Fig. 2a, b). We further confirmed that the S2P<sub>50</sub> NPs could efficiently protect the siRNA from enzymatic degradation when incubated in serum solution (Supplementary Fig. 2c). Importantly, our siRNA NPs were able to escape from the endosomes after 4 h incubation with HeLa-Luc cells as demonstrated by the confocal microscopy images (Supplementary Fig. 3).

We then evaluated the gene silencing effect of siRNA NPs in HeLa-Luc cells. Our siLuc-loaded S2P<sub>50</sub> NPs were able to effectively silence the Luc expression of HeLa-Luc cells at a siLuc concentration as low as 12.5 nM as demonstrated by the bioluminescence imaging (Fig. 4e, f). In addition, the *in vitro* biocompatibility of the siRNA NPs was assessed by



different assays. The HeLa-Luc, RAW 264.7 and HEK-293 cells treated with S2P<sub>50</sub> NPs containing varying concentrations of siLuc (6.25-50 nM) displayed high cell viabilities of more than 80% (Supplementary Fig. 4a). The treatment of RAW 264.7 cells with S2P<sub>0</sub> and S2P<sub>50</sub> NPs (siLuc 50 nM) induced minimal cell apoptosis comparable to that induced by PBS (Supplementary Fig. 4b). Both of the S2P<sub>0</sub> and S2P<sub>50</sub> NPs (siLuc 50 nM) did not show any significant inhibitions to the growth of RAW 264.7 cells after monitoring their proliferation as long as 7 days (Supplementary Fig. 4c).

Next, we assessed the *in vivo* applications of our polymer-lipid hybrid NPs using the treatment of atherosclerosis in a mouse model as an example. To this end, we first established a mouse model of atherosclerosis by feeding the *Ldlr* gene-deficient (*Ldlr*<sup>-/-</sup>) mice with a high-fat Western diet (WD) for 8 weeks. The resulting WD-fed *Ldlr*<sup>-/-</sup> mice were intravenously injected with either PBS, S2P<sub>50</sub> siCtrl NPs, or S2P<sub>50</sub> siCamk2g NPs. After the treatment, the sections of the aortic roots of the mice were isolated for analysis of the therapeutic effect of NPs. The immunofluorescence imaging revealed that the S2P<sub>50</sub> siCamk2g NPs could efficiently silence the expression of the CaMKII $\gamma$  in the lesional macrophages, resulting in the activation of the MerTK pathway (Fig. 5a-d). TUNEL assay showed that the S2P<sub>50</sub> siCamk2g NPs-induced activation of the MerTK pathway could significantly promote the lesional efferocytosis (Fig. 5e-f). Consequently, plaque stability was improved in the S2P<sub>50</sub> siCamk2g NP-treated groups as evidenced by the smaller necrotic area and thicker fibrous cap compared to that of the PBS or S2P<sub>50</sub> siCtrl NP-treated groups (Fig. 5g-j).

We also investigated the *in vivo* toxicity of our polymer-lipid hybrid NPs under the same conditions (*e.g.* dose, frequency) as the assay for therapeutic effect evaluation. To this end, BALB/c mice were injected intravenously with PBS or S2P<sub>50</sub> siCamk2g NPs for 4 weeks, two doses per week. After the treatment, sections of the major organs and blood of the mice were isolated and collected for analysis of the toxicity of NPs. No obvious damage to the major organs of NP-treated mice was observed as indicated by the H&E staining images (Fig. 6a). In addition, blood from both PBS- and NP-treated mice showed no significant differences in a wide range of hematological parameters, including alanine aminotransferase (ALT), aspartate aminotransferase (AST), blood urea nitrogen (BUN), hematocrit (HCT), hemoglobin (Hb), lymphocyte count (LY), mean corpuscular volume (MCV), mean corpuscular hemoglobin (MCH), mean corpuscular hemoglobin concentration (MCHC), neutrophil count (NPC), red blood cell count (RBC), and white blood cell count (WBC) (Fig. 6b).

## Supplementary Material

Refer to Web version on PubMed Central for supplementary material.

## Acknowledgements

This study is supported by the American Heart Association (AHA) Collaborative Science Award (No. 2018A004190; to W. Tao), Harvard Medical School/Brigham and Women's Hospital Khoury Innovation Award (No. 2020A003219; to W. Tao), Anesthesiology Department Basic Scientist Grant (No. 2420 BPA075; to W. Tao), Center for Nanomedicine Research Fund No.2019A014810; to W. Tao), National Institutes of Health (HL127464, HL145228, and HL087123; to I. Tabas), and National Institutes of Health (R00 HL145131; to A. Yurdagul). W.

Tao also received a start-up package (for 3 years) from the Department of Anesthesiology, Perioperative and Pain Medicine to establish his independent research laboratory at Harvard Medical School and Brigham and Women's Hospital.

## Data availability

The authors declare that the main data discussed in this protocol are available in the supporting primary research papers (<https://doi.org/10.1126/scitranslmed.aay1063>) and <https://doi.org/10.1073/pnas.1303958110>). The raw datasets are available for research purposes from the corresponding authors upon reasonable request.

## Appendix

### Related links

Key references using this protocol

- Tao W et al. *Sci. Transl. Med* 12, eaay1063 (2020) (DOI:10.1126/scitranslmed.aay1063): <https://stm.sciencemag.org/content/12/553/eaay1063> [PubMed: 32718990]
- Zhu X et al. *Proc. Natl. Acad. Sci* 112, 7779–7784 (2015) (DOI:10.1073/pnas.1505629112): <https://www.pnas.org/content/112/25/7779> [PubMed: 26056316]
- Islam MA et al. *Nat. Biomed. Eng* 2, 850–864 (2018) (DOI:10.1038/s41551-018-0284-0): <https://www.nature.com/articles/s41551-018-0284-0> [PubMed: 31015614]
- Kong N et al. *Sci. Transl. Med* 11, eaaw1565 (2019) (DOI:10.1126/scitranslmed.aaw1565): <https://stm.sciencemag.org/content/11/523/eaaw1565> [PubMed: 31852795]


## Reference

- Schultze JL, Schmieder A & Goerdts S Macrophage activation in human diseases. *Semin. Immunol* 27, 249–256 (2015). [PubMed: 26303100]
- Moore Kathryn J. & Tabas I Macrophages in the Pathogenesis of Atherosclerosis. *Cell* 145, 341–355 (2011). [PubMed: 21529710]
- Ross R Atherosclerosis — An Inflammatory Disease. *N. Engl. J. Med* 340, 115–126 (1999). [PubMed: 9887164]
- Kuznetsova T, Prange KHM, Glass CK & de Winther MPJ Transcriptional and epigenetic regulation of macrophages in atherosclerosis. *Nat. Rev. Cardiol* 17, 216–228 (2020). [PubMed: 31578516]
- Shi J, Kantoff PW, Wooster R & Farokhzad OC Cancer nanomedicine: progress, challenges and opportunities. *Nat. Rev. Cancer* 17, 20–37 (2017). [PubMed: 27834398]
- Zhu X et al. Long-circulating siRNA nanoparticles for validating Prohibitin1-targeted non-small cell lung cancer treatment. *Proc. Natl. Acad. Sci* 112, 7779–7784 (2015). [PubMed: 26056316]
- Davidson BL & McCray PB Current prospects for RNA interference-based therapies. *Nat. Rev. Genet* 12, 329–340 (2011). [PubMed: 21499294]
- Whitehead KA, Langer R & Anderson DG Knocking down barriers: advances in siRNA delivery. *Nat. Rev. Drug Discovery* 8, 129–138 (2009). [PubMed: 19180106]
- Castanotto D & Rossi JJ The promises and pitfalls of RNA-interference-based therapeutics. *Nature* 457, 426–433 (2009). [PubMed: 19158789]
- Chen J et al. Synthesis and use of an amphiphilic dendrimer for siRNA delivery into primary immune cells. *Nat. Protoc* 16, 327–351 (2021). [PubMed: 33277630]
- Loynachan CN et al. Renal clearable catalytic gold nanoclusters for in vivo disease monitoring. *Nat. Nanotechnol* 14, 883–890 (2019). [PubMed: 31477801]
- Terrazas M & Kool ET RNA major groove modifications improve siRNA stability and biological activity. *Nucleic Acids Res.* 37, 346–353 (2008). [PubMed: 19042976]

13. Sioud M, Furset G & Cekaite L Suppression of immunostimulatory siRNA-driven innate immune activation by 2'-modified RNAs. *Biochem. Biophys. Res. Commun* 361, 122–126 (2007). [PubMed: 17658482]
14. Choung S, Kim YJ, Kim S, Park H-O & Choi Y-C Chemical modification of siRNAs to improve serum stability without loss of efficacy. *Biochem. Biophys. Res. Commun* 342, 919–927 (2006). [PubMed: 16598842]
15. Kanasty R, Dorkin JR, Vegas A & Anderson D Delivery materials for siRNA therapeutics. *Nat. Mater* 12, 967–977 (2013). [PubMed: 24150415]
16. Yin H et al. Non-viral vectors for gene-based therapy. *Nat. Rev. Genet* 15, 541–555 (2014). [PubMed: 25022906]
17. Mullard A FDA approves landmark RNAi drug. *Nat. Rev. Drug Discovery* 17, 613–614 (2018).
18. Frank-Kamenetsky M et al. Therapeutic RNAi targeting PCSK9 acutely lowers plasma cholesterol in rodents and LDL cholesterol in nonhuman primates. *Proc. Natl. Acad. Sci* 105, 11915–11920 (2008). [PubMed: 18695239]
19. Zimmermann TS et al. RNAi-mediated gene silencing in non-human primates. *Nature* 441, 111–114 (2006). [PubMed: 16565705]
20. Tadin-Strapps M et al. siRNA-induced liver ApoB knockdown lowers serum LDL-cholesterol in a mouse model with human-like serum lipids. *J. Lipid Res* 52, 1084–1097 (2011). [PubMed: 21398511]
21. Pan H et al. Anti-JNK2 peptide–siRNA nanostructures improve plaque endothelium and reduce thrombotic risk in atherosclerotic mice. *Int. J. Nanomed* 13, 5187 (2018).
22. Flores AM et al. Nanoparticle therapy for vascular diseases. *Arterioscler., Thromb., Vasc. Biol* 39, 635–646 (2019). [PubMed: 30786744]
23. Sager HB et al. RNAi targeting multiple cell adhesion molecules reduces immune cell recruitment and vascular inflammation after myocardial infarction. *Sci. Transl. Med* 8, 342ra380–342ra380 (2016).
24. Dahlman JE et al. In vivo endothelial siRNA delivery using polymeric nanoparticles with low molecular weight. *Nat. Nanotechnol* 9, 648 (2014). [PubMed: 24813696]
25. Leuschner F et al. Therapeutic siRNA silencing in inflammatory monocytes in mice. *Nat. Biotechnol* 29, 1005–1010 (2011). [PubMed: 21983520]
26. Tao W et al. siRNA nanoparticles targeting CaMKII $\gamma$  in lesional macrophages improve atherosclerotic plaque stability in mice. *Sci. Transl. Med* 12, eaay1063 (2020). [PubMed: 32718990]
27. Doran AC et al. CAMKII $\gamma$  suppresses an efferocytosis pathway in macrophages and promotes atherosclerotic plaque necrosis. *J. Clin. Investig* 127, 4075–4089 (2017). [PubMed: 28972541]
28. Doran AC, Yurdagul A & Tabas I Efferocytosis in health and disease. *Nat. Rev. Immunol* 20, 254–267 (2020). [PubMed: 31822793]
29. Finn AV, Nakano M, Narula J, Kolodgie FD & Virmani R Concept of vulnerable/unstable plaque. *Arterioscler., Thromb., Vasc. Biol* 30, 1282–1292 (2010). [PubMed: 20554950]
30. Islam MA et al. Restoration of tumour-growth suppression in vivo via systemic nanoparticle-mediated delivery of PTEN mRNA. *Nat. Biomed. Eng* 2, 850–864 (2018). [PubMed: 31015614]
31. Kong N et al. Synthetic mRNA nanoparticle-mediated restoration of p53 tumor suppressor sensitizes p53-deficient cancers to mTOR inhibition. *Sci. Transl. Med* 11, eaaw1565 (2019). [PubMed: 31852795]
32. Baden LR et al. Efficacy and safety of the mRNA-1273 SARS-CoV-2 vaccine. *N. Engl. J. Med* 384, 403–416 (2021). [PubMed: 33378609]
33. Polack FP et al. Safety and efficacy of the BNT162b2 mRNA Covid-19 vaccine. *N. Engl. J. Med* 383, 2603–2615 (2020). [PubMed: 33301246]
34. Tang Z et al. A materials-science perspective on tackling COVID-19. *Nat. Rev. Mater*, 1–14 (2020).
35. Wang Z et al. mRNA vaccine-elicited antibodies to SARS-CoV-2 and circulating variants. *Nature* 592, 616–622 (2021). [PubMed: 33567448]


36. Bradley T et al. Antibody responses after a single dose of SARS-CoV-2 mRNA vaccine. *N. Engl. J. Med* 384, 1959–1961 (2021). [PubMed: 33755375]
37. Tang Z et al. Insights from nanotechnology in COVID-19 treatment. *Nano today* 36, 101019 (2021). [PubMed: 33178330]
38. Song MS, Salmena L & Pandolfi PP The functions and regulation of the PTEN tumour suppressor. *Nat. Rev. Mol. Cell Biol* 13, 283–296 (2012). [PubMed: 22473468]
39. Kowalski PS, Rudra A, Miao L & Anderson DG Delivering the messenger: Advances in technologies for therapeutic mRNA delivery. *Mol. Ther* 27, 710–728 (2019). [PubMed: 30846391]
40. Ahmadi S et al. Stimulus-responsive sequential release systems for drug and gene delivery. *Nano Today* 34, 100914 (2020). [PubMed: 32788923]
41. Liu Y & Chen X Macrophage-targeted gene therapy to improve the atherosclerotic plaque stability. *Matter* 3, 621–622 (2020).
42. Zhao W, Hou X, Vick OG & Dong Y RNA delivery biomaterials for the treatment of genetic and rare diseases. *Biomaterials* 217, 119291 (2019). [PubMed: 31255978]
43. Wang Y & Yu C Emerging concepts of nanobiotechnology in mRNA delivery. *Angew. Chem. Int. Ed* 59, 23374–23385 (2020).
44. Vaughan HJ, Green JJ & Tzeng SY Cancer-targeting nanoparticles for combinatorial nucleic acid delivery. *Adv. Mater* 32, 1901081 (2020).
45. Zhou Y et al. Blood-brain barrier-penetrating siRNA nanomedicine for Alzheimer's disease therapy. *Sci. Adv* 6, eabc7031 (2020). [PubMed: 33036977]
46. Truong NP et al. An influenza virus-inspired polymer system for the timed release of siRNA. *Nat. Commun* 4, 1902 (2013). [PubMed: 23695696]
47. Hu Q, Li H, Wang L, Gu H & Fan C DNA nanotechnology-enabled drug delivery systems. *Chem. Rev* 119, 6459–6506 (2019). [PubMed: 29465222]
48. Hu Q, Wang S, Wang L, Gu H & Fan C DNA nanostructure-based systems for intelligent delivery of therapeutic oligonucleotides. *Adv. Healthcare Mater* 7, 1701153–1701172 (2018).
49. Lee H et al. Molecularly self-assembled nucleic acid nanoparticles for targeted in vivo siRNA delivery. *Nat. Nanotechnol* 7, 389–393 (2012). [PubMed: 22659608]
50. Thai HBD et al. Kidney-targeted cytosolic delivery of siRNA using a small-sized mirror DNA tetrahedron for enhanced potency. *ACS Cent. Sci* 6, 2250–2258 (2020). [PubMed: 33376785]
51. Bujold KE, Hsu JCC & Sleiman HF Optimized DNA “nanosuitcases” for encapsulation and conditional release of siRNA. *J. Am. Chem. Soc* 138, 14030–14038 (2016). [PubMed: 27700075]
52. Rahman MA et al. Systemic delivery of bc12-targeting siRNA by DNA nanoparticles suppresses cancer cell growth. *Angew. Chem. Int. Ed* 56, 16023–16027 (2017).
53. Semple SC et al. Rational design of cationic lipids for siRNA delivery. *Nat. Biotechnol* 28, 172–176 (2010). [PubMed: 20081866]
54. Kulkarni JA, Witzigmann D, Chen S, Cullis PR & van der Meel R Lipid nanoparticle technology for clinical translation of siRNA therapeutics. *Acc. Chem. Res* 52, 2435–2444 (2019). [PubMed: 31397996]
55. Whitehead KA et al. Degradable lipid nanoparticles with predictable in vivo siRNA delivery activity. *Nat. Commun* 5, 4277 (2014). [PubMed: 24969323]
56. Wittrup A et al. Visualizing lipid-formulated siRNA release from endosomes and target gene knockdown. *Nat. Biotechnol* 33, 870–876 (2015). [PubMed: 26192320]
57. Sahay G et al. Efficiency of siRNA delivery by lipid nanoparticles is limited by endocytic recycling. *Nat. Biotechnol* 31, 653–658 (2013). [PubMed: 23792629]
58. Liu S et al. Membrane-destabilizing ionizable phospholipids for organ-selective mRNA delivery and CRISPR–Cas gene editing. *Nat. Mater* 20, 701–710 (2021). [PubMed: 33542471]
59. Cheng Q et al. Selective organ targeting (SORT) nanoparticles for tissue-specific mRNA delivery and CRISPR–Cas gene editing. *Nat. Nanotechnol* 15, 313–320 (2020). [PubMed: 32251383]
60. Xu X et al. Enhancing tumor cell response to chemotherapy through nanoparticle-mediated codelivery of siRNA and cisplatin prodrug. *Proc. Natl. Acad. Sci* 110, 18638–18643 (2013). [PubMed: 24167294]

61. Wilson DS et al. Orally delivered thioketal nanoparticles loaded with TNF- $\alpha$ -siRNA target inflammation and inhibit gene expression in the intestines. *Nat. Mater* 9, 923–928 (2010). [PubMed: 20935658]
62. Díez S, Miguéliz I & De Iarduya CT Targeted cationic poly (D, L-lactic-co-glycolic acid) nanoparticles for gene delivery to cultured cells. *Cell. Mol. Biol. Lett* 14, 347–362 (2009). [PubMed: 19194666]
63. Shi J, Xiao Z, Votruba AR, Vilos C & Farokhzad OC Differentially charged hollow core/shell lipid-polymer-lipid hybrid nanoparticles for small interfering RNA delivery. *Angew. Chem. Int. Ed* 50, 7027–7031 (2011).
64. Mui BL et al. Influence of polyethylene glycol lipid desorption rates on pharmacokinetics and pharmacodynamics of siRNA lipid nanoparticles. *Mol. Ther.--Nucleic Acids* 2, e139 (2013). [PubMed: 24345865]
65. Love KT et al. Lipid-like materials for low-dose, in vivo gene silencing. *Proc. Natl. Acad. Sci* 107, 1864–1869 (2010). [PubMed: 20080679]
66. Schöttler S et al. Protein adsorption is required for stealth effect of poly(ethylene glycol)- and poly(phosphoester)-coated nanocarriers. *Nat. Nanotechnol* 11, 372–377 (2016). [PubMed: 26878141]
67. Petros RA & DeSimone JM Strategies in the design of nanoparticles for therapeutic applications. *Nat. Rev. Drug Discovery* 9, 615–627 (2010). [PubMed: 20616808]
68. Moghimi SM, Hunter AC & Murray JC Long-circulating and target-specific nanoparticles: theory to practice. *Pharmacol. Rev* 53, 283–318 (2001). [PubMed: 11356986]
69. Gilleron J et al. Image-based analysis of lipid nanoparticle-mediated siRNA delivery, intracellular trafficking and endosomal escape. *Nat. Biotechnol* 31, 638 (2013). [PubMed: 23792630]
70. Ahmed SA, Goyal RM Jr & Walsh JE A new rapid and simple non-radioactive assay to monitor and determine the proliferation of lymphocytes: an alternative to [<sup>3</sup>H] thymidine incorporation assay. *J. Immunol. Methods* 170, 211–224 (1994). [PubMed: 8157999]
71. Miller LR et al. Considering sex as a biological variable in preclinical research. *FASEB J.* 31, 29–34 (2017). [PubMed: 27682203]
72. Fredman G et al. Targeted nanoparticles containing the proresolving peptide Ac2-26 protect against advanced atherosclerosis in hypercholesterolemic mice. *Sci. Transl. Med* 7, 275ra220–275ra220 (2015).
73. Zhu X et al. Surface de-PEGylation controls nanoparticle-mediated sirna delivery in vitro and in vivo. *Theranostics* 7, 1990–2002 (2017). [PubMed: 28638484]
74. Liu J et al. Fast and efficient CRISPR/Cas9 genome editing in vivo enabled by bioreducible lipid and messenger RNA nanoparticles. *Adv. Mater* 31, 1902575 (2019).
75. Ensuring accurate resource identification. *Nat. Protoc* 15, 1879–1880 (2020). [PubMed: 32494032]




**Fig. 1. Schematic of CaMKII $\gamma$ -mediated plaque stability.**

The introduction of siCamk2g into the macrophages silences CaMKII $\gamma$ , resulting in the activation of the MerTK pathway which further enhances efferocytosis. The increased efferocytosis leads to a decrease in plaque necrosis and fibrous cap thinning, which synergistically improve the plaque stability.<sup>26,27</sup>



**Fig. 2. Workflow of the protocol.**


The protocol comprises five main stages, including: **a**, synthesis of components (steps 1-24); **b**, preparation of the targeted siRNA NPs (steps 25-52); **c**, physicochemical characterization of the siRNA NPs (steps 53-108); **d**, *in vitro* characterization of the siRNA NPs (steps 109-164); and **e**, *in vivo* therapeutic effect evaluation of the siCamk2g NPs (steps 165-189).



**Fig. 3. Synthesis of cationic lipid and ligand.**


**a**, G0-C14. Adapted with permission from ref.<sup>60</sup>, National Academy of Sciences; **b**, DSPE-PEG-S2P. Adapted with permission from ref.<sup>26</sup>, American Association for the Advancement of Science.






**Fig. 4. Physicochemical and in vitro characterization of siRNA NPs.**

**a-b**, Size and zeta potential of siRNA NPs (S2P<sub>0</sub> and S2P<sub>50</sub>) were determined by DLS. Polydispersity index (PDI).  $n = 3$ . **c**, siRNA encapsulation efficiency (EE) of S2P<sub>0</sub> and S2P<sub>50</sub> NPs was measured by plate reader. Dy647-siRNA was loaded into the NPs (n.s., not significant).  $n = 6$ . **d**, Transmission electron microscopy (TEM) images of S2P<sub>0</sub> and S2P<sub>50</sub> NPs. **e-f**, *In vitro* gene silencing of siRNA NPs. HeLa (HeLa-Luc) cells were incubated with S2P<sub>50</sub> NPs containing varying concentrations of siLuc for 24 h, and their luciferase-expressing was characterized by bioluminescence imaging. The average radiance of the cells was quantified.  $n = 3$ . Data are presented as mean  $\pm$  SEM. \*\*\*\* $P < 0.0001$  in one-way ANOVA. **a-f**, Adapted with permission from ref.<sup>26</sup>, American Association for the Advancement of Science.



**Fig. 5. In vivo therapeutic effect of siCamk2g NPs.**

Eight-week WD-fed *Ldlr*<sup>-/-</sup> mice were injected intravenously with two doses of PBS, S2P<sub>50</sub> siCtrl NPs (1 nmol siRNA per injection), or S2P<sub>50</sub> siCamk2g NPs (1 nmol siRNA per injection) per week for 4 weeks. *n* = 7 to 9 mice per group. **a-d**, Sections of the aortic root were immunostained for CaMKII $\gamma$ , MerTK, and Mac2. Nuclei were stained with DAPI. Scale bars, 200  $\mu$ m. The mean fluorescence intensity (MFI) of each group was normalized to the PBS group. **e-f**, TUNEL assay was used to characterize the apoptotic cells (ACs) of sections of the aortic root. The lesional efferocytosis was quantified as the ratio of TUNEL<sup>+</sup> cells (red) in contact with Mac2<sup>+</sup> macrophages (green) to isolated TUNEL<sup>+</sup> cells. Arrows indicate isolated or free TUNEL<sup>+</sup> cells, and arrowheads indicate TUNEL<sup>+</sup> cells in contact with Mac2<sup>+</sup> macrophages. Scale bars, 50  $\mu$ m. **g-j**, Sections of the aortic root were stained with picrosirius red to determine the lesion size and necrotic core area. The fibrous (collagen) cap thickness is quantified from one midpoint and two shoulder regions of the plaque and normalized to the lesion size. Scale bars, 200  $\mu$ m. \**P* < 0.05 in one-way ANOVA. Data are presented as mean  $\pm$  SEM. **a-j**, Adapted with permission from ref.<sup>26</sup>, American Association for the Advancement of Science. All the animal procedures were conducted in accordance with the protocols reviewed and approved by IACUC at Harvard Medical School and Columbia University Irving Medical Center.



**Fig. 6. In vivo toxicity assessment of siCamk2g NPs.**

Six-week-old female BALB/c mice were injected intravenously with two doses of PBS or S2P<sub>50</sub> siCamk2g NPs (1 nmol siRNA per injection) per week for 4 weeks.  $n = 3$  mice per group. The mice were euthanized by CO<sub>2</sub> inhalation three days after the last administration. **a**, the major organs were sectioned followed by H&E staining. Scale bars, 100  $\mu$ m. **b**, Blood was collected and various hematological parameters were evaluated, including alanine aminotransferase (ALT), aspartate aminotransferase (AST), blood urea nitrogen (BUN), hematocrit (HCT), hemoglobin (Hb), lymphocyte count (LY), mean corpuscular volume (MCV), mean corpuscular hemoglobin (MCH), mean corpuscular hemoglobin concentration (MCHC), neutrophil count (NPC), red blood cell count (RBC), and white blood cell count (WBC). Data are presented as mean  $\pm$  SEM. The differences between PBS and S2P<sub>50</sub> siCamk2g NPs groups are not statistically significant for all parameters (Student's  $t$  test). **a-b**, Adapted with permission from ref.<sup>26</sup>, American Association for the Advancement of Science. All the animal procedures were conducted in accordance with the protocols reviewed and approved by IACUC at Harvard Medical School and Columbia University Irving Medical Center.

**Table 1 |**

The siRNA sequences used in this protocol.

NP formulations	siRNA	Strands	Sequences
S2P <sub>50</sub> siLuc NPs	siLuc	Sense	5'-CUU ACG CUG AGU ACU UCG AdTdT-3'
		Antisense	5'-UCG AAG UAC UCA GCG UAA GdTdT-3'
S2P <sub>50</sub> siCtrl NPs	siLuc used as a control	Sense	5'-CUU ACG CUG AGU ACU UCG AdTdT-3'
		Antisense	5'-UCG AAG UAC UCA GCG UAA GdTdT-3'
S2P <sub>50</sub> Dy647-siLuc NPs	Dy647-siLuc	Sense	5'-Dy647-CUU ACG CUG AGU ACU UCG AdTdT-3'
		Antisense	5'-UCG AAG UAC UCA GCG UAA GdTdT-3'
S2P <sub>50</sub> siCamk2g NPs	siCamK2g	Sense	5'-AAC GUG GUA CAU AAU GCU ACA-3'
		Antisense	5'-UGU AGC AUU AUG UAC CAC GUU-3'

**Table 2 |**

The composition of different siRNA NPs.

NP formulations	Weight ( $\mu\text{g}$ )			
	siRNA/G0-C14 (ratio)	PLGA	DSPE-PEG	DSPE-PEG-S2P (from Step 24)
S2P <sub>0</sub> NPs	25/250	2500	1000	0
S2P <sub>10</sub> NPs	25/250	2500	900	100
S2P <sub>30</sub> NPs	25/250	2500	700	300
S2P <sub>50</sub> NPs	25/250	2500	500	500

Author Manuscript

Author Manuscript

Author Manuscript

Author Manuscript

**Table 3 |**

Troubleshooting table.

Steps	Problem	Possible reason	Solution
10	No product	G0-C14 is stuck on the silica gel	Increase the percentage of MeOH to increase the polarity of the eluent
14	Low yield of G0-C14 product	The amines of G0 are partially oxidized, resulting in a yellow G0 solution	Always store the G0 solution in a 4 °C fridge and flush the bottle with argon after each use to protect the amines from oxidization.
24	Low yield of DSPE-PEG-S2P product	The peptide S2P is not completely dissolved.	Warm up the solution at 37 °C in a water bath or add more HEPES buffer.
42	Precipitation of siRNA	The vigorous mixing process changes the homogenous state of the siRNA/G0-C14 solution	Mixing the siRNA/G0-C14 solution by gentle pipetting. Do not vortex the solution
57	Poor reproducibility of DLS results	The siLuc NPs are too diluted and below the DLS detection limit or the extreme dilution destabilizes the NPs.	Adjust the dilution factor according to the detection limit of the DLS instrument and do not dilute too much.
68	Poor contrast of TEM images	The volume of staining solution is not large enough or incubation time is too short	Increase the volume of the staining solution or increase the incubation time
	Very few samples in TEM images	The sample is too diluted	Increase the concentration of the sample, or repeat the sample loading procedure several times
	The background is too dark in TEM images	The incubation time for the staining solution is too long	Decrease the incubation time of the staining solution
91	The gel bands are too weak	The siRNA concentration is too low	Before gel electrophoresis, concentrate the siRNA samples using Amicon tubes (MWCO, 10 kDa)
101	Low siRNA EE	Fail to transfer all the Dy647-siLuc NPs from the Amicon tube	Pipette the Dy647-siLuc NPs solution vigorously before transferring; Also, wash the tube one or two times with a small volume of PBS.
121	Low fluorescence signals of Dy647	The concentration of Dy647-siLuc NPs is too low or the cell density is too high	Raise the concentration of Dy647-siLuc NPs or decrease the cell density.
131	Low bioluminescence signal	The bioluminescence signal decays	Prepare and warm up the imaging system before adding the substrate. Upon adding the substrate, finish the measurement within 10 min.
	Low gene silencing efficiency	G0-C14 is degraded	The G0-C14 can be easily degraded after dissolution in the solvent. Aliquot the G0-C14 solution and store aliquots at -20 °C. Do not use the same batch more than three times.
152	All the samples are PI positive	The PI staining process is too long	Shorten the PI staining process; Finish the measurement within 1 h after the PI staining
165	The centrifugation takes too much time	The concentration of siRNA NPs is too high, making the ultrafiltration difficult	Split the samples into more Amicon tubes for washing, combine and transfer all the samples to one Amicon tube to adjust the concentration of NPs for <i>in vivo</i> injection
166	The volume of the NPs is too large for injection	The siRNA concentration is too low	Add less PBS to redisperse the NPs during the preparation of NPs
176	Weak fluorescence signal of the sections	Low activity of the primary antibody or the primary antibody concentration is too low	Handle the primary antibodies on ice to retain their optimal activity; Decrease the dilution factor of the primary antibodies.