



# Accurate, Sensitive, and Rapid Detection of *Pseudomonas aeruginosa* Based on CRISPR/Cas12b with One Fluid-Handling Step

## <sup>©</sup>Xiaotong Qiu,<sup>a</sup> Xueping Liu,<sup>b</sup> Ruixue Wang,<sup>c</sup> Xiao Ma,<sup>d</sup> Lichao Han,<sup>a</sup> Jiang Yao,<sup>b</sup> Zhenjun Li<sup>a</sup>

<sup>a</sup>State Key Laboratory for Infectious Disease Prevention and Control, National Institute for Communicable Disease Control and Prevention, Chinese Center for Disease Control and Prevention, Beijing, China

<sup>b</sup>School of Laboratory Medicine and Life Sciences, Wenzhou Medical University, Wenzhou, China

<sup>c</sup>Shanxi Bethune Hospital, Shanxi Academy of Medical Sciences, Tongji Shanxi Hospital, Third Hospital of Shanxi Medical University, Taiyuan, China

<sup>d</sup>Department of Clinical Laboratory, Beijing Friendship Hospital, Capital Medical University, Beijing, China

Xiaotong Qiu and Xueping Liu contributed equally to this work. The author order of equal contribution was determined by both authors.

**ABSTRACT** *Pseudomonas aeruginosa* is a major bacterial pathogen causing nosocomial infections and accounts for morbidity and mortality among patients with cystic fibrosis. An accurate, sensitive, and rapid method to detect P. aeruginosa is critical for the early control of infection and patient management. In this study, we established a P. aeruginosa clustered regularly interspaced short palindromic repeats testing in one pot (CRISPR-top) assay which detected P. aeruginosa with one fluid-handling step in one tube. The reaction was performed isothermally within 1 h; thus, specific instruments were not required. The optimal reaction conditions of this assay were determined to be a temperature of 55°C; working concentrations of 1  $\mu$ M for the forward inner primer and backward inner primer, 0.5  $\mu$ M for the loop forward primer and loop backward primer, and 0.25  $\mu$ M for the forward outer primer and backward outer primer; as well as a 2  $\mu$ M concentration single-stranded DNA reporter molecules. In terms of specificity, our assay showed 100% inclusivity and exclusivity among 48 strains, including 15 P. aeruginosa clinical isolates and 33 non-P. aeruginosa strains. The limit of detection of our method was 10 copies per reaction mixture. Forty-six human sputum specimens from patients with respiratory symptoms were tested. Using the results of quantitative realtime PCR as the gold standard, our method showed 85.3% (29/34) sensitivity, 100% (12/ 12) specificity, a positive predictive value of 100% (29/29), and a negative predictive value of 70.6% (12/17). In summary, the P. aeruginosa CRISPR-top assay developed in the present study is a high-efficiency alternative tool for the accurate and rapid detection of P. aeruginosa, especially in resource-limited settings.

**IMPORTANCE** This study reports a *P. aeruginosa* CRISPR-top assay which can precisely identify *P. aeruginosa* using nucleic acids from pure cultures or clinical samples in one pot with one fluid-handling step. The *P. aeruginosa* CRISPR-top reaction is suitable for on-site testing, and its diagnostic performance can be compared with that of qPCR.

**KEYWORDS** CRISPR, Cas12b, LAMP, *Pseudomonas aeruginosa*, nucleic acid detection, diagnosis, accurate diagnosis

P seudomonas aeruginosa, a Gram-negative aerobe, is a major bacterial pathogen causing nosocomial infections. In the United States, *P. aeruginosa* is the leading cause of ventilator-associated pneumonia in long-term care hospitals and the second most common cause in intensive care units (1). In China, *P. aeruginosa* is the third most common infection in adults with pneumonia (2). In addition, chronic lung colonization

Editor Silvia T. Cardona, University of Manitoba Copyright © 2023 Qiu et al. This is an openaccess article distributed under the terms of the Creative Commons Attribution 4.0 International license.

Address correspondence to Zhenjun Li, lizhenjun@icdc.cn. The authors declare no conflict of interest. Received 2 September 2022 Accepted 13 December 2022 Published 9 January 2023 by *P. aeruginosa* is the major cause of morbidity and mortality among patients with cystic fibrosis (CF) (3). *P. aeruginosa* also causes urinary tract infection (4), sepsis (5), and skin and subcutaneous tissue infections (6). Recently, the emergence of multidrug-resistant *P. aeruginosa*, especially carbapenem-resistant *P. aeruginosa*, has become a challenge for global health (7). Therefore, accurate and rapid detection of *P. aeruginosa* is critical for the early control of infection and patient management.

To date, the routine detection and identification of *P. aeruginosa* has relied on microbiological culture in most laboratories. However, the conventional culture-based identification methods, for example, biochemical identification and matrix-assisted laser desorption ionization–time of flight mass spectrometry (MALDI-TOF MS), are time-consuming (>2 days) and dependent on professional operation and laboratory apparatus (8). Thus, the traditional detection method is not suitable for rapid diagnosis. Although quantitative real-time PCR (qPCR) technology is not culture dependent, has been proven to be a valuable tool to detect *P. aeruginosa*, and is superior to culture-based methods (9), it still requires expensive and sophisticated thermal cyclers.

Recently, clustered regularly interspaced short palindromic repeats (CRISPR)/ CRISPR-associated (CRISPR/Cas) systems have paved the way for next-generation nucleic acid detection. Several CRISPR-based diagnostic platforms, for example, CRISPR/Cas12a-based DNA endonuclease-targeted CRISPR trans reporter (DETECTR) and CRISPR/Cas13a-based specific high-sensitivity enzymatic reporter unLOCKing (SHERLOCK) have been used to detect multiple pathogens and single-nucleotide polymorphisms, showing single-base mismatch specificity and attomolar sensitivity (10, 11). These Cas effectors possess unique collateral cleavage activities that transiently cleave fluorophore quencher-labeled single-stranded DNA (ssDNA) or RNA reporter molecules once activated, releasing a fluorescent signal. The use of guide RNA (gRNA) makes CRISPR/Cas-based detection extremely specific. However, most CRISPR-based detection platforms require two separate fluid-handling steps, namely, preamplification and CRISPR-based detection steps, resulting in multiple operations with the potential for contamination (12-14). More recently, a one-pot, one-step CRISPR-based nucleic acid detection technology, named CRISPR-top, has been successfully devised and applied to the diagnosis of coronavirus disease 2019 (15). AapCas12b, a thermophilic RNA-guided endonuclease from the type V-B CRISPR/Cas system of Alicyclobacillu acidophilus, works at temperatures from 37°C to 60°C (16), which overlaps with the thermal profile of the Bst 2.0 enzyme in loop-mediated isothermal amplification (LAMP) at 55°C to 60°C (17). The CRISPR-top technology employs AapCas12b and couples LAMP and the Cas12b reaction in one pot to achieve accurate detection of pathogens, thus avoiding the crosscontamination caused by multiple steps.

In this study, we applied the CRISPR-top diagnosis platform to establish an accurate, sensitive, and rapid *P. aeruginosa* detection assay, termed the *P. aeruginosa* CRISPR-top assay (Fig. 1). *oprL*, a *P. aeruginosa*-specific and conserved gene encoding a peptidogly-can-associated lipoprotein (18), was selected as the target gene. The optimal reaction conditions of the *P. aeruginosa* CRISPR-top assay were determined, and its diagnostic performance was evaluated using clinical specimens and compared with that of qPCR.

#### RESULTS

**Confirmation of LAMP primers.** To confirm the validity of the primers (Table 1), LAMP reactions were performed at 60°C for 1 h. A total of four reaction mixtures, including one with 10 ng/ $\mu$ L ATCC 27853 genomic DNA, one with 10 pg/ $\mu$ L ATCC 27853 genomic DNA, one with 10 pg/ $\mu$ L ATCC 27853 genomic DNA, and two with deionized water (DW), were tested. The turbidities of the reaction mixtures containing 10 ng/ $\mu$ L or 10 pg/ $\mu$ L *P. aeruginosa* genomic DNA began to increase significantly at 21 min and 38 min, respectively, whereas the turbidities of the two negative controls did not increase within 1 h (see Fig. S1 in the supplemental material). Thus, the validity of the LAMP primer set was confirmed.

**Optimal reaction conditions.** First, the optimal reaction temperature was determined. Genomic DNA from *P. aeruginosa* ATCC 27853 was used in each positive reaction and DW comprised the negative control. The reactions were conducted from 53°C



FIG 1 Schematic diagram of *P. aeruginosa* CRISPR-top workflow (created using Figdraw).

to 58°C with 1°C increments. The experiments at each temperature were repeated three times. As shown in Fig. 2a, 55°C was the optimum reaction temperature of the *P. aeruginosa* CRISPR-top assay, as it produced the highest fluorescence values and the shortest take-off time.

Second, the optimal primer concentration was determined. Six LAMP primers were diluted to 100  $\mu$ M using DW and premixed as 40  $\mu$ L each of forward inner primer (FIP) and backward inner primer (BIP), 20  $\mu$ L each of loop forward (LF) and loop backward (LB) primers, and 10  $\mu$ L each of forward outer primer (F3) and backward outer primer (B3). Serial volumes of the LAMP primer premixture (0.3 to 2.1  $\mu$ L, in 0.3- $\mu$ L intervals) were used to determine the optimal primer concentration for the *P. aeruginosa* CRISPR-top assay. *P. aeruginosa* ATCC 27853 and DW served as the positive control and the negative control, respectively. As shown in Fig. 2b, the reaction mixtures containing 0.6 to 1.5  $\mu$ L of the primer premixture resulted in effective detection. Considering both the cost and effectiveness, 0.9  $\mu$ L of the primer premixture per reaction was determined as the optimum primer concentration, corresponding to working concentrations of 1  $\mu$ M FIP and BIP, 0.5  $\mu$ M LB and LF, and 0.25  $\mu$ M F3 and B3.

Finally, the optimal ssDNA reporter molecule concentration was determined. The ssDNA reporter molecule was diluted to 100  $\mu$ M, and serial volumes (0.2 to 2.0  $\mu$ L, in 0.3- $\mu$ L

| rimer Type              |                         | Sequence (5′–3′) <sup>a</sup>                  |        |
|-------------------------|-------------------------|--|--------|
| Pa-F3                   | Forward outer primer    | CGCGTAGTGCTGGAAGG                              | 17 nt  |
| Pa-B3                   | Backward outer primer   | GGTTCTGAGCCCAGGACTG                            | 19 nt  |
| Pa-FIP                  | Forward inner primer    | TAGCGCTGAACGGCCTTG <b>TTC</b> GGCACCCGCGAGTACA | 37 nt  |
| Pa-BIP                  | Backward inner primer   | GCAGGGTGTTTCGCCGTCGTCGTGGCCGGTAG               | 32 nt  |
| Pa-LF                   | Loop forward primer     | GCACGACGCTCGCCCAGAG                            | 19 nt  |
| Pa-LB                   | Loop backward primer    | GGTAAAGAGCGTCCGGTCG                            | 19 nt  |
| gRNA                    | Guide RNA               | GUCUAGAGGACAGAAUUUUUCAACGGGUGUGCCA             | 111 nt |
| -                       |                         | AUGGCCACUUUCCAGGUGGCAAAGCCCGUUGAGCUUCU         |        |
|                         |                         | CAAAUCUGAGAAGUGGCACGGCACCCGCGAGUACAAUAU        |        |
| ssDNA reporter molecule | ssDNA reporter molecule | 6-FAM-TTATTATTAT-BHQ1                          | 10 nt  |

<sup>a</sup>The inserted PAM site (TTC) is in bold.



**FIG 2** Optimal reaction conditions of the *P. aeruginosa* CRISPR-top assay. Error bars represent the means  $\pm$  standard error of means (SEM); n = 3 technical replicates. (a) Optimal reaction temperature. The reactions were performed under temperatures ranging from 53°C to 58°C at 1°C intervals. (b) Optimal primer concentration. Serial volumes of LAMP primer premixture (0.3 to 2.1  $\mu$ L at 0.3- $\mu$ L intervals) were used to prepare the reaction mixtures. (c) Optimal ssDNA reporter molecule concentration. Serial volumes of 100  $\mu$ M ssDNA probe (0.2 to 2.0  $\mu$ L at 0.3- $\mu$ L intervals) were used to prepare the reaction mixtures.

intervals) were examined. The positive reactions and the negative control were prepared as described above. As shown in Fig. 2c, 0.5  $\mu$ L of the ssDNA reporter molecule per reaction mixture, i.e., a working concentration of 2  $\mu$ M, showed the optimum detection results.

**Specificity and sensitivity of the** *P. aeruginosa* **CRISPR-top assay.** Forty-eight strains, including 15 *P. aeruginosa* clinical isolates and 33 non-*P. aeruginosa* strains (Table S1), were used to verify the specificity of this method. *P. aeruginosa* ATCC 27853 and DW were used as templates in the positive control and the negative control reactions, respectively. The results revealed that all the *P. aeruginosa* strains showed positive results, and all the non-*P. aeruginosa* strains and DW showed negative results (Fig. 3).

In this study, 10-fold serial dilutions of recombinant plasmids containing the targeted fragment of the *P. aeruginosa oprL* gene (from  $10^8$  to 1 copies/ $\mu$ L) and DW (negative control) were used to determine the analytical sensitivity of the method. Finally, we determined that this method could detect at least 10 copies per reaction mixture (Fig. 4).

**Application of the** *P. aeruginosa* **CRISPR-top assay to clinical samples.** A total of 46 human sputum specimens from patients with respiratory symptoms were tested by *P. aeruginosa* CRISPR-top assay, qPCR, and culture. As a result, 76.1% (35/46) of the samples showed coincident results using all three methods; 13% (6/46) showed positive results by CRISPR-top and qPCR but were negative by culture; and 10.9% (5/46) showed negative results by CRISPR-top and culture but were positive by qPCR. Taking the result of qPCR as the gold standard, the sensitivity and specificity of the *P. aeruginosa* CRISPR-top assay were 85.3% (29/34) and 100% (12/12), respectively; the positive predictive value and the negative predictive value were 100% (29/29) and 70.6% (12/17), respectively (Table 2). More details of these results are shown in Fig. 5.



FIG 3 Specificity of the *P. aeruginosa* CRISPR-top assay. PC, positive control (*P. aeruginosa* ATCC 27853); NC, negative control (DW). 1, Klebsiella pneumoniae; 2, Acinetobacter baumannii; 3, Streptococcus pneumoniae; 4, Moraxella catarrhalis; 5, Staphylococcus aureus; 6, Escherichia coli; 7, Mycobacterium tuberculosis; 8, Klebsiella aerogenes; 9, Klebsiella oxytoca; 10, Serratia marcescens; 11, Acinetobacter pizzerialum; 12, Acinetobacter junii; 13, Streptococcus mitis; 14, Streptococcus agalactiae; 15, Steptococcus progenes; 16, Streptococcus salivarius; 17, Streptococcus epidermidis; 22, Stenotophomonas maltophilia; 23, Nocardia cyriacigeorgica; 24, Nocardia farcinica; 25, Corynebacterium striatum; 26, Corynebacterium simulans; 27, Corynebacterium propinquum; 28, Corynebacterium aurimucosum; 29, Enterococcus faecalis; 30, Aeromonas caviae; 31, Elizabethkingia anophelis; 32, Ralstonia mannitolilytica; 33, Rothia kristinae; 34 to 48, *P. aeruginosa* clinical isolates. Error bars represent means ± SEM (*n* = 3 technical replicates).

### DISCUSSION

Several studies have suggested that early detection and adequate treatment are important to prevent or postpone chronic colonization of *P. aeruginosa* in patients with CF, and molecular detection methods, e.g., qPCR, are superior to culture to detect early colonization (9, 19, 20). However, PCR-based molecular diagnostic techniques require sophisticated instruments, and thus their application is restricted in some resource-limited regions and in the field. CRISPR-based diagnostics were reported as one of seven technologies to watch in 2022 (21). Isothermal amplification coupled with CRISPR-based detection has shown higher sensitivity and specificity than isothermal amplification alone in pathogen detection (22). In this report, a CRISPR-assisted detection method, a *P. aeruginosa* CRISPR-top assay, was developed. The whole assay process can be completed within 1 h at a constant temperature and without sophisticated instruments, making this method suitable for point-of-care tests.

Optimization of reaction conditions is critical for a detection system to achieve good efficiency. The *P. aeruginosa* CRISPR-top assay is performed at a constant temperature; thus, a suitable reaction temperature is key to acquiring more effective detection. In this study, the optimal reaction temperature was selected based on the final fluorescence values and the take-off time. Although a positive result can be observed within 40 min, 1 h is still recommended for our assay to achieve maximum sensitivity. Additionally, appropriate working concentrations of the primers and the ssDNA reporter molecule are important. Too much or too little of the oligonucleotides might decrease the effectiveness of detection (Fig. 2b and c), as reported in a previous study (22).



**FIG 4** Sensitivity of the *P. aeruginosa* CRISPR-top assay. Bars represent means  $\pm$  SEM (*n* = 3 technical replicates). NC, negative control. \*, *P* < 0.05; \*\*, *P* < 0.01; \*\*\*, *P* < 0.001; ns, no significance (two-tailed Student's *t* test).

|                        | qPCR         |              |             |             |                  |                  |
|------------------------|--------------|--------------|-------------|-------------|------------------|------------------|
| CRISPR-top result      | No. positive | No. negative | Sensitivity | Specificity | PPV <sup>a</sup> | NPV <sup>a</sup> |
| Positive               | 29           | 0            |             |             |                  |                  |
| Negative               | 5            | 12           | 85.3%       | 100%        | 100%             | 70.6%            |
| Total ( <i>n</i> = 46) | 34           | 12           |             |             |                  |                  |

TABLE 2 Comparison between performance of the P. aeruginosa CRISPR-top assay and qPCR

<sup>a</sup>PPV, positive predictive value; NPV, negative predictive value.

To validate the feasibility of our method in clinical practice, 46 human sputum specimens from patients with respiratory symptoms were tested. qPCR is the gold standard for conventional nucleic acid detection; therefore, we evaluated the consistency of CRISPR-top and qPCR for diagnosis. Our results indicated that the CRISPR-top technique was a potential alternative to gPCR, although the sensitivity of our method was slightly low, which probably resulted from LAMP amplification and CRISPR-Cas12 collateral cleavage occurring simultaneously. Six positive samples tested by our assay were not validated by culture. There may be three explanations. (i) The reaction mixtures were cross-contaminated. This hypothesis, however, can probably be excluded, as all necessary precautions to avoid contamination were taken, and the six "false-positive" samples were also positive when tested using qPCR. (ii) These "false-positive" samples corresponded to a very low inoculum and/or transient colonization (19). The ultrasensitivity of our method means that infection with a low dose of bacteria or transient colonization would be identified by our assay, but not by culture. (iii) The *P. aeruginosa* strains in the six samples were in a viable but noncultivable (VBNC) state. qPCR has been validated to detect VBNC P. aeruginosa (23); thus, the fact that all six samples were positive by qPCR tends to support this hypothesis. Furthermore, five samples showed positive results by qPCR, whereas they were negative by our CRISPR-based approach and culture. Although we cannot exclude the possibility that the limit of detection (LOD) of qPCR is lower than our method and the five positive qPCR samples corresponded to much lower inoculums, in our previous study (13) qPCR showed a similar LOD to a CRISPR-based assay but generated some false-positive results.

Some isothermal amplification technologies, such as recombinase polymerase amplification (24) and multiple cross-displacement amplification (25), have been applied for rapid and sensitive *P. aeruginosa* detection; however, previous research indicated that a relatively high rate of false-positive results might occur in the isothermal amplification assay, resulting in misdiagnosis (26, 27). With LAMP, the false-positive results are usually caused by nonspecific primer interactions (26) and aerosol contamination (28). Careful primer design is the key to avoiding false positives. Moreover, combined with CRISPR-based detection, the false-positive results can be alleviated by the complementary base pairing of the gRNA and the target, which is an advantage of CRISPR-based diagnostics (29). Although a simple, low-cost, and ultrasensitive DNA probe based on a lateral flow biosensor (LFB) with CRISPR/Cas12a and LAMP has been



**FIG 5** Validation of the *P. aeruginosa* CRISPR-top assay using clinical samples. The horizontal dashed line indicates the cutoff value set based on the mean + 3 (standard deviations). Bars represent means  $\pm$  SEM (n = 3 technical replicates). PC, positive control; NC, negative control.

developed with proven effectiveness in *P. aeruginosa* detection, this method has three operation steps, namely, LAMP amplification, Cas12 reaction, and LFB readout (30). After amplification, the LAMP products are transferred to a Cas reaction system, and then the reaction products are tested using LFB. These multiple manual operations might increase the risk of environmental contamination and cross-contamination. In contrast to the existing CRISPR-based diagnosis platform (30), the *P. aeruginosa* CRISPR-top assay combines LAMP amplification and Cas collateral cleavage in one reaction tube, thereby decreasing the risk of environmental contamination and cross-contamination to the maximum extent, achieving accurate, sensitive, and rapid detection of *P. aeruginosa*. Additionally, reaction mixture preparation should be strictly separated from sampling and reaction, which is another way to avoid contamination.

Although CRISPR-mediated detection technologies have become powerful molecular diagnosis tools, the preamplification step is controversial. Some scholars have studied amplification-free CRISPR/Cas detection systems to simplify the process and decrease the risk of cross-contamination (31, 32). However, integrating preamplification with detection is another valuable way to solve the issue. Although it has relatively low sensitivity and a long reaction time compared with other two-step CRISPR-Cas detection methods (14, 33, 34), the *P. aeruginosa* CRISPR-top assay, with only one fluid-handing step, significantly simplifies the manual operations and saves human resources. According to the principle of this method, a fluorescence reader and a lateral flow biosensor can be used to obtain the results (15). Thus, a temperature-controlled small device and portable strips are sufficient for the *P. aeruginosa* CRISPR-top detection, making it suitable for point-of-care testing. Moreover, *P. aeruginosa* CRISPR-top detection is low cost. Each reaction costs approximately 3.50 USD; therefore, this assay is suitable for resource-limited settings.

The *P. aeruginosa* CRISPR-top assay is simple, specific, and of low cost; however, there are a few aspects that could be further refined in the future. More clinical sample types, especially samples from sterile parts, such as blood, pleural effusion, and urine, should be analyzed to provide more clinical diagnostic significance. Additionally, a nucleic acid extraction-free CRISPR-based detection platform could be developed to further simplify operations.

In conclusion, the developed *P. aeruginosa* CRISPR-top assay is a simple, rapid, and reliable method that can detect *P. aeruginosa* in one reaction tube and at a constant temperature. This method is comparable with qPCR; thus, it could be a valuable alternative tool for the rapid and accurate detection of *P. aeruginosa*, especially for on-site testing and in resource-limited settings.

#### **MATERIALS AND METHODS**

**Reagents and instruments.** Wizard Genomic DNA purification kits (catalog number A1125; Promega, Madison, WI, USA) were used for DNA extraction from pure cultures and clinical specimens. DNA isothermal amplification kits (catalog number HT0600; HuiDeXin, Tianjin, China) were used for the LAMP reaction and the *P. aeruginosa* CRISPR-top assay. AapCas12b enzyme with Cas12b buffer (catalog number HT100008; HuiDeXin) was used for gRNA/Cas12b complex preparation. Premix *Ex Taq* 2× (probe qPCR, catalog number RR390A; TaKaRa, Dalian, China) was used for qPCR detection. A spectrophotometer with fixed wavelengths of 260 and 280 nm (Nanodrop LITE; Thermo Fisher Scientific, Shanghai, China) was used to determine the concentration of double-stranded DNA. A real-time turbidimeter (Loopamp LA-320c; Eiken, Tokyo, Japan) was used to monitor the LAMP products in primer screening. An isothermal metal bath (catalog number 88870004; Thermo Fisher Scientific) was used for gRNA/Cas12b complex preparation. A real-time fluorescence qPCR instrument (QuantStudio 6 Flex; Applied Biosystems, Foster City, CA, USA) was used for qPCR and as the fluorescence reader in the CRISPR-top reaction. MALDI-TOF MS (microflex LRF; Bruker, Bremen, Germany) was used to identify the strains isolated from clinical samples.

**Bacterial strains and DNA extraction.** *P. aeruginosa* reference strain ATCC 27853 was used to establish the *P. aeruginosa* CRISPR-top assay. Fifteen *P. aeruginosa* clinical isolates and 33 non-*P. aeruginosa* strains were used to verify the specificity of the method. The strains used in this study are listed in Table S1. All strains were stored in 20% (wt/vol) glycerol broth at  $-70^{\circ}$ C and were identified by MALDI-TOF MS before DNA extraction. The concentrations of extracted genomic DNA from these strains were greater than 10 ng/µL, as determined using the Nanodrop LITE spectrophotometer.

**Primer and gRNA design.** A LAMP primer set targeting the *P. aeruginosa oprL* gene (GenBank accession number 882991) was designed using Primer Explorer version 5 software (http://primerexplorer.jp/lampv5e/index.html), and a TTC sequence (the standard protospacer adjacent motif [PAM] sequence of



**FIG 6** LAMP primers and gRNA design regions. A part of the *oprL* sequence (bp 277 to 490) is shown. Right-pointing arrows and left-pointing arrows indicate sense and complementary strands, respectively.

AapCas12b) was artificially inserted between the F1c and F2 regions of the forward inner primer (FIP). Oligo Analyzer version 3.1 (Integrated DNA Technologies, Coralville, IA, USA) was used to assess the secondary structure of the primers. The gRNA and ssDNA reporter molecule were designed based on the principle of CRISPR-top (15). The sequences of the LAMP primers, gRNA, and ssDNA reporter molecule are listed in Table 1 and the design sites of the primers and gRNA are shown in Fig. 6. The LAMP primers and the ssDNA reporter molecule (6-fluorescein amidite [6-FAM]/Black Hole Quencher 1 [BHQ1] labeled) were synthesized and purified to high-performance liquid chromatography (HPLC) grade by Sangon Biotech (Shanghai, China). The gRNA was synthesized and purified to RNase-free HPLC grade by GenScript Biotech (Nanjing, China).

**LAMP reaction.** To confirm the feasibility of the LAMP primer set, LAMP amplification was performed at 60°C for 1 h using genomic DNA from *P. aeruginosa* ATCC 27853 (10 ng/ $\mu$ L and 10 pg/ $\mu$ L) and nuclease-free DW. The 25- $\mu$ L LAMP reaction mixture contained 12.5  $\mu$ L of 2× isothermal reaction buffer, 8 units of *Bst* 2.0 DNA polymerase, 1.6  $\mu$ M (each) FIP and BIP, 0.8  $\mu$ M (each) LB and LF, 0.4  $\mu$ M (each) F3 and B3, 1  $\mu$ L of DNA or DW (negative control), and DW up to 25  $\mu$ L. The products of LAMP amplification were monitored using the Loopamp LA-320c real-time turbidimeter.

**Standard P. aeruginosa CRISPR-top assay.** First, the gRNA/Cas12b complex was prepared in advance, as described for a previous study (22). In brief, the mixture consisted of 15 pmol AapCas12b and 2  $\mu$ L of gRNA (10  $\mu$ M) in Cas12b buffer, which was incubated at 37°C for 10 min. The prepared complex was used immediately or within 12 h (stored at 4°C). Then, the *P. aeruginosa* CRISPR-top reaction mixture was prepared. The mixture, at a final volume of 25  $\mu$ L, contained 12.5  $\mu$ L of 2× isothermal reaction buffer, 8 units of *Bst* 2.0 DNA polymerase, 1  $\mu$ M (each) FIP and BIP, 0.5  $\mu$ M (each) LB and LF, 0.25  $\mu$ M (each) F3 and B3, 2  $\mu$ M ssDNA reporter molecule, 3.5  $\mu$ L gRNA/Cas12b complex, 1  $\mu$ L DNA template or DW (negative control), and DW up to 25  $\mu$ L. The reaction was performed at 55°C for 1 h and monitored using a fluorescence reader.

**Quantitative real-time PCR.** A qPCR assay was performed using the primers and probe targeting the *oprL* gene, as described previously (19, 35). The qPCR mixture comprised 12.5  $\mu$ L of 2× Premix *Ex Taq* (probe qPCR), 0.3  $\mu$ M oprL-F primer (5'-AACAGCGGTGCCGTTGAC-3'), 0.3  $\mu$ M oprL-R primer (5'-GTCGGAGCTGTCGTACTCGAA-3'), 0.2  $\mu$ M hydrolysis probe (5'-FAM-TGAGCGACGAAGCC-BHQ1-3'), and 1  $\mu$ L of DNA template or DW (negative control) and was made up to a final volume of 25  $\mu$ L with DW. Cycling was performed on the QuantStudio 6 Flex real-time qPCR instrument with an initial hold at 95°C for 10 min, followed by 40 cycles at 95°C for 15 s and 60°C for 1 min (19).

**Specificity and sensitivity of the** *P. aeruginosa* **CRISPR-top assay.** Genomic DNAs from 15 *P. aeruginosa* clinical isolates and 33 non-*P. aeruginosa* strains (Table S1) were extracted and used to determine the specificity of the assay. Among the 33 non-*P. aeruginosa* strains, pneumonia-related pathogenic bacteria (for example, *Klebsiella pneumoniae, Streptococcus pneumoniae, Staphylococcus aureus,* and *Moraxella catarrhalis*) and other common bacteria (for example, *Escherichia coli, Staphylococcus epidermidis, Enterococcus faecalis,* etc.) were included. *P. aeruginosa* ATCC 27853 and DW served as the positive control and the negative control, respectively.

To determine the sensitivity of the *P. aeruginosa* CRISPR-top assay, the targeted 200-bp fragment of the *oprL* gene was cloned into a pUC57 vector to construct a standard plasmid (33). The copy number of DNA molecules was calculated according to the following equation: DNA copy number (in copies per microliter) =  $[6.02 \times 10^{23} \times \text{DNA} \text{ concentration}$  (in nanograms per microliter)  $\times 10^{-9}$ ]/[plasmid length (in base pairs)  $\times 660$ ].

Ten-fold serial dilutions of standard plasmids (from  $10^8$  to 1 copies/ $\mu$ L) and DW (negative control) were used as the templates to determine the LOD of this assay. The experiments were repeated three times.

**Evaluation of effectiveness of the** *P. aeruginosa* **CRISPR-top assay using clinical samples.** To evaluate the effectiveness of the *P. aeruginosa* CRISPR-top assay toward clinical samples, 46 human sputum specimens from patients with respiratory symptoms were collected from two hospitals in China. Each sputum specimen was inoculated on a Columbia blood agar plate, a chocolate agar plate, and an eosin-methylene blue agar plate, and plates were incubated at 35°C for 24 to 72 h. The isolated colonies from the sputum samples were identified using MALDI-TOF MS. To prepare DNA templates, 1 mL each of sputum sample was first digested using 4% sodium hydroxide solution, and then the total DNA was extracted using the DNA purification kit according to the technical manual. The sputum samples were

stored at  $-70^{\circ}$ C until nucleic acid extraction. The same DNA templates were used for the *P. aeruginosa* CRISPR-top testing and qPCR.

**Ethics statement.** This study was approved by the Research Ethics Committee of National Institute for Communicable Disease Control and Prevention, Chinese Center for Disease Control and Prevention (permit number ICDC-2019015). All experiments were performed according to relevant regulations.

Data availability. The data that support the findings of this study are available from the corresponding author upon reasonable request.

#### SUPPLEMENTAL MATERIAL

Supplemental material is available online only. **SUPPLEMENTAL FILE 1**, PDF file, 0.1 MB.

#### ACKNOWLEDGMENTS

We thank Figdraw for their assistance with creating Fig. 1. This work was supported by grants from the National Key R&D Program of China (grant numbers 2021YFC2301105, 2021YFC2301001) and the Strategic Research and Consulting Project of Chinese Academy of Engineering (grant number 2022-XBXD-18).

We declare that the research was conducted in the absence of any commercial or financial relationships that could be construed as a potential conflict of interest.

#### REFERENCES

- Weiner-Lastinger LM, Abner S, Edwards JR, Kallen AJ, Karlsson M, Magill SS, Pollock D, See I, Soe MM, Walters MS, Dudeck MA. 2020. Antimicrobial-resistant pathogens associated with adult healthcare-associated infections: summary of data reported to the National Healthcare Safety Network, 2015–2017. Infect Control Hosp Epidemiol 41:1–18. https://doi .org/10.1017/ice.2019.296.
- 2. Li ZJ, Zhang HY, Ren LL, Lu QB, Ren X, Zhang CH, Wang YF, Lin SH, Zhang XA, Li J, Zhao SW, Yi ZG, Chen X, Yang ZS, Meng L, Wang XH, Liu YL, Wang X, Cui AL, Lai SJ, Jiang T, Yuan Y, Shi LS, Liu MY, Zhu YL, Zhang AR, Zhang ZJ, Yang Y, Ward MP, Feng LZ, Jing HQ, Huang LY, Xu WB, Chen Y, Wu JG, Yuan ZH, Li MF, Wang Y, Wang LP, Fang LQ, Liu W, Hay SI, Gao GF, Yang WZ, Chinese Centers for Disease Control and Prevention Etiology of Respiratory Infection Surveillance Study Team. 2021. Etiological and epidemiological features of acute respiratory infections in China. Nat Commun 12:5026. https://doi.org/10.1038/s41467-021-25120-6.
- Deschaght P, Van Daele S, De Baets F, Vaneechoutte M. 2011. PCR and the detection of Pseudomonas aeruginosa in respiratory samples of CF patients. A literature review. J Cyst Fibros 10:293–297. https://doi.org/10 .1016/j.jcf.2011.05.004.
- Mittal R, Aggarwal S, Sharma S, Chhibber S, Harjai K. 2009. Urinary tract infections caused by Pseudomonas aeruginosa: a minireview. J Infect Public Health 2:101–111. https://doi.org/10.1016/j.jiph.2009.08.003.
- Fabre V, Amoah J, Cosgrove SE, Tamma PD. 2019. Antibiotic therapy for Pseudomonas aeruginosa bloodstream infections: how long is long enough? Clin Infect Dis 69:2011–2014. https://doi.org/10.1093/cid/ciz223.
- Spernovasilis N, Psichogiou M, Poulakou G. 2021. Skin manifestations of Pseudomonas aeruginosa infections. Curr Opin Infect Dis 34:72–79. https:// doi.org/10.1097/QCO.00000000000717.
- Tenover FC, Nicolau DP, Gill CM. 2022. Carbapenemase-producing Pseudomonas aeruginosa: an emerging challenge. Emerg Microbes Infect 11: 811–814. https://doi.org/10.1080/22221751.2022.2048972.
- Tramper-Stranders GA, van der Ent CK, Wolfs TF. 2005. Detection of Pseudomonas aeruginosa in patients with cystic fibrosis. J Cyst Fibros 4(Suppl 2):37–43. https://doi.org/10.1016/j.jcf.2005.05.009.
- Boutin S, Weitnauer M, Hassel S, Graeber SY, Stahl M, Dittrich AS, Mall MA, Dalpke AH. 2018. One time quantitative PCR detection of Pseudomonas aeruginosa to discriminate intermittent from chronic infection in cystic fibrosis. J Cyst Fibros 17:348–355. https://doi.org/10.1016/j.jcf.2017.12.007.
- Chen JS, Ma E, Harrington LB, Da Costa M, Tian X, Palefsky JM, Doudna JA. 2018. CRISPR-Cas12a target binding unleashes indiscriminate singlestranded DNase activity. Science 360:436–439. https://doi.org/10.1126/ science.aar6245.
- Gootenberg JS, Abudayyeh OO, Lee JW, Essletzbichler P, Dy AJ, Joung J, Verdine V, Donghia N, Daringer NM, Freije CA, Myhrvold C, Bhattacharyya RP, Livny J, Regev A, Koonin EV, Hung DT, Sabeti PC, Collins JJ, Zhang F. 2017. Nucleic acid detection with CRISPR-Cas13a/C2c2. Science 356: 438–442. https://doi.org/10.1126/science.aam9321.

- Chen W, Luo H, Zeng L, Pan Y, Parr JB, Jiang Y, Cunningham CH, Hawley KL, Radolf JD, Ke W, Ou J, Yang J, Yang B, Zheng H. 2022. A suite of PCR-LwCas13a assays for detection and genotyping of Treponema pallidum in clinical samples. Nat Commun 13:4671. https://doi.org/10.1038/s41467 -022-32250-y.
- Qiu X, Xu S, Liu X, Ren H, Han L, Li Z. 2022. CRISPR/Cas12a-based diagnostic platform accurately detects Nocardia farcinica targeting a novel species-specific gene. Front Cell Infect Microbiol 12:884411. https://doi.org/ 10.3389/fcimb.2022.884411.
- Zhu X, Wang X, Li S, Luo W, Zhang X, Wang C, Chen Q, Yu S, Tai J, Wang Y. 2021. Rapid, ultrasensitive, and highly specific diagnosis of COVID-19 by CRISPR-based detection. ACS Sens 6:881–888. https://doi.org/10.1021/ acssensors.0c01984.
- Li S, Huang J, Ren L, Jiang W, Wang M, Zhuang L, Zheng Q, Yang R, Zeng Y, Luu LDW, Wang Y, Tai J. 2021. A one-step, one-pot CRISPR nucleic acid detection platform (CRISPR-top): application for the diagnosis of COVID-19. Talanta 233:122591. https://doi.org/10.1016/j.talanta.2021.122591.
- Yang H, Gao P, Rajashankar KR, Patel DJ. 2016. PAM-dependent target DNA recognition and cleavage by C2c1 CRISPR-Cas endonuclease. Cell 167:1814–1828.e12. https://doi.org/10.1016/j.cell.2016.11.053.
- Notomi T, Okayama H, Masubuchi H, Yonekawa T, Watanabe K, Amino N, Hase T. 2000. Loop-mediated isothermal amplification of DNA. Nucleic Acids Res 28:E63. https://doi.org/10.1093/nar/28.12.e63.
- De Vos D, Lim A, Jr, Pirnay JP, Struelens M, Vandenvelde C, Duinslaeger L, Vanderkelen A, Cornelis P. 1997. Direct detection and identification of Pseudomonas aeruginosa in clinical samples such as skin biopsy specimens and expectorations by multiplex PCR based on two outer membrane lipoprotein genes, oprl and oprL. J Clin Microbiol 35:1295–1299. https://doi.org/10.1128/jcm.35.6.1295-1299.1997.
- 19. Le Gall F, Le Berre R, Rosec S, Hardy J, Gouriou S, Boisrame-Gastrin S, Vallet S, Rault G, Payan C, Hery-Arnaud G. 2013. Proposal of a quantitative PCR-based protocol for an optimal Pseudomonas aeruginosa detection in patients with cystic fibrosis. BMC Microbiol 13:143. https://doi.org/10.1186/1471-2180-13-143.
- Hery-Arnaud G, Nowak E, Caillon J, David V, Dirou A, Revert K, Munck MR, Frachon I, Haloun A, Horeau-Langlard D, Le Bihan J, Danner-Boucher I, Ramel S, Pelletier MP, Rosec S, Gouriou S, Poulhazan E, Payan C, Ferec C, Rault G, Le Gal G, Le Berre R. 2017. Evaluation of quantitative PCR for early diagnosis of Pseudomonas aeruginosa infection in cystic fibrosis: a prospective cohort study. Clin Microbiol Infect 23:203–207. https://doi.org/ 10.1016/j.cmi.2016.11.016.
- 21. Eisenstein M. 2022. Seven technologies to watch in 2022. Nature 601: 658–661. https://doi.org/10.1038/d41586-022-00163-x.
- Qiu X, Liu X, Ma X, Wang R, Chen S, Li F, Li Z. 2022. One-pot isothermal LAMP-CRISPR-based assay for Klebsiella pneumoniae detection. Microbiol Spectr https://doi.org/10.1128/spectrum.01545-22.

- Mangiaterra G, Amiri M, Di Cesare A, Pasquaroli S, Manso E, Cirilli N, Citterio B, Vignaroli C, Biavasco F. 2018. Detection of viable but non-culturable Pseudomonas aeruginosa in cystic fibrosis by qPCR: a validation study. BMC Infect Dis 18:701. https://doi.org/10.1186/s12879-018-3612-9.
- 24. Yang H, Wang Y, Yang Q, Fan H, Wang L, Zhang T, Li Z, Liu G, Zhao P, Wu H, Dong J, Liang W. 2021. A rapid and sensitive detection method for Pseudomonas aeruginosa using visualized recombinase polymerase amplification and lateral flow strip technology. Front Cell Infect Microbiol 11:698929. https://doi.org/10.3389/fcimb.2021.698929.
- Zhao F, Niu L, Nong J, Wang C, Wang J, Liu Y, Gao N, Zhu X, Wu L, Hu S. 2018. Rapid and sensitive detection of Pseudomonas aeruginosa using multiple cross displacement amplification and gold nanoparticle-based lateral flow biosensor visualization. FEMS Microbiol Lett 365 . https://doi .org/10.1093/femsle/fny147.
- Hardinge P, Murray JAH. 2019. Reduced false positives and improved reporting of loop-mediated isothermal amplification using quenched fluorescent primers. Sci Rep 9:7400. https://doi.org/10.1038/s41598-019-43817-z.
- Wu H, Zhao P, Yang X, Li J, Zhang J, Zhang X, Zeng Z, Dong J, Gao S, Lu C. 2020. A recombinase polymerase amplification and lateral flow strip combined method that detects Salmonella enterica serotype Typhimurium with no worry of primer-dependent artifacts. Front Microbiol 11:1015. https://doi.org/10.3389/fmicb.2020.01015.
- Wong YP, Othman S, Lau YL, Radu S, Chee HY. 2018. Loop-mediated isothermal amplification (LAMP): a versatile technique for detection of micro-organisms. J Appl Microbiol 124:626–643. https://doi.org/10.1111/ jam.13647.
- Trung NT, Son LHP, Hien TX, Quyen DT, Bang MH, Song LH. 2022. CRISPR-Cas12a combination to alleviate the false-positive in loop-mediated isothermal amplification-based diagnosis of Neisseria meningitidis. BMC Infect Dis 22:429. https://doi.org/10.1186/s12879-022-07363-w.
- Mukama O, Wu J, Li Z, Liang Q, Yi Z, Lu X, Liu Y, Liu Y, Hussain M, Makafe GG, Liu J, Xu N, Zeng L. 2020. An ultrasensitive and specific point-of-care

CRISPR/Cas12 based lateral flow biosensor for the rapid detection of nucleic acids. Biosens Bioelectron 159:112143. https://doi.org/10.1016/j .bios.2020.112143.

- 31. Liu TY, Knott GJ, Smock DCJ, Desmarais JJ, Son S, Bhuiya A, Jakhanwal S, Prywes N, Agrawal S, Diaz de Leon Derby M, Switz NA, Armstrong M, Harris AR, Charles EJ, Thornton BW, Fozouni P, Shu J, Stephens SI, Kumar GR, Zhao C, Mok A, Iavarone AT, Escajeda AM, McIntosh R, Kim S, Dugan EJ, Pollard KS, Tan MX, Ott M, Fletcher DA, Lareau LF, Hsu PD, Savage DF, Doudna JA, IGI Testing Consortium. 2021. Accelerated RNA detection using tandem CRISPR nucleases. Nat Chem Biol 17:982–988. https://doi .org/10.1038/s41589-021-00842-2.
- 32. Fozouni P, Son S, Diaz de Leon Derby M, Knott GJ, Gray CN, D'Ambrosio MV, Zhao C, Switz NA, Kumar GR, Stephens SI, Boehm D, Tsou CL, Shu J, Bhuiya A, Armstrong M, Harris AR, Chen PY, Osterloh JM, Meyer-Franke A, Joehnk B, Walcott K, Sil A, Langelier C, Pollard KS, Crawford ED, Puschnik AS, Phelps M, Kistler A, DeRisi JL, Doudna JA, Fletcher DA, Ott M. 2021. Amplification-free detection of SARS-CoV-2 with CRISPR-Cas13a and mobile phone microscopy. Cell 184:323–333.e9. https://doi.org/10.1016/j.cell .2020.12.001.
- Deng Z, Hu H, Tang D, Liang J, Su X, Jiang T, Hu X, Ying W, Zhen D, Xiao X, He J. 2022. Ultrasensitive, specific, and rapid detection of Mycoplasma pneumoniae using the ERA/CRISPR-Cas12a dual system. Front Microbiol 13:811768. https://doi.org/10.3389/fmicb.2022.811768.
- Li F, Xiao J, Yang H, Yao Y, Li J, Zheng H, Guo Q, Wang X, Chen Y, Guo Y, Wang Y, Shen C. 2022. Development of a rapid and efficient RPA-CRISPR/ Cas12a assay for Mycoplasma pneumoniae detection. Front Microbiol 13: 858806. https://doi.org/10.3389/fmicb.2022.858806.
- 35. Joly B, Pierre M, Auvin S, Colin F, Gottrand F, Guery B, Husson MO. 2005. Relative expression of Pseudomonas aeruginosa virulence genes analyzed by a real time RT-PCR method during lung infection in rats. FEMS Microbiol Lett 243:271–278. https://doi.org/10.1016/j.femsle.2004.12.012.