



A structural basis for prion strain diversity

In the format provided by the
authors and unedited

Table of Contents:

Supplementary Table 1. Determination of prion infectivity titre in cell culture

Supplementary Table 2. Cryo-EM data collection, refinement and validation statistics

Supplementary Figure 1. 2D classification of ME7 prion fibrils

Supplementary Figure 2. Quality of the ME7 map and model

Supplementary references

Supplementary Table 1. Determination of prion infectivity titre in cell culture

Dilution	Positive Wells (w/w ₀)			Percentage Positivity		
	RML BH I6200*	ME7 BH I21487	Purified ME7 (1X)†	RML BH I6200*	ME7 BH I21487	Purified ME7 (1X)†
1.10 ⁻⁴	24 / 24	24 / 24	24 / 24	100.0%	100.0%	100.0%
3.10 ⁻⁵	24 / 24	24 / 24	24 / 24	100.0%	100.0%	100.0%
1.10 ⁻⁵	20 / 24	24 / 24	15 / 24	83.3%	100.0%	62.5%
3.10 ⁻⁶	5 / 24	16 / 24	5 / 24	20.8%	66.7%	20.8%
1.10 ⁻⁶	2 / 24	9 / 24	1 / 24	8.3%	37.5%	4.2%
3.10 ⁻⁷	0 / 24	2 / 24	0 / 24	0.0%	8.3%	0.0%
1.10 ⁻⁷	0 / 24	1 / 24	0 / 24	0.0%	4.2%	0.0%
NI	0 / 24	0 / 24	0 / 25	0.0%	0.0%	0.0%

log TCIU ml⁻¹‡	5.8	6.1	5.5
----------------------------------	------------	------------	------------

Measurement of total protein and PrP content:

Total protein (mg ml⁻¹)	10.0	10.0	0.0027
PrP content (µg ml⁻¹)	5.4	11.0	2.7

Calculation of specific infectivity:

log TCIU mg⁻¹ protein	4.8	5.1	8.1
log TCIU mg⁻¹ PrP	8.1	8.1	8.1

(w/w₀), the number of prion infected wells (w) versus the total number of inoculated wells (w₀); BH, brain homogenate; TCIU, tissue culture infectious units; NI, non-prion inoculated


* Prion titre of 10^{7.3±0.5} (mean + s.d., n=6) intracerebral LD₅₀ units/ml in Tg20 mice¹

† Purified material was resuspended to the same volume as the 10% (w/v) BH from which it was derived


‡ TCIU calculated according to the Poisson distribution²

Supplementary Table 2. Cryo-EM data collection, refinement and validation statistics

	ME7 (EMD-15043) (PDB 8A00)
Data collection and processing	
Magnification	105,000 x
Voltage (kV)	300
Electron exposure (e ⁻ /Å ²)	49.4
Defocus range (μm)	from -2.4 to -0.9
Pixel size (Å)	0.828
Symmetry imposed	C1
Initial particle images (no.)	223,065
Final particle images (no.)	40,239
Map resolution (Å)	2.6
FSC threshold	0.143
Map resolution range (Å)	2.6-3
Refinement	
Initial model used (PDB code)	7QIG
Model resolution (Å)	2.6
FSC threshold	0.143
Model resolution range (Å)	121.7 - 2.6
Map sharpening <i>B</i> factor (Å ²)	-26.75
Model composition	
Non-hydrogen atoms	6456
Protein residues	408
Ligands	none
<i>B</i> factors (Å²)	
Protein	24.18-67.46
r.m.s. deviations	
Bond lengths (Å)	0.002
Bond angles (°)	0.580
Validation	
MolProbity score	1.49
Clashscore	2.79
Poor rotamers (%)	0
Ramachandran plot	
Favored (%)	93.53
Allowed (%)	6.47
Disallowed (%)	0



Supplementary Figure 1. 2D classification of ME7 prion fibrils. Gallery of representative 2D classes (box size: 397.5 x 397.5 Å).



Supplementary Figure 2. Quality of the ME7 map and model. **a** Top, Fourier Shell Correlation (FSC) plots for independent half-reconstructions, determined in Relion 4.0-beta. The final plot (black line) is corrected for overfitting with high-resolution noise substitution. Bottom, FSC plot for the model:map fit. **b** Local resolution of the protein density calculated with Relion 4.0-beta LocRes. **c** Atomic B-factor values colour-coded on the solvent-excluded model surface.

Supplementary references

- 1 Wenborn, A. et al. A novel and rapid method for obtaining high titre intact prion strains from mammalian brain. *Sci Rep* **5**, 10062 (2015).
- 2 Edgeworth, J. A., Jackson, G. S., Clarke, A. R., Weissmann, C., & Collinge, J. Highly sensitive, quantitative cell-based assay for prions adsorbed to solid surfaces. *Proc Natl Acad Sci USA* **106**, 3479-3483 (2009).