

Table 1 Antibody dilutions for western blot analysis

Antibody name	Dilution
TRAF6	1/ 5000
MEK1/2	1/ 3000
p-p65, p-I κ B α , I κ B α , p-MEK1/2	1/ 2500
β -actin, MKK7, p-ASK1, p-ERK1, Bad, Bax	1/ 2000
p-MKK7, JNK1, IKK β , IL-10	1/ 1500
DUSP9, p65, p-IKK β , IKK, p-IKK, IL-1 β , p-TAK1, TAK1, ASK1, p-p38, ERK1, Bcl2, TNF- α , c-Caspase3	1/ 1000
p-JNK1, p38	1/ 500

TRAF6, p-MKK7, p-JNK1, p-MEK1/2, IL-10, IL-1 β , c-Caspase-7 were from Abcam, Cambridge, UK; β -actin, p-IKK, IKK, IKK β , p-I κ B α , I κ B α , p-p65, Bad, p65, MKK7, TAK1, Bcl2, p-TAK1, TNF- α , MEK1/2, c-Caspase-3 were from Beyotime, Shanghai, China; DUSP9, p-IKK β (Invitrogen, Carlsbad, USA), p-ERK1 were from Merck, Shanghai, China; and JNK1, p-p38, p38, ERK1, ASK1, p-ASK1, p-MKK4, and MKK4 were from Zenbio, Chengdu, China.