

Supplementary Table S3. Vectra Polaris marker combinations used to identify immune cell types (supporting methods and Fig. 5).

Cell type	Vectra Polaris marker combination
Total T cells	CD45+/CD3+
CD8 T cells	CD45+/CD3+/CD8+
CD4 T cells	CD45+/CD3+/CD4+
T regulatory cells	CD45+/CD4+/FoxP3+
Macrophages	CD45+/F480+
B cells	CD45+/B220+
Total Leukocytes	CD45+
Total non-leukocytes	CD45-