

## SUPPLEMENTARY TABLE

**Supplementary Table 1. The primer sequences used in RT-PCR.**

---

miR-670-3p	Forward, 5'-CTGATCGTGAGGAGAGTGT-3'; Reverse, 5'-GGTCTTCGACATCGGGGCGG-3'
NLRP3	Forward, 5'-CAACCTCACGTCACACTGCT-3'; Reverse, 5'-TTTCAGACAACCCCAGGTTC-3'
ASC	Forward, 5'-CTGACGGATGAGCAGTACCA-3'; Reverse, 5'-CAGGATGATTTGGTGGGATT-3'
Caspase-1	Forward, 5'-TGTAATGAAGACTGCTACCTGGC-3'; Reverse, 5'-CCTTGTTTCTGAGGATGAAGGATGT-3'
GSDMD	Forward, 5'-ATGGATGGGCAGATACAGGG-3'; Reverse, 5'-TGCTGCAGGACTTTGTGTTC-3'
NOX2	Forward, 5'-ACAAGGTTTATGACGATGAGCC-3'; Reverse, 5'-TTGAGCAACACGCACTGGAA-3'
U6	Forward, 5'-CGCTTCGGCAGCACATATAC-3'; Reverse, 5'-TTCACGAATTTGCGTGTCAT-3'
GAPDH	Forward, 5'-GGAGCGAGATCCCTCCAAAAT-3'; Reverse, 5'-GGCTGTTGTCATACTTCTCATGG-3'

---