

Supplementary Table 1. Antibodies and staining conditions used in this study

Antibody	Animal	Clonality	Clone	Dilution	Antigen retrieval	Vendor
CD68	Mouse	Monoclonal	KP-1	1:500	CB (pH 6); heating in an autoclave for 20 min	DAKO, Santa Clara, CA, USA
CD163	Mouse	Monoclonal	10D6	1:100	CB (pH 6); heating in an autoclave for 20 min	Novocastra, Leica Biosystems, Nussloch, Germany
CD3	Rabbit	Polyclonal		1:500	CB (pH 6); heating in an autoclave for 20 min	DAKO, Santa Clara, CA, USA
CD20	Mouse	Monoclonal	L26	1:1000	CB (pH 6); heating in an autoclave for 20 min	DAKO, Santa Clara, CA, USA