

Supplementary Table 2 - Antibody reagents for FACS sorting.

| <u>Antibody (Clone)</u> | <u>Fluorochrome</u> | <u>Isotype</u> | <u>Source</u> | <u>Dilution</u> |
|-------------------------|---------------------|-----------------------|---------------|-----------------|
| Anti-human CD19 (HIB19) | PE | Mouse IgG1, κ | BioLegend | 1:50 |
| Anti-human CD3 (HIT3a) | FITC | Mouse IgG2a, κ | BioLegend | 1:50 |
| Anti-human CD14 (63D3) | APC | Mouse IgG1, κ | BioLegend | 1:50 |
| DAPI | | | Sigma-Aldrich | 1:10000 |