

S1 Table: Antibodies used for Flow Cytometry Analysis of Contaminating Cell Populations

Staining step	Fluorophore	Specificity	Supplier	Clone	Catalogue Number
Surface Staining	BUV496	CD16	BD	3G8	612945
	V450	CD66b	BD	G10F5	561649
Staining After Perm/Wash Incubation	BV650	CD15	BioLegend	W6D3	323034
	Spark Blue 550	CD45	BioLegend	2D1	368549
	PE-eFluor-610	CD14	Thermo Fisher (Invitrogen)	61D3	61-0149-42
	PE-Cy5	CD3	BD	UCHT1	561007