

Supplementary Table S1 Flow cytometry antibodies.

Antigen	Fluorophore	Clone	Company	Dilution
CD3	BV785	UCHL1	BL	1:100
	APC-Fire750	SK7	BD	1:100
CD4	BUV395	SK3	BD	1:100
CD8	PE-Cy5	HIT8a	BL	1:100
CD11c	PE-Cy7	B-ly6	BD	1:100
CD14	V450	MomegaP9	BD	1:100
CD16	AF700	3G8	BL	1:100
CD18	PE-Cy7	TS1/18	BL	1:100
CD19	BV605	HIB19	BL	1:100
	APC-Fire750	SJ25C1	BD	1:100
CD39	BV510	AI	BL	1:100
CD45	PerCP-Cy5.5	HI30	BL	1:100
CD56	BV650	5.1H11	BL	1:100
	PE	B159	BD	1:100
CD64	APC-Fire750	10.1	BL	1:100
CD103	PE-CF594	Ber-ACT8	BL	1:100
CD158a	FITC	HP-MA4	BL	1:100
CD158b	PE	DX27	BL	1:100
HLA-DR	V500	G46-6	BD	1:100
Perforin	AF647	deltaG9	BD	1:100

Antibodies used for flow cytometry and cell sorting. To stain intracellularly with anti-Perforin, the BD Biosciences Cytofix/Cytoperm kit was used following manufacturer instructions exactly. BD, BD Biosciences; BL, Biolegend.