

**Supplemental Table 3.** Modified carbapenem inactivation method using cefiderocol disks.<sup>a</sup>

<b>Isolate</b>	<b>TSB alone</b>	<b>TSB + EDTA</b>	<b>IDMH alone</b>	<b>IDMH + EDTA</b>
No isolate added to broth	25	25	25	27
PA01	26	25	22	25
Positive control <sup>b</sup>	19	27	20	30
BSL-9	25	26	25	25
BSL-15	23	25	25	26

<sup>a</sup>All disk diameters are reported in millimeters.

<sup>b</sup>*P. aeruginosa* ST773 isolate producing NDM-1, VIM-5, and OXA-21