

Supporting information

Lead optimization of the 5-phenylpyrazolopyrimidinone NPD-2975 towards compounds with improved antitrypanosomal efficacy

Yang Zheng,^{a,1} Magali van den Kerkhof,^{b,1} Mohamed Ibrahim,^a Iwan J.P. de Esch,^a Louis Maes,^b Geert Jan Sterk,^a Guy Caljon^{b,*} and Rob Leurs^{a,*}

^a Amsterdam Institute of Molecular and Life Sciences, Division of Medicinal Chemistry, Faculty of Science, Vrije Universiteit Amsterdam, De Boelelaan 1108, 1081 HZ Amsterdam, the Netherlands.

^b Laboratory of Microbiology, Parasitology and Hygiene (LMPH), University of Antwerp, Universiteitsplein 1, 2610 Wilrijk, Belgium.

¹ Authors contributed equally to this manuscript.

* Corresponding author: guy.caljon@uantwerpen.be (Guy Caljon); r.leurs@vu.nl (Rob Leurs).

Table of contents

Table S1. <i>In vitro</i> metabolism: intrinsic clearance of 11 , 31a-c and 35 in comparison with 2 .	Page S2
Table S2. Measured brain concentrations of 31c in the pharmacokinetic experiments at 24 h after treatment.	Page S3
Figure S1. Detection of parasite infection by SL-RNA detection by qPCR	Page S4
Figure S2-S141. LCMS, ¹ H NMR, ¹³ C NMR and 1D NOESY spectroscopy data.	Page S5-101

Table S1. *In vitro* metabolism: intrinsic clearance of **11**, **31a-c** and **35** in comparison with **2**.

Species	Phase I/II	Time	2 (%)		11 (%)		31a (%)		31b (%)		31c (%)		35 (%)	
			Mean	STEDV	Mean	STEDV	Mean	STEDV	Mean	STEDV	Mean	STEDV	Mean	STEDV
Mouse	I	0	100	-	100	-	100	-	100	-	100	-	100	-
		15	73	3.9	69	3	2.7	1	43	11	108	11.8	113	7.5
		30	55	2.6	54	1.5	0.8	0.02	22	6.4	105	3	110	5.9
		60	33	4.8	31	0.4	0.6	0.2	7	1.8	105	8.8	110	3.6
	II	0	100	-	100	-	100	-	100	-	100	-	100	-
		15	93	6	102	2.4	85	13.4	107	0.7	94	1.2	56	2
		30	81	3.9	108	15	83	15.2	99	11.3	95	8.7	25	3.4
		60	71	7.6	97	8.9	70	4.2	83	0.1	90	11	10	1.5
Human	I	0	100	-	100	-	100	-	100	-	100	-	100	-
		15	100	12	97	2.8	96	4.5	89	4.4	112	10.1	113	2.5
		30	101	5.6	93	2.8	103	13.6	90	3.9	117	13.7	107	4
		60	93	8.5	91	5.5	87	3.5	86	5.8	114	7	109	7
	II	0	100	-	100	-	100	-	100	-	100	-	100	-
		15	74	0.8	100	5.1	81	0.7	99	12.9	100	3.5	82	8.4
		30	75	6.2	98	1.7	79	3.3	118	3	96	8.9	67	8
		60	69	2	99	3.8	82	10.8	95	0.6	99	7.7	46	0.1

Table S2. Measured brain concentrations of **31c** in the pharmacokinetic experiments at 24 h after treatment.

Compound	Administration route	Mouse	Brain levels (ng/mL)	Stdev	Brain levels (μ M)	Stdev
31c	PO	m1	24	23	0.088	0.083
		m2	188	80	0.690	0.294
		m3	30	28	0.111	0.103
	IP	m1	11	4	0.040	0.016
		m2	8	-	0.029	-
		m3	91	101	0.335	0.370







Figure S1. Ct-values of the SL-RNA qPCR of various organs collected from treated animals after infection with *T. brucei*. Results are expressed as the Ct-value for each individual measurement. Results were obtained from three independent repeats. Und./Uns.: Undetected or unspecific.

Acquired by : Admin
 Date Acquired : 20/8/2018 4:45:22 PM
 Sample Name : YAZH01-240-1
 Sample ID :
 Tray# : 1
 Vial# : 17
 Injection Volume : 2
 Data File : C:\LabSolutions\Data\2018\2018-wk34\YAZH01-240-1.lcd
 Background File : blanco 200818.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 21/8/2018 9:01:55 AM



1 - PDA Multi 1 / 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		2.768	833110	99.270
2		3.994	6128	0.730





MS Spectrum Table


#1 Ret.Time:
BG Mode: Calc 2.710-<->3.000(272-<->301)
Mass Peaks: 10 Base Peak: 213.00(333159) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	213.00	333159	100.00				6	258.15	12526	3.76			
2	214.05	35617	10.69				7	276.10	19140	5.75			
3	215.10	3891	1.17				8	277.15	3742	1.12			
4	235.10	28364	8.51				9	447.25	35732	10.73			
5	236.00	3950	1.19				10	448.15	6737	2.02			

Figure S2. LCMS spectrum of compound 4.





Acquired by : Admin
 Date Acquired : 4/6/2019 12:23:28 PM
 Sample Name : YAZH01-246
 Sample ID :
 Tray# : 1
 Vial# : 14
 Injection Volume : 2
 Data File : C:\LabSolutions\Data\2019\2019-wk23\YAZH01-246.lcd
 Background File :
 Method File : Method SCAN BASE standard.lcm
 Report Format : DefaultLCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 4/6/2019 1:37:50 PM



PDA Ch1 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		2.289	4043163	97.221
2		2.779	28032	0.674
3		2.974	87555	2.105





MS Spectrum Table




#1 Ret Time:
 BG Mode: Calc 2.090-<->2.640(210-<->265)
 Mass Peaks: 6 Base Peak: 183.00(630062) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	151.95	17274	2.74				4	184.05	75227	11.94			
2	170.00	40497	6.43				5	204.95	7757	1.23			
3	183.00	630062	100.00				6	256.15	16405	2.60			

Figure S5. LCMS spectrum of compound 5.



Acquired by : Admin
 Date Acquired : 13/3/2019 4:10:07 PM
 Sample Name : YAZH_28_20190313
 Sample ID :
 Tray# : 1
 Vial# : 18
 Injection Volume : 2
 Data File : C:\LabSolutions\Data\2019\2019-wk10\YAZH_28_20190313.lcd
 Background File : blanco_13032019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 13/3/2019 4:23:24 PM



MS Spectrum Table

#1 Ret Time:
BG Mode: Calc 3.920-<->4.290(393-<->430)
Mass Peaks: 8 Base Peak: 287.10(825619) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	287.10	825619	100.00				5	310.00	9548	1.16			
2	288.10	162597	19.69				6	328.10	9457	1.15			
3	289.10	14314	1.73				7	350.10	35884	4.35			
4	309.05	36108	4.37				8	595.30	14294	1.73			

Figure S8. LCMS spectrum of compound 7.




Figure S9. $^1\text{H NMR}$ spectrum of compound 7.






Figure S10. $^{13}\text{C NMR}$ spectrum of compound 7.

Acquired by : Admin
 Date Acquired : 3/9/2018 2:19:49 PM
 Sample Name : YAZH01-248-1
 Sample ID :
 Tray# : 1
 Vial# : 13
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2018\2018-wk36\YAZH01-248-1.lcd
 Background File : blanco 03092018.lcd
 Method File : Method SCAN ACID standard MW50.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 3/9/2018 2:52:36 PM




1 | PDA Multi 1 / 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		4.455	1443459	100.000



#1 Ret.Time:Averaged 4.500-4.520(Scan#:451-453)
Mass Peaks:11 Base Peak:228.05(343753) Polarity:Pos Segment1 - Event1



MS Spectrum Table

#1 Ret.Time:
BG Mode:Calc 4.400-<->4.630(441<->464)
Mass Peaks:11 Base Peak:228.05(343753) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	59.05	9624	2.80				7	229.05	31162	9.07			
2	83.05	95735	27.85				8	269.15	106700	31.04			
3	84.05	3885	1.13				9	270.15	15272	4.44			
4	105.05	5189	1.51				10	288.30	3630	1.06			
5	151.00	3535	1.03				11	338.30	5872	1.71			
6	228.05	343753	100.00										

Figure S11. LCMS spectrum of compound 8.

YAZH01-248-1/7

^1H NMR (600 MHz, Chloroform-*d*) δ 3.98 (s, 3H), 3.96 (s, 3H), 3.43 (hept, $J = 7.1$ Hz, 1H), 1.30 (d, $J = 6.9$ Hz, 6H).




Figure S12. ^1H NMR spectrum of compound 8.

YAZH01-248-1/8

^{13}C NMR (151 MHz, CDCl_3) δ 159.3, 153.3, 132.1, 131.9, 53.7, 39.2, 26.5, 21.5.




Figure S13. ^{13}C NMR spectrum of compound 8.

YAZH01-248-1 1D NOE/2
YAZH01-248-1 1D NOE







Figure S14. 1D NOESY spectrum of compound **8**.

Acquired by : Admin
 Date Acquired : 3/9/2018 2:38:23 PM
 Sample Name : YAZH01-248-2
 Sample ID :
 Tray# : 1
 Vial# : 14
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2018\2018-wk36\YAZH01-248-2.lcd
 Background File : blanco 03092018.lcd
 Method File : Method SCAN ACID standard MW50.lcm
 Report Format : Default.LCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-01072015a.lct
 Processed by : Admin
 Modified Date : 3/9/2018 2:53:35 PM



1: PDA Multi 1 / 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		3.955	1168339	100.000



MS Spectrum Table

#1 Ret.Time:
BG Mode:Calc 3.920<>4.150(393<>416)
Mass Peaks:17 Base Peak:228.05(435416) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	59.05	11139	2.56				10	228.05	435416	100.00			
2	83.05	30697	7.05				11	229.10	43077	9.89			
3	105.05	10233	2.35				12	230.05	7305	1.68			
4	196.00	47038	10.80				13	269.15	17962	4.13			
5	203.00	93195	21.40				14	275.50	4931	1.13			
6	204.20	8468	1.94				15	276.10	10672	2.45			
7	214.10	23006	5.28				16	296.10	15696	3.60			
8	223.50	18131	4.16				17	361.10	4417	1.01			
9	223.60	24965	5.73										

Figure S15. LCMS spectrum of compound 12.

YAZH01-248-2/1

^1H NMR (600 MHz, Chloroform-*d*) δ 3.94 (s, 3H), 3.94 (s, 3H), 3.48 (hept., $J = 7.2$ Hz, 1H), 1.40 (d, $J = 7.2$ Hz, 6H).




Figure S16. ^1H NMR spectrum of compound **12**.

YAZH01-248-2/2

^{13}C NMR (151 MHz, CDCl_3) δ 160.8, 146.3, 137.4, 132.1, 53.1, 39.1, 25.8, 19.4.




Figure S17. ^{13}C NMR spectrum of compound **12**.






Figure S18. 1D NOESY spectrum of compound **12**.

Acquired by : Admin
 Date Acquired : 13/9/2018 4:39:27 PM
 Sample Name : YAZH01-257
 Sample ID :
 Tray# : 1
 Vial# : 27
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2018\2018-wk37\YAZH01-257.lcd
 Background File : blanco 13092018.lcd
 Method File : Method SCAN ACID standard MW100.lcm
 Report Format : Default.LCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 13/9/2018 5:21:56 PM




#1 PDA Multi 1 / 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		1.994	10665	1.228
2		2.137	838126	96.509
3		4.762	19653	2.263



#1 Ret.Time:Averaged 2.190-2.210(Scan#:220-222)

Mass Peaks:7 Base Peak:183.05(315118) Polarity:Pos Segment1 - Event1



MS Spectrum Table

#1 Ret.Time:

BG Mode:Calc 2.100<->2.350(211<->236)

Mass Peaks:7 Base Peak:183.05(315118) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	166.00	92920	29.49				5	221.10	3883	1.23			
2	167.10	9564	3.04				6	224.05	9755	3.10			
3	183.05	315118	100.00				7	385.15	6722	2.13			
4	184.05	44037	13.97										

Figure S19. LCMS spectrum of compound 10.




Figure S20. ¹H NMR spectrum of compound 10.







Figure S21. ¹³C NMR spectrum of compound 10.

Acquired by : Admin
 Date Acquired : 9/10/2018 10:39:03 AM
 Sample Name : YAZH01-262
 Sample ID :
 Tray# : 1
 Vial# : 6
 Injection Volume : 3
 Data File : C:\LabSolutions\Data\2018\2018-wk41\YAZH01-262.lcd
 Background File : blanco 09102018.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 9/10/2018 11:37:02 AM



1 PDA Multi 1 / 254nm 4nm
PDA Ch1 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		4.486	351625	100.000



MS Spectrum Table
#1 Ret Time:
BG Mode Calc 4.440->4.730(445->474)
Mass Peaks: 5 Base Peak: 287.05(398536) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	287.05	398536	100.00				4	328.15	5201	1.31			
2	288.15	60500	15.18				5	350.10	7664	1.92			
3	289.30	5307	1.33										

Figure S22. LCMS spectrum of compound 11.




Figure S23. ¹H NMR spectrum of compound **11**.







Figure S24. ¹³C NMR spectrum of compound **11**.

Acquired by : Admin
 Date Acquired : 13/9/2018 4:48:05 PM
 Sample Name : YAZH01-258
 Sample ID :
 Tray# : 1
 Vial# : 28
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2018\2018-wk37\YAZH01-258.lcd
 Background File : blanco 13092018.lcd
 Method File : Method SCAN ACID standard MW100.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 13/9/2018 5:22:57 PM



#1 PDA Multi 1 / 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		1.779	453050	99.404
2		2.469	2715	0.596



MS Spectrum Table

#1 Ret Time:
BG Mode:Calc 1.740<->2.070(175<->208)
Mass Peaks:10 Base Peak:183.05(279952) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	166.10	35487	12.68				6	197.00	5190	1.85			
2	167.00	3201	1.14				7	205.00	31743	11.34			
3	183.05	279952	100.00				8	206.00	3173	1.13			
4	184.05	65044	23.23				9	246.10	6568	2.35			
5	185.05	5591	2.00				10	387.15	13120	4.69			

Figure S25. LCMS spectrum of compound 14.

YAZH01-258.1.fid




Figure S26. ¹H NMR spectrum of compound 14.

YAZH01-258/2







Figure S27. ¹³C NMR spectrum of compound 14.

Acquired by : Admin
 Date Acquired : 15/10/2018 10:05:15 AM
 Sample Name : YAZH01-263
 Sample ID :
 Tray# : 1
 Vial# : 7
 Injection Volume : 3
 Data File : C:\LabSolutions\Data\2018\2018-wk42\YAZH01-263.lcd
 Background File : blanco_15102018.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 15/10/2018 10:36:12 AM



Peak#	Name	Ret. Time	Area	Area %
1		4.115	1342813	99.184
2		4.580	11050	0.816



MS Spectrum Table
#1 Ret Time: BG Mode:Calc 4.060-->4.400(407-->441)
Mass Peaks:11 Base Peak:287.10(814033) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	287.10	814033	100.00				7	332.15	13825	1.70			
2	288.15	127144	15.62				8	350.10	106657	13.10			
3	289.10	14865	1.83				9	351.10	17400	2.14			
4	309.10	35012	4.30				10	595.25	50713	6.23			
5	328.10	73982	9.09				11	596.25	23883	2.93			
6	329.15	14796	1.82										

Figure S28. LCMS spectrum of compound 15.




Figure S29. ¹H NMR spectrum of compound 15.





Figure S30. ¹³C NMR spectrum of compound 15.



Acquired by : Admin
 Date Acquired : 17/4/2019 8:59:50 AM
 Sample Name : LABI-044
 Sample ID :
 Tray# : 1
 Vial# : 1
 Injection Volume : 5
 Data File : C:\LabSolutions\Data\2019\2019-wk16\LABI-044.lcd
 Background File : blanco 17042019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultLCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 19/6/2019 3:30:24 PM



PDA Ch1 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		3.400	1205913	100.000

Peak Table



MS Spectrum Table

#1 Ret.Time:
 BG Mode: Calc 3.360<->3.670(337<->368)
 Mass Peaks: 5 Base Peak: 272.10(595466) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	272.10	595466	100.00				4	313.10	47323	7.95			
2	273.10	116312	19.53				5	314.10	8585	1.44			
3	274.10	10013	1.68										

Figure S31. LCMS spectrum of compound 17.




Figure S32. ¹H NMR spectrum of compound 17.






Figure S33. ¹³C NMR spectrum of compound 17.

Acquired by : Admin
 Date Acquired : 7/3/2019 4:10:39 PM
 Sample Name : YAZH_23-1
 Sample ID :
 Tray# : 1
 Vial# : 17
 Injection Volume : 4
 Data File : C:\LabSolutions\Data\2019\2019-wk09\YAZH_23-1.lcd
 Background File : blanco 07032019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 7/3/2019 4:29:18 PM



PDA Ch1 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		5.077	21334096	100.000



MS Spectrum Table

#1 Ret.Time:
 BG Mode:Calc 5.020<->5.570(503<->558)
 Mass Peaks:5 Base Peak:256.05(1386430) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	256.05	1386430	100.00				4	297.10	296569	21.39			
2	257.10	198226	14.30				5	298.10	45252	3.26			
3	258.10	25936	1.87										

Figure S34. LCMS spectrum of compound **18a**.

YAZH_23-1/1




Figure S35. ¹H NMR spectrum of compound **18a**.

YAZH_23-1/3




Figure S36. ¹³C NMR spectrum of compound **18a**.







Figure S37. 1D NOESY spectrum of compound **18a** and its regioisomer.

Acquired by : Admin
 Date Acquired : 18/6/2019 9:45:35 AM
 Sample Name : MOIB_23_1
 Sample ID :
 Tray# : 1
 Vial# : 3
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk25\MOIB_23_1.lcd
 Background File : blanco 18062019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 18/6/2019 10:59:57 AM



MS Spectrum Table

#1 Ret. Time:
 BG Mode: Calc 5.580-<->5.880(559-<->589)
 Mass Peaks: 10 Base Peak: 284.10(196283) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	280.25	3989	2.03				6	296.15	4447	2.27			
2	282.10	5645	2.88				7	323.15	2573	1.31			
3	284.10	196283	100.00				8	325.15	57615	29.35			
4	285.05	35064	17.86				9	326.10	14474	7.37			
5	286.20	3600	1.83				10	418.30	2919	1.49			

Figure S38. LCMS spectrum of compound **18b**.




Figure S39. ¹H NMR spectrum of compound **18b**.




Figure S40. ¹³C NMR spectrum of compound **18b**.







Figure S41. 1D NOESY spectrum of compound **18b** and its regioisomer.

Acquired by : Admin
 Date Acquired : 19/6/2019 10:34:53 AM
 Sample Name : YAZH_43-1
 Sample ID :
 Tray# : 1
 Vial# : 8
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk25\YAZH_43-1.lcd
 Background File : blanco 19062019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 19/6/2019 4:22:46 PM



PDA Ch1 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		5.699	4966643	100.000



MS Spectrum Table

#1 Ret.Time:
 BG Mode: Calc 5.650-<->6.020(566-<->603)
 Mass Peaks: 15 Base Peak: 284.10(224105) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	181.95	3035	1.35				9	285.10	41155	18.36			
2	200.05	4316	1.93				10	286.35	3826	1.71			
3	240.90	4325	1.93				11	296.20	2329	1.04			
4	242.05	4400	1.96				12	323.05	3143	1.40			
5	269.25	2361	1.05				13	325.10	21110	9.42			
6	280.15	3036	1.35				14	326.15	3795	1.69			
7	282.30	5961	2.66				15	418.00	2245	1.00			
8	284.10	224105	100.00										

Figure S42. LCMS spectrum of compound **18c**.




Figure S43. ¹H NMR spectrum of compound **18c**.




Figure S44. ¹³C NMR spectrum of compound **18c**.







Figure S45. 1D NOESY spectrum of compound **18c**.

Acquired by : Admin
 Date Acquired : 21/3/2019 12:05:23 PM
 Sample Name : MOIB_1_pure_biggest_peak
 Sample ID :
 Tray# : 1
 Vial# : 9
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk11\MOIB_1_pure_biggest_peak.lcd
 Background File : blanco 21032019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 30/3/2021 4:44:59 PM



Peak#	Name	Ret. Time	Area	Area %
1		4.125	71087	6.087
2		4.606	1096754	93.913



MS Spectrum Table

#1 Ret. Time:
 BG Mode: Calc 4.570<->4.950 (458<->496)
 Mass Peaks: 5 Base Peak: 316.10 (579289) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	240.15	7836	1.35				4	317.10	82577	14.25			
2	247.05	8419	1.45				5	318.15	11594	2.00			
3	316.10	579289	100.00										

Figure S46. LCMS spectrum of compound **18d**.




Figure S47. ¹H NMR spectrum of compound **18d**.




Figure S48. ¹³C NMR spectrum of compound **18d**.







Figure S49. 1D NOESY spectrum of compound **18d** and its regioisomer.

Acquired by : Admin
 Date Acquired : 26/3/2019 4:14:25 PM
 Sample Name : MOIB_2_pure_big
 Sample ID :
 Tray# : 1
 Vial# : 32
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk13\MOIB_2_pure_big.lcd
 Background File : blanco 26032019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 6/4/2021 10:55:03 AM



Peak#	Name	Ret. Time	Area	Area %
1		4.022	714	0.098
2		5.014	1235	0.169
3		5.752	728241	99.733



MS Spectrum Table

#2 Ret. Time:
 BG Mode: Calc 5.730<->5.960(574<->597)
 Mass Peaks: 13 Base Peak: 468.20(50927) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	217.90	1105	2.17				8	468.20	50927	100.00			
2	282.10	845	1.66				9	469.15	11685	22.94			
3	322.75	1627	3.19				10	470.20	3443	6.76			
4	329.10	873	1.71				11	506.25	4792	9.41			
5	348.05	675	1.33				12	507.25	1404	2.76			
6	457.40	5176	10.16				13	508.45	1305	2.56			
7	458.25	1603	3.15										

Figure S50. LCMS spectrum of compound **18e**.




Figure S51. ¹H NMR spectrum of compound **18e**.




Figure S52. ¹³C NMR spectrum of compound **18e**.







Figure S53. 1D NOESY spectrum of compound **18e** and its regioisomer.

Acquired by : Admin
 Date Acquired : 2/7/2019 11:38:22 AM
 Sample Name : YAZH_44
 Sample ID :
 Tray# : 1
 Vial# : 3
 Injection Volume : 3
 Data File : C:\LabSolutions\Data\2019\2019-wk27\YAZH_44-20190702.lcd
 Background File : 19.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultLCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg\01072015a.lct
 Processed by : Admin
 Modified Date : 2/7/2019 12:25:47 PM



PDA Ch1 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		4.047	6035513	95.733
2		4.223	269043	4.267



MS Spectrum Table

#1 Ret.Time:
 BG Mode:Calc 4.000<->4.220(401<->423)
 Mass Peaks:10 Base Peak:241.05(208535) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	198.90	3110	1.49				6	241.05	208535	100.00			
2	209.05	7479	3.59				7	242.00	22572	10.82			
3	209.95	2828	1.36				8	282.05	29765	14.27			
4	224.00	26262	12.59				9	283.10	4444	2.13			
5	225.05	4312	2.07				10	354.10	2765	1.33			

Figure S54. LCMS spectrum of compound 19c.




Figure S55. ¹H NMR spectrum of compound **19c**.







Figure S56. ¹³C NMR spectrum of compound **19c**.

Acquired by : Admin
 Date Acquired : 1/5/2019 4:50:28 PM
 Sample Name : MOIB_9-6
 Sample ID :
 Tray# : 1
 Vial# : 54
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk18\MOIB_9-6.lcd
 Background File : blanco 01052019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 6/4/2021 10:57:17 AM



Peak#	Name	Ret. Time	Area	Area %
1		4.374	1419030	99.676
2		4.691	4610	0.324



MS Spectrum Table

#1 Ret Time:
 BG Mode: Calc 4.310<->4.690(432<->470)
 Mass Peaks: 8 Base Peak: 333.10(505858) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	255.90	6300	1.25				5	336.20	5295	1.05			
2	333.10	505858	100.00				6	360.10	23540	4.65			
3	334.15	92961	18.38				7	380.85	15211	3.01			
4	335.15	11248	2.22				8	526.30	5702	1.13			

Figure S57. LCMS spectrum of compound 19e.




Figure S58. $^1\text{H NMR}$ spectrum of compound **19e**.







Figure S59. $^{13}\text{C NMR}$ spectrum of compound **19e**.

Acquired by : Admin
 Date Acquired : 26/6/2019 10:12:06 AM
 Sample Name : MOIB_20_pure
 Sample ID :
 Tray# : 1
 Vial# : 3
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk26\MOIB_20_pure.lcd
 Background File : Blanco 26062019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultLCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 26/6/2019 10:41:04 AM



Peak#	Name	Ret. Time	Area	Area %
1		2.844	1429628	100.000



MS Spectrum Table

#1 Ret.Time:
 BG Mode:Calc 2.790<->3.120(280<->313)
 Mass Peaks:14 Base Peak:270.05(544210) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	199.05	12940	2.38				8	307.95	8087	1.49			
2	253.05	28192	5.18				9	315.10	12831	2.36			
3	270.05	544210	100.00				10	333.10	83124	15.27			
4	271.10	63109	11.60				11	334.15	11363	2.09			
5	272.10	8015	1.47				12	431.50	10826	1.99			
6	292.05	119659	21.99				13	431.85	7007	1.79			
7	293.05	16508	3.03				14	561.25	10716	1.97			

Figure S60. LCMS spectrum of compound 19f.




Figure S61. ^1H NMR spectrum of compound **19f**.







Figure S62. ^{13}C NMR spectrum of compound **19f**.

Acquired by : Admin
 Date Acquired : 21/6/2019 1:13:59 PM
 Sample Name : MOIB_27-af
 Sample ID :
 Tray# : 1
 Vial# : 23
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk25\MOIB_27-af.lcd
 Background File : blanco_an 21062019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 21/6/2019 1:27:50 PM



Peak#	Name	Ret. Time	Area	Area %
1		2.904	1902079	100.000



MS Spectrum Table

#1 Ret.Time:
 BG Mode: Calc 2.860<->3.300(287<->331)
 Mass Peaks: 7 Base Peak: 211.05(1030108) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	194.00	309501	30.05				5	213.10	13033	1.27			
2	195.05	36760	3.57				6	252.10	37148	3.61			
3	211.05	1030108	100.00				7	253.05	14027	1.36			
4	212.05	159664	15.50										

Figure S63. LCMS spectrum of compound **20b**.




Figure S64. ¹H NMR spectrum of compound **20b**.







Figure S65. ¹³C NMR spectrum of compound **20b**.

Acquired by : Admin
 Date Acquired : 10/4/2019 9:48:23 AM
 Sample Name : MOIB_6_pure
 Sample ID :
 Tray# : 1
 Vial# : 44
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk15\MOIB_6_pure.lcd
 Background File : blanco 10042019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 30/3/2021 5:22:22 PM



Peak#	Name	Ret. Time	Area	Area %
1		2.426	794189	98.234
2		2.956	14276	1.766



MS Spectrum Table

#1 Ret Time: 2.488
 BG Mode: Calc 2.380<->2.740(239<->275)
 Mass Peaks: 8 Base Peak: 227.10(974979) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	178.05	16005	1.64				5	228.05	121523	12.46			
2	210.05	73336	7.52				6	229.10	11674	1.20			
3	211.10	12632	1.30				7	249.05	60850	6.24			
4	227.10	974979	100.00				8	290.10	32523	3.34			

Figure S66. LCMS spectrum of compound **20d**.




Figure S67. ¹H NMR spectrum of compound **20d**.







Figure S68. ¹³C NMR spectrum of compound **20d**.

Acquired by : Admin
 Date Acquired : 3/5/2019 4:04:24 PM
 Sample Name : MOIB_13-5
 Sample ID :
 Tray# : 1
 Vial# : 42
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk18\MOIB_13-5.lcd
 Background File : blanco 03052019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultLCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 3/5/2019 4:13:44 PM



Peak#	Name	Ret. Time	Area	Area %
1		3.586	1786961	100.000



MS Spectrum Table

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	286.05	521756	16.33				4	304.15	606797	18.99			
2	287.10	88123	2.76				5	305.15	63853	2.00			
3	303.15	3195339	100.00										

Figure S69. LCMS spectrum of compound 20e.




Figure S70. ¹H NMR spectrum of compound 20e.







Figure S71. ¹³C NMR spectrum of compound 20e.

Acquired by : Admin
 Date Acquired : 17/4/2019 1:15:29 PM
 Sample Name : YAZH_40_20190417
 Sample ID :
 Tray# : 1
 Vial# : 7
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk16\YAZH_40_20190417.lcd
 Background File : blanco_17042019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 17/4/2019 2:32:24 PM



Peak#	Name	Ret. Time	Area	Area %
1		4.813	187701	100.000



MS Spectrum Table

#1 Ret Time:
 BG Mode: Calc 4.770-<->5.070(478-<->508)
 Mass Peaks: 5 Base Peak: 301.10(387117) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	301.10	387117	100.00				4	336.10	3948	1.02			
2	302.10	72822	18.81				5	366.25	9420	2.43			
3	303.15	7328	1.89										

Figure S72. LCMS spectrum of compound **21a**.




Figure S73. ¹H NMR spectrum of compound 21a.







Figure S74. ¹³C NMR spectrum of compound 21a.

Acquired by : Admin
 Date Acquired : 27/6/2019 2:18:26 PM
 Sample Name : MOIB_28_rc_pure
 Sample ID :
 Tray# : 1
 Vial# : 43
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk26\MOIB_28_rc_pure.lcd
 Background File : blanco 27062019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultLCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 27/6/2019 2:31:01 PM



Peak Table

Peak#	Name	Ret. Time	Area	Area %
1		5.170	1889998	100.000



MS Spectrum Table

#1 Ret. Time:
 BG Mode:None
 Mass Peaks:8 Base Peak:315.10(917553) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	315.10	917553	100.00				5	337.05	18540	2.02			
2	316.10	194142	21.16				6	378.05	41095	4.48			
3	317.00	20939	2.28				7	379.10	9874	1.08			
4	318.20	11618	1.27				8	446.00	9436	1.03			

Figure S75. LCMS spectrum of compound **21b**.




Figure S76. ¹H NMR spectrum of compound **21b**.







Figure S77. ¹³C NMR spectrum of compound **21b**.

Acquired by : Admin
 Date Acquired : 10/7/2019 10:42:25 AM
 Sample Name : YAZH_46
 Sample ID :
 Tray# : 1
 Vial# : 8
 Injection Volume : 2
 Data File : C:\LabSolutions\Data\2019\2019-wk28\YAZH_46.lcd
 Background File : blanco10072019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : Default.LCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 10/7/2019 11:01:10 AM



Peak#	Name	Ret. Time	Area	Area %
1		5.169	1664	0.175
2		5.338	945899	99.746
3		5.725	748	0.079



MS Spectrum Table

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	273.00	21147	2.03				4	317.15	21843	2.09			
2	315.10	1043830	100.00				5	378.15	18923	1.81			
3	316.10	207896	19.92										

Figure S78. LCMS spectrum of compound **21c**.




Figure S79. ¹H NMR spectrum of compound **21c**.






Figure S80. ¹³C NMR spectrum of compound **21c**.


Acquired by : Admin
 Date Acquired : 26/4/2019 12:01:34 PM
 Sample Name : MOIB_10_final_pure
 Sample ID :
 Tray# : 1
 Vial# : 16
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk17\MOIB_10_final_pure.lcd
 Background File : blanco 26042019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 17/6/2019 4:58:25 PM



Peak#	Name	Ret. Time	Area	Area %
1		4.469	5959	0.482
2		4.581	1231102	99.518



#:5 Ret.Time:Averaged 4.620-4.640(Scan#:-463-465)
 BG Mode:Calc 4.510<->4.860(452<->487)



MS Spectrum Table

#:5 Ret.Time:
 BG Mode:Calc 4.510<->4.860(452<->487)

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	331.10	1435605	100.00				4	353.10	20172	1.41			
2	332.10	336947	23.47				5	394.15	17736	1.24			
3	333.10	36941	2.57										

Figure S81. LCMS spectrum of compound **21d**.




Figure S82. ¹H NMR spectrum of compound 21d.







Figure S83. ¹³C NMR spectrum of compound 21d.

Acquired by : Admin
 Date Acquired : 7/5/2019 5:45:32 PM
 Sample Name : MOIB_14_RC_pure
 Sample ID :
 Tray# : 1
 Vial# : 64
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk19\MOIB_14_RC_pure.lcd
 Background File : blanco 07052019 lcd
 Method File : Method SCAN.ACID standard.lcm
 Report Format : Default1.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 17/6/2019 5:01:55 PM



PeakTable

Peak#	Name	Ret. Time	Area	Area %
1		5.348	268843	98.593
2		5.574	3836	1.407



MS Spectrum Table

#6 Ret.Time:
 BG Mode: Calc 5.300<->5.580(531<->559)
 Mass Peaks: 7 Base Peak: 407.15(132273) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	341.30	1343	1.02				5	429.10	11031	8.34			
2	407.15	132273	100.00				6	430.20	3538	2.67			
3	408.15	33503	25.33				7	470.20	3787	2.86			
4	409.35	4694	3.55										

Figure S84. LCMS spectrum of compound **21e**.




Figure S85. ^1H NMR spectrum of compound **21e**.






Figure S86. ^{13}C NMR spectrum of compound **21e**.


Acquired by : Admin
 Date Acquired : 10/7/2019 5:27:37 PM
 Sample Name : MOIB_33_rc
 Sample ID :
 Tray# : 1
 Vial# : 51
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk28\MOIB_33_rc.lcd
 Background File : blanco10072019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : Default.LCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 11/7/2019 2:39:12 PM



Peak#	Name	Ret. Time	Area	Area %
1		3.610	3354	1.806
2		3.839	178768	96.262
3		3.936	3589	1.932



#2 Ret. Time: Averaged 3.890-3.910(Scan#:390-392)
Mass Peaks:22 Base Peak:344.05(236385) Polarity:Pos Segment1 - Event1



#2 Ret. Time:
BG Mode:Calc 3.790<->4.090(380<->410)

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	240.80	4303	1.82				12	407.15	20315	8.59			
2	273.05	8266	3.50				13	408.15	6037	2.55			
3	327.05	4572	1.93				14	444.10	2386	1.01			
4	344.05	236385	100.00				15	542.75	9356	3.96			
5	345.15	50190	21.23				16	543.75	3479	1.47			
6	346.15	6162	2.61				17	709.35	12682	5.36			
7	366.10	61037	25.82				18	710.35	2952	1.25			
8	367.05	11640	4.92				19	710.65	2410	1.02			
9	371.30	3347	1.42				20	714.30	4332	1.83			
10	371.80	3246	1.37				21	714.95	5069	2.14			
11	382.15	2841	1.20				22	787.25	3681	1.56			

Figure S87. LCMS spectrum of compound 21f.




Figure S88. ¹H NMR spectrum of compound **21f**.







Figure S89. ¹³C NMR spectrum of compound **21f**.

Acquired by : Admin
 Date Acquired : 11/7/2019 11:19:41 AM
 Sample Name : MOIB_36_final
 Sample ID :
 Tray# : 1
 Vial# : 12
 Injection Volume : 3
 Data File : C:\LabSolutions\Data\2019\2019-wk28\MOIB_36_final.lcd
 Background File : blanco 11072019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 31/3/2021 10:51:15 AM



PDA Ch1 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		4.048	1436731	97.987
2		5.080	29518	2.013



MS Spectrum Table

#1 Ret. Time:
 BG Mode Calc 3.980<->4.370 (399<->438)
 Mass Peaks: 7 Base Peak: 317.10 (916094) Polarity: Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	317.10	916094	100.00				5	380.20	17058	1.86			
2	318.10	162880	17.78				6	448.00	13229	1.44			
3	319.05	17836	1.95				7	458.00	9755	1.06			
4	339.10	16003	1.75										

Figure S90. LCMS spectrum of compound 22.




Figure S91. ¹H NMR spectrum of compound 22.







Figure S92. ¹³C NMR spectrum of compound 22.

Acquired by : Admin
 Date Acquired : 29/5/2019 12:37:27 PM
 Sample Name : MOIB_18_pure
 Sample ID :
 Tray# : 1
 Vial# : 21
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk22\MOIB_18_pure.lcd
 Background File : blanco 29052019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 29/5/2019 1:03:53 PM



PeakTabl

Peak#	Name	Ret. Time	Area	Area %
1		3.578	13616793	100.000



MS Spectrum Table

#1 RetTime:
BG Mode:Calc 3.560<->3.830(357<->384)
Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	181.95	10119	5.30				6	215.00	19639	10.29			
2	184.10	2766	1.45				7	215.90	2171	1.14			
3	199.95	16007	8.39				8	241.00	6945	3.64			
4	201.15	2044	1.07				9	255.15	4946	2.59			
5	214.00	190889	100.00										

Figure S93. LCMS spectrum of compound **23**.




Figure S94. ¹H NMR spectrum of compound **23**.






Figure S95. ¹³C NMR spectrum of compound **23**.

Acquired by : Admin
 Date Acquired : 16/7/2018 10:51:14 AM
 Sample Name : LABA-022
 Sample ID :
 Tray# : 1
 Vial# : 9
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2018\2018-wk29\LABA-022.lcd
 Background File : blanco 16072018.lcd
 Method File : Method SCAN ACID standard neg.lcm
 Report Format : Default.LCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 16/7/2018 11:08:48 AM




Peak#	Name	Ret. Time	Area	Area %
1		3.189	765106	100.000



#1 Ret.Time: Averaged 3.230-3.250(Scan#: 324-326)

Mass Peaks:24 Base Peak:179.00(141779) Polarity:Neg Segment1 - Event1



MS Spectrum Table



#1 Ret.Time:

BG Mode:Calc 3.140<->3.460(315<->347)


Mass Peaks:24 Base Peak:179.00(141779) Polarity:Neg Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	179.00	141779	100.00				13	381.25	9413	6.64			
2	180.25	11417	8.05				14	382.20	3611	2.55			
3	181.20	1755	1.24				15	413.10	2066	1.46			
4	220.15	2520	1.78				16	439.15	1523	1.07			
5	224.15	3579	2.52				17	443.25	1420	1.00			
6	225.05	34564	24.38				18	486.00	1824	1.29			
7	247.05	4057	2.86				19	561.95	2607	1.84			
8	279.10	2230	1.57				20	585.05	2845	2.01			
9	280.00	2476	1.75				21	645.70	2233	1.57			
10	358.10	1454	1.03				22	795.40	2152	1.52			
11	359.15	67854	47.86				23	858.30	2564	1.81			
12	360.15	11941	8.42				24	1051.00	2181	1.54			





Figure S96. LCMS spectrum of compound 27a.



Acquired by : Admin
 Date Acquired : 23/10/2018 2:30:32 PM
 Sample Name : ODVI01_018 crude 231018
 Sample ID :
 Tray# : 1
 Vial# : 17
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2018\2018-wk43\ODVI01_018 crude 231018.lcd
 Background File : 23102018 blanco.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 23/10/2018 2:40:02 PM



Peak#	Name	Ret. Time	Area	Area %
1		3.465	452342	97.867
2		4.012	1072	0.232
3		4.319	8786	1.901



MS Spectrum Table

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	177.05	47070	8.04				5	236.05	317612	54.23			
2	178.05	6493	1.11				6	237.10	42704	7.29			
3	195.05	585623	100.00				7	389.25	35464	6.06			
4	196.05	79866	13.64				8	390.30	8125	1.39			

Figure S99. LCMS spectrum of compound **27b**.







Figure S100. ¹H NMR spectrum of compound **27b**.

Acquired by : Admin
 Date Acquired : 24/7/2018 10:42:48 AM
 Sample Name : LABA-026-CRUDE-PR
 Sample ID :
 Tray# : 1
 Vial# : 15
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2018\2018-wk30\LABA-026-CRUDE-PR.lcd
 Background File : blanco 24072018.lcd
 Method File : Method SCAN ACID standard neg.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 24/7/2018 12:11:18 PM



11 PDA Multi 1 / 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		3.252	1116635	100.000



MS Spectrum Table
#1 Ret. Time:
BG Mode:Calc 3.210<->3.480(322<->349)
Mass Peaks:14 Base Peak:224.05(105117) Polarity:Neg Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	180.05	89072	84.74				8	301.85	1934	1.84			
2	181.05	11197	10.65				9	321.95	3422	3.26			
3	197.95	3988	3.79				10	439.75	1971	1.88			
4	223.10	1357	1.29				11	449.70	1347	1.28			
5	224.05	105117	100.00				12	470.70	1270	1.21			
6	225.05	23822	22.66				13	471.70	3498	3.33			
7	225.95	2303	2.19				14	508.60	1327	1.26			

Figure S101. LCMS spectrum of compound **28a**.




Figure S102. ^1H NMR spectrum of compound **28a**.







Figure S103. ^{13}C NMR spectrum of compound **28a**.

Acquired by : Admin
 Date Acquired : 24/5/2018 12:30:40 PM
 Sample Name : YAZH01-126
 Sample ID :
 Tray# : 1
 Vial# : 12
 Injection Volume : 3
 Data File : C:\LabSolutions\Data\2018\2018-wk21\YAZH01-126.lcd
 Background File : blanco_24052018.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 24/5/2018 2:23:18 PM



Peak#	Name	Ret. Time	Area	Area %
1		3.259	3089630	96.427
2		4.166	87842	2.742
3		4.777	26639	0.831



MS Spectrum Table

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	179.95	1133	12.40				8	261.85	1089	11.92			
2	195.95	2092	22.90				9	262.85	1727	18.91			
3	214.00	9135	100.00				10	263.85	1142	12.50			
4	215.05	1374	15.04				11	273.85	1949	21.34			
5	222.00	3478	38.07				12	298.80	1197	13.10			
6	254.85	2347	25.69				13	339.75	1002	10.97			
7	257.90	1253	13.72				14	493.50	1137	12.45			

Figure S104. LCMS spectrum of compound **28c**.




Figure S105. ¹H NMR spectrum of compound **28c**.






Figure S106. ¹³C NMR spectrum of compound **28c**.

Acquired by : Admin
 Date Acquired : 17/5/2018 4:33:53 PM
 Sample Name : YAZH01-215
 Sample ID :
 Tray# : 1
 Vial# : 12
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2018\2018-wk20\YAZH01-215-1.lcd
 Background File : blanco 170518.lcd
 Method File : Method SCAN ACID standard neg.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 17/5/2018 6:33:28 PM



PDA Ch1 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		2.231	783773	100.000



MS Spectrum Table

#1 Ret.Time:
 BG Mode:None
 Mass Peaks:757 Base Peak:183.85(23845) Polarity:Neg Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	151.00	478	2.00				15	169.00	336	1.41			
2	153.00	288	1.21				16	170.00	401	1.68			
3	154.00	360	1.51				17	172.00	389	1.63			
4	155.00	258	1.08				18	174.00	423	1.77			
5	156.00	383	1.61				19	176.00	249	1.04			
6	158.00	460	1.93				20	177.00	363	1.52			
7	159.00	395	1.66				21	179.00	368	1.54			
8	160.00	305	1.28				22	181.00	404	1.69			
9	161.00	352	1.48				23	183.00	439	1.84			
10	163.00	487	2.04				24	183.85	23845	100.00			
11	164.00	470	1.97				25	184.90	3538	14.84			
12	165.00	368	1.54				26	185.90	1748	7.33			
13	167.00	520	2.18				27	187.90	313	1.31			
14	168.00	501	2.10				28	188.90	460	1.93			

Figure S107. LCMS spectrum of compound **28d**.




Figure S108. ¹H NMR spectrum of compound **28d**.









Figure S109. ¹³C NMR spectrum of compound **28d**.



Acquired by : Admin
 Date Acquired : 8/21/2017 3:12:01 PM
 Sample Name : YAZH01-134
 Sample ID :
 Tray# : 1
 Vial# : 40
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2017\2017-wk34\YAZH01-134.lcd
 Background File : BLANCO21082017.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 8/21/2017 3:40:37 PM





Peak#	Name	Ret. Time	Area	Area %
1		1.596	36034	4.373
2		1.687	787999	95.627




MS Spectrum Table

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	171.00	20339	100.00				7	216.10	7347	36.12			
2	172.20	2238	11.00				8	229.05	1522	7.48			
3	193.05	9397	46.20				9	234.00	4866	23.92			
4	194.00	1006	4.95				10	359.85	1012	4.98			
5	209.00	2641	12.98				11	363.05	11123	54.69			
6	211.00	1085	5.33				12	364.10	2179	10.71			

Figure S113. LCMS spectrum of compound 29e.





Acquired by : Admin
 Date Acquired : 8/25/2017 11:08:22 AM
 Sample Name : YAZH01-137-1
 Sample ID :
 Tray# : 1
 Vial# : 12
 Injection Volume : 10
 Data File : C:\LabSolutions\Data\2017\2017-wk34\YAZH01-137-1.lcd
 Background File : BLANCO25082017.lcd
 Method File : Method SCAN ACID standard MW50.lcm
 Report Format : DefaultLCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 8/25/2017 12:14:03 PM



Peak Table

Peak#	Name	Ret. Time	Area	Area %
1		0.711	3474924	96.648
2		1.002	120517	3.352



MS Spectrum Table

#1 Ret.Time:
BG Mode:Calc 0.730<->0.990(74<->100)
Mass Peaks:32 Base Peak:141.00(571800) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	124.00	33633	5.88				9	169.05	94928	16.60			
2	141.00	571800	100.00				10	170.10	12503	2.19			
3	142.00	137717	24.08				11	183.10	25925	4.53			
4	143.00	10576	1.85				12	191.05	17087	2.99			
5	152.10	9772	1.71				13	195.15	6916	1.21			
6	156.00	19768	3.46				14	204.05	14362	2.51			
7	163.00	127066	22.22				15	205.05	94767	16.57			
8	163.95	11210	1.96				16	206.05	17858	3.12			

Figure S116. LCMS spectrum of compound 30e.







Figure S117. ¹H NMR spectrum of compound 30e.

Acquired by : Admin
 Date Acquired : 17/7/2018 3:20:53 PM
 Sample Name : LABA-025-ENDPR
 Sample ID :
 Tray# : 1
 Vial# : 32
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2018\2018-wk29\LABA-025-ENDPR.lcd
 Background File : blanco 17072018.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : Default.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 17/7/2018 3:56:31 PM



MS Spectrum Table
#1 Ret Time:
BG Mode: Calc 4.270-<->4.490 (428-<->450)

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	299.15	141864	100.00				7	340.20	14881	10.49			
2	300.15	27404	19.32				8	341.15	3574	2.52			
3	301.35	2464	1.74				9	362.15	20335	14.33			
4	321.10	6833	4.82				10	363.15	4944	3.49			
5	322.15	2042	1.44				11	619.25	4795	3.38			
6	337.15	2872	2.02				12	620.20	2208	1.56			

Figure S118. LCMS spectrum of compound **31a**.




Figure S119. ^1H NMR spectrum of compound 31a.








Figure S120. ^{13}C NMR spectrum of compound 31a.

Acquired by : Admin
 Date Acquired : 10/12/2018 9:07:44 AM
 Sample Name : ODVI01_044 3-6 10122018
 Sample ID :
 Tray# : 1
 Vial# : 2
 Injection Volume : 2
 Data File : C:\LabSolutions\Data\2018\2018-wk50\ODVI01_044 3-6 10122018.lcd
 Background File : blanco 10122018.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 10/12/2018 9:17:37 AM



1 PDA Multi 1 / 254nm 4nm
 PDA Ch1 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		3.761	942	0.351
2		4.530	266986	99.649



MS Spectrum Table
 #1 Ret Time:
 BG Mode:Calc 4.470<->4.780(448<->479)
 Mass Peaks:12 Base Peak:313.10(356819) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	313.10	356819	100.00				7	354.15	43530	12.20			
2	314.10	75772	21.24				8	355.15	8991	2.52			
3	315.05	8093	2.27				9	376.10	44498	12.47			
4	335.15	21537	6.04				10	377.05	11658	3.27			
5	336.15	5160	1.45				11	647.20	6549	1.84			
6	351.15	6088	1.71				12	652.50	5918	1.66			

MS Spectrum Table
 #2 Ret Time:
 BG Mode:None
 Mass Peaks:7 Base Peak:258.15(21219) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	258.15	21219	100.00				5	332.10	1029	4.85			
2	259.35	2860	13.48				6	515.80	1603	7.55			
3	321.25	1058	4.99				7	537.75	1177	5.55			
4	331.10	2095	9.87										

Figure S121. LCMS spectrum of compound **31b**.




Figure S122. ¹H NMR spectrum of compound **31b**.






Figure S123. ¹³C NMR spectrum of compound **31b**.

Acquired by : Admin
 Date Acquired : 7/8/2018 10:24:01 AM
 Sample Name : YAZH01-233
 Sample ID :
 Tray# : 1
 Vial# : 6
 Injection Volume : 3
 Data File : C:\LabSolutions\Data\2018\2018-wk32\YAZH01-233.lcd
 Background File : blanco 070818.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultL.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 7/8/2018 10:33:54 AM




Peak#	Name	Ret. Time	Area	Area %
1		4.342	1344980	100.000



MS Spectrum Graph

#1 Ret. Time: Averaged 4.400-4.420(Scan#:441-443)
 Mass Peaks:15 Base Peak:287.15(323082) Polarity:Pos Segment1 - Event1



MS Spectrum Table

#1 Ret. Time:
 RG Mode:Calc 4 290<->4 600(430<->461)

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	287.15	323082	100.00				9	350.10	85019	26.31			
2	288.15	68146	21.09				10	351.20	15835	4.90			
3	289.25	5304	1.64				11	456.95	4410	1.36			
4	309.10	40327	12.48				12	595.40	10554	3.27			
5	310.10	7211	2.23				13	596.45	3293	1.02			
6	325.05	5165	1.60				14	600.45	7979	2.47			
7	328.15	70804	21.92				15	600.70	10022	3.10			
8	329.15	14802	4.58										

Figure S124. LCMS spectrum of compound **31c**.




Figure S125. $^1\text{H NMR}$ spectrum of compound **31c**.







Figure S126. $^{13}\text{C NMR}$ spectrum of compound **31c**.

Acquired by : Admin
 Date Acquired : 27/6/2018 12:16:59 PM
 Sample Name : YAZH01-221
 Sample ID :
 Tray# : 1
 Vial# : 46
 Injection Volume : 3
 Data File : C:\LabSolutions\Data\2018\2018-wk26\YAZH01-221.lcd
 Background File : blanco 27062018.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : Default.CMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 27/6/2018 12:35:32 PM



PeakTable

Peak#	Name	Ret. Time	Area	Area %
1		3.520	963705	99.477
2		4.018	5062	0.523



MS Spectrum Table

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	259.10	451571	100.00				6	300.15	43860	9.71			
2	260.10	73084	16.18				7	301.15	8388	1.86			
3	261.05	7254	1.61				8	322.15	67636	14.98			
4	281.10	27348	6.06				9	323.05	17785	3.94			
5	297.00	6123	1.36				10	539.30	18577	4.11			

Figure S127. LCMS spectrum of compound **31d**.




Figure S128. ¹H NMR spectrum of compound **31d**.







Figure S129. ¹³C NMR spectrum of compound **31d**.

Acquired by : Admin
 Date Acquired : 20/10/2017 4:13:41 PM
 Sample Name : YAZH01-153-4
 Sample ID :
 Tray# : 1
 Vial# : 14
 Injection Volume : 10
 Data File : C:\LabSolutions\Data\2017\2017-wk42\YAZH01-153-4.lcd
 Background File : blanco 20102017.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultLCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 7/6/2019 5:51:55 PM



PDA Ch1 230nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		2.755	39494	0.905
2		3.140	4178158	95.689
3		3.296	148718	3.406



MS Spectrum Table

#:1 Ret.Time:
 BG Mode:Calc 3.080<->3.620(309<->363)
 Mass Peaks:14 Base Peak:313.10(616217) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	245.05	450100	73.04				8	315.15	13252	2.15			
2	246.05	69475	11.27				9	331.10	20953	3.40			
3	267.10	6364	1.03				10	335.10	10290	1.67			
4	271.10	11398	1.85				11	376.20	6628	1.08			
5	286.05	24835	4.03				12	647.35	93552	15.18			
6	313.10	616217	100.00				13	648.35	31493	5.11			
7	314.10	112209	18.21				14	649.45	6417	1.04			

Figure S130. LCMS spectrum of compound **31e**.




Figure S131. ¹H NMR spectrum of compound **31e**.






Figure S132. ¹³C NMR spectrum of compound **31e**.

Acquired by : Admin
 Date Acquired : 18/6/2019 3:39:49 PM
 Sample Name : VUF16375
 Sample ID :
 Tray# : 1
 Vial# : 35
 Injection Volume : 5
 Data File : C:\LabSolutions\Data\2019\2019-wk25\VUF16375.lcd
 Background File : blanco 18062019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : DefaultLCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 18/6/2019 3:51:25 PM




PDA Ch1 254nm 4nm

Peak#	Name	Ret. Time	Area	Area %
1		2.906	7092524	98.714
2		3.170	46741	0.651
3		3.517	45662	0.636



#:1 Ret.Time:Averaged 2.960-2.980(Scan#:297-299)
 BG Mode:Calc 2.870<->3.180(288<->319)
 Mass Peaks:17 Base Peak:231.00(441377) Polarity:Pos Segment1 - Event1



MS Spectrum Table

#:1 Ret.Time:
BG Mode:Calc 2.870<->3.180(288<->319)
Mass Peaks:17 Base Peak:231.00(441377) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	231.00	441377	100.00				10	294.00	107749	24.41			
2	231.95	58109	13.17				11	294.95	19940	4.52			
3	232.90	4546	1.03				12	309.00	12625	2.86			
4	253.05	69571	15.76				13	330.95	14689	3.33			
5	254.05	9430	2.14				14	405.10	5598	1.27			
6	269.00	14447	3.27				15	480.45	4520	1.02			
7	272.05	22163	5.02				16	483.15	12558	2.85			
8	273.00	4580	1.04				17	484.15	4810	1.09			
9	276.10	6628	1.50										

Figure S133. LCMS spectrum of compound **31f**.




Figure S134. ¹H NMR spectrum of compound **31f**.






Figure S135. ¹³C NMR spectrum of compound **31f**.

Acquired by : Admin
 Date Acquired : 21/2/2019 2:41:40 PM
 Sample Name : LABI-009-NANP
 Sample ID :
 Tray# : 1
 Vial# : 11
 Injection Volume : 1
 Data File : C:\LabSolutions\Data\2019\2019-wk07\LABI-009-NANP.lcd
 Background File : blanco 21022019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : Default.LCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 19/7/2019 12:06:34 PM



PeakTab1

Peak#	Name	Ret. Time	Area	Area %
1		3.095	2251	0.291
2		3.229	8753	1.130
3		3.509	763389	98.579



MS Spectrum Table

#1 Ret.Time:
BG Mode:Calc 3.460<->3.790(347<->380)
Mass Peaks:13 Base Peak:232.05(477051) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	185.00	31347	6.37				8	254.05	26742	5.61			
2	204.05	59882	12.55				9	273.10	113864	23.87			
3	205.05	7954	1.67				10	274.10	16207	3.40			
4	232.05	477051	100.00				11	295.10	62053	13.01			
5	233.05	70683	14.82				12	296.05	10037	2.10			
6	234.05	5072	1.06				13	485.25	5272	1.11			
7	245.05	22245	4.66										

Figure S136. LCMS spectrum of compound 33.




Figure S137. ¹H NMR spectrum of compound **33**.






Figure S138. ¹³C NMR spectrum of compound **33**.

Acquired by : Admin
 Date Acquired : 18/3/2019 2:12:47 PM
 Sample Name : LABI-22-FINAL-2
 Sample ID :
 Tray# : 1
 Vial# : 34
 Injection Volume : 3
 Data File : C:\LabSolutions\Data\2019\2019-wk11\LABI-22-FINAL-2.lcd
 Background File : blanco 18032019.lcd
 Method File : Method SCAN ACID standard.lcm
 Report Format : Default.LCMS.lcr
 Tuning File : C:\LabSolutions\Tuning\Tuning-ESI-pos-neg01072015a.lct
 Processed by : Admin
 Modified Date : 18/3/2019 2:23:04 PM




PeakTable

Peak#	Name	Ret. Time	Area	Area %
1		4.274	2829474	100.000



MS Spectrum Graph



MS Spectrum Table

#1 Ret.Time:
BG Mode:Calc 4.210<->4.520(422<->453)
Mass Peaks:12 Base Peak:307.05(475327) Polarity:Pos Segment1 - Event1

#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic	#	m/z	Abs.Inten.	Rel.Inten.	Charge	Polarity	Monoisotopic
1	307.05	475327	100.00				7	348.10	79010	16.62			
2	308.05	105556	22.21				8	349.10	18277	3.85			
3	309.40	8327	1.75				9	370.10	47873	10.07			
4	328.95	28045	5.90				10	371.10	11539	2.43			
5	330.00	5891	1.24				11	635.20	9222	1.94			
6	344.95	5676	1.19				12	640.50	20656	4.35			

Figure S139. LCMS spectrum of compound 35.




Figure S140. ¹H NMR spectrum of compound 35.




Figure S141. ¹³C NMR spectrum of compound 35.