

Figure S1




Figure S1. Binding Profile of PEG-reactive Mabs against Fluorophores used in B cell Sorting. PEG Specific mAbs (DH1511-DH1516), anti-MPER mAb (DH1318.2), and anti-influenza HA (CH65), were tested by ELISA for binding to commercial Streptavidin-fluorochrome conjugates SAV-Vio Bright 515 and SAV-Alexa Fluor647. **(A)** All six PEG-specific mAbs bound to SAV-Vio Bright 515 whereas DH1318.2 and CH65 exhibited minimal binding. **(B)** PEG-specific mAbs, anti-MPER mAb DH1318.2, and CH65 exhibited minimal binding to SAV-Alexa Fluor647.