

*Supplementary Table 1: MRI acquisition protocols for Centers A and B*

	Center A			Center B		
	T2w	DWI	DCE	T2w	DWI	DCE
<b>Field of view (mm)</b>	160 x 160	160 x 160	200 x 200	180 x 180	180 x 180	180 x 180
<b>In plane resolution (mm)</b>	0.36 x 0.36	2.23 x 2.23	0.71 x 0.71	0.38 x 0.38	1.88 x 1.88	0.94 x 0.94
<b>section thickness/gap (mm)</b>	3/0	3/0	3/0	3/0	3/0	3/0
<b>Repetition time (ms)</b>	3020	7000	3.6	5403	2919	4.7
<b>Echo time (ms)</b>	85	101	1.3	100	61	2.2
<b>Flip angle (degrees)</b>	160	90	20	90	90	10
<b>b-values (s/mm<sup>2</sup>)</b>	-	0, 800; 1000 (calculated)	-	-	50, 500, 1000; 1700 (calculated)	-

*T2w: T2weighted image, DWI: Diffusion-weighted imaging, DCE: Dynamic Contrast Enhanced*