

630 **Supplementary files**

631

632 **Supplementary Figure 1. qPCR primers for iPSCs.**

| iPSCs-QPCR primer | primer (5' to 3')           |
|-------------------|-----------------------------|
| <i>Foxg1-F</i>    | CCTGCCCTGTGAGTCTTTAAG       |
| <i>Foxg1-R</i>    | GTTCACTTACAGTCTGGTCCC       |
| <i>Neun-F</i>     | GTAGAGGGACGGAAAATTGAGG      |
| <i>Neun-R</i>     | CATAGAATTTCAGGCCCGTAGAC     |
| <i>Brn2-F</i>     | AAAGTAACTGTCAAATGCGCG       |
| <i>Brn2-R</i>     | GCTGTAGTGGTTAGACGCTG        |
| <i>Gfap-F</i>     | CCTCCAGCGATTCAACCTTT        |
| <i>Gfap-R</i>     | GAAGCTCCAAGATGAAACCAAC      |
| <i>Co2-F</i>      | CCCCACATTAGGCTTAAAAACAGAT   |
| <i>Co2-R</i>      | TATACCCCCCGGTCTGTAGCGGT     |
| <i>Dloop-F</i>    | TATCTTTGGCGGTATGCACTTAACAGT |

633

|                 |                            |
|-----------------|----------------------------|
| <i>Dloop-R</i>  | TGATGAGATTAGTAGTATGGGAGTGG |
| <i>Shh-F</i>    | CGGAGCGAGGAAGGGAAAG        |
| <i>Shh-R</i>    | TTGGGGATAAACTGCTTAGGC      |
| <i>Gli1-F</i>   | TCTGCCCCATTGCCCACTTG       |
| <i>Gli1-R</i>   | TACATAGCCCCAGGCCATACCTC    |
| <i>Patch1-F</i> | CGGCGTTCTCAATGGGCTGGTTT    |
| <i>Patch1-R</i> | GTGGGGCTGCTGTTCGGGTCG      |
| <i>Actin-F</i>  | TCACCCACACTGTGCCATCTACGA   |
| <i>Actin-R</i>  | CAGCGGAACCGCTCATTGCCAATGG  |

635 **Supplementary Figure 2. qPCR primers for wild-type and IOB mouse.**

636

| WT and IOB mouse -QPCR primer | primer (5' to 3')          |
|-------------------------------|----------------------------|
| <i>neun-F</i>                 | GTTGCCTACCGGGGTGCACAC      |
| <i>neun-R</i>                 | TGCTCCAGTGCCGCTCCATAAG     |
| <i>pax6-F</i>                 | GTTCCCTGTCCTGTGGACTC       |
| <i>pax6-R</i>                 | ACCGCCCTGGTTAAAGTCT        |
| <i>brn2-F</i>                 | AGAGCCCAGGGCAGAAAAGT       |
| <i>brn2-R</i>                 | GGCGCTCTGGTTAAAGGAG        |
| <i>gfap-F</i>                 | GGACTGAACCATGTCCCTTGTC     |
| <i>gfap-R</i>                 | AGGCTAGCTCTATCGGTATAACCTAA |
| <i>Mito1-F</i>                | CTAGAAACCCGAAACCAAA        |
| <i>Mito1-R</i>                | CCAGCTATCACCAAGCTCGT       |
| <i>Cox1-F</i>                 | TGCTAGCCGCAGGCATTACT       |
| <i>Cox1-R</i>                 | CGGGATCAAAGAAAGTTGTGT      |
| <i>gli1-F</i>                 | TTATGGAGCAGCCAGAGAGA       |
| <i>gli1-R</i>                 | GAGCCCGCTCTTTGTTAAT        |
| <i>gli2-F</i>                 | TGAAGGATTCCCTGCTCGTG       |
| <i>gli2-R</i>                 | GAAGTTTCCAGGACAGAACCA      |
| <i>gli3-F</i>                 | AAGCGGTCCAAGATCAAGC        |
| <i>gli3-R</i>                 | TGTTCCCTCCGGCTGTTTC        |
| <i>patch1-F</i>               | TGACAAAGCCGACTACATGC       |
| <i>patch1-R</i>               | AGCGTACTCGATGGGCTCT        |
| <i>gapdh-F</i>                | CATCACTGCCACCCAGAAGACTG    |
| <i>gapdh-R</i>                | ATGCCAGTGAGCTTCCGTTAG      |

637

638 **Supplementary Figure 3. Genotyping processing for wild-type and IOB  
639 mouse.**

640 (a) Genotyping primers for wild-type and IOB mouse. (b) Genotyping  
641 processing for wild-type and IOB mouse.

642

643

644 **(a) Genotyping primers**

| OCRL Wild-type (250 bp) | primer (5' to 3')         |
|-------------------------|---------------------------|
| GT-F                    | CCCTTTCATCTGTTAGGAGAAAT   |
| GT-R                    | GCATGGTAAACGCACTATGTGG    |
| OCRL KO (480 bp)        | primer (5' to 3')         |
| OCRL KO -F              | GCCCTTGATTCTAATCCCTTTCATC |
| OCRL KO -R              | TCTGAGCCCAGAAAGCGAAG      |

645

| Inpp5b (wt 150 bp/ mut 220 bp) | primer (5' to 3')     |
|--------------------------------|-----------------------|
| E25-F                          | TAAAGTCTGAAAATCCAAGGC |
| E34-R                          | CTCATTCTCCTTGATTCCAA  |

646

| INPP5B HumanBAC (240 bp) | primer (5' to 3')   |
|--------------------------|---------------------|
| in Exon 1F               | CCACCCCCACGATTGACTC |
| in Exon 1R               | GGTGTCCCAGCCCTCAG   |


647

648 **(b) Genotyping processing**

649 Step 1 - OCRL WT



Step 2- OCRL KO




650


651 Step 3- Inpp5b (wt 150 bp/ mut 220 bp)  
652 bp)

Step 4- INPP5B HumanBAC (240

653



654



655