

Table 3. Results from the Taq-Man assay on for leukemias; two positive for the *TCF3/PBX1* and two negative

Sample	Detector	CT	Mean	DCT	DDCT	Ratio sample/ref.
<i>TCF3/PBX1</i>	FGF9	26.45	26.38	11.06	-11.525	2946.9
	18S	15.3	15.32			
<i>TCF3/PBX1</i>	FGF9	27.04	27.106667	10.51667	-12.0683	4294.7
	18S	16.61	16.59			
Other, ref.	FGF9	38.3	38.275	22.585	0	1
	18S	15.7	15.69			
Other	FGF9	34.42	34.53	18.04	-4.545	23.3
	18S	16.45	16.49			