

**Supplemental Table 2: LC/MS gradient for the separation of isoprenoid phosphates.**

	<b>Gradient composition</b>			
<b>Time [min]</b>	<b>% A</b>	<b>% B</b>	<b>% C</b>	<b>% D</b>
0.0	95	2	3	0
0.3	95	2	3	0
0.5	93	2	0	5
1.4	55	40	0	5
2.2	40	55	0	5
2.5	0	95	0	5
3.0	0	95	0	5
3.1	0	99	0	1
3.6	0	99	0	1
3.7	95	2	3	0
7.0	95	2	3	0

See Methods for details.