

Additional file 8 –Primers used in the PCR experiments.

Number		Primer		Products
ATCC12228 RP62A				
SE1029	SERP0918	Sense	5'-GGCGTTAATGCTGAACAATTT-3'	2192
		Anti	5'-TCGTTGATTTCGCTTGAAAAA-3'	
SE0231	SERP2349	Sense	5'-TCCTGTGATCATTTCGTTTGG-3'	760
		Anti	5'-AGGTTGGGGTCATTTTGATG-3'	
SE0800	SERP0689	Sense	5'-ACCCACTCAAATGAGAAGTTA-3'	981
		Anti	5'-ACCATTTGCTATTTTAAAGAC-3'	
SE1546	SERP1399	Sense	5'-TGAGCACTTCCTCATCAACAA-3'	846
		Anti	5'-TTTGGCAAGGAGAATATGGG-3'	
SE2085	SERP2099	Sense	5'-CTACTCTTAGCACCACCACCG-3'	1224
		Anti	5'-AGGGTGAACGATGGTGAAAG-3'	
SE1042	SERP0931	Sense	5'-TCAATTGTTGTTCCAGGCAT-3'	424
		Anti	5'-ATCGATCCTTTGGATGATGG-3'	
SE0977	SERP0866	Sense	5'-ACGTAGGGGTCTCCCTGTTT-3'	421
		Anti	5'-TTTAATGGCTGGTGCTGTTG-3'	
SE1170	SERP1049	Sense	5'-CACTGCTTAATTGGGTCTACA-3'	1422
		Anti	5'-GGAATGGTTTGGTTTGTGATGA-3'	
SE0502	SERP0385	Sense	5'-CCCCAAATGAAAAAGTGGAA-3'	394
		Anti	5'-AGTCACTTTTAGGTTCGCGTTC-3'	
SE0546	SERP0431	Sense	5'-TAAACGAGCCAAACAAGGC-3'	1083
		Anti	5'-GCCTCTTGGTAAGTGGTTGG-3'	
SE0265		Sense	5'-TGTTGATTTGAGGTTGCTCAT-3'	422
		Anti	5'-CGAACGGAACCTTGAACAAT-3'	
	SERP2313	Sense	5'-CGAACGGAACCTTCAACAAT-3'	
		Anti	5'-TGTTGATTTGAGGTTGCTCAT-3'	
SE0378		Sense	5'-TTTGAAGATGGGCATTTCCCT-3'	470
		Anti	5'-GCATAAACTGCTATTGCTTG-3'	
	SERP0259	Sense	5'-TGAAAAAGACAATGTTTCATCA-3'	
		Anti	5'-CGAATTTCAAAAGGCGAGAG-3'	
SE0039		Sense	5'-AGATGGCAGTTTTTGTGTCG-3'	392
		Anti	5'-GGACGTTGGTTCGCATATTT-3'	
	SERP2495	Sense	5'-GGACGTTGGTTCGCATATTT-3'	
		Anti	5'-AAGATGGCAGTTTTTGTGTTG-3'	