

DryingStage Active Tray 989

Post Dry Information

Oven Temperature: 92.0 °C

Placed in Oven: 08/23/2006 1pm 0
mm/dd/yyyy Hour Min

Drying Period: 20.0 hours

Digestion Acid: Conc nitric acid
0.7 mL

Tube Information

Click the button below to edit tube level information such as dry weights, standard reference materials, and digestion acid values.

TUBE EDITOR

SaveStage

DryingStage EXIT TUBE EDITOR
Active Tray 989

Tube Information

The following list displays the weighed tubes only by default. To edit the non-weighed tubes of this experiment click the view non-weighed tubes link at the bottom.

Weighed Tubes

Tube ID	Std Mtrl	Dig Acid	Dry Weight
2	N/A	Conc nitric acid 0.7 mL	0.00212
11	N/A	Conc nitric acid 0.7 mL	0.0021
20	N/A	Conc nitric acid 0.7 mL	0.00169
29	N/A	Conc nitric acid 0.7 mL	0.00284
38	N/A	Conc nitric acid 0.7 mL	0.0022
47	N/A	Conc nitric acid 0.7 mL	0.00302
56	N/A	Conc nitric acid 0.7 mL	0.00274