

Cells	MethoCult GF plus (E-, GEMM-, GM-, M-CFCs) #cells/1 ml/dish	MegaCult + SCF, EPO (E-, Mk-CFCs) #cells/0.75 ml/chamber
hESC/OP9 coculture		
- days 4-5	200000	NA
- day 6	100000	NA
- days 7-8	50000	NA
- day 9	20000	NA
- CD34 ⁻ CD43 ⁻ (days 6, 9)	250000	NA
CD34⁺ cell subsets		
- CD43 ⁺	1000	NA
- CD43 ⁻ KDR ⁺	20000	NA
- CD43 ⁻ KDR ⁻	20000	NA
CD43⁺ cell subsets		
- CD41a/CD235a ⁺ CD45 ⁻	5000	500
- CD41a/CD235a ⁻ CD45 ⁻	400	500
- CD41a/CD235a ⁻ CD45 ⁺	200	500
- CD41a/CD235a ⁺ CD45 ⁺	5000	500

NA indicates not applicable.