


**B**

	G2A Lysates	Pool 2	
		Pre-sort	Post-sort



**Silver staining**—Silver staining of proteins separated by NuPAGE 4-12 % Bis-Tris gel from MDCK cells stably expressing G2A Naked2-EGFP. Whole cell lysates (10  $\mu$ g), pre-sort Pool 2 (10  $\mu$ g) and post-sort Pool 2 (15  $\mu$ g) were loaded in each lane. Note enrichment of bands in different regions of the gel in post-sort compared to pre-sort. Protein concentrations were determined using Micro BCATM Protein Assay Kit (Pierce).