

Table 2. Primers utilized for real-time PCR analysis

Amplicon	Forward Primer	Reverse Primer
SFRP1	5'-GCCCCGAGATGCTTAAGTGTGACAA-3'	5'- ACTCGCTGGCACAGAGATGTTCAA-3'
Cyclin D1	5'-AACTACCTGGACCGCTTCCT -3'	5'- CCACTTGAGCTTGTTTACCA-3'
E-cadherin	5'-ACACCATCCTCAGCCAAGA-3'	5'-CGTAGGGAAACTCTCTCGGT-3'
vimentin	5'-AAAGTGTGGCTGCCAAGAAC-3'	5'-AGCCTCAGAGAGGTCAGCAA-3'
Snail	5'-GCTGCAGGACTCTAATCCAGA-3'	5'-ATCTCCGGAGGTGGGATC-3'
Slug	5'-CATGCCTGTCATACCACAAC-3'	5'-GGTGTCTAGATGGAGGAGGG-3'
ZEB1	5'- GGGAGGAGCAGTGAAAGAGA-3'	5'-TTTCTTGCCCTTCCTTTCTG-3'
ZEB2	5'- AAGCCAGGGACAGATCAGC-3'	5'-CCACACTCTGTGCATTTGAACT-3'
GAPDH	5'-GAAGATGGTGATGGGATTTC-3'	5'-GAAGGTGAAGGTCGGAGT -3'