

Supplemental Table 1. PCR primer pairs used to detect progestin binding proteins in GL and GL5 cells			
Target Gene	Forward Primer	Reverse Primer	Product Size
<i>PGR</i>	GTGCAAGGTTGGAGACAGCT	TTTGCCCTTCAGAAGCGGAC	213 bp
<i>PGRMC1</i>	GCCTGGATAAGGAAGCACTG	TAATCATTTCCTCCGGGCACT	202 bp
<i>Paqr7</i>	CGCTCTTCTGGAAGCCGTAC ATCTATG	CAGACGGTGGGTCCAGACAT TCAC	100 bp
<i>Paqr8</i>	GTCAACGTCTGGACCCATT	GTGGGAGAGCTCTGACTTGG	192 bp
<i>Paqr5</i>	CAGCTGTTTCACGTGTGTGT GATCCTG	GGACAGAAGTATGGCTCCAG CTATCTGAG	120 bp