

Assay	Site	L primer	PCR forward (L) primers	H primer	PCR reverse (H) primers	Size
K1	497/498del	L00460	cttttggcggatgcacttttaac	H00536	gtgtgtgtgtgctgggtagg	77
K1	1189	L01142	actacgagccacagcttaaaactc	H01220	tatcgattacagaacaggctcctc	67
K1	5913	L05846	cctgtctttagatttacagtccaatgc	H05950	ccaatgtctttgtggttttagag	118
K1	8137	L08105	gcaattcccggacgtctaa	H08205	atgaaactgtggtttgctccac	94
K1	11299	L11266	caccctaggctcactaaacattcta	H11335	gttggctcaggagtttgatagttct	70
K1	12308	L12257	ctttctcaacttttaaggataacagc	H12336	ttacttttatttggagttgcacca	80
K1	12705	L12663	ctcagacccaacattaatcagttc	H12766	agccgatgaacagttggaa	104
K1	16224	L16186	ccctcccatgcttacia	H16291	gtgggtaggtttgttggtatccta	106
K1	16311	L16257	caccctcaccactagga	H16363	gacgagaagggttgactgtaat	107
K1	16362	L16320	ccattaccgtacatagcacattaca	H16424	atattgatttcacggaggatggt	105