

Table S5.**Sequences of newly-designed primers used for qPCR analyses.**

Ero1lb	<i>Forward</i>	5'-AGCTGGTCCCAGTGACAGAAA-3'
	<i>Reverse</i>	5'-ACAGCCCTGCATTACAGAGGA-3'
Nnat	<i>Forward</i>	5'-TCTGGACCAAGTCGGAACAGT-3'
	<i>Reverse</i>	5'-TGTAGGGCTTCAGCTGGTGAT-3'
Pdia4	<i>Forward</i>	5'-GTAGCGAAGATCGATGCAACC-3'
	<i>Reverse</i>	5'-CAGGATCTTGATGGTGGGGTA-3'
Pdia6	<i>Forward</i>	5'-AAACATCAGTCCCTGGGAGGT-3'
	<i>Reverse</i>	5'-CGGCATCTACAATGGCTTCTC-3'
Edem3	<i>Forward</i>	5'-CACCTTGATCCTCGAGTTTGC-3'
	<i>Reverse</i>	5'-TCGCTGTCTTTTCTCCAGAG-3'
Dnajc10	<i>Forward</i>	5'-ATCGGAGCTGTCAACTGTGGT-3'
	<i>Reverse</i>	5'-ACTGCAGCCATTCCAGACCTA-3'
Atf4	<i>Forward</i>	5'-GCAGTGTTGCTGTAACGGACA-3'
	<i>Reverse</i>	5'-CGCTGTTCAGGAAGCTCATCT-3'
Eif4ebp1	<i>Forward</i>	5'-TGCTACTGCAGGGGAATGACT-3'
	<i>Reverse</i>	5'-ACCCAGCCTAAGGAAAGATGG-3'
Trib3	<i>Forward</i>	5'-CTGTGCCAAAGACCAAAGTCC-3'
	<i>Reverse</i>	5'-CTTGAATCCCACAGCCAAGTC-3'
HPRT	<i>Forward</i>	5'-TGCTCGAGATGTCATGAAGGA-3'
	<i>Reverse</i>	5'-CCAGCAGGTCAGCAAAGAACT-3'