

An On-Line, High-Pressure Digestion System for Protein Characterization by Hydrogen Deuterium Exchange

Lisa M. Jones, Hao Zhang, Ilan Vidavsky, and Michael L. Gross

Department of Chemistry, Washington University, St. Louis, MO 63130

Supplementary material

EXPERIMENTAL SECTION

Digestion Apparatus. The apparatus was plumbed with silica tubing (O.D. 360 μm , I.D. 150 μm), and a 20 μL sample loop was used for protein injection. For digestion, a constant flow of solvent A (99.9% H_2O , 0.1% formic acid) was kept through the system. An 8 min gradient from 0-100% solvent B (99.9% acetonitrile, 0.1% formic acid) with a flow rate of 40 $\mu\text{L}/\text{min}$ was run prior to MS analysis.

Back Exchange Correction. To account for back exchange during analysis, the deuterium content was calculated from:

$$D_{\text{H}\rightarrow\text{D}} = \frac{m - m_{0\%}}{m_{100\%} - m_{0\%}} \times N$$

where, m , $m_{0\%}$, and $m_{100\%}$ are the centroids for each peptide peak after a specified HDX period, zero-time control, and the 100% control. The 100% exchange control was the result for 48 hr exchange. The average back exchange for both the high and low pressure digestion was ~20%.