

Supplementary Table I. Antibodies used in flow cytometry.

Antigen-Fluorochrome	Manufacturer	Clone	Isotype	Dilution
CD3-Pacific Blue	eBioscience	OKT3	IgG2a	1:40
CD4-PE-Cy7	eBioscience	RPA-T4	IgG1	1:200
CD8-PerCP-Cy5.5	BD Pharmingen	SK1	IgG1	1:200
IL-17 Alexa 488	eBioscience	eBio64 DEC17	IgG1	1:20
IL-4-PE	BD Pharmingen	3010.211	IgG1	1:20
IL-22-APC	R&D	142928	IgG1	1:20
IFN γ -Alexa 700	BD Pharmingen	B27	IgG1	1:500

19 **Supplementary Figure 1. Similar composition of cells isolated from atopic**
20 **dermatitis and psoriasis lesions.** Percentages of isolated T cells, CD4⁺ T
21 helper cells and CD8⁺ T cytotoxic cells are statistically similar between atopic
22 dermatitis and psoriasis skin lesions.

23

ACCEPTED MANUSCRIPT

Supplementary Figure 1

Composition of cells isolated from skin lesions

