

## Supplementary Table 2

Details of external (ex) and internal (in) primer pairs used for single cell multiplex RT-PCR.

Template	ex/in	Sense primer	Anti-sense primer
CD3 $\epsilon$	ex	ATAGGCGGTGATGAGGATGA	TAGTCTGGGTTGGGAACAGG
CD3 $\epsilon$	in	ATGTCTGCTACCCCAGAGGA	TGGCCTCTCCTTGTTTTGTC
IL-10	ex	GAGAACAGCTGCACCCACTT	GGCCTTGCTCTTGTTTTCAC
IL-10	in	TGCCTTCAGCAGAGTGAAGA	GGGAAGAAATCGATGACAGC
granzyme A	ex	GAACAAAAGGTCCCAGGTCA	GCACGAGTCTCTCCACCTC
granzyme A	in	TTGGGGCTCACTCAATAACC	GGCTTCCAGCACAAACCATA
granzyme B	ex	CCCTGGGAAAACATTCACAC	ATTACAGCGGGGGCTTAGTT
granzyme B	in	TGAAGATGACAGTGCAGGAAG	CGTTTCATGGTTTTCTTTATCCA
perforin	ex	CATGTAACCAGGGCCAAAGT	CCGAACAGCAGGTCGTTAAT
perforin	in	GCTGGACGTGACTCCTAAGC	AGGGCCCTCTTGAAGTCAG