

Table S1. Primers used for PCR DIG labeling and sequencing,

PCR Reaction	Predict size (bp)	Primers	Primer sequences (5'-3')	Reference
DIG probe netB	384	AKP78 AKP79	GCTGGTGCTGGAATAAATGC TCGCCATTGAGTAGTTTCCC	Keyburn et al., 2008
DIG probe hdhA	387	HDHA-F HDHA-R	CAATTGTTACTGCAGCTTCAAGAGG TAATTTTTACCGCCACCTTCATAATA	This study
Long Range PCR-5' NEL_1	8503	SP-F DCM-R	GTGCAGTTACTATTATAGGACC TTGAAGCACTAAGGATTATAGA	This study
Long Range PCR-3' NEL_1	2677	BLAC-F LEXA-F	AACTACTTAATAGACACAGGAA TACAGGATCAGTATCATATACC	This study
CP1 pcr product 4214 bp	3000	WBRA-F LKI-R	GCAGATTCACCTCCTGTAACAA CAGTTCATATGTATGTGTTGAC	This study