

**Table S1. Primers used for RT-PCR**

Primer	Sequence (5' to 3')
<i>rp49</i> -forward	ATCGGTTACGGATCGAACAAG
<i>rp49</i> -reverse	GTAAACGCGGTTCTGCAT
<i>Barley</i> -forward	GGAGGAAGAGCGTGTTGTAGG
<i>Barley</i> -reverse	ACAACACCGACATGCTCTACAG
<i>White</i> -forward	GCAACCAAGTAAATCAACTGC
<i>White</i> -reverse	GTTTTGGCACAGCACTTTGTG
<i>bFtz-F1</i> -forward	TGGCGTACTTTTAGCGTCCT
<i>bFtz-F1</i> -reverse	AATACAAGAATCGATCTTCAAGTGG
<i>E74</i> -forward	AACGTGCAATCTCTTGAAAGG
<i>E74</i> -reverse	AGATAAAGGTACTIONTAGAGAACGCAAC
<i>E93</i> -forward	AGAACGCGTTGCTGAAGAAT
<i>E93</i> -reverse	GATTGCTCTGGCTGATCTCC
<i>BR-C</i> -forward	GCCCCTCAGTTCATCATGTC
<i>BR-C</i> -reverse	CTGCAACTGTGTGTTGAGACC