

Supplemental Table S2 Cell-Type Specific Marker Probes

Gene	Primer Fwd 5'->3'-OH	Taqman Probe 5'->3'-OH	Primer Rev 5'->3'-OH
18S (internal standard)	CGCCTACCACATCCAAGGAA	CAGCAGGCCGCGCAAATTACCCA	GGGCTCGAAAGAGTCCTGT
NF-H	ACACAGATGCTATGCGCTCG	CCCAAGAGGAGATAACTGAGTACCGGCG	CTCCTTGGTGCTTTTCAGGG
Nf-L	GAGGAAGAGATCGCCGAGCT	AGGTGTCCTCCAAGCCCGACCTC	TCCTTGAGAGCGGCGG
Nf-M	GATTTAGCACATTTTCAGGAAGCA	TGTACACACACCCGACAGCCCTCAGTCAC	GGCTCGACTTTGGTCTTCTG
PLP	GCAGTCTCTGCGCTGATGC	AGAATGTATGGTGTCTCCCATGGAATGCTTTT	GGCTGTTTTCAGATGGACA
MBP	AGCAAGTACCATGGACCATGC	CATGGCTTCTCCCAAGGCACAGAG	CACCCCTGTACCAGCTAAAG
GFAP	GATGGAGCTCAATGACCGCT	CTACATCGAGAAGTCCGCTTCCTGGAAC	TGGCTCGAAGCTGGTTCAG
Nestin	ACCAGCGACATGGAGGGT T	TTTTCGATGTGGGAGCTCAATCGACG	CTGCTCCTCCAGCGTCTTG
beta2/beta3 tubulin	CAGTTCACAGCCATGTTCCG	AAGGCCTTCTCCACTGGTACACGC	TCATTCATGTTGCTCTCGGC