

Supplemental Figure 1




Figure S1: Characterization of A β O_s preparation. Representative western immunoblot of an A β O preparation probed with rabbit polyclonal anti-oligomer antibody. The left panel shows a lightly exposed film to facilitate evaluation of the size range of higher molecular weight oligomers, and the right panel shows a stronger film exposure of the same membrane to allow detection of low-n oligomers (dimers, trimers, tetramers). The size of aggregates varied from dimers (~ 9 kDa) to higher molecular weight oligomers (~50-100 kDa).