

Table S6. Taqman probes for ChIP-PCR measurements

Name	Sequence 5'->3'
DbpI2 fwd	TGGGACGCCTGGGTACAC
DbpI2 rev	GGGAATGTGCAGCACTGGTT
DbpI2 probe	FAM-CCCTGATCCGACCTCCAGGCAA-TAMRA
Per1 fwd	AGCCCTCTCAGCCTATGAGAAAGT
Per1 rev	CCCGCCCTGCCTAAATCA
Per1 probe	FAM-TTTAGGGCAGGGCTGGCATTGCG-TAMRA
Gapdh fwd	CATGGCCTTCCGTGTTCTTA
Gapdh rev	CCTGCTTCACCACCTTCTTGA
Gapdh probe	FAM-CCGCCTGGAGAAACCTGCCAAGTATG-TAMRA