



Supplemental Fig. S2. Physical interaction between AtNug2 and various AtL10a homologs. Pull-down assay of GST-AtNug2 with [³⁵S]methionine-labeled AtL10aA, AtL10aB or AtL10aC. Radiolabeled AtL10a proteins were incubated with GST alone (lanes 2, 5 and 8) or GST-AtNug2 (lanes 3, 6 and 9). Reaction mixtures were incubated with glutathione beads, centrifuged and then the supernatant was discarded. Residual precipitates were boiled and subjected to SDS-PAGE. Five percent of total ³⁵S-labeled AtL10a was loaded onto lanes 1, 4 and 7, and 50% of each sample was used for pull-down assays (lanes 3, 6 and 9).