

Supplemental Table. DNA probes used for EMSA in Supplemental Fig. 1. Only the top strand sequence is shown. These were annealed with their complementary bottom strands and used in EMSA.

Name	Top strand sequence
4rA	5' - TTTTGGGAAGGAAGGAAGGAATTTTTTTTTTTTTTTT - 3'
4rB	5' - TTTTTTGGGAAGGAAGGAAGGAATTTTTTTTTTTTTTTT - 3'
4rC	5' - TTTTTTTTTGGGAAGGAAGGAAGGAATTTTTTTTTTTTTTTT - 3'
4rD	5' - TTTTTTTTTTTGGGAAGGAAGGAAGGAATTTTTTTTTTTTTTTT - 3'
4rE	5' - TTTTTTTTTTTTTTGGGAAGGAAGGAAGGAATTTTTTTTTTTTTTTT - 3'
4rF	5' - TTTTTTTTTTTTTTTTTGGGAAGGAAGGAAGGAATTTTTTTTTTTTTTTT - 3'
4rG	5' - TTTTTTTTTTTTTTTTTTTGGGAAGGAAGGAAGGAATTTTTTTTTTTTTTTT - 3'