

Table S2. The primers used for qPCR experiments

Gene	Forward	Reverse
<i>nvd</i>	5'-GGAAGCGTTGCTGACGACTGTG-3'	5'-TAAAGCCGTCCACTTCCTGCGA-3'
<i>dib</i>	5'-TGCCCTCAATCCCTATCTGGTC-3'	5'-ACAGGGTCTTCACACCCATCTC-3'
<i>spok</i>	5'-GCCGAGCTAAATTTCTCCGCTT-3'	5'-TATCTCTTGGGCACACTCGTG-3'
<i>phm</i>	5'-GGATTTCTTTCGGCGCGATGTG-3'	5'-TGCCTCAGTATCGAAAAGCCGT-3'
<i>sad</i>	5'-CCGCATTACAGCAGTCAGTGG-3'	5'-ACCTGCCGTGTACAAGGAGAG-3'
<i>rpl23</i>	5'-GACAACACCGGAGCCAAGAACC-3'	5'-GTTTGCCTGCGCAATAACCAC-3'
<i>ptth</i>	5'-CCACGAGCTCATTTCGTAACCTTGC-3'	5'-CACTCCACATCCCACAGAGATGGCGATGG-3'