

Table S2. Primer sets for ChIP analysis

Gene	Primers	
<i>Per1</i> CRE	forward	5'- AACACGATGTTCCCTAGTGC-3'
	reverse	5'- ATCAGTCCCCCACAGTAACC-3'
<i>Per2</i> CRE	forward	5'- GAGACTCTGCCAGGTGGATG-3'
	reverse	5'- CACCCTGCACACACCCAAAG-3'
<i>Cry2</i> CRE	forward	5'- TGTA ACTCCAACCCCCAGAAA-3'
	reverse	5'- AGTCTCTGCCCCGACTCTCTC-3'
<i>Bmal1</i> CRE	forward	5'- TCAAATCCCAGCAACCACATG-3'
	reverse	5'- TTTTACTTGTACCCTGATGCTTAAT-3'
<i>Per2</i> E-box	forward	5'- GGTTCCGCCCCGCCAGTATGC-3'
	reverse	5'- CCGTCACTTGGTGCGCTCGGC-3'