

Supplementary Table S3: Dilutions used and manufacturers' information for antibodies used in immunohistochemistry (IHC) staining

<i>Tumour markers</i>	Dilution for IHC	Clone	Manufacturer
CD44	1:1000	SFF-304	Bender MedSystems, San Diego, CA
Musashi-1	1:600	AB5977	Millipore, Billerica, MA
CD133	1:20	AC133	MiltenyiBiotec, Auburn, CA
Ki67	1:1000	ab15580	Abcam, Cambridge, MA
PCNA	1:1000	NCL-PCNA	Novocasta, UK