



sFig. 3. Analysis of recombinant Ad-FK and CD46ex-Fc proteins. (A, B) Individual FK proteins were produced using the Baculovirus expression system. Purified FK proteins from Ad3, Ad5, Ad7, Ad11, and Ad35 were analyzed by either 12.5% reducing PAGE and sypro ruby red staining (A), or by 10% native PAGE and Coomassie Blue staining (B). (C) Analysis of recombinant CD46ex-huFc. Three μ g of purified CD46ex-huFc were analyzed by 10% native PAGE and Coomassie Blue staining.